B-1000 Series

INSTRUCTION MANUAL

| |
|--------------|
| Model |
| B-1000POL-I |

Version: 1
Issued: 17, 06, 2015

Table of Contents

Warning

Symbols and conventions

Safety Information

Intended use

List of accessories and spare parts

Overview

Unpacking

Assembling

Using the microscope

Maintenance

Troubleshooting

Equipment disposal

Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.


Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

List of accessories and spare parts

| CAT. NO. | DESCRIPTION |
|----------|---|
| M-1001 | WF10x/22mm eyepiece (pair) |
| M-1004.N | Centering telescope, 30mm diameter |
| M-781 | WF10x/22mm micrometer eyepiece (10mm, 0.1mm div.) |
| M-1011 | Trinocular Head (3 positions) |
| M-1012 | Binocular ERGO head |
| M-1033 | Bertrand Lens with Analyzer and slot for slides (with Lambda, 1/4 Lambda and Quartz Edge) |
| M-1034 | Incident Polarizing Attachment |
| M-1044 | Quintuple nosepiece with centrable positions for POL objectives |
| M-1090 | 5x IOS LWD POL PLAN objective |
| M-1091 | 10x IOS LWD POL PLAN objective |
| M-1092 | 20x IOS LWD POL PLAN objective |
| M-1093 | 50x IOS LWD POL PLAN objective |
| M-1145 | Rotating Stage, centrable |
| M-1146 | Attachable mechanical stage for rotating Stage |
| M-005 | Micrometric slide, 26x76mm, range 1mm, div. 0,01mm |
| M-690 | Eyecup (pair) |
| M-619 | Photo tube adapter for full frame SLR camera |
| M-173 | Photo tube adapter for APS-C SLR camera |
| M-699 | Photo tube adapter for DIGI digital camera series |
| M-620 | CCD camera adapter for 1/3" sensor |
| M-620.1 | CCD camera adapter for 1/2" sensor |
| M-114 | CCD camera adapter 0,45x |
| M-113.1 | Ring adapter, 30mm (for monocular and binocular microscopes) |
| M-617.1N | Phase contrast set with IOS PLAN objectives 40x |

Overview

Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Assembling

Once you open the box, these are the microscope's components:

- ① Condenser
- ② Optical head
- ③ Main body
- ④ Eyepieces
- ⑤ Power supply
- ⑥ Stage
- ⑦ Objectives
- ⑧ Bertrand lens
- ⑨ Retardation plates
- ⑩ Polarizer and analyzer
- ⑪ Polarizing attachment

Put the microscope stand on a solid table.

First insert the incident polarization attachment, using the 2,5mm Allen wrench to tighten the screw. All screws are already inserted into each threaded hole.

(Fig.1)

Fig.1

Insert the Bertrand lens attachment, using the 3mm Allen wrench to tighten the screw.

(Fig.2)

Fig.2


Insert the optical head above the Bertrand lens, using the other **3mm** Allen wrench to tighten the screw. (Fig.3)

Insert both eyepieces into the tubes of the optical head. (Fig.4)

Insert the condenser under the stage: push until it is well inserted into its holder (under the condenser there is a pin that must fully enter the guide of the holder). (Fig. 5)

Mount the rotating stage: at the bottom of the stage there is a spring, push this spring toward the stage support ①, then push the stage downward ②. (Fig.6)

Screw each objective into the thread of the nosepiece, in order of magnification. (Fig.7)

Fig.7

Insert the retardation plates by sliding them into the slot under the optical head. (Fig.8)

Fig.8

Center each objective using the provided screwdrivers (see chapter using the microscope). (Fig.9)

Fig.9

Insert the jack for incident illumination. (Fig.10)


Fig.10

Insert the power supply jack on the connector at the rear.
(Fig.11)

The rotating stage is provided with:
a) Stage centering knobs;
b) Stage lock knob;
c) 45° rotation clicks;
d) Optional X-Y translator.
(Fig.12)

The Bertrand lens attachment is provided with
a) Analyzer filter;
b) Retardation plates;
c) Lens focus adjustment;
d) Lens X-Y centering screws.
(Fig.13)

Using the microscope

Control keyboard

B-1000 illumination can be managed through the keyboard placed on the left of the stand:

1) ON-OFF: press this button to turn on or off the LED illuminator.

2) BOOST: press this button in order to increase the brightness (useful for high-magnification objectives or very opaque specimens).

2) AUTO-OFF: if you want the illuminator to switch off automatically, press this button until 15, 30 or 60 minutes delay is set. After this period of time, the light will turn off. You have to press the ON-OFF button to turn it on again. (Fig.14)


Fig.14

Do not enable boost mode while observing with low magnification objectives (4x, 10x) with fully open diaphragm: the high brightness may hurt user's eyes.

LED settings

In order to activate the illuminator, insert the plug of the external power supply into the socket and turn on the main switch on the side of the main body. Press the ON-OFF button on the control keyboard and turn the brightness adjustment knob to a brightness suitable for observation. (Fig.15)

Position I: transmitted illumination.
Position II: incident illumination.

Fig.15

Preliminary settings

Before trying to focus the sample, ensure that:

- The field diaphragm is fully open. (Fig.16)
- The aperture diaphragm is fully open. (Fig.17)
- No filter is inserted under the observation head.

Fig.16

Fig.17

- The analyzer filter is in OUT position (lever fully pulled out). (Fig.18)
- Photo port closed (lever on observation head fully in).

- Bertrand lens in OUT position (you read "0" on the Bertrand lens disc under the observation head).
- (If you use 5X objective) Swing out lens of the condenser in OUT position. (Fig.19)

Adjust the observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw. (Fig.20)

Adjust interpupillary distance

Hold the right and left parts of the observation head using both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. (Fig.21)

Place the specimen on the stage

Fix the specimen slide on the mechanical stage using the slide-clamp. Ensure that the specimen is centred over the stage opening. (Fig.22)


Fig.22

Focus tension adjustment

Turn the tension-adjust knob to get a suitable tension for the focus system. (Fig.23)
NOTE: if the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the friction.


Fig.23

Focus-stop knob

Loosen the focus-stop knob, rotate the coarse focusing knob to bring the slide into focus with a 5X or 10X objective. Then rotate the focus-stop in order to block the height of the stage. This simplifies the next focusing operations. The focus-stop knob is also useful to avoid accidental contacts between objective and specimen. (Fig.24)

Fig.24

Diopter adjustment

Adjust the fine focusing knob to get the image sharp and clear while observing with your right eye, then turn the left diopter ring to a sharp and clear image also with the other eye. The highpoint eye-pieces allow the user to wear glasses. (Fig.25)

NOTE: For the optimal parafocality of the image, it's suggested to wear your glasses during the normal use of the microscope.

Fig.25


Centering the condenser

Insert the swing-out lens of the condenser ① and fully close the field diaphragm ②. Rotate the condenser adjustment knob ③ until you see a sharp image of the closed field diaphragm (a bright spot of light). Act on the condenser centering screws ④ until you move the bright spot in the center of the image field. Then re-open the field diaphragm ②. (Fig.26)

Fig.26

CENTERING THE CONDENSER

Adjust the aperture of the iris diaphragm under the condenser to set the numerical aperture of the illuminator, thus controlling image contrast and resolution. It is suggested to set the aperture diaphragm between 70% and 80% of objective's N.A. (Fig.27)


Fig.27

Centering the stage

With the Bertrand lens out, focus on your slide. In this procedure we'll center the optical axis of each objective with the rotation axis of the stage. Insert the 10x objective, rotating the nosepiece. Look at your sample while continuously rotating the stage clockwise and then counter-clockwise by a little angle (e.g. 30° or 45°). During these oscillations you should spot a point on the sample that doesn't rotate on a circumference but only revolving about itself (it is the center of the rotations). (Fig.28)



Fig.28

Using the **stage centering screws**, bring this point in the center of the field of view. In this way the mechanical center of rotation of the stage coincides with the system's optical axis. (Fig.29)

Insert another objective, rotating the nosepiece. Repeat the above operation (stage rotations), and bring the point into the center using the **objective centering screws** (not the stage centering screws). (Fig.30)

NOTE :The stage can be locked in its position using the stage lock knob.

Bertrand lens

In order to have a conoscopic view of your sample, you can insert a Bertrand lens in the optical path by rotating the disc under the observation head.

The Bertrand lens can be fine aligned along Z axis by rotating this disc. (Fig.31)

The lens can also be centered in X-Y axis by means of two screws located at the back of the lens. (Fig.32)

Polarizer and analyzer filters

Insert or remove the lower polarizer as required. It can also be rotated at any angle. (Fig.33)

Fig.33

The analyzer filter can be inserted by pushing the slider inward. By rotating the front disc, it can also be rotated at any angle. (Fig.34)

Fig.34

Retardation filter plates

The microscope is provided with these retardation plates: lambda, lambda/4 and quartz wedge. You can insert them in the optical path, sliding the plate holder in the slot under the observation head. (Fig.35)


Fig.35

Incident illuminator

Pull the lever in R position, in order to insert the prism for reflected light. (Fig.36)


Push the main switch to position II and use the brightness adjustment knob to set the desired level. (Fig.37)


Open/close the field and aperture diaphragms as needed. (Fig.38)

- ① Field diaphragm
- ② Aperture diaphragm


You can center both diaphragms with the screws (one on each side). (Fig.39)

In order to perform an incident polarization analysis, insert the polarization filter into the attachment. You can use the analyzer and Bertrand lens as described for transmitted light observation. (Fig.40)

- ① Polarizer filter

Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

| PROBLEM | CAUSE | SOLUTION |
|--|--|--|
| 1. Optical System | | |
| LED does not light. | Power cord is unplugged. | Plug power cord into the power outlet. |
| LED operates, but field of view remains dark. | Aperture and field iris diaphragms are not opened wide enough. | Adjust them to proper sizes. |
| | Condenser is lowered too much. | Adjust the condenser height position. |
| | Light path selector knob is set to the camera position. | Move the knob to the eye position. |
| Field of view is obscured or not evenly illuminated | Light path selector knob is in an intermediate position. | Set the knob according to the observation method. |
| | Revolving nosepiece is not correctly engaged. | Make sure that the revolving nosepiece clicks properly into place. |
| | Condenser is not attached properly. | Re-attach it. |
| | Revolving nosepiece is not attached properly. | Push the side dovetail all the way until it is stopped. |
| | An objective that falls outside of the condenser's illumination range is used. | Use a condenser to match the purpose. |
| | Condenser is not properly centered. | Center the condenser. |
| | Field iris diaphragm is stopped down too far. | Open the field iris diaphragm until it circumscribes the field. |
| Dirt or dust is visible in the field of view. | Dirt/dust on the eyepieces | Clean thoroughly. |
| | Dirt on the top surface of the condenser | |
| | Dirt/dust on the specimen | |
| Visibility is poor. · Image is not poor. · Contrast is poor. · Details are indistinct. · Image glares. | Condenser is lowered too far. | Adjust the condenser height position. |
| | Aperture iris diaphragm is stopped down too far. | Open aperture iris diaphragm. |
| | Revolving nosepiece is not mounted properly. | Push the slide dovetail all the way until it is stopped. |
| | Front lens of objective is dirty. | Clean objective. |
| | Immersion oil is not being used with an oil immersion objective. | Use immersion oil. |
| | Immersion oil contains bubbles. | Remove the bubbles. |
| | Recommended immersion oil is not used. | Use the provided immersion oil. |
| | Dirt/dust on specimen. | Clean it. |
| | Dirt/dust on condenser | |
| | Inappropriate object side or cover glass thickness. | Replace with glass of recommended thickness. |

| PROBLEM | CAUSE | SOLUTION |
|--|--|---|
| One side of image is blurred. | Objective is not correctly engaged in light path. | Make sure that revolving nosepiece clicks into place correctly. |
| | Revolving nosepiece is not correctly mounted. | Push slide dovetail all the way until it is stopped. |
| | Stage is not correctly mounted. | Re-attach it. |
| | Specimen is not correctly mounted on stage. | Place specimen correctly on top of stage and secure it with slide holder. |
| Image appears to waver. | Revolving nosepiece is not correctly mounted. | Push slide dovetail all the way until it is stopped. |
| | Objective is not correctly engaged in light path. | Make sure that revolving nosepiece clicks into place correctly. |
| | Condenser is not properly centered. | Center the condenser. |
| Field of view becomes only slightly brighter when the voltage is raised. | Condenser is not properly centered. | Center the condenser. |
| | Condenser is lowered too far. | Adjust the condenser height position. |
| 2. Coarse/Fine Adjustment | | |
| Coarse adjustment knob is hard to turn. | Tension adjustment ring is tightened excessively. | Loose ring. |
| | You are trying to raise stage while focus-lock lever is kept locked. | Unlock focus-lock lever |
| Stage drifts down by itself or focus is lost during observation. | Tension adjustment ring is too loose. | Tighten ring. |
| Coarse adjustment will not go all the way up. | Focus-lock lever is locked at a too low height. | Unlock focus-lock lever. |
| Coarse adjustment will not go all the way down. | Condenser holder is too low. | Raise condenser holder. |
| Objective makes contact with specimen before focus is obtained. | Specimen is mounted upside down. | Mount specimen correctly. |
| 3. Observation Tube | | |
| Field of view of one eye does not match that of the other. | Interpupillary distance is incorrect. | Adjust interpupillary distance. |
| | Incorrect diopter adjustment. | Adjust diopter. |
| | Different eyepieces are used on left and right. | Change on eyepiece to match the other so that both sides are the same type. |
| | Your view is not accustomed to microscope observation. | Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope. |
| 4. Stage | | |
| Image shifts when you touch stage. | Stage is not properly mounted. | Clamp stage. |
| Specimen stops midway on the X-axis traverse. | Specimen is not correctly positioned. | Place specimen correctly. |
| X- and Y-axis knobs are too tight or too loose. | Tension of X- and Y-axis knobs is too high or too low. | Adjust tension. |

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.