B-380 Series

INSTRUCTION MANUAL

Model
B-383FL

Version: 2
Issued: 19, 03, 2015

Table of Contents

Warning

Symbols and conventions

Safety Information

Intended use

Package Contents

Unpacking

Product specifications

Overview

Using the microscope

Using the fluorescence

Maintenance

Troubleshooting

User replaceable accessories and spare parts

Equipment disposal

Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely.

Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

Package Contents

DESCRIPTION	QUANTITY
Microscope stand with nosepiece, stage, condenser	1
Optical head (trinocular)	1
Objective E-PLAN IOS 4x	1
Objective E-PLAN IOS 10x	1
Objective E-PLAN IOS 20x	1
Objective E-PLAN IOS 40x	1
Objective E-PLAN IOS 100x (oil)	1
Eyepiece WF10x/20mm	2
Fluorescence attachment	1
Fluorescence external power supply	1
Mercury lamp house	1
Dust cover	1
Power supply output 6Vdc	1

Unpacking


The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.


Place the observation head onto the top of the arm and tighten the lock-screw. Insert the eyepieces into the eye tubes and fix them with the tiny screws on the side of the eye tubes. Remove the plastic protection film from the microscope stage.

Product specifications

Head:	Trinocular, 30° inclined, 360° rotating. Interpupillary adjustment 48-75 mm.
Eyepieces:	WF10X/20mm.
Nosepiece:	5-position reversed revolving nosepiece. Ball bearing rotation.
Objectives:	IOS E-PLAN 4x/0.10, 10x/0.25, 20x/0.40, 40x/0.65, 100x/1.25 (oil).
Focusing system:	Coaxial coarse and fine.
Stage:	Double layer mechanical sliding stage, dimensions 216x150mm, moving range 78x54mm. Belt-drive in X direction..
Condenser:	Abbe condenser, sliding-in, N.A. 1.25 with centering system.
Illumination:	Transmitted light: X-LED ³ , with manual brightness control. Epi-fluorescence: HBO 100W high pressure mercury bulb.
Power supply:	External power supply: Input 100-240Vac 50-60Hz / Output 6Vdc 1A Fluorescence Power Supply 100W. Timer and current display.

Overview

Using the microscope

Adjust the observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw.

Place the specimen on the stage

Lock the specimen slide on the mechanical stage using the slide clamp. Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage

Lamp settings

The microscope is fitted with a white LED illuminator. Before turning on the light, read the section about electrical safety precautions. Connect the external power supply into the socket at the base. Turn the brightness adjustment knob to a brightness suitable for observation.

Adjust interpupillary distance

Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen.

Focus and diopter adjustment


Turn the dioptic adjustment ring on the left eyepiece to the zero position. Turn the coarse focus knob in order to focus the slide with an objective with low magnification. Adjust the fine focus knob until you obtain a clear and defined picture observing with the right eye, and then act on the left dioptic compensation ring observing with the left eye. When the image appears in focus, choose the necessary objective with the revolving nosepiece. Turn the tension-adjust-knob to get a suitable tension for the focus system.

Condenser adjustment

Turn the condenser adjustment knob to move the condenser up or down for a clear and uniform sample illumination. The two condenser centring screws can be used to centre the condenser.

Set the numerical aperture

Adjust the aperture of the iris diaphragm to set the numerical aperture of the illuminator, thus controlling image contrast and resolution.

Using the fluorescence

Assembling the epi-fluorescence attachment

Extract the attachment from the special packaging and place it on top of the stand of the microscope tightening the locking screw. Then place the optical head above the attachment, setting it with the appropriate locking screw. In the rear part of the attachment slide the lamp housing and connect the cable to the external power supply. Open the lamp housing door and insert the mercury lamp (included in the package) between the two terminals (note that the two poles of the lamp have different sizes so that the insertion in the terminals is unique).

Use and alignment of the lamp

Press the ON/OFF button on the external power supply. Check that the indicator of current after a few seconds shows a value equal to approximately 4A. If the indicator has a value of 0A, press the TRIGGER button. Place on the stage your fluorescence sample and focus. Act on the lamp centering screws until you reach the maximum intensity of emitted fluorescence from the sample. Finally turn the lamp focusing knob in order to achieve uniformity of lighting on the entire field of view.

Warning: for optimal use of the fluorescence wait 10-15 minutes since power-up in order to allow the correct thermal equilibrium of the lamp.

Insertion and replacement of the mercury lamp

If the current drops below 3.5A, replace the lamp. Lifetime is approximately 400 hours.

Take care to cool the lamp for at least 30 minutes before opening its door.

After replacing the lamp, press the reset button to reset the time usage counter.

Unplug any cable from the epi-illuminator, then unscrew the lamp door and open it.

The lamp holder will come out together with the door. Remove the old lamp or the plastic piece (in case this is the first insertion).

It is recommended to use protective glasses because, if not handled with due caution, lamp could crash and injure the operator with glass splinters. Also do not touch the lamp glass during the extraction out of its box, taking care to not exert pressure on any part of it. Plug one lamp end into the flexible metal holder (it fits only in one way), and then insert the other end in the fixed part.

Screw it just enough to keep the lamp in position. If it is too tight, the lamp could break. Replace the lamp door in its position and connect the cable on the side of the housing to the external power supply.

Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
LIGHT DOESN'T TURN ON	Power supply not connected	Check that the 6Vdc power supply jack is well inserted on the rear of the microscope
		Rotate the brightness adjustment control and check if there is an increase in the light output
IMAGE CANNOT BE SEEN OR IS DARK.	The iris diaphragm aperture is not completely opened	Completely open the iris diaphragm aperture
	Brightness level is low	Rotate the brightness adjustment potentiometer
	Objective is not aligned with the optical axis	Rotate the nosepiece until an objective is well inserted in the optical path (it "clicks")
IMAGE IS UNCLEAR, BLURRED OR HAS INSUFFICIENT CONTRAST.	Objectives or filters are dirty	Wipe them clean
	The iris diaphragm aperture is not opened correctly	Open the iris diaphragm aperture completely
	Condenser at wrong height	Rotate the condenser knob until you see a uniform illumination

User replaceable accessories and spare parts

CAT. NO.	DESCRIPTION
M-160	Eyepiece WF10x/20mm.
M-161	Eyepiece WF15x.
M-162	Eyepiece WF20x.
M-163	Micrometer eyepiece WF10x/20mm.
M-005	26x76 mm micrometric slide. Range 1 mm, div. 0,01 mm.
M-144	Objective 4x/0,10 E-PLAN IOS.
M-145	Objective 10x/0,25 E-PLAN IOS.
M-146	Objective 20x/0,40 E-PLAN IOS.
M-147	Objective 40x/0,65 E-PLAN IOS
M-335	Objective PLAN IOS MET 50x/0.75
M-149	Objective 60x/0,80 E-PLAN IOS.
M-148	Objective 100x/1,25 E-PLAN IOS (Oil).
M-181	Complete Phase Contrast Set with IOS PLAN obj. 10x, 20x, 40x, 100x , with Darkfield position
M-174.1	Polarizing set, filters only (for B-380 series).
M-175	Rotating table for polarizing set.
M-185	Darkfield condenser for dry objectives.
M-173	Photo adapter for APS-C and Full Frame Reflex cameras
M-114	C-Mount adapter for 1/2" sensor.
M-116	C-Mount adapter for 2/3" sensor.
M-069	Solar battery pack.
M-151	HBO100W high-pressure mercury bulb for fluoresc. (for B-383FL only).

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.