B-380 Series

INSTRUCTION MANUAL

Model
B-383LD1
B-383LD2

Version: 2
Issued: 18, 03, 2015

Table of Contents

Warning

Symbols and conventions

Safety Information

Intended use

Package Contents

Unpacking

Product specifications

Overview

Using the microscope

Using the fluorescence

Maintenance

Troubleshooting

User replaceable accessories and spare parts

Equipment disposal

Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely.

Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

Package Contents

DESCRIPTION	QUANTITY
Microscope stand with nosepiece, stage, condenser	1
Optical head (trinocular)	1
Objective E-PLAN IOS 4x	1
Objective E-PLAN IOS 10x	1
Objective E-PLAN IOS 20x	1
Objective E-PLAN IOS 40x	1
Objective PLAN IOS MET 50x	1
Eyepiece WF10x/20mm	2
Fluorescence attachment	1
Dust cover	1
Power supply output 6Vdc	1

Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Place the observation head onto the top of the arm and tighten the lock-screw. Insert the eyepieces into the eye tubes and fix them with the tiny screws on the side of the eye tubes. Remove the plastic protection film from the microscope stage.

Product specifications


B-383LD1


Head:	Trinocular, 30° inclined, 360° rotating. Interpupillary adjustment 48-75 mm.
Eyepieces:	WF10X/20mm.
Nosepiece:	5 position on reversed revolving nosepiece. Ball bearing rotation.
Objectives:	IOS E-PLAN 4x/0.10, 10x/0.25, 20x/0.40, 40x/0.65, and PLAN 50x/0.75 (no cover slide).
Focusing system:	Coaxial coarse and fine..
Stage:	Double layer mechanical sliding stage, dimensions 216x150mm, moving range 78x54mm. Belt-drive in X direction.
Condenser:	Abbe condenser, sliding-in, N.A. 1.25 with centering system.
Illumination:	Transmitted light: X-LED ³ , with manual brightness control. Epi-fluorescence: high-power blue LED.
Power supply:	External power supply: Input 100-240Vac 50-60Hz / Output 6Vdc 1A.

B-383LD2

Head:	Trinocular, 30° inclined, 360° rotating. Interpupillary adjustment 48-75 mm.
Eyepieces:	WF10X/20mm.
Nosepiece:	5-position reversed revolving nosepiece. Ball bearing rotation.
Objectives:	IOS E-PLAN 4x/0.10, 10x/0.25, 20x/0.40, 40x/0.65, and PLAN 50x/0.75 (no cover slide).
Focusing system:	Coaxial coarse and fine.
Stage:	Double layer mechanical sliding stage, dimensions 216x150mm, moving range 78x54mm. Belt-drive in X direction.
Condenser:	Abbe condenser, sliding-in, N.A. 1.25 with centering system.
Illumination:	Transmitted light: X-LED ³ , with manual brightness control. Epi-fluorescence: high-power white LED (for B and G filtersets).
Power supply:	External power supply: Input 100-240Vac 50-60Hz / Output 6Vdc 1A.

Overview

Using the microscope

Adjust the observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw.

Place the specimen on the stage

Lock the specimen slide on the mechanical stage using the slide clamp. Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage

LED lamp settings – Transmitted illumination

The microscope is fitted with high brightness white LED illuminator. Press the selector switch on the back of the main body to I position. Turn the brightness adjustment knob to a brightness suitable for observation.

LED lamp settings – Epillumination

The microscope is fitted with high brightness LED illuminator. Press the selector switch on the back of the main body to II position. Turn the brightness adjustment knob to a brightness suitable for observation.

Adjust interpupillary distance

Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen.

Focus and diopter adjustment


Turn the dioptic adjustment ring on the left eyepiece to the zero position. Turn the coarse focus knob in order to focus the slide with an objective with low magnification. Adjust the fine focus knob until you obtain a clear and defined picture observing with the right eye, and then act on the left dioptic compensation ring observing with the left eye. When the image appears in focus, choose the necessary objective with the revolving nosepiece. Turn the tension-adjust-knob to get a suitable tension for the focus system.

Condenser adjustment

Turn the condenser adjustment knob to move the condenser up or down for a clear and uniform sample illumination. The two condenser centring screws can be used to centre the condenser.

Set the numerical aperture

Adjust the aperture of the iris diaphragm to set the numerical aperture of the illuminator, thus controlling image contrast and resolution.

Using the fluorescence

Assembling the epi-fluorescence attachment

Take the attachment out of its packaging and put it onto the microscope stand, then tighten the attachment locking screw.

Place the optical head onto the top of the attachment and tighten the optical head locking screw. Connect the cable from the epi-fluorescence attachment to the output jack on the rear of the microscope.

Lighting up the fluorescence LED light

Connect the external 6Vdc power supply to a wall socket, and plug its output to the jack on the rear of the microscope:

In order to light up the epi-fluorescence LED, press the main switch on the rear in position II. Then rotate the brightness adjustment knob on the left to the desired value. Pull the filter selector (1) to right or left in order to insert the fluorescence filter into the optical path. Put the filter selector (1) in the center if you want to use the brightfield transmitted illumination. Unlike a mercury lamp system, B-383LD LED illumination doesn't need any power-up time for heating, and can be used immediately after switching on. Also, the LED source is pre-aligned in factory and doesn't need any alignment operation.

Looking at the fluorescent specimen

Focus on your sample, and adjust the light intensity as needed through the brightness adjustment knob. In order to improve the darkness of the background (thus improving contrast), it is strongly suggested to put the provided black plate under the stage:

Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

Troubleshooting

Review the information in the table below to troubleshoot operating problems.


PROBLEM	CAUSE	SOLUTION
LIGHT DOESN'T TURN ON	Power supply not connected	Check that the 6Vdc power supply jack is well inserted on the rear of the microscope
		Rotate the brightness adjustment control and check if there is an increase in the light output
IMAGE CANNOT BE SEEN OR IS DARK.	The iris diaphragm aperture is not completely opened	Completely open the iris diaphragm aperture
	Brightness level is low	Rotate the brightness adjustment potentiometer
	Objective is not aligned with the optical axis	Rotate the nosepiece until an objective is well inserted in the optical path (it "clicks")
IMAGE IS UNCLEAR, BLURRED OR HAS INSUFFICIENT CONTRAST.	Objectives or filters are dirty	Wipe them clean
	The iris diaphragm aperture is not opened correctly	Open the iris diaphragm aperture completely
	Condenser at wrong height	Rotate the condenser knob until you see a uniform illumination

User replaceable accessories and spare parts

CAT. NO.	DESCRIPTION
M-160	Eyepiece WF10x/20mm.
M-161	Eyepiece WF15x.
M-162	Eyepiece WF20x.
M-163	Micrometer eyepiece WF10x/20mm.
M-005	26x76 mm micrometric slide. Range 1 mm, div. 0,01 mm.
M-144	Objective 4x/0,10 E-PLAN IOS.
M-145	Objective 10x/0,25 E-PLAN IOS.
M-146	Objective 20x/0,40 E-PLAN IOS.
M-147	Objective 40x/0,65 E-PLAN IOS
M-149	Objective 60x/0,80 E-PLAN IOS.
M-148	Objective 100x/1,25 E-PLAN IOS (Oil).
M-335	Objective PLAN IOS MET 50x/0.75.
M-173	Photo adapter for APS-C and Full Frame Reflex cameras.
M-114	C-Mount adapter for 1/2" sensor.
M-069	Solar battery pack.
M-181	Complete Phase Contrast Set with IOS PLAN obj. 10x, 20x, 40x, 100x , with Darkfield position.
M-174.1	Polarising set, filters only (for B-380 series).
M-175	Rotating table for polarising set.
M-185	Darkfield condenser for Dry Objectives.
M-116	C-Mount adapter for 2/3" sensor.

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.