**B-383 Series**

# INSTRUCTION MANUAL

<b>Model</b>
B-383MET

Ver. 2.4    2024

---

## Table of Contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
4.	Unpacking	5
5.	Intended use	5
6.	Symbols and conventions	5
7.	Instrument description	6
8.	Assembling	8
9.	Summary of transmitted light observation procedures	10
10.	Use of the microscope in transmitted light	11
10.1	Turning on the instrument	11
10.2	Light intensity adjustment	11
10.3	Adjusting the interpupillary distance	11
10.4	Diopter adjustment	11
10.5	Coarse focus tension adjustment	12
10.6	Coarse upper limit lever	12
10.7	Stage	12
10.7.1	Use of the glass plate	12
10.8	Centering the condenser	13
10.9	Aperture diaphragm	13
11.	Summary of reflected light observation procedures	14
12.	Use of the microscope in reflected light	15
12.1	Turning on the instrument	15
12.2	Light intensity adjustment	15
12.3	Use of field diaphragm (FS)	15
12.4	Effects of field diaphragm	15
12.5	Use of aperture diaphragm (AS)	16
12.6	Use of color filters	16
12.7	Use of polarized light	16
13.	Microphotography	18
13.1	Use of C-mount cameras	18
13.2	Use of reflex cameras	18
14.	Maintenance	19
15.	Troubleshooting	20
	Equipment disposal	22

---

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

## 2. Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

### 3. Package content

- ① Microscope frame
- ② Reflected light Illuminator
- ③ Observation head
- ④ Eyepieces
- ⑤ Objectives
- ⑥ Tension adjustment tool
- ⑦ Allen wrench
- ⑧ Color filters
- ⑨ Analyzer + polarizer
- ⑩ Glass plate
- ⑪ Dust cover
- ⑫ Power supply

---

## 4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 5. Intended use

### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

## 6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION


This symbol indicates a potential risk and alerts you to proceed with caution.


### ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

## 7. Instrument description

Opposite side

## 8. Assembling

1. Insert the reflected light illuminator ① on the frame and lock the locking knob ② with the provided Allen wrench. (Fig. 1)

2. Connect the illuminator cable ③ to the socket ④ placed on the back side of the frame. (Fig. 2)

3. Insert the optical head above the epi-illuminator and tighten the locking screw. (Fig. 3)
- **Hold the head with one hand during the locking in order to avoid that the head falls.**

4. Insert both eyepieces into the tubes of the optical head. (Fig. 4)

5. Insert the polarizer. (Fig. 5)


6. Remove the protection cap and insert analyzer. (Fig. 6-7)

7. Insert the power supply jack in the connector placed on the back side of the frame. (Fig. 8)

## 9. Summary of transmitted light observation procedures

## 10. Use of the microscope in transmitted light

### 10.1 Turning on the instrument

Turn the main switch ① on the rear of the stand to "I". (Fig. 9)

- There is a three-position switch on the rear of the stand: position "I" turns on the transmitted light, position "II" turns on the reflected light and position "O" turns off the microscope.

### 10.2 Light intensity adjustment

Operate on the light intensity adjustment knob ② to increase / decrease the illumination voltage. (Fig. 10)

### 10.3 Adjusting the interpupillary distance

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

- **The graduation on the interpupillary distance indicator ③, pointed by the spot "." on the eyepiece holder, shows the distance between the operator's eyes. (Fig. 11)**

The range of the interpupillary distance is 48-75 mm.

### 10.4 Diopter adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
  2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the diopter adjustment ring ④ to compensate. (Fig. 12)
- **The adjustment range is  $\pm 5$  diopter. The number indicated on the adjustment ring graduation should correspond to the operator's diopter correction.**

## 10.5 Coarse focus tension adjustment

- **Adjust the tension using the provided tool.**

The tension of the coarse focusing knob is factory preset.


1. To modify the tension according to personal's needs, rotate the ring ①. (Fig. 13)
- Clockwise rotation increases the tension.
  - If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

## 10.6 Coarse upper limit lever

The upper limit knob has two functions: prevent the contact between slide and objective and acts as "focus memory".

1. After focusing the specimen, pull the lever ② toward the front of the microscope and lock it. (Fig. 14)
- In this way the focus upper limit is set.
2. Now one can lower the stage with coarse focus knob, replace the specimen and raise again the stage up to the upper limit: specimen will be in approximate focus and will need a fine adjustment to get the proper focus.
- **Fine focus movement is not affected by the coarse focus lock.**
  - **To unlock, move the lever in the opposite direction to the one used for the locking.**

## 10.7 Stage

Stage accepts metallurgical samples with thickness max 15 mm, placed on a standard support 26 x 76 mm.

1. Open the spring arm of the slide holder ③ and place frontally the slides on the stage. (Fig. 15)
  2. Gently release the spring arm of the slide holder.
- **A sudden release of the spring arm could cause the falling of the slide.**

### 10.7.1 Use of the glass plate

When using embedded metallographic specimens it may be necessary to use the square glass plate.

1. Open the spring arm of the slide holder □ and place frontally the glass plate by aligning the positioning notch ④ so that the movable arm fits directly into the notch. (Fig. 16)
  2. Gently release the spring arm of the slide holder.
- **The glass plate has a positioning notch to facilitate insertion of the movable arm.**
  - **This notch is provided and does not indicate a glass breakage.**

## 10.8 Centering the condenser

- The condenser is mounted and pre-centered before shipment from the factory.
- To remove the condenser, use a 1.5 mm Allen key and turn the fixing screw on the right-hand side of the condenser holder.

If re-centering is necessary, proceed as follow:


1. Insert the 5x objective into the optical path (if the 5x objective is not available, use the lower magnification objective).
2. Focus the specimen.
3. Close the aperture diaphragm using the ① dial, moving the dial towards the "4" value for the 5X objective. (Fig. 17)
4. Raise the condenser to the end of its travel by turning the condenser height adjustment screw ② on the left side of the condenser holder.
5. Centre the condenser using the centering screws ③ until the field of view is evenly lit (no lighter or darker areas within the field of view).
6. Then open the diaphragm completely.


## 10.9 Aperture diaphragm

The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.

1. Move the aperture ring ① (Fig. 18) to the value corresponding to the objective in use. You can move the dial to lower or higher values to adjust the observation to your preference.
- With low contrast specimens set the numerical aperture value (printed on the condenser ring) to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's ring in order to obtain an image like the one in fig. 19.

## 11. Summary of reflected light observation procedures

## 12. Use of the microscope in reflected light

### 12.1 Turning on the instrument

Turn the main switch ① on the rear of the stand to "II". (Fig. 20)

### 12.2 Light intensity adjustment


Operate on the light intensity adjustment knob ② to increase / decrease the illumination voltage. (Fig. 21)

### 12.3 Use of field diaphragm (FS)


• **The field and aperture diaphragms on the illuminator for reflected light are pre-centered before shipment from the factory and do not require further centering.**

1. Move the selector switch ③ on the reflected light illuminator toward left to position . (Fig. 22)
2. Place the sample on the stage, insert the 10x objective into the optical path and focus.
3. In normal use, open the aperture until the image circumscribes the field of view.

### 12.4 Effects of field diaphragm

The field diaphragm adjusts the illuminated area to obtain an image with high contrast. Adjust the field diaphragm according to the lens in use until the iris diaphragm circumscribes the field of view to eliminate unnecessary light to the eyepieces. (Fig. 23)

### 12.5 Use of aperture diaphragm (AS)

The numerical aperture value (N.A.) of the aperture diaphragm affects image contrast. Increasing or decreasing this value according to the numerical aperture of the lens changes the resolution, contrast and depth of field of the image.

1. Move the aperture ring ④ (Fig. 24) to open or close the diaphragm to the desired position.
- For low-contrast specimens, set the numerical aperture value to approximately 70%-80% of the objective's N.A. If necessary, remove an eyepiece and, looking into the empty eyepiece holder, adjust the condenser ring until an image like that in Fig. 19 is obtained.


Fig. 24

### 12.6 Use of color filters

1. Insert the filter ① in one of the two slots ② placed in the upper part of the reflected light illuminator. (Fig. 25).

Filter	Application
Blue	Converts the source color temperature to daylight
Yellow	Contrast filter for the observation of wafers and semiconductors
Frosted white	Reduces lighting irregularities by creating a homogeneous field of view

Fig. 25

### 12.7 Use of polarized light

1. Insert the polarizer all the way into the optical path. (Fig. 26)

Fig. 26


2. Insert the analyzer all the way into the optical path. (Fig. 27)

Fig. 27


3. Put a flat mirror on the stage and get into approximate focus the mirror surface.
4. While observing into the eyepieces, rotate the scale of the analyzer ① until you achieve the darkest image. (Fig. 28)
  - Now the light extinction (so-called "crossed Nicol position") is achieved.
5. Put a specimen on the stage, focus the specimen and begin observation.

## 13. Microphotography

### 13.1 Use of C-mount cameras

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 29)

2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 30)

### 13.2 Use of reflex cameras

1. Insert the Reflex adapter ① into the relay tube to the microscope ②.
2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
3. Connect the Reflex camera ④ to the "T2" ring just installed (Fig. 31).
4. Mount the other end of the relay tube ② into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 29)
  - "T2" ring is not provided along with the microscope, but is commercially available.
  - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
  - To calculate the magnification of the camera: objective magnification \* camera magnification \* lens magnification.
  - **If using an SLR camera, mirror movement may cause the camera to vibrate.**
  - **We suggest lifting the mirror, using long exposure times and a remote cord.**

---

## 14. Maintenance

### Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

### To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

### Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

**For the best results, use the OPTIKA cleaning kit (see catalogue).**

If you need to send the microscope to Optika for maintenance, please use the original packaging.

## 15. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
<b>I. Optical Section:</b>		
LED operates, but field of view remains dark	Power supply is unplugged	Connect
	Brightness is too low	Set brightness to a proper level
	Color filter or oblique illumination slider are in a wrong position	Move to a click stop
Field of view is obscured or not evenly illuminated	Revolving nosepiece is not correctly engaged	Make sure that the revolving nosepiece clicks properly into place
	Oblique illumination slider is in a wrong position	Move the slider to a fully open position
Dirt or dust is visible in the field of view	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture diaphragm is stopped down too far	Open aperture diaphragm
	The field diaphragm (both for reflected light and transmitted light) is not well centered	Set the diaphragm according to Koehler settings
Visibility is poor <ul style="list-style-type: none"> <li>• Image is not good</li> <li>• Contrast is poor</li> <li>• Details are indistinct</li> <li>• Image glares</li> </ul>	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture diaphragm is too closed or too open	Adjust aperture diaphragm
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly
	For reflected light observation, the specimen must not have any coverglass	Use specimens without coverglass
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position
One side of the image is unfocused	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Slide is mounted not in a flat position (tilted)	Place the specimen in a flat position on the stage
<b>II. Mechanical Section:</b>		
Coarse focus knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
Focus is unstable	Tension adjustment ring is too loose	Tighten tension adjustment ring
<b>III. Electrical Section</b>		
LED doesn't turn on	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection
<b>IV. Observation tube</b>		
Field of view of one eye does not match that of the other	Interpupillary distance is incorrect	Adjust interpupillary distance
	Incorrect diopter adjustment	Adjust diopter
	Your view is not accustomed to microscope observation	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope

---

<b>V. Microphotography</b>		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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**OPTIKA® S.r.l.**