B-380 Series

INSTRUCTION MANUAL

Model
B-383POL

Version: 3
Issued: 18, 03, 2015

Table of Contents

Warning

Symbols and conventions

Safety Information

Intended use

Package Contents

Unpacking

Product Specifications

Overview

Using the microscope

Maintenance

Troubleshooting

User replaceable accessories and spare parts

Equipment disposal

Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position.

Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely.

Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

Package Contents

DESCRIPTION	QUANTITY
Microscope stand with head, focus and rotating stage	1
WF10x/20 eyepiece	2
Bertrand lens attachment	1
Objective E-PLAN IOS POL 4x	1
Objective E-PLAN IOS POL 10x	1
Objective E-PLAN IOS POL 40x	1
Objective E-PLAN IOS POL 60x	1
Condenser with rotating polarizer	1
Retardation plates (λ , $\lambda/4$, quartz wedge)	3
Allen wrench for objective centering	2
Power supply 6Vdc	1
Dust cover	1

Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

1. First of all remove all the components from the original package and store the box in a safe and dry place.
2. Place the main body of the microscope on a flat and sturdy surface.

Don't connect the power supply cable at this stage.

3. Place the attachment with analyzer /Bertrand Lens onto the main body and turn the head locking screw (1).

4. Insert the optical head and lock the screw.

5. After removing the protection caps, fully insert both eyepieces.

6. Connect the external power supply cable to the socket on the back of the microscope.


If there is no need to install additional accessories, the instrument is now ready for use.


Product specifications

Head	Trinocular, 30° inclined, 360° rotating.
Eyepieces	WF10X/20mm with cross hair.
Bertrand lens	Swing-out type; centrable.
Polarizing attachment	0°-90° rotating analyzer filter. Tint plates included: 1° order red (λ), $\lambda/4$, quartz wedge.
Nosepiece	4-positions with centering mechanism for all objectives.
Objectives	E-PLAN IOS POL (Strain-free) 4x/0.10, 10x/0.25, 40x/0.65, 60x/0.80.
Focusing system	Coaxial coarse and fine.
Stage	160mm dia.; 360° rotating with stop knob and 0.1° Vernier.
Condenser	N.A. 1.25, with iris diaphragm, focusable and centrable. With rotating polarizer filter.
Illumination	X-LED ³ , with manual brightness control.
Power supply	External power supply: Input 100-240Vac 50-60Hz / Output 6Vdc 1A

Overview

Using the microscope

Observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw.

Place the specimen on the stage

Lock the specimen slide on the mechanical stage using the slide clamp. Ensure that the specimen is centred over the stage opening by adjusting the coaxial knobs of the stage.

Illumination system settings

The microscope is fitted with a white LED illuminator. Turn the brightness adjustment knob to a brightness suitable for observation.

Adjust interpupillary distance

Hold the right and left parts of the observation head with both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen.

The white dot (°) placed on the right eyepiece shows the set interpupillary distance. Just remember this value to help on later settings.

Focus tension adjustment

The tension of the coarse focusing knob is preset by factory.

To change the tension according to your preference, just rotate the knob (3) clockwise in order to increase it.

Excessive tension could damage the mechanism of focus. A too loosed tension causes the descent of the stage by gravity or a sudden loss of focus. In this case, rotate the knob (3) to increase the tension.

TENSION
KNOB (3)

COARSE FOCUS

FINE FOCUS

Dioptric adjustment

Turn the dioptric adjustment ring on the left eyepiece halfway. Turn the coarse focus knob in order to focus the slide with an objective with low magnification. Adjust the fine focus knob until you obtain a clear and defined picture observing with the right eye, and then repeat the operation with the left dioptric compensation ring and the left eye. When the image appears in focus, choose the necessary objective with the revolving nosepiece.

Condenser

The condenser, placed below the stage, is equipped with the following features:

APERTURE DIAPHRAGM
ADJUSTMENT RING (4)

ROTATING POLARIZER WITH
ANGULAR SCALE (9)

Raise or lower the condenser through the knob (5) to obtain a clear and uniform illumination of the sample. To center the condenser: completely close the iris diaphragm (4). Using the condenser centering screws (2), move the diaphragm in the center of the field of view. Then gradually expand the diaphragm until it is tangent to the edges of the field of view. If necessary, you can perform an additional adjustment. The condenser is centered when the edges of iris diaphragm are tangent to the field of view.

In the bottom part of the condenser there is a polarizer filter, which can be rotated through 360°, with a reference notch for a precise angle measurement.

Numerical aperture setting

The value of the numerical aperture (N.A.) of the diaphragm is an indication of the contrast of the illumination system. Matching the value of illumination system's N.A. with that of the objective ensures the best results in terms of contrast and image quality.

To set the numerical aperture of the illuminator, adjust the opening of the iris diaphragm (4). In this way you control contrast and image resolution.

For samples with low contrast set the iris to about 75% of the value of the objective's numerical aperture.

HEIGHT ADJUSTMENT FOR
CONDENSER (5)


System centering

In order to correctly center the whole optical system, first remove the Bertrand lens out of the path by rotating its disc. Remove also the lower polarizer filter.

The rotating stage is factory pre-centered and so it doesn't need any additional operation.

It's necessary to center the objectives following this procedure:

1. Place a specimen on the stage.
2. Insert 10x objective.
3. Insert the two little wrenches into the holes for nosepiece centering (7).
4. Focus your sample while continuously rotating the stage clockwise and then counter-clockwise by a little angle (e.g. 30° or 45°). During these oscillations you should spot a point on the sample that doesn't rotate on a circumference but only revolving about itself (fixed point). With the two objective's centering wrenches, bring this point in the center of the crossline.
5. Insert the 20x objective.
6. Insert the two little wrenches into the holes for nosepiece centering (7), above the 20x objective.
7. Rotate the two wrenches until the same point of the sample is brought in the center of the crossline.
8. Repeat the operation for the other objectives.

Polarized light observation

In order to correctly use the polarized light, it's necessary to perform a setting for "light extinction":

1. Remove from the optical path the sample, retardation plates, etc...
2. Insert the 10x objective.
3. Insert in the optical path the analyzer filter and rotate the rotation scale (8) in correspondence to position "0".
4. Put the angular scale (9) of the polarizer filter to "0" position.
5. Loosen the lock screw of the polarizer and, looking into the eyepieces, rotate the scale of the polarizer until you obtain the complete extinction of the light (totally dark).
6. Then lock the polarizer screw.

Bertrand lens

Bertrand lens is employed in conoscopic observation mode (lens inserted) or orthoscopic mode (lens not inserted).

Orthoscopic observation

The path of the polarized light is parallel to the optical axis, thus allowing the observation of the optical characteristics of the specimen. You can use objectives from 4x to 100x.

1. Remove the Bertrand lens out from optical path.
 2. Insert the analyzer filter using the proper lever (12) and look at the sample.
 3. Insert a retardation plate (11) pushing it completely inside the slot under the head.
- To remove the plate from the optical path, pull it outward till the first click.

Conoscopic observation

1. Use objectives from 20X to 100X.
2. Insert the analyzer using the proper lever (12) and put it in the extinction position.
3. Insert the Bertrand lens.
4. Insert the objective of choice (from 20X to 100X).
5. Open the aperture diaphragm.
6. It's possible to center the conoscopic image adjusting the screws (10) placed under of the head.


IMAGE CONOSCOPIC

In this manner it's possible to look at the conoscopic images that form on the back focal plane of the objective. If the peripheral part of the image is too dark, you can just move the condenser vertically to find the optimal position for the illumination.

Retardation plates

Retardation plates (1 lambda, corresponding to 530nm or Red 1° order, 1/4 lambda corresponding to 137nm, and quartz wedge) are used to make the polarization analysis more precise and refined.

Insert or remove the plate in the optical path of the microscope by pushing it into the slot (11) under the head. Retardation plates have no visible effect when used in normal brightfield light.

Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
LIGHT DOESN'T TURN ON	Power supply not connected	Check that the 6Vdc power supply jack is well inserted on the rear of the microscope
		Rotate the brightness adjustment control and check if there is an increase in the light output
IMAGE CANNOT BE SEEN OR IS DARK.	The iris diaphragm aperture is not completely opened	Completely open the iris diaphragm aperture
	Brightness level is low	Rotate the brightness adjustment potentiometer
	Objective is not aligned with the optical axis	Rotate the nosepiece until an objective is well inserted in the optical path (it "clicks")
IMAGE IS UNCLEAR, BLURRED OR HAS INSUFFICIENT CONTRAST.	Objectives or filters are dirty	Wipe them clean
	The iris diaphragm aperture is not opened correctly	Open the iris diaphragm aperture completely
	Condenser at wrong height	Rotate the condenser knob until you see a uniform illumination

User replaceable accessories and spare parts

CAT. NO.	DESCRIPTION
M-160	Eyepiece WF10x/20mm.
M-161	Eyepiece WF15x.
M-162	Eyepiece WF20x.
M-163	Micrometer eyepiece WF10x/20mm.
M-005	26x76 mm micrometric slide. Range 1 mm, div. 0,01 mm.
M-144P	Objective 4x/0,10 E-PLAN IOS POL.
M-145P	Objective 10x/0,25 E-PLAN IOS POL.
M-146P	Objective 20x/0,40 E-PLAN IOS POL.
M-147P	Objective 40x/0,65 E-PLAN IOS POL.
M-149P	Objective 60x/0,80 E-PLAN IOS POL.
M-148P	Objective 100x/1,25 E-PLAN IOS (Oil).
M-185	Darkfield condenser for dry objectives.
M-173	Photo tube adapter for reflex camera.
M-114	CCD camera adapter.
M-069	Solar battery pack.

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.