EDUCAM STUDENT

INSTRUCTION MANUAL

	Model
	4083.2
	4083.3
	4083.5

Version: 2
Issued: 04, 09, 2014

Description

The **EDUCAM STUDENT** video camera is especially designed to meet the various requirements in the educational field.

When it is connected to a professional monitor or simply to a TV set, **EDUCAM** is able to carry out many different functions. It can be used in the following ways:

- as an episcopo, for the reproduction of texts, documents, photographs, etc.
- as an enlarger, to enlarge small objects, insects, minerals, etc.
- for video-microscopy, connectable both to biological or stereo microscopes using the adapters included.
- as a overhead projector, for the projection of tracings
- as a camera for teleconferences, that can be connected to a computer and to the Internet
- as a camera for video-conferences, used in conferences, assemblies, meetings, etc.
- as a camera for filming, with the help of a video recorder.

This video camera enables you to record clearly even when there is little light available, thanks to its extremely sensitive electronic circuits. The special lens enables you to focus from 0,76 mm, up to an infinite distance.

The special flexible arm and the heavy weight of the base, gives the **EDUCAM** its versatile, sturdy and stabile characteristics. The camera head, placed on the end of the flexible arm, can be placed in any position, even projecting from the base, without jeopardising the system's stability. A special joint allows the head to rotate avoiding connector breakages inside the flexible arm.

Overview 4083.5

4083.5

Overview

Technical specifications

Sensor	1/3" colour CCD element
Horizontal resolution	420 TV lines
Total resolution	244 k pixel
Objective	8 mm, F:1.2
Sensitivity	min. 0,08 lux at F:1.2
Picture controls	electronic shutter
Gain control	automatic
White balance	automatic
Video signal	PAL compatible - VHS
Working distance	from 0,76 mm to infinite
Magnifications	90x at 0,76 mm (with 21" monitor)
Power supply	12 V DC, 150 mA, with external power adapter
Switching on	via 2-position switch
PC link	standard USB 2.0
Length of flexible arm	50 cm (Student model), 65 cm (Student Pro model)
Base diameter	17 cm
Weight	3.4 Kg

Assembly

Power supply

The Educam Student video camera requires 12 V DC working voltage and comes with an external power supply. Attach the video camera output cable of the Educam Student video camera to the external power supply socket, and plug the power supply directly into a 220 V mains socket.

Video signal

The video signal from the Educam video camera is transmitted to a TV monitor. The video output signal cable can be found on the back of the camera base. It is fitted with a yellow plug that should be plugged into the monitor/colour TV directly or via the SCART adapter provided (in the yellow video socket).

Finally, if using the SCART adapter, it has to be inserted into the SCART socket at the back of the TV set. See the TV manufacturer's instruction manual for details of how to switch on the TV and use its SCART circuit.

Operation

Switching on and off the video camera

The Educam video camera has a switch at the rear of its base. Switch to "1" to turn on the video camera and to "0" to turn it off.

Using the video camera


Position the camera head with the lens and the electronic sensor directly in front of the object that you wish to observe using its flexible arm. The lens provided with the camera has a manual focusing and must therefore be adjusted every time the distance between the object and the camera is changed. This can be done by simply turning the lens on the camera head. The object must be at least 0.76 mm away from the lens. There are no limits to the maximum distance.

Capturing video from a microscope


The Educam video camera has been designed for effective use with microscopes.

The following special adapters (provided) are required for this application:

Adapter for biological microscopes: To be used directly in the optical tube of the microscope after having extracted the eyepiece. This adapter contains a WF10x eyepiece that needs to be replaced with the one belonging to the microscope.


Adapter for stereo microscopes: To be used directly on one of the optical cylinders on a stereo microscope after having extracted the eyepiece. The biological microscope adapter needs to be inserted inside this adapter.

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.