LAB Series

INSTRUCTION MANUAL

	Model
	LAB-10
	LAB-20
	LAB-30

Version: 2
Issued: 10, 01, 2017

Table of Contents

Warning

Symbols and conventions

Safety Information

Intended use

List of accessories and spare parts

Overview

Unpacking

Assembling

Maintenance

Troubleshooting

Equipment disposal

Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.


Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

List of accessories and spare parts

CAT. NO.	DESCRIPTION
ST-001.1	Eyepieces (pair) wf5x/22mm.
ST-401	Eyepieces (pair) WF10x/20 mm.
ST-402	Eyepieces (pair) WF15x/15 mm.
ST-403	Eyepieces (pair) WF20x/10 mm.
ST-405	Eyepiece micrometer WF10x/20mm.
ST-085	Additional lens 0.5x (w.d. 165mm) (Only for LAB-20 and LAB-30.)
ST-086	Additional lens 1.5x (w.d. 47mm) (Only for LAB-20 and LAB-30.)
ST-415	Eyecups (pair).
DC-002	Plastic dust cover, medium
ST-092	Protective glass for stereohead.
ST-417	Round plastic stage
M-114	C-Mount adapter for 1/2" sensor.
M-115	C-Mount adapter for 1/3" sensor.
M-116	C-Mount adapter for 2/3" sensor.
M-173	Photo adapter for APS-C and Full Frame Reflex cameras.
M-113.1	Ring adapter, 30mm (for monocular and binocular microscopes).
ST-418	Focusable C-Mount adapter for 1/3" sensor.
ST-419	Focusable C-Mount adapter for 1/2" sensor.
15104	Lens cleaner, 50ml.
VP-LAB	IQ/OQ/PQ Validation Protocols.

Overview

Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Assembling

Once you open the box, these are the microscope's components:

1. Insert the stereo head into the holder on the stand, and tighten the side screw. (Fig.1)


2. Insert the eyepieces into the tubes of the head. (Fig.2)

Using the microscope

1. Place the specimen on the stage

Place the specimen on the microscope stage and lock the specimen using the slide-clamp if it is necessary. Ensure that the specimen is centred over the stage opening. (Fig.3)

2. Turn on the light

The microscope comes with an electrical illuminator. Insert the plug of the cable into the power socket and turn on the switch on the main body and select your light source. (Fig.4) Turn the light intensity control to a brightness suitable for observation. (Fig.5)

3. Adjust interpupillary distance

Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. If two circles appear, the interpupillary distance is too big, and if two overlapped circles appear, the interpupillary distance is too small. (Fig.6)

4. Focus

Rotate the focusing knob to bring the sample into focus with the highest magnification. (Fig.7)

5. Dioptic compensation

This compensation makes it possible for people with glasses to adjust the microscope to their eyes and use the microscope without glasses. Turn the zoom down to the lowest magnification. Adjust the diopter compensation ring of the right eyepiece until the image of the right eyepiece is clear and sharp. Repeat the procedure for the left eyepiece. Then, check the focus of the image for the whole zoom range. It should now be perfectly parfocal (focus is always maintained during the change of magnification). (Fig.8)

6. Magnification

Select the desired magnification by adjusting the zoom knob. Change the eyepieces and/or add an appropriate magnifying objective lens if necessary. (Fig.9)

Total magnification used can be calculated as:


Eyeiece magnification x Zoom magnification x Objective lens magnification.

7. Contrast disc

You can use the black/white disc in order to enhance the contrast of the image when using the incident illumination. (Fig.10).

If you want to look at bright opaque objects, place the disc with the black side up. (Fig.11)


8. Working distance

Normal working distance for the standard configuration is 80 mm (LAB-10) and 85mm (LAB-20 & LAB-30). (Fig.12)

9. Video capturing

Two kinds of observation heads, binocular or trinocular, are available for the LAB series. They can all be connected to cameras via an adaptor.

Please refer to the adaptor and camera manuals for further details.

Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
The illumination is ON, but the field of view is dark.	The plug is not connected to the illumination	Connect the cable
	The brightness is too low	Adjust to a proper setting
The edge of the field of view is vignetted or the brightness is asymmetric.	The incident illuminator is not correctly oriented	Change the angle of the incident illuminator
Dust and stains can be seen in the field of view.	There are stains and dust on the specimen	Clean the specimen
	There are stains and dust on the eyepiece	Clean the eyepiece
Poor image quality: The image is not sharp The contrast is not high The details are not clear	The lenses (condenser, objective, eyepieces) is dirty	Thoroughly clean all the optical system
One side of the image is out of focus.	The specimen is out of place (tilted)	Place the specimen flat on the stage.
II. Mechanical Section:		
The coarse focus knob is hard to turn.	The tension adjustment collar is too tight	Loosen the tension adjustment collar
The focus is unstable.	The tension adjustment collar is too loose	Tighten the tension adjustment collar
III. Electric section:		
The LED doesn't turn on.	No power supply	Check the power cord connection
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord
IV. Viewing tube assembly:		
The field of view of the two eyes is different	The interpupillar distance is not correct	Adjust the interpupillar distance
	The dioptic correction is not right	Adjust the dioptic correction
	The viewing technique is not correct, and the operator is straining the eyesight	When look into the objective, do not stare at the specimen but look at the whole field of view. Periodically, move the eyes away to look at a distant object, then back into the objective
V. Microphotography and video:		
The image is unfocused	Incorrect focussing	Adjusting the focus system as in the present manual
The edge of the image is unfocussed	To some degree, it is inherent to the nature of achromatic objectives	The problem can be minimized by a correct setting of the aperture diaphragm
Bright patches appear on the image	Stray light is entering the microscope through the eyepieces and through the camera viewfinder	Cover the eyepieces and the viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA S.r.l.