B-1000 Series

INSTRUCTION MANUAL

Model
B-1000
B-1000BF
B-1000PH
B-1000TI-2
B-1000TI-3
B-1000TI-5
B-1000TI-10

Ver. 2.1 2019

Summary

1. Warning	3
2. Symbols and conventions	3
3. Safety Information	3
4. Intended use	3
5. Overview	4
5.1 Manual version	4
5.2 Motorised version	6
6. Unpacking	8
7. Assembling	8
7.1 Manual version	8
7.2 Motorised version	9
7.3 Assembling the microscope	10
7.4 Only for motorised version	12
8. Summary of brightfield observation procedures	13
9. Use of the microscope	14
9.1 General switch on	14
9.2 Control keyboard	14
9.3 Brightness adjustment	14
9.4 Adjust the observation head	15
9.5 Adjust the interpupillary distance	15
9.6 Diopter adjustment	15
9.7 Use of eyeshields	16
9.8 Focus tension adjustment	16
9.9 Focus stop lever	16
9.10 Stage	17
9.11 Condenser centering	17
9.12 Use of oil immersion objective	18
9.13 Use of ALC system	19
9.14 Nosepiece rotation	19
9.15 Focusing	19
9.16 Stage	19
10. Use of universal condenser for brightfield / darkfield / phase contrast	20
10.1 Brightfield observation (BF)	20
10.2 Darkfield observation (DF)	20
10.3 Phase contrast observation (PH)	21
10.4 Use of green filter	22
11. DIC observation	23
11.1 Koehler DIC transmitted light	23
1.2 Nomarski DIC transmitted light	24
12. Microphotography	26
12.1 Installing the "C" mount adapter	26
12.2 Use of reflex camera	26
13. Maintenance	27
14. Troubleshooting	28
Equipment disposal	30

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

3. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

4. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.


IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

5. Overview

5.1 Manual version

5.2 Motorised version

Only the parts related to the motor are highlighted; all the other components of the microscope remain unchanged compared to the manual version.

X-AXIS MOVEMENT
KNOB (MANUAL
MOVEMENT)

STAGE
CONNECTION
CABLE

NOSEPIECE
ROTATION
BUTTONS

6. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

7. Assembling

7.1 Manual version

Once opened the box, the microscope parts (manual version) are the following:

- ① Frame
- ② Objectives
- ③ Stage
- ④ Condenser

- ⑤ Observation head
- ⑥ Eyepieces
- ⑦ ALC system (M-1030) (Optional)
- ⑧ Power supply

7.2 Motorised version

Once opened the box, the microscope parts (motorised version) are the following:

- ① Frame
- ② Objectives
- ③ Stage
- ④ Condenser
- ⑤ Observation head
- ⑥ Eyepieces

- ⑦ ALC system (M-1030) (Optional)
- ⑧ Microscope power supply
- ⑨ Motorised parts power supply
- ⑩ Serial cable
- ⑪ PS/2 mouse

7.3 Assembling the microscope

1. Put the microscope stand on a solid table. Insert M-1030 attachment (if provided) above the stand, using the 2mm Allen wrench to tighten the screw. (Fig.1)

2. Connect the cable of the ALC (Automatic Light Control) system to the socket placed on the right side of the frame. (Fig.2)


3. Insert the optical head above the attachment, using the 2mm Allen wrench to tighten the screw. (Fig.3)


4. Insert eyepieces into the empty eyepiece sleeves. (Fig. 4)

5. Insert the condenser under the stage: position until it is well inserted into its holder (under the condenser there is a pin that must fully enter the guide of the holder). (Fig. 5)
6. Lock the condenser fixing knob ①.

7. Mount the stage: lower the support using the coarse focus knob, then place the stage and firmly tighten the lock screw ②. (Fig. 6)

8. Screw the objectives into the revolving nose-piece, in order of magnification. (Fig. 7)

9. Insert the power supply jack in the socket placed in the rear side of the frame. (Fig. 8)

7.4 Only for motorised version

10. Assemble the stage in the same way as the manual version. Check the perfect alignment of the rear part of the stage with the rear arm of the frame. An imperfect alignment could lead to an incorrect functioning of the system. (Fig. 9)


11. Connect the cable ① from the stage to the frame and tight the locking screws of the connectors ②. (Fig. 10)

12. Connect the provided cables: ③ 12V power supply for the motorised parts; ④ 6V microscope power supply; ⑤ serial cable; ⑥ PS/2 mouse. Connect power cables as the last step. (Fig. 11)


8. Summary of brightfield observation procedures

9. Use of the microscope

9.1 General switch on


To activate the transmitted light illuminator, insert the plug of the external power supply into the mains socket and turn the main switch ①, located on the left side of the stand, to the position “1”. (Fig. 12)

9.2 Control keyboard

B-1000 illumination can be managed through the keyboard placed on the left of the stand. (Fig. 13)

- **ON-OFF** (②): press this key (after switching the main switch on “1”) to turn ON or OFF the microscope LED.
- **BOOST** (③): press this button in order to increase the brightness (useful for high-magnification objectives or very opaque specimens).
⚠ Do not enable boost mode while observing with low magnification objectives (4x, 10x) with fully open diaphragm: the high brightness may hurt user’s eyes.
- **AUTO OFF** (④): if you want the illuminator to switch off automatically, press this button until 15, 30 or 60 minutes delay is set. After this period of time, the light will turn off. You have to press the ON-OFF button to turn it on again.

9.3 Brightness adjustment

Use the brightness adjustment dial ⑤ on the left side of the microscope to increase or decrease the light intensity on the specimen. (Fig. 14)


9.4 Adjust the observation head

Loosen the lock-screw ①, turn the observation head to a comfortable position for observation, and then lock the lock-screw. (Fig. 15)

9.5 Adjust the interpupillary distance

Hold the right and left parts of the observation head using both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. (Fig. 16)

9.6 Diopter adjustment

Adjust the fine focusing knob to get the image sharp and clear while observing with your right eye, then turn the left diopter ring ② to obtain a sharp and clear image also with the left eye. (Fig. 17)
The highpoint eyepieces allow the user to wear glasses

- NOTE: For the optimal parfocality of the image, it's suggested to wear your glasses during the normal use of the microscope.

9.7 Use of eyeshields

- **Use with eyeglasses**

Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses. (Fig. 18)

- **Use without eyeglasses**

Raise eyeshields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation. (Fig. 19)

9.8 Focus tension adjustment

Turn the tension adjustment knob to get a suitable tension for the focus system. (Fig. 20)

- **NOTE:** if the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the friction.

9.9 Focus stop lever

Using the focusing knobs focus the specimen with a 4x or 10x objective. Then rotate the focus stop ① in order to block the height of the stage. This simplifies the next focusing operations. (Fig. 21)

The focus stop lever is also useful to avoid accidental contacts between objective and specimen.

9.10 Stage

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverslide 0,17mm. (Fig. 22)
It is possible to place two slides side by side on the stage.

1. Open the spring arm of the slide holder ① and place from the front the slide on the stage.
2. Gently release the spring arm of the slide holder.
 - **A sudden release of the spring arm could cause the falling of the slide.**

9.11 Condenser centering


1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Insert the front lens of the swing-out condenser ①. (Fig. 23)
3. Rotate the field diaphragm ring ② in the direction showed by the arrow, to fully close the diaphragm.
4. Rotate the condenser height adjustment knob ③ to focus the edges of the diaphragm.
5. Rotate the two centering screws ④ to bring the bright spot in the center of the field of view.
6. Gradually open the diaphragm. The condenser is centered when the diaphragm image is symmetrical to the field of view.
7. In normal use, open the diaphragm until it circumscribes the field of view.

Effect of field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image.

Set the diaphragm according to the objective in use until it circumscribes the field of view, in order to eliminate unnecessary light to eyepieces. (Fig. 24)


Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image
- With low contrast specimens set the numerical aperture value ⑤ (printed on the condenser ring) to about 70%-80% of the objective's N.A. (Fig. 25). If necessary, remove one eyepiece and, looking into empty sleeve, adjust the condenser's ring in order to obtain an image like the one in fig. 26.

Example: with objective PLAN 40x / 0.65 set the scale to $0.65 \times 0.8 = 0,52$

9.12 Use of oil immersion objective

1. Focus the specimen with a low power objective.
2. Lower the stage (remembering to lock the coarse upper limit knob).
3. Put a drop of oil (provided) on the area of the specimen to be observed. (Fig. 27)
 - **Make sure that there are no oil bubbles. Air bubbles in the oil damage the image quality.**
 - To check for bubbles: remove an eyepiece, fully open the aperture diaphragm and observe the objective exit pupil. (The pupil must be round and bright).
 - To remove the bubbles, gently move the nose-piece to the right and left to move the immersion objective a few times and allow the air bubbles to move.
4. Insert immersion objective.
5. Return the stage to the upper focusing point and obtain an optimal focus using the fine focus knob.
6. After use, gently remove the oil with a soft paper towel or a lightly moistened optic paper with a mixture of ethyl ether (70%) and absolute alcohol (30%).
 - **The immersion oil, if not immediately cleaned, could crystallize creating a glass-like layer. In this situation the observation of the specimen would be difficult if not impossible due to the presence of an additional thickness on the objective.**

9.13 Use of ALC system

1. Adjust the desired brightness through the eyepieces using the light intensity dial (chapter 9.3).
2. Press the MEM key ① to store this setting (Fig. 28). The light on the microscope will turn off for some seconds, then will turn on again.
 - **The settings could not be working when the light intensity is too low or too high. This is not a defect.**
3. The LED ② of the AUTO key ③ will light up to show that the system is now active.
4. Now the system will automatically adapt the brightness to the eyepieces when an objective is changed, when the aperture diaphragm is used or when another specimen is placed on the stage.
5. Pressing the AUTO key, the ALC system will be disabled, but keeping in memory the setting of the illumination previously done.
6. Pressing the AUTO key again the setting is recalled.
 - **When ALC system is active the light intensity dial is not active.**
 - To store a new setting, repeat step 1 and 2. This will overwrite the old setting and will store a new one.


Fig. 28

Only for motorised version

9.14 Nosepiece rotation

1. To change magnification it is possible operate on the nosepiece movement buttons located on the right side of the frame. (Fig. 29). Orange button ④ rotates the nosepiece clockwise, while the blue button ⑤ rotates the nosepiece counter-clockwise.
2. As an alternative it is possible operate on the right and left mouse buttons.

Fig. 29

9.15 Focusing

Focus motor is activated through the mouse wheel. Front or rear rotation of the mouse wheel raises or lowers the stage. (Fig. 30)

Fig. 30

9.16 Stage

1. Stage is moved through the mouse movement. A mouse movement to the front or to the back ① causes a stage movement of the stage along the Y axis, while a left or right movement ② causes a stage movement of the stage along the X axis. (Fig. 31)
2. It is always possible, however, operate on the translation knobs of the stage for a manual movement.


Fig. 31

10. Use of universal condenser for brightfield / darkfield / phase contrast

Universal condenser provided with B-1000 allows observation in brightfield, darkfield and phase contrast.

Observation mode	Condenser turret position
Brightfield	BF (Fig. 32)
Darkfield	DF (Fig. 33)
Phase contrast 10x	10/20 (Fig. 34)
Phase contrast 20x	10/20 (Fig. 34)
Phase contrast 40x	40 (Fig. 35)
Phase contrast 100x	100 (Fig. 36)

10.1 Brightfield observation (BF)





1. Rotate the condenser turret to insert the "BF" position.
2. Now repeat the steps described in the procedure "*Summary of brightfield observation procedures*" at page 10.

10.2 Darkfield observation (DF)

1. Rotate the condenser turret to insert the "DF" position.
 2. Open the aperture diaphragm.
 3. Place a specimen on the stage and focus.
 4. Observing into eyepieces raise or lower the condenser until a homogeneous illumination of the specimen can be achieved, thus obtaining a proper darkfield effect.
- **Darkfield requires a huge amount of light. Switching from darkfield to brightfield, one could be dazzled. Do not keep your eyes on the eyepieces when moving the condenser turret from DF to BF.**
 - **"Dry" darkfield observation, that is, without the use of oil, is only possible with objectives with N.A. lower than 0.7.**
 - **Observing in darkfield, it may be necessary to raise the condenser from the normal position to obtain a more homogeneous illumination. This is not a defect.**

10.3 Phase contrast observation (PH)

1. Center the condenser as already described at page 14.
2. Rotate the condenser turret to insert the "10/20" position
3. Insert 10x objective into the light path.
4. Place a specimen on the stage and focus.
5. Remove one eyepiece and insert the centering telescope. (Fig. 37)
6. Rotate the upper part of the centering telescope until the two phase rings (one dark and one bright) visible in the telescope are in focus. (Fig. 37)
7. Using centering screws on the condenser ① (Fig. 38), center the phase rings to make the bright ring ② be concentric to the dark ring ③. (Fig- 39-40)
8. Insert 20x objective (do not rotate the condenser turret) and check the centering of the two rings. (Fig. 40)
9. Repeat the same operation with other objectives to check the ring centering: 40x objective – turret position "40", 100x objective – turret position "100".
10. At the end remove the centering telescope, reinstall the eyepiece and begin observation.
 - **With 40x and 100x objectives it may be useful to slightly raise the condenser, to obtain a better projection of the phase rings. This is not a defect.**
 - **With the 4X objective, the condenser could have a dark halo at the periphery of the field of view. This is not to be considered a defect.**

10.4 Use of green filter

- Green filter is used to increase the contrast of the image during phase contrast observation.
- Place the filter on the field lens of the microscope and begin the observation. (Fig. 41)
- For brightfield or darkfield observation it advisable to remove the green filter from the light path.

11. DIC observation

The microscope allows the observation in Differential Interferential Contrast (DIC) with two different methods: Koehler DIC and Nomarski DIC.

The Koehler DIC method is the simplest both from the point of view of installation and from the point of view of use, while the Nomarski DIC method provides for a more complex setup.

Both methods work in transmitted light but can be used in combination with fluorescence observation reflected light and therefore the transmitted light only set-up is different from the combined fluorescence observation set-up.



11.1 Koehler DIC transmitted light

The observation in Koehler DIC in transmitted light requires the kit consisting of the following accessories: Polarizer ①, Analyzer for transmitted light ②, Interferential green filter ③ and DIC slider ④. (Fig. 42)

1. Place the polarizer on the lens at the base of the microscope.
2. Remove the dummy slider from the nosepiece and insert the analyzer into the dummy slider, then insert the ⑤ assembly into the slot ⑥. (Fig. 43)
3. Remove the slide from the stage.
4. Rotate the polarizer at the base of the microscope to achieve maximum darkening of the eyepieces.

5. Once the maximum darkening is found, remove the slider from the nosepiece, remove the analyzer from the dummy slider and insert it into the DIC prism. Now insert the DIC slider ⑦ into the slot ⑥. (Fig. 44)
6. Close the condenser aperture diaphragm a little.

7. Put a specimen on the stage and focus.
8. Begin the observation by rotating the DIC ⑧ slider knob to obtain a three-dimensional sample effect. (Fig. 44)
 - For a better effect on the image you can use the green filter IF550 which must be placed on top of the polarizer.


1.2 Nomarski DIC transmitted light


The observation in Nomarski DIC in transmitted light requires the kit consisting of the following accessories: Universal condenser ① (containing the dedicated DIC prisms according to the objectives in use), Analyzer for transmitted light ②, DIC slider ③. (Fig. 46)

Universal Condenser Controls


- Ⓐ Optical inserts markers
- Ⓑ Aperture diaphragm scale
- Ⓒ Aperture diaphragm lever
- Ⓓ Top lens
- Ⓔ Top lens lever
- Ⓕ Optical inserts centering screws


- Ⓖ Polarizer rotation fixing screw
- Ⓗ Polarizer rotation knob
- Ⓘ Polarizer in/out knob
- ⓵ Polarizer slider locking screw
- Ⓚ Polarizer
- Ⓛ Indicator markers

1. Using the knob ①, insert the polarizer ② embedded in the condenser and loosen the polarizer rotation fixing screw ③. (Fig. 47)

2. Remove the dummy slider from the nosepiece and insert the analyzer into the dummy slider, then insert the ④ assembly into the slot ⑤. (Fig. 48)

3. Remove the slide from the stage.
4. Turn the polarizer knob (H) under the condenser to achieve maximum darkening of the eyepieces, and then tighten the polarizer locking screw (G). (Fig. 49)

Fig. 49

5. Once the maximum darkening is found, remove the slider from the nosepiece, remove the analyzer from the dummy slider and insert it into the DIC prism. Now insert the DIC slider (6) into the slot (5). (Fig. 50)


Fig. 50

6. Rotate the condenser turret (7) to insert the DIC prism matching the objective in use. (Fig. 51)
- **The condenser is supplied with magnetic markers. Each marker is specific to the type of insert mounted in the condenser (DIC, PH, DF, etc.).**

Fig. 51

7. Put a specimen on the stage and focus.
8. Begin the observation by turning the knob on the DIC slider (8) to obtain a three-dimensional effect of the sample. (Fig. 52)

Fig. 52

12. Microphotography

12.1 Installing the “C” mount adapter

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 45)
2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 46)

12.2 Use of reflex camera

1. Insert the Reflex adapter ② into the relay tube to the microscope ①.
2. Screw the “T2” ring ③ (not provided) to the reflex adapter.
3. Connect the Reflex camera ④ to the “T2” just installed. (Fig. 47)
 - “T2” ring is not provided along with the microscope, but is commercially available.
 - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
 - To calculate the magnification of the camera: objective magnification * camera magnification * lens magnification.
 - **If using an SLR camera, mirror movement may cause the camera to vibrate. We suggest lifting the mirror, using long exposure times and a remote cord.**

13. Maintenance

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

14. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark.	Power supply is unplugged.	Connect
	Brightness is too low	Set brightness to a proper level
The edge of the field of view is vignetted or the brightness is asymmetric.	The incident illuminator is not correctly oriented.	Change the angle of the incident illuminator
	The turret of the phase contrast condenser is in a wrong position	Move the turret to a click stop
Dirt or dust is visible in the field of view.	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture diaphragm is stopped down too far	Open aperture diaphragm
	The condenser is not well centered or it is in a wrong height	Set the condenser according to Koehler settings
Visibility is poor. <ul style="list-style-type: none"> • Image is not good. • Contrast is poor. • Details are indistinct. • Image glares 	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture diaphragm is too closed or too open	Adjust aperture diaphragm
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly
	For transmitted light observation the coverglass thickness must not exceed 0.17 mm	Use a coverglass with thickness 0,17 mm
	For phase contrast observation a brightfield objective is used instead of a phase contrast objective	Use a phase contrast objective
	Phase rings of objective and condenser are not well centered	Operate on centering screws to obtain a proper centering
	Objective in use is not compatible with condenser phase ring	Use a compatible objective
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position
One side of the image is out of focus.	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The specimen is out of place (tilted)	Place the specimen flat on the stage.
	The optical performance of the sample cover glass is poor	Use a cover glass of better quality
II. Mechanical Section:		
The coarse focus knob is hard to turn.	The tension adjustment collar is too tight	Loosen the tension adjustment collar
The focus is unstable.	The tension adjustment collar is too loose	Tighten the tension adjustment collar
III. Electric section		
The LED doesn't turn on.	No power supply	Check the power cord connection
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord

IV. Viewing tube assembly		
The field of view of the two eyes is different	The interpupillar distance is not correct	Adjust the interpupillar distance
	The dioptic correction is not right	Adjust the dioptic correction
	The viewing technique is not correct, and the operator is straining the eyesight	When look into the objective, do not stare at the specimen but look at the whole field of view. Periodically, move the eyes away to look at a distant object, then back into the objective
V. Microphotography and video		
The image is unfocused	Incorrect focussing	Adjusting the focus system as in the present manual
The edge of the image is unfocussed	To some degree, it is inherent to the nature of achromatic objectives	The problem can be minimized by a correct setting of the aperture diaphragm
Bright patches appear on the image	Stray light is entering the microscope through the eyepieces and through the camera viewfinder	Cover the eyepieces and the viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.