B-1000 Series

INSTRUCTION MANUAL

Model
B-1000FL-LED

Ver. 2.1 2020

Table of Contents

1. Warning	3
2. Symbols and conventions	3
3. Safety Information	3
4. Intended use	3
5. Instrument description	4
6. Unpacking	6
7. Assembling	6
7.1 Assembling the microscope	7
7.1.1 Manual version	7
7.1.2 Motorised version	9
8. Brightfield (transmitted light) observation procedures	10
9. Use of the microscope in Brightfield (transmitted light)	11
9.1 General switch on	11
9.2 Control keyboard	11
9.3 Brightness adjustment	11
9.4 Adjust the observation head	12
9.5 Adjust the interpupillary distance	12
9.6 Diopter adjustment	12
9.7 Use of eyeshields	12
9.8 Light path selection	13
9.9 Coarse focus tension adjustment	13
9.10 Focus stop lever	14
9.11 Stage	14
9.12 Centering the condenser	15
9.13 Effect of field diaphragm	15
9.14 Aperture diaphragm	15
9.15 Use of oil immersion objective	16
9.16 Only for motorised version	17
9.16.1 Nosepiece rotation	17
9.16.2 Focusing	17
9.16.3 Stage	17
10. Fluorescence (reflected light) observation procedures	18
11. Use of the microscope in Fluorescence (reflected light)	19
11.1 LED switch on	19
11.2 Use of fluorescence	19
11.3 Use of light excluding plate	20
12. Microphotography	21
12.1 Use of "C"mount camera	21
12.2 Use of Reflex camera	21
13. Maintenance	22
14. Troubleshooting	23
Equipment disposal	25

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

2. Symbols and conventions

The following chart is an illustrated glossary of the symbols used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed carefully.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

3. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

4. Intended use


Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.


IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

5. Instrument description

Opposite side

6. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

7. Assembling

Once opened the box, the microscope parts are the following:

- | | |
|-------------------------------|---|
| ① Frame | ⑧ Light excluding plate |
| ② Objectives | ⑨ Dust cover |
| ③ Stage | ⑩ Allen wrench |
| ④ Condenser | ⑪ Immersion oil (if 100x obj. is included in the configuration) |
| ⑤ Observation head | ⑫ Power supply (2 units) |
| ⑥ Eyepieces | |
| ⑦ LED fluorescence attachment | |

7.1 Assembling the microscope

7.1.1 Manual version

1. Put the microscope stand on a solid table. Insert fluorescence attachment above the stand, using the Allen wrench to tighten the screw. (Fig. 1)


2. Insert the optical head above the attachment, using the Allen wrench to tighten the screw. (Fig. 2)


3. Insert eyepieces into the empty eyepiece sleeves. (Fig. 3)

4. Insert the condenser under the stage: position until it is well inserted into its holder (under the condenser there is a pin that must fully enter the holder guide). (Fig. 4)
5. Lock the condenser fixing knob ①.

6. Mount the stage: lower the support using the coarse focus knob, then place the stage and tighten the lock screw ①. (Fig. 5)

7. Screw each objective into the thread of the nose-piece, clockwise with increasing magnification. (Fig. 6)


8. Insert the power supply jack in the socket placed in the rear side of the of the microscope: one for transmitted light and one for fluorescence. (Fig. 7-8)

7.1.2 Motorised version

1. Assemble the stage in the same way as the manual version. Check the perfect alignment of the rear part of the stage with the rear arm of the frame. An imperfect alignment could lead to an incorrect functioning of the system. (Fig. 9)


2. Connect the cable ① from the stage to the frame and tight the locking screws of the connectors ②. (Fig. 10)

3. Connect the provided cables: ③ 12V power supply for the motorised parts; ④ 6V microscope power supply; ⑤ serial cable; ⑥ PS/2 mouse. (Fig. 11)
- **Connect power cables as the last step.**


8. Brightfield (transmitted light) observation procedures

9. Use of the microscope in Brightfield (transmitted light)

9.1 General switch on


1. To activate the transmitted light illuminator, put the main switch ①, located on the left side of the stand, to the position "1". (Fig. 12)

9.2 Control keyboard

B-1000 illumination can be managed through the keyboard placed on the left of the stand. (Fig. 13)

- **ON-OFF** (②): press this key (after switching the main switch on "1") to turn ON or OFF the microscope LED.
 - **BOOST** (③): press this button in order to increase the brightness (useful for high-magnification objectives or very opaque specimens).
- ⚠ Do not enable boost mode while observing with low magnification objectives (4x, 10x) with fully open diaphragm: the high brightness may hurt user's eyes.**
- **AUTO OFF** (④): if you want the illuminator to switch off automatically, press this button until 15, 30 or 60 minutes delay is set. After this period of time, the light will turn off. You have to press the ON-OFF button to turn it on again.

9.3 Brightness adjustment

1. Use the brightness adjustment dial ⑤ on the left side of the microscope to increase or decrease the light intensity on the specimen. (Fig. 14)

9.4 Adjust the observation head


1. Loosen the locking screw ①, turn the observation head to a comfortable position for observation, and then lock the locking screw again. (Fig. 15)

9.5 Adjust the interpupillary distance

1. Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.
- **The graduation on the interpupillary distance indicator ②, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator’s eyes. (Fig. 16)**

The range of the interpupillary distance is 48-75 mm.

9.6 Diopter adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
 2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptic adjustment ring ③ to compensate. (Fig. 17)
- **The adjustment range is ± 5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator’s dioptic correction.**

9.7 Use of eyeshields

• Use with eyeglasses

1. Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses. (Fig. 18)

- **Use without eyeglasses**

1. Raise eyeshields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation. (Fig. 19)

Fig. 19

9.8 Light path selection

- The observation head is equipped with an optical path selector that allows the light to be splitted to the eyepieces and to the photo / TV port.
1. Move the selector ① to one of the three possible positions to split the light. (Fig. 20)

POSITION	LIGHT
IN	100% EYEPIECES
MIDDLE	50% EYEPIECES / 50% TV
OUT	100% TV

Fig. 20

9.9 Coarse focus tension adjustment

The tension of the coarse focusing knob is factory preset.

1. To modify the tension according to personal's needs, rotate the ring ②. (Fig. 21)
- Clockwise rotation increases the tension.
 - If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.



Fig. 21

9.10 Focus stop lever

The upper limit knob has two functions: prevent the contact between slide and objective and acts as “focus memory”.

1. After focussing the specimen, pull the lever ① toward the front of the microscope and lock it. (Fig. 22).
- In this way the focus upper limit is set.
2. Now one can lower the stage with coarse focus knob, replace the specimen and raise again the stage up to the upper limit: specimen will be in approximate focus and will need a fine adjustment to get the proper focus..
- **Fine focus movement is not affected by the coarse focus lock.**
 - **To unlock, move the lever in the opposite direction to the one used for the locking.**
- **Two stoppers ② are inserted on the frame. DO NOT REMOVE THE TWO STOPPERS.**

9.11 Stage

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverslide 0,17 mm. (Fig. 23)
It is possible to place two slides side by side on the stage.

1. Open the spring arm of the slide holder ① and place from the front the slide on the stage.
 2. Gently release the spring arm of the slide holder.
- **A sudden release of the spring arm could cause the falling of the slide.**

9.12 Centering the condenser


1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Insert the front lens of the swing-out condenser ①. (Fig. 24)
3. Rotate the field diaphragm ring ② clockwise, to fully close the diaphragm.
4. Rotate the condenser height adjustment knob ③ to focus the edges of the diaphragm.
5. Rotate the two centering screws ④ to bring the bright spot in the centre of the field of view.
6. Gradually open the diaphragm. The condenser is centred when the diaphragm image is symmetrical to the field of view.
7. In normal use, open the diaphragm until it circumscribes the field of view.

9.13 Effect of field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image.



Set the diaphragm according to the objective in use until it circumscribes the field of view, in order to eliminate unnecessary light to eyepieces. (Fig. 25)

9.14 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image
- With low contrast specimens set the numerical aperture value ⑤ (printed on the condenser ring) to about 70%-80% of the objective's N.A. (Fig. 26). If necessary, remove one eyepiece and, looking into empty eyepiece sleeve, adjust the condenser's ring in order to obtain an image like the one in Fig. 27.

Example: with objective PLAN 40x / 0.65 set the scale to $0.65 \times 0.8 = 0,52$

9.15 Use of oil immersion objective

1. Focus the specimen with a low power objective.
2. Lower the stage (remember to lock the coarse upper limit knob).
3. Put a drop of oil (provided) on the area of the specimen to be observed. (Fig. 28)
 - **Make sure that there are no oil bubbles. Air bubbles in the oil damage the image quality.**
 - To check for bubbles: remove an eyepiece, fully open the aperture diaphragm and observe the objective exit pupil. (The pupil must be round and bright).
 - To remove the bubbles, gently move the nose-piece to the right and to the left to move the immersion objective a few times and allow the air bubbles to move away.
4. Insert immersion objective in the light path.
5. Return the stage to the upper focusing point and obtain an optimal focus using the fine focus knob.
6. After use, gently remove the oil with a soft paper towel or a lightly moistened optic paper with a mixture of ethyl ether (70%) and absolute alcohol (30%).
 - **Immersion oil, if not immediately cleaned, could crystallize creating a glass-like layer. In this situation the observation of the specimen would be difficult if not impossible due to the presence of an additional thickness on the objective.**

Fig. 28

9.16 Only for motorised version

9.16.1 Nosepiece rotation

1. To change magnification it is possible operate on the nosepiece movement buttons located on the right side of the frame. (Fig. 29). Orange button ① rotates the nosepiece clockwise, while the blue button ② rotates the nosepiece counter-clockwise.
2. As an alternative it is possible operate on the right and left mouse buttons.


9.16.2 Focusing

- Focus motor is activated through the mouse wheel.
1. Front or rear rotation of the mouse wheel raises or lowers the stage. (Fig. 30)


9.16.3 Stage

1. Stage is moved through the mouse movement. A mouse movement to the front or to the back ③ causes a stage movement of the stage along the Y axis, while a left or right movement ④ causes a stage movement of the stage along the X axis. (Fig. 31)
- It is always possible, however, operate on the translation knobs of the stage for a manual movement.

10. Fluorescence (reflected light) observation procedures

11. Use of the microscope in Fluorescence (reflected light)

11.1 LED switch on

1. Turn the main switch ①. (Fig. 32)
2. Adjust to the desired brightness by rotating the dial ①.

11.2 Use of fluorescence

The filter turret is provided with 4 positions.

- Position "1" is free and is used for Brightfield
 - Position "2" contains B (Blue) filter
 - Position "3" contains G (Green) filter
 - Position "4" is free and is used for Brightfield
 - **It is not possible to install additional fluorescence filters.**
1. Move the filter selector ② in the desired position. (Fig. 33)

FILTER NAME	EXCITATION FILTER	DICHROIC MIRROR	BARRIER FILTER	APPLICATIONS
B	450-490 nm	495 nm	520LP nm	<ul style="list-style-type: none">• FITC: fluorescent antibodies• Achridine orange: DNA, RNA• Auramine
G	500-540 nm	560 nm	575LP nm	<ul style="list-style-type: none">• Rhodamine, TRITC: fluorescent antibodies• Propidium iodide: DNA, RNA• RFP

11.3 Use of light excluding plate

- **Microscope is provided with a light excluding plate that can be placed on the stage and prevents flare and reflections coming from the condenser front lens.**

The plate can be used in two different ways.

1. Mode n° 1: place the plate on the stage (under the slide holder) and place the slide directly over the plate. (Fig. 34)
 2. Mode n° 2: lower the condenser and insert the plate between the two layers of the stage. (Fig. 35).
- **In both cases it is possible to move the sample using the stage X-Y translation knobs.**


Fig. 34

Fig. 35

12. Microphotography

12.1 Use of “C”mount camera

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 36)

2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 37)

12.2 Use of Reflex camera

1. Insert the Reflex adapter ① into the relay tube to the microscope ②.
2. Screw the “T2” ring ③ (not provided) to the reflex adapter.
3. Connect the Reflex camera ④ to the “T2” ring just installed (Fig. 38).
4. Mount the other end of the relay tube ② into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 36)
 - “T2” ring is not provided along with the microscope, but is commercially available.
 - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
 - To calculate the magnification of the camera: objective magnification * camera magnification * lens magnification.
 - **If using an SLR camera, mirror movement may cause the camera to vibrate.**
 - **We suggest lifting the mirror, using long exposure times and a remote cord.**

13. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

14. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical System:		
LED operates, but field of view remains dark	Brightness is too low	Set brightness to a proper level
	Aperture and field diaphragms are not opened wide enough	Adjust them to proper sizes
	Condenser is lowered too much.	Adjust the condenser height position
	Fluorescence filter selector is not in a click stop	Move the selector to a click stop
	Fluorescence filter is not suitable for the specimen	Use a suitable filter
	Light path selector knob is set to the camera position	Move the knob to the eye position
Field of view is obscured or not evenly illuminated	Light path selector knob is in an intermediate position	Set the knob according to the observation method
	Revolving nosepiece is not correctly engaged	Make sure that the revolving nosepiece clicks properly into place
	Condenser is not attached properly	Re-attach it
	Revolving nosepiece is not attached properly	Push the side dovetail all the way until it is stopped
	An objective that falls outside of the condenser's illumination range is used	Use a condenser to match the purpose
	Condenser is not properly centered	Center the condenser
	Field iris diaphragm is stopped down too far	Open the field iris diaphragm until it circumscribes the field
	Fluorescence filter selector is not in a click stop	Move the selector to a click stop
Dirt or dust is visible in the field of view	Dirt/dust on the eyepieces	Clean thoroughly
	Dirt/dust on the condenser surface	
	Dirt/dust on the specimen	
Image looks double	Aperture iris diaphragm is stopped down too far	Open aperture iris diaphragm
	The condenser is not well centered or it is in a wrong height	Set the condenser according to Koehler settings
Low image quality <ul style="list-style-type: none"> • Image is not sharp • Contrast is poor • Details are indistinct • Image glares 	Condenser is lowered too much	Adjust the condenser height position
	Aperture diaphragm is stopped down too far	Open aperture diaphragm
	Revolving nosepiece is not mounted properly	Push the slide dovetail all the way until it is stopped
	Front lens of objective is dirty	Clean objective
	Immersion oil is not being used with an oil immersion objective	Use immersion oil
	Immersion oil contains bubbles	Remove the bubbles
	Recommended immersion oil is not used	Use the provided immersion oil
	For transmitted light observation, coverglass thickness must not exceed 0.17 mm	Use a coverglass with thickness 0.17 mm

One side of the image is unfocused	Revolving nosepiece is not correctly mounted	Push slide dovetail all the way until it is stopped
	Stage is not correctly mounted	Re-attach it
	Specimen is not correctly mounted on stage	Place specimen correctly on top of stage and secure it with slide holder
Image appears to waver	Revolving nosepiece is not corrected mounted	Push slide dovetail all the way until it is stopped
	Objective is not correctly engaged in light path	Make sure that revolving nosepiece clicks into place correctly
	Condenser is not properly centered	Center the condenser
Field of view becomes only slightly brighter when the voltage is raised.	Condenser is not properly centered	Center the condenser
	Condenser is lowered too far	Adjust the condenser height position
II. Mechanical Section:		
Coarse adjustment knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
	You are trying to raise stage while focus-lock lever is kept locked	Unlock focus-lock lever
Stage drifts down by itself or focus is lost during observation	Tension adjustment ring is too loose	Tighten ring
Coarse adjustment will not go all the way up	Focus-lock lever is locked at a too low height	Unlock focus-lock lever
Coarse adjustment will not go all the way down	Condenser holder is too low	Raise condenser holder
Objective makes contact with specimen before focus is obtained	Specimen is mounted upside down	Mount specimen correctly
III. Electrical Section:		
LED doesn't turn on	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection
IV. Observation tube:		
The field of view of one eye does not match the field of view of the other eye	Interpupillary distance is incorrect	Adjust interpupillary distance
	Incorrect diopter adjustment	Adjust diopter
	Your view is not accustomed to microscope observation	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope
V. Microphotography:		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.