B-1000 Series

INSTRUCTION MANUAL

Model
B-1000FL-LED

Ver. 2.0 2019

Table of Contents

Warning

Symbols and conventions

Safety Information

Intended use

List of accessories and spare parts

Overview

Unpacking

Assembling

Using the microscope

Maintenance

Troubleshooting

Equipment disposal

Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.


IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

List of accessories and spare parts

CAT. NO.	DESCRIPTION
M-005	Micrometric slide, 26x76mm, range 1mm, div. 0,01mm
M-613	Polarizing set (filters only).
M-615	Lambda filter for polarizing set.
M-617.1N	Phase contrast set with IOS PLAN objective 40x.
M-690	Eyecup (pair).
M-619	Photo tube adapter for full frame SLR camera.
M-173	Photo adapter for APS-C and Full Frame Reflex cameras.
M-620	C-Mount adapter for 1/3" sensor.
M-620.1	C-Mount adapter for 1/2" sensor.
M-699	Universal adapter for OPTIKAM & DIGI cameras (only models with eyepiece adapter, 23mm).
M-114	C-Mount adapter for 1/2" sensor.
M-116	C-Mount adapter for 2/3" sensor.
M-113.1	Ring adapter, 30mm (for monocular and binocular microscopes).
15008	OPTIKA immersion oil, 10ml.

Overview


Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Assembling

Once you open the box, these are the microscope's components:

Put the microscope stand on a solid table.
First insert the LED attachment, use the 2mm Allen wrench to tighten the screw. All screws are already inserted into each threaded hole. (Fig.1)


Insert the optical head above the attachment, using the **2mm Allen wrench** to tighten the screw. (Fig.2)


Insert both eyepieces into the tubes of the optical head. (Fig.3)

Insert the condenser under the stage: position until it is well inserted into its holder (under the condenser there is a pin that must fully enter the guide of the holder). (Fig. 4)

Mount the stage: lower the support using the coarse focus knob, then place the stage and firmly tighten the lock screw. (Fig.5)

Screw each objective into the thread of the nosepiece, in order of magnification. (Fig.6)

Insert the power supply jack on the connector at the rear: one for brightfield illumination, one for illumination. (Fig.7)

Using the microscope

Control keyboard

B-1000 illumination can be managed through the keyboard placed on the left of the stand:

1) ON-OFF: press this button to turn on or off the LED illuminator.

2) BOOST: press this button in order to increase the brightness (useful for high-magnification objectives or very opaque specimens).

3) AUTO-OFF: if you want the illuminator to switch off automatically, press this button until 15, 30 or 60 minutes delay is set. After this period of time, the light will turn off. You have to press the ON-OFF button to turn it on again. (Fig.8)


Fig.8

Do not enable boost mode while observing with low magnification objectives (4x, 10x) with fully open diaphragm: the high brightness may hurt user's eyes.

LED settings

In order to activate the transmitted light illuminator, insert the plug of the external power supply into the socket and turn on the main switch on the side of the main body. Press the ON-OFF button on the control keyboard and turn the brightness adjustment knob to a brightness suitable for observation. (Fig.9)

Fig.9

Adjust the observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw. (Fig.10)

Fig.10

Adjust interpupillary distance

Hold the right and left parts of the observation head using both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. (Fig.11)

Fig.11

Place the specimen on the stage

Fix the specimen slide on the mechanical stage using the slide-clamp. Ensure that the specimen is centred over the stage opening. (Fig.12)

Fig.12

Focus tension adjustment

Turn the tension-adjust knob to get a suitable tension for the focus system. (Fig.13)
NOTE: if the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the friction.


Fig.13

Focus-stop knob

Loosen the focus-stop knob, rotate the coarse focusing knob to bring the slide into focus with a 4X or 10X objective. Then rotate the focus-stop in order to block the height of the stage. This simplifies the next focusing operations. The focus-stop knob is also useful to avoid accidental contacts between objective and specimen. (Fig.14)

Fig.14

Diopter adjustment

Adjust the fine focusing knob to get the image sharp and clear while observing with your right eye, then turn the left diopter ring to a sharp and clear image also with the other eye. The highpoint eye-pieces allow the user to wear glasses. (Fig.15)

NOTE :For the optimal parafocality of the image, it's suggested to wear your glasses during the normal use of the microscope.

Fig.15

Centering the condenser

Insert the swing-out lens of the condenser ① and fully close the field diaphragm ②. Rotate the condenser adjustment knob ③ until you see a sharp image of the closed field diaphragm (a bright spot of light). Act on the condenser centering screws ④ until you move the bright spot in the center of the image field. Then re-open the field diaphragm ②. (Fig.16)



Fig.16

CENTERING THE CONDENSER

Adjust the aperture of the iris diaphragm under the condenser to set the numerical aperture of the illuminator, thus controlling image contrast and resolution. It is suggested to set the aperture diaphragm between 70% and 80% of objective's N.A. (Fig.17)


Fig.17

LED

Insert the plug of the power supply (on the rear).

Adjust the optical head and focus on your sample using the brightfield illumination, as described above.


Switch off the brightfield transmitted illumination and place the black slide under the stage (Fig.18)

Move the filter selector lever ① to the desired position: B (blue excitation) or G (green excitation) (Fig.19)

Rotate the light intensity knob ① on the right until you reach a sufficient brightness. (Fig.20)

When you move the filter selector lever, the corresponding LED switches on automatically.

The LEDs are factory-aligned and don't need any set-up operation.

Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
1. Optical System		
LED does not light.	Power cord is unplugged.	Plug power cord into the power outlet.
LED operates, but field of view remains dark.	Aperture and field iris diaphragms are not opened wide enough.	Adjust them to proper sizes.
	Condenser is lowered too much.	Adjust the condenser height position.
	Light path selector knob is set to the camera position.	Move the knob to the eye position.
Field of view is obscured or not evenly illuminated	Light path selector knob is in an intermediate position.	Set the knob according to the observation method.
	Revolving nosepiece is not correctly engaged.	Make sure that the revolving nosepiece clicks properly into place.
	Condenser is not attached properly.	Re-attach it.
	Revolving nosepiece is not attached properly.	Push the side dovetail all the way until it is stopped.
	An objective that falls outside of the condenser's illumination range is used.	Use a condenser to match the purpose.
	Condenser is not properly centered.	Center the condenser.
	Field iris diaphragm is stopped down too far.	Open the field iris diaphragm until it circumscribes the field.
Dirt or dust is visible in the field of view.	Dirt/dust on the eyepieces	Clean thoroughly.
	Dirt on the top surface of the condenser	
	Dirt/dust on the specimen	
Visibility is poor. · Image is not poor. · Contrast is poor. · Details are indistinct. · Image glares.	Condenser is lowered too far.	Adjust the condenser height position.
	Aperture iris diaphragm is stopped down too far.	Open aperture iris diaphragm.
	Revolving nosepiece is not mounted properly.	Push the slide dovetail all the way until it is stopped.
	Front lens of objective is dirty.	Clean objective.
	Immersion oil is not being used with an oil immersion objective.	Use immersion oil.
	Immersion oil contains bubbles.	Remove the bubbles.
	Recommended immersion oil is not used.	Use the provided immersion oil.
	Dirt/dust on specimen.	Clean it.
	Dirt/dust on condenser	
	Inappropriate object slide or cover glass thickness.	Replace with glass of recommended thickness.

PROBLEM	CAUSE	SOLUTION
One side of image is blurred.	Objective is not correctly engaged in light path.	Make sure that revolving nosepiece clicks into place correctly.
	Revolving nosepiece is not correctly mounted.	Push slide dovetail all the way until it is stopped.
	Stage is not correctly mounted.	Re-attach it.
	Specimen is not correctly mounted on stage.	Place specimen correctly on top of stage and secure it with slide holder.
Image appears to waver.	Revolving nosepiece is not correctly mounted.	Push slide dovetail all the way until it is stopped.
	Objective is not correctly engaged in light path.	Make sure that revolving nosepiece clicks into place correctly.
	Condenser is not properly centered.	Center the condenser.
Field of view becomes only slightly brighter when the voltage is raised.	Condenser is not properly centered.	Center the condenser.
	Condenser is lowered too far.	Adjust the condenser height position.
2. Coarse/Fine Adjustment		
Coarse adjustment knob is hard to turn.	Tension adjustment ring is tightened excessively.	Loose ring.
	You are trying to raise stage while focus-lock lever is kept locked.	Unlock focus-lock lever
Stage drifts down by itself or focus is lost during observation.	Tension adjustment ring is too loose.	Tighten ring.
Coarse adjustment will not go all the way up.	Focus-lock lever is locked at a too low height.	Unlock focus-lock lever.
Coarse adjustment will not go all the way down.	Condenser holder is too low.	Raise condenser holder.
Objective makes contact with specimen before focus is obtained.	Specimen is mounted upside down.	Mount specimen correctly.
3. Observation Tube		
Field of view of one eye does not match that of the other.	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Different eyepieces are used on left and right.	Change on eyepiece to match the other so that both sides are the same type.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope.
4. Stage		
Image shifts when you touch stage.	Stage is not properly mounted.	Clamp stage.
Specimen stops midway on the X-axis traverse.	Specimen is not correctly positioned.	Place specimen correctly.
X- and Y-axis knobs are too tight or too loose.	Tension of X- and Y-axis knobs is too high or too low.	Adjust tension.

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.