**B-1000 Series**

# INSTRUCTION MANUAL

Model
B-1000POL
B-1000POL-I

Ver. 2.4    2023

---

## Table of Contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
3.1	B-1000POL	4
3.2	B-1000POL-I	5
4.	Unpacking	6
5.	Intended use	6
6.	Symbols and conventions	6
7.	Instrument description	7
7.1	B-1000POL	7
7.2	B-1000POL-I	9
8.	Assembling	11
8.1	Microscope assembling	11
8.1.1	B-1000POL	11
8.1.2	B-1000POL-I	13
9.	Transmitted light brightfield observation procedures	17
9.1	B-1000POL	17
9.2	B-1000POL-I	18
10.	Use of the microscope (transmitted light brightfield)	19
10.1	General switch on	19
10.2	Control keyboard	19
10.3	Brightness adjustment	19
10.4	Adjust the observation head	19
10.5	Adjust the interpupillary distance	20
10.6	Diopter adjustment	20
10.7	Use of eye shields	20
10.8	Light path selection	21
10.9	Coarse focus tension adjustment	21
10.10	Focus stop lever	21
10.11	Centering the condenser	22
10.12	Effects of the field diaphragm	22
10.13	Aperture diaphragm	22
10.14	Stage	23
10.14.1	45° click stop lever	23
11.	Polarized transmitted light observation procedures	24
11.1	B-1000POL	24
11.2	B-1000POL-I	25
12.	Polarized reflected light observation procedures	26
12.1	B-1000POL-I	26
13.	Use of the microscope in polarized light	27
13.1	Centering of rotatable stage	27
13.2	Nosepiece centering	28
13.3	Checking the extinction of the light	28
13.3.1	B-1000POL	28
13.3.2	B-1000POL-I	29
13.4	Centering of reflected light diaphragms	30
13.4.1	Field diaphragm (FS)	30
13.4.2	Aperture diaphragm (AS)	30
13.5	Use of tint plates	31
13.6	Use of the Bertrand lens	31
13.7	Use of diffusing filter (B-1000POL-I)	31
14.	Microphotography	32
14.1	Use of C-mount cameras	32
14.2	Use of Reflex cameras	32
15.	Maintenance	33
16.	Troubleshooting	34
	Equipment disposal	36

---

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

## 2. Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

### 3. Package content

#### 3.1 B-1000POL

① Microscope frame

② Observation head

③ Condenser

④ Eyepieces

⑤ Objectives

⑥ Rotatable stage + mechanical stage

⑦ Bertrand lens

⑧ Analyzer

⑨ Tint plates

⑩ Dummy slider

⑪ Nosepiece centering screws

⑫ Allen wrench

⑬ Power supply

⑭ Dust cover

### 3.2 B-1000POL-I

- ① Microscope frame
- ② Eyepieces
- ③ Objectives
- ④ Observation head
- ⑤ Bertrand lens
- ⑥ Analyzer
- ⑦ Condenser
- ⑧ Rotatable stage + mechanical stage

- ⑨ Tint plates
- ⑩ Diffusing filter
- ⑪ Reflected light illuminator
- ⑫ Reflected light polarizer
- ⑬ Nosepiece centering screws
- ⑭ Allen wrench
- ⑮ Power supply
- ⑯ Dust cover

---

## 4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 5. Intended use

### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

## 6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


### ELECTRICAL SHOCK


This symbol indicates a risk of electrical shock.

## 7. Instrument description

### 7.1 B-1000POL


**B-1000POL (opposite side)**

7.2 B-1000POL-I

## 8. Assembling

### 8.1 Microscope assembling

#### 8.1.1 B-1000POL

1. Insert the Bertrand lens ① on the frame and lock the locking screw ② with the provided Allen wrench. (Fig. 1)

2. Insert the optical head on the Bertrand lens and lock the locking screw with the provided Allen wrench. (Fig. 2)

- **Hold the head with one hand during the locking in order to avoid that the head falls.**

3. Insert both eyepieces into the tubes of the optical head. (Fig. 3)

- **One of the two eyepieces is equipped with a crosshair for centering the entire optical system. It is advisable to insert the eyepiece with crosshair in the right eyepiece sleeve.**

4. Insert the condenser under the stage: position until it is well inserted into its holder (under the condenser there is a pin that must fully enter the holder guide). (Fig. 4)
5. Lock the condenser fixing knob ①.

6. Mount the rotating stage: at the bottom of the stage there is a spring, push this spring toward the stage support ①, then push the stage downward ②. (Fig. 5)

7. Screw each objective into the thread of the nosepiece, clockwise with increasing magnification. (Fig. 6)

- The nosepiece is centrable and has small holes that are used for centering each position.
- One position has no centering holes. This position is “pre-centered” and serves as a reference for all other positions. It is recommended that the 10X be placed in this position.

8. Remove the dummy slider from Bertrand lens ③ and insert the analyser ④. (Fig. 7 - 8)

9. Insert the power supply jack in the connector placed at the rear side of the microscope. (Fig. 9)

### 8.1.2 B-1000POL-I

1. Insert the reflected light illuminator ① on the frame and lock the locking screw ② with the provided Allen wrench. (Fig. 10)

2. Connect the illuminator cable jack to the socket ③ placed in the back side of the frame. (Fig. 11)

3. Insert the Bertrand lens ④ on the reflected light illuminator and lock the locking screw ⑤ with the provided Allen wrench. (Fig. 12)

4. Insert the optical head on the Bertrand lens and lock the locking screw with the provided Allen wrench. (Fig. 13)

  - **Hold the head with one hand during the locking in order to avoid that the head falls.**

5. Insert both eyepieces into the tubes of the optical head. (Fig. 14)

  - **One of the two eyepieces is equipped with a crosshair for centering the entire optical system. It is advisable to insert the eyepiece with crosshair in the right eyepiece sleeve.**

6. Insert the condenser under the stage: position until it is well inserted into its holder (under the condenser there is a pin that must fully enter the holder guide). (Fig. 15)
7. Lock the condenser fixing knob ①.

8. Mount the rotating stage: at the bottom of the stage there is a spring, push this spring toward the stage support ①, then push the stage downward ②. (Fig. 16)


9. Screw each objective into the thread of the nosepiece, clockwise with increasing magnification. (Fig. 17)

- The nosepiece is centrable and has small holes that are used for centering each position.
- One position has no centering holes. This position is “pre-centered” and serves as a reference for all other positions. It is recommended that the 10X be placed in this position.

10. Insert reflected light polarizer ①. (Fig. 18)

11. Remove the dummy slider from Bertrand lens ② and insert the analyser ③. (Fig. 19 - 20)


- 
12. Insert the power supply jack in the connector placed at the rear side of the microscope. (Fig. 21).


## 9. Transmitted light brightfield observation procedures

### 9.1 B-10000POL

9.2 B-1000POL-I


## 10. Use of the microscope (transmitted light brightfield)

### 10.1 General switch on

To activate the transmitted light illuminator, put the main switch ①, located on the left side of the stand, to the position “1”. (Fig. 22)


- **For B-1000POL-I model only.** There is a three-position switch on the left side of the stand: position “I” turns on the transmitted light, position “II” turns on the reflected light and position “O” turns off the microscope.

### 10.2 Control keyboard

B-1000 illumination can be managed through the keyboard placed on the left of the stand. (Fig. 23)

- **ON-OFF (②):** press this key (after switching the main switch on “1”) to turn ON or OFF the microscope LED.
- **BOOST (③):** press this button in order to increase the brightness (useful for high-magnification objectives or very opaque specimens).  
**Do not enable boost mode while observing with low magnification objectives (4x, 10x) with fully open diaphragm: the high brightness may hurt user’s eyes.**
- **AUTO OFF (④):** if you want the illuminator to switch off automatically, press this button until 15, 30 or 60 minutes delay is set. After this period of time, the light will turn off. You have to press the ON-OFF button to turn it on again.

### 10.3 Brightness adjustment

Use the brightness adjustment dial ⑤ on the left side of the microscope to increase or decrease the light intensity on the specimen. (Fig. 24)

### 10.4 Adjust the observation head

Loosen the locking screw ①, turn the observation head to a comfortable position for observation, and then lock the locking screw again. (Fig. 25)

### 10.5 Adjust the interpupillary distance

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

- **The graduation on the interpupillary distance indicator ②, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator’s eyes. (Fig. 26)**

The range of the interpupillary distance is 48-75 mm.

### 10.6 Diopter adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
  2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptic adjustment ring ③ to compensate. (Fig. 27)
- **The adjustment range is  $\pm 5$  diopter. The number indicated on the adjustment ring graduation should correspond to the operator’s dioptic correction.**

### 10.7 Use of eye shields

- **Use with eyeglasses**

Fold rubber eye shields with both hands. Folded eye shields avoid scratching the lenses of eyeglasses. (Fig. 28)

- **Use without eyeglasses**

Raise eye shields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation. (Fig. 29)

### 10.8 Light path selection

- The observation head is equipped with an optical path selector that allows the light to be distributed to the eyepieces and to the photo / TV port.
- Move the selector ① to one of the three possible positions to distribute the light. (Fig. 30)

POSITION	LIGHT
IN	100% EYEPIECES
MIDDLE	50% EYEPIECES / 50% TV
OUT	100% TV

Fig. 30

### 10.9 Coarse focus tension adjustment

The tension of the coarse focusing knob is factory preset.

- To modify the tension according to personal's needs, rotate the ring ①, placed on the right side of the frame just beside the coarse focus knob. (Fig. 31)
- Clockwise rotation increases the tension.
  - If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.


Fig. 31

### 10.10 Focus stop lever

The upper limit knob has two functions: prevent the contact between slide and objective and acts as "focus memory".

- After focussing the specimen, pull the lever ① toward the front of the microscope and lock it. (Fig. 32).
- In this way the focus upper limit is set.
- Now one can lower the stage with coarse focus knob, replace the specimen and raise again the stage up to the upper limit: specimen will be in approximate focus and will need a fine adjustment to get the proper focus.
- Fine focus movement is not affected by the coarse focus lock.**
  - To unlock, move the lever in the opposite direction to the one used for the locking.**
- Two stoppers ② are inserted on the frame. DO NOT REMOVE THE TWO STOPPERS.**


Fig. 32

### 10.11 Centering the condenser

1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Insert the front lens of the swing-out condenser ①. (Fig. 33)
3. Rotate the field diaphragm ring ② in counterclockwise direction, to fully close the diaphragm.
4. Rotate the condenser height adjustment knob ③ to focus the edges of the diaphragm.
5. Rotate the two centering screws ④ to bring the bright spot in the centre of the field of view.
6. Gradually open the diaphragm. The condenser is centred when the diaphragm image is symmetrical to the field of view. (Fig. 34)
7. In normal use, open the diaphragm until it circumscribes the field of view.

Fig. 33

### 10.12 Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image.

Set the diaphragm according to the objective in use until it circumscribes the field of view, in order to eliminate unnecessary light to eyepieces.


Fig. 34

### 10.13 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
- With low contrast specimens set the numerical aperture value ① (printed on the condenser ring) to about 70%-80% of the objective's N.A. (Fig. 35) If necessary, remove the eyepiece and, looking into the empty sleeve, adjust the condenser's ring in order to obtain an image like the one in fig. 36.

**Example:** with objective PLAN 40x / 0.65 set the scale to 0.65 x 0.8 = 0.52

Fig. 35



Fig. 36


## 10.14 Stage

- The rotatable stage is equipped with an X-Y mechanical stage applied directly on the stage surface.
  - The movement is made by means of the two knobs ① (Fig. 37) that allow to move the sample without affecting the rotation on the 360° of the stage itself.
1. Open the spring arm of the slide holder ② and place from the front the slide on the stage.
- **Gently release the spring arm of the slide holder.**
  - **A sudden release of the spring arm could cause the falling of the slide.**


### 10.14.1 45° click stop lever

- This lever is used to enable or disable the 45° click stop function during stage rotation.
  - When the function is enabled, the operator will hear a “click” every 45° by rotating the stage.
1. Pull the lever ③ towards the front of the microscope to enable the function. (Fig. 38)
  2. Push the lever towards the rear to disable the function.


# 11. Polarized transmitted light observation procedures

## 11.1 B-1000POL


11.2 B-1000POL-I

## 12. Polarized reflected light observation procedures

### 12.1 B-1000POL-I

### 13. Use of the microscope in polarized light

- The system allows observation in Orthoscopy (crossed Nicol) or in Conoscopy (crossed Nicol with the use of the Bertrand lens).
- For optimal performance in polarized light microscopy, accurate optical adjustments are essential before beginning the observation.

#### 13.1 Centering of rotatable stage

1. Loosen the stage rotation locking screw ① and rotate the stage until the graduated scale of the stage ② and the vernier scale are aligned ③ in the "0" position. (Fig. 39)
- This operation serves to ensure a standard reference position for centering the rotating stage.


Fig. 39

2. Put the 10x objective in the optical path.
3. Move a small recognizable detail in the field of view ④ using the X-Y movement knobs on the mechanical stage and place it at the centre of the eyepiece crosshair. (Fig. 40)
4. Focus on this detail.
5. By turning the stage, the detail in focus will draw a circle ⑤. (Fig. 40)
6. Move again the stage on the "0" position and lock the locking screw ①. (Fig. 39)


Fig. 40

7. Using the centering screws of the stage ⑥ move the detail in diametrically opposite direction to the described circle. The movement must be about half the diameter of the circle drawn. (Fig. 41)
8. Using the X-Y movement knobs on the mechanical stage place the detail in the centre of the eyepiece crosshair. Loosen the stage locking screw again and rotate the stage again.
9. If the centering has been carried out correctly, by turning the stage, the image of the focused detail does not move with respect to the centre of the crosshair. If not, repeat the operations described by 1. to 8. until you reached the perfect coincidence of the centre of rotation of the stage with the centre of the crosshair, for which the specimen remains in the centre of the crosshair when turning the stage.


Fig. 41

## 13.2 Nosepiece centering

1. After centering the stage with the 10x objective, return the recognizable detail used for centering in the centre of the crosshair.
2. Rotate the nosepiece by inserting all the other objectives in the optical path and check that the detail is always in the centre of the crosshair.
3. If this is not the case, use the centering screws ① on the nosepiece to ensure that all the objectives are perfectly centred with respect to the optical axis. (Fig. 42)

## 13.3 Checking the extinction of the light

### 13.3.1 B-1000POL

1. Insert the swing-out polarizer ① built in the condenser. (Fig. 43)
2. Remove the specimen from the light path and insert 10x.
3. Loosen the locking screw of the polarizer ② and check that the scale is in the "0" position ③. (Fig. 44)
4. Insert rotatable analyser in the light path, loosen the locking screw ④ and put the vibration scale on 0° ⑤, then lock the locking screw ④. (Fig. 45)
5. Rotate the polarizer scale ③ to obtain total extinction (total dark in the eyepieces). Then tighten the screw ②. (Fig. 44)

- It may happen that the scale of the polarizer is not perfectly aligned on the reference mark but is shifted by one or two notches. This is not a defect but it is due to the mechanical alignment of the polarizers during assembly.



### 13.3.2 B-1000POL-I

#### Extinction in reflected light

1. Move the selector ① on the reflected light illuminator in the fully inserted position, corresponding to letter “R”. (Fig. 46)

2. Insert the reflected light polarizer ②. (Fig. 47)
3. Place on the stage a flat mirror and insert 10x objective.

4. Insert rotatable analyser in the light path, loosen the locking screw ③ and turn the vibration direction scale ④ until the eyepieces are completely dark, then lock the locking screw ③. (Fig. 48)
- It may happen that the scale of the polarizer is not perfectly aligned on the reference mark but is shifted by one or two notches. This is not a defect but it is due to the mechanical alignment of the polarizers during assembly.



#### Extinction in transmitted light

1. Move the selector ① on the reflected light illuminator in the fully removed position, corresponding to letter “T”. (Fig. 46)
2. Repeat procedure described in steps 1. to 5. for the B-1000POL.

## 13.4 Centering of reflected light diaphragms


### 13.4.1 Field diaphragm (FS)

1. Move the selector ① on the reflected light illuminator in the fully inserted position, corresponding to letter "R". (Fig. 46)
2. Place a specimen on the stage, insert 10x objective and focus the specimen.
3. Rotate the knurled ring of the field diaphragm ② in the direction shown by the arrow to fully close the diaphragm. (Fig. 49)
4. Using the provided Allen wrench use the two centering screws ③ to bring the bright spot in the centre of the field of view.
5. Gradually open the diaphragm. The illuminator is centred when the diaphragm image is symmetrical to the field of view. (Fig. 50)
6. In normal use, open the diaphragm until it circumscribes the field of view.

### 13.4.2 Aperture diaphragm (AS)

1. Rotate the knurled ring of the aperture diaphragm ④ in the direction shown by the arrow to fully close the diaphragm.
2. Remove one eyepiece.
3. While observing in the empty eyepiece holder, use the provided Allen wrench and use the centering screws ⑤ to bring the bright spot in the centre of the field of view. (Fig. 51)
4. The illuminator is centred when the diaphragm image is symmetrical to the field of view.
  - The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
  - With low contrast specimens rotate the ring on the illuminator to set the numerical aperture value to about 70%-80% of the objective's N.A. If necessary, remove one eyepiece and, looking into empty sleeve, adjust the ring of the illuminator in order to obtain an image like the one in Fig. 28.

### 13.5 Use of tint plates

Three tint plates are supplied with the microscope:


- Lambda ( $\lambda$ ) plate (1<sup>st</sup> order Red)
  - Lambda/4 ( $\lambda/4$ ) plate
  - “Quartz wedge” (Q) plate
1. Insert in the right slot of the Bertrand lens ① one of the tint plates ②. (Fig. 52)
  2. In polarized light, inserting one of the plates will have chromatic effects on the specimen.
    - Using the  $\lambda$  plate (also called 1<sup>st</sup> order Red) the specimen will take a magenta tinge.
    - Using the  $\lambda/4$  plate the specimen will take a colour tending to pale yellow.
    - Using the Q plate the specimen will show a series of coloured bands that will fade as the plate is inserted.

### 13.6 Use of the Bertrand lens

Bertrand lens allows observation in Orthoscopy and Conoscopy. In the disengaged position (“O”) the lens allows observation in Orthoscopy, while in the engaged position (“B”) it is possible to make observations in Conoscopy.

1. Rotate the upper knurled ring of the Bertrand lens ① to engage the “B” position. (Fig. 53)
  2. Using one objective from 20x to 60x, focus the conoscopic image using the focus ring ②.
  3. If the conoscopic image is not perfectly centred with respect to the optical axis, centre the image using the centering screws.
- By turning the stage you will see black fringes that will appear and disappear depending on the rotation of the stage. These fringes are the crystallization axes of that specific crystal. (Fig. 54)

### 13.7 Use of diffusing filter (B-1000POL-I)

Depending on the type of sample to be observed, it may be useful to remove or insert the diffusing filter ④ that is available in the back of the illuminator.


1. Insert the slide into the illuminator until it reaches the end of its stroke to insert the diffuser filter in the optical path. (Fig. 55)
2. Pull out one click (until the first “click”) the slide to remove the filter from the optical path, but always leaving the slide in place.
3. If you intend to completely remove the slide from the illuminator, remove it completely from its housing.

## 14. Microphotography

### 14.1 Use of C-mount cameras

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 56)

2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 57)

### 14.2 Use of Reflex cameras

1. Insert the Reflex adapter ① into the relay tube to the microscope ②.
  2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
  3. Connect the Reflex camera ④ to the "T2" just installed (Fig. 58).
  4. Mount the other end of the relay tube ② into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 56)
- "T2" ring is not provided with the microscope, but is commercially available.
  - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
  - To calculate the magnification of the camera: objective magnification \* camera magnification \* lens magnification.
  - **If using an SLR camera, mirror movement may cause the camera to vibrate.**
  - **We suggest lifting the mirror, using long exposure times and a remote cord.**

---

## 15. Maintenance

### Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

### To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

### Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

**For the best results, use the OPTIKA cleaning kit (see catalogue).**

If you need to send the microscope to Optika for maintenance, please use the original packaging.

## 16. Troubleshooting

Review the information in the table below to troubleshoot operating problems.


PROBLEM	CAUSE	SOLUTION
<b>I. Optical Section:</b>		
LED operates, but field of view remains dark.	Power supply is unplugged.	Connect
	Brightness is too low	Set brightness to a proper level
	Bertrand lens is in	Remove Bertrand lens from light path
	You are in a position of extinction	Disengage analyzer from the light path
Field of view is obscured or not evenly illuminated	Revolving nosepiece is not correctly engaged.	Make sure that the revolving nosepiece clicks properly into place.
	Tint plate, a filter or Bertrand lens are in an intermediate position.	Move to a click stop
Dirt or dust is visible in the field of view.	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture iris diaphragm is stopped down too far.	Open aperture iris diaphragm.
	The condenser is not well centered or it is in a wrong height	Set the condenser according to Kohler settings.
Visibility is poor. <ul style="list-style-type: none"> <li>• Image is not sharp.</li> <li>• Contrast is poor.</li> <li>• Details are indistinct.</li> </ul>	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture iris diaphragm is too closed or too open.	Adjust aperture iris diaphragm.
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly.
	For transmitted light observation, the coverglass thickness must not exceed 0.17mm	Use a coverglass with thickness 0.17mm
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position.
One side of the image is unfocused	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Slide is mounted not in a flat position (tilted)	Place the specimen in a flat position on the stage
	Poor quality of the glass slide	Use a glass slide with higher quality
Conoscopic image cannot be seen	Swing-out lens of the condenser is not in the light path	Insert the lens in the light path
	Bertrand lens is not in the light path.	Insert the lens in the light path.
Total extinction cannot be obtained	Analyzer is not in the light path	Insert analyzer in the light path.
<b>II. Mechanical Section:</b>		
Coarse focus knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
Focus is unstable	Tension adjustment ring is too loose	Tighten tension adjustment ring
<b>III. Electrical Section:</b>		
LED doesn't turn on.	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection

<b>IV. Observation tube:</b>		
Field of view of one eye does not match that of the other.	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope.
<b>V. Microphotography:</b>		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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**OPTIKA® S.r.l.**