B-1000 Series

INSTRUCTION MANUAL

| Model |
|-------------|
| B-1000Ti-2 |
| B-1000Ti-3 |
| B-1000Ti-5 |
| B-1000Ti-10 |

Version: 2
Issued: 16, 02, 2015

Description

The following assembly procedure can be used for the multidiscussion systems B-1000Ti-2, B-1000Ti-3, B-1000Ti-5 and B-1000Ti-10. The assembly sequence is identical for all systems.

The microscope is a precision scientific instrument designed for long lasting life with a minimal maintenance, being built following the best optical and mechanical standards and developed for a daily use. Optika reminds that the present manual contains important informations for a safe use and a correct maintenance of the instrument. It has to be available for each user.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

Unpacking and assembly

The multidiscussion system is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the upper part of the container. Take care not to let optical items fall out and get damaged. Remove each part from the housing and place them on a sturdy table.

- 2.1 Install the microscope body as indicated in the specific instruction manual.
Do not install the observation head.
- 2.2 Place the splitter attachment of the multidiscussion system and tighten the lock screw on the right side of the microscope (Fig. 1).

POINTER POSITION CONTROLLER (1)

POINTER INTENSITY CONTROLLER (2)

Fig. 1

- 2.3 In order to use the arrow which will be used as a pointer, connect the 5Vdc power supply to the rear Jack of the splitter attachment. (Fig. 2). Now use the pointer position controller (1) or the pointer intensity controller (2) to adjust the position and intensity of the pointer.

Fig. 2

2.4 Begin to place the connection tubes on both sides of the splitter attachment. Firmly screw the black ring nut. (Fig. 3-4).

- **Every connection is identified by a letter printed on both sides of the connection. Make sure to match the letters in order to correctly assemble the microscope.**

Fig. 3

2.5 Repeat the same procedure on the other side of the splitter.


Fig. 4

2.6 Insert the second part of the extension tube (Fig. 6).

- **At the end of the first extension tube there is a lens (Fig. 5). Make sure it is free from dirt, dust or other contaminants before to proceed with the assembling of the second extension tube.**


Fig. 5

2.7 Fully insert the second extension tube in the right position. Using the provided Allen wrench (small one) lock the fixing screws (1) to block the extension tube (Fig. 6).


Fig. 6

- 2.8** Adjust the height of the multi-head holder. Loosen the base fixing knob (1), unscrew the base (2) in order to reach the desired height, then lock the knob. (Fig. 7).
Make sure that each extension tube is perfectly horizontal.

Fig. 7

- 2.9** Insert the optical head (always matching the reference letters) and screw the black fixing ring. (Fig. 8).


Fig. 8

- 2.10** Insert the binocular heads, matching the reference letters (Fig. 9).


Fig. 9

- 2.11** Insert the provided eyepieces (WF10X/20) into binocular heads. (Fig. 10).

Fig. 10

- 2.12 Repeat all these operations for all the other observation points.


Fig. 11

- 2.13 (only B-1000Ti-10) mount the second splitter attachment over the pointer attachment, then fix the extension tube rotating the black ring nut (always following the coupling catters): Fig. 11 and Fig. 12


Fig. 12

- 2.14 (only B-1000Ti-10) repeat for the 5 back heads the same procedure of the front heads, always following the coupling letters and numbers.

Fig. 13

2.15 Install the trinocular head over the splitter, directly onto the microscope body. (Fig. 14)

Fig. 14

2.16 Once completely assembled, the B-1000Ti-5 looks like this.

Fig. 15

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OPTIKA® S.r.l.