**B-150V Series**

# INSTRUCTION MANUAL

Model
B-151R-PLV
B-152R-PLV
B-159R-PLV
B-151V
B-153V

Ver. 1.3    2023

## Table of contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
3.1	B-151V / B-151R-PLV	4
3.2	B-152R-PLV / B-153V	4
3.3	B-159R-PLV	5
4.	Unpacking	6
5.	Intended use	6
6.	Symbols and conventions	6
7.	Instrument description	7
7.1	B-151V / B-151R-PLV	7
7.2	B-152R-PLV / B-153V / B-159R-PLV	8
8.	Assembling	10
8.1	Assembling the microscope with observation head	10
8.2	Assembling the microscope with digital head	10
8.3	Polarizing set (optional)	11
9.	Use of the microscope	12
9.1	Light intensity adjustment	12
9.2	Coarse focus tension adjustment	12
9.3	Stage	12
9.4	Use of diffusing filter	12
9.5	Use of oil immersion objective	13
9.6	Aperture diaphragm	13
9.7	Use with rechargeable batteries	14
9.8	Use of the polarizer (optional)	14
10.	Use of the digital head	15
10.1	Operation procedure	15
10.1.1	Connect the power adapter to the camera	15
10.1.2	Use of the rear buttons	15
10.1.3	Snap	15
10.2	Menu and Functions	16
10.3	Operating Instructions	17
10.3.1	White Balance	17
10.3.2	Exposure	17
10.3.3	Line Measurement	18
a)	Cross Line	18
b)	Calibration	18
c)	Measurement	19
10.3.4	Color Adjustment	19
10.3.5	Record Video	19
10.3.6	Other Functions	20
a)	Monochrome	20
b)	Flip	20
c)	Field of View (FOV)	20
10.3.7	Other Settings	20
a)	Power Frequency	20
b)	SD Management	21
c)	Timed Snap	21
d)	Language	21
e)	Factory Reset	22
f)	Software Upgrade	22
f.1)	System normally upgraded	22
f.2)	System upgrade failed	22
g)	Version	23
10.4	Failure analysis and troubleshooting	23
11.	Maintenance	24
12.	Troubleshooting	25
	Equipment disposal	26

---

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

## 2. Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

### 3. Package content

#### 3.1 B-151V / B-151R-PLV

- |                               |                               |
|-------------------------------|-------------------------------|
| ① Frame                       | ⑦ Green filter                |
| ② Monocular observation head  | ⑧ Power supply for microscope |
| ③ Eyepiece                    | ⑨ Camera + display            |
| ④ Tension adjustment tool     | ⑩ Power supply for display    |
| ⑤ Objectives (4X / 10X / 40X) | ⑪ Diffusing filter            |
| ⑥ Dust cover                  | ⑫ SD card + SD card reader    |

#### 3.2 B-152R-PLV / B-153V

- |                                |                            |
|--------------------------------|----------------------------|
| ① Frame                        | ⑥ Dust cover               |
| ② Monocular observation head   | ⑦ Green filter             |
| ③ Eyepiece                     | ⑧ Power supply             |
| ④ Tension adjustment tool      | ⑨ Camera + display         |
| ⑤ Objectives                   | ⑩ Power supply for display |
| • B-152 (4X / 10X / 40X)       | ⑪ Diffusing filter         |
| • B-153 (4X / 10X / 40X / 60X) | ⑫ SD card + SD card reader |

### 3.3 B-159R-PLV

- ① Frame
- ② Binocular observation head
- ③ Eyepieces
- ④ Tension adjustment tool
- ⑤ Objectives (4X / 10X / 40X / 100X)
- ⑥ Dust cover
- ⑦ Green filter
- ⑧ Power supply
- ⑨ Immersion oil
- ⑩ Camera + display
- ⑪ Power supply for display
- ⑫ Diffusing filter
- ⑬ SD card + SD card reader

---

## 4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 5. Intended use

### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

## 6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


### ELECTRICAL SHOCK


This symbol indicates a risk of electrical shock.

## 7. Instrument description


### 7.1 B-151V / B-151R-PLV

7.2 B-152R-PLV / B-153V / B-159R-PLV

SD CARD SLOT

SNAP BUTTON

ON/OFF BUTTON

MENU CONTROL  
BUTTONS

POWER SUPPLY  
CONNECTOR

## 8. Assembling

### 8.1 Assembling the microscope with observation head

1. Remove the dust cap from the microscope frame and from the bottom of the observation head.
  2. Insert the optical head above the stand and tighten the screw ①. (Fig. 1)
- **Hold the head with one hand during the locking in order to avoid that the head falls.**

3. Insert the eyepiece into the empty tube of the optical head. (Fig. 2)

4. Insert the microscope power supply jack in the socket placed at the rear side of the microscope. (Fig. 3)

### 8.2 Assembling the microscope with digital head

1. Remove the dust cap from the microscope frame and from the bottom of the digital head.
  2. Insert the digital head above the stand and tighten the screw ①. (Fig. 4)
- **Hold the head with one hand during the locking in order to avoid that the head falls.**

### 8.3 Polarizing set (optional)


1. Place the polarizer on the light exit ① at the base of the microscope. (Fig. 5)

2. Loosen the head fixing knob ② and remove the head or the digital head from the microscope frame. (Fig. 6)

3. Insert the analyzer ③ into the hole inside the frame. (Fig. 7)
4. Put back the head into its original position and lock the fixing knob.

## 9. Use of the microscope

### 9.1 Light intensity adjustment

Operate on the light intensity adjustment dial ① to turn ON / OFF the microscope and to increase / decrease the illumination voltage. (Fig. 8)

Fig. 8

### 9.2 Coarse focus tension adjustment

- **Adjust the tension using the provided tool.**

The coarse knob tension is pre-set in the factory.

1. To modify the tension according to personal's needs, rotate the ring ② using the provided tool (Fig. 9).
- Clockwise rotation increases the tension.
- If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

Fig. 9

### 9.3 Stage

#### All models except B-151 series

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverslide 0,17mm. (Fig. 10)

1. Open the spring arm of the slide holder ③ and place the slide from the front on the stage.
2. Gently release the spring arm of the slide holder.
- **A sudden release of the spring arm could cause the falling of the slide.**

Fig. 10

### 9.4 Use of diffusing filter

#### All models with Digital head

In order to ensure an even illumination with low magnification objectives, it is strongly suggested to use a diffusing filter to be placed in the filter holder placed under the condenser.

1. Rotate the filter holder ④ placed below the condenser.
2. Insert the diffusing filter ⑤ in its place. (Fig. 11)
3. Put the filter holder into the light path.

Fig. 11

## 9.5 Use of oil immersion objective


### B-159R-PLV only

1. Focus the specimen with a low power objective.
  2. Lower the stage.
  3. Put a drop of oil (provided) on the area of the specimen to be observed. (Fig. 12)
- **Make sure that there are no oil bubbles. Air bubbles in the oil damage the image quality.**
  - To check for bubbles: remove an eyepiece, fully open the aperture diaphragm and observe the objective exit pupil. (The pupil must be circular and bright).
  - To remove the bubbles, gently move the nosepiece to the right and left to move the immersion objective a few times and allow the air bubbles to move.
4. Insert immersion objective.
  5. Return the stage to the upper focusing point and obtain an optimal focus using the fine focus knob.
  6. After use, gently remove the oil with a soft paper towel or a lightly moistened optic paper with a mixture of ethyl ether (70%) and absolute ethyl alcohol (30%).
- **The immersion oil, if not immediately cleaned, could crystallize creating a glass-like layer. In this situation the observation of the specimen would be difficult (even not impossible) due to the presence of an additional thickness on the objective.**

## 9.6 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image. Move the diaphragm lever ① (Fig. 13) toward left or right to decrease or increase the N.A. value.
- With low contrast specimens set the numerical aperture to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's diaphragm in order to obtain an image like the one in Fig. 14.

## 9.7 Use with rechargeable batteries

This operation can be done only with “R” models.

When the microscope is plugged with the power supply, the LED indicator for the battery recharge ② is lit. (Fig. 15)

- Power supply connected, battery low: LED is RED. When charging is completed, it turns GREEN.
- Power supply not connected, microscope off: LED is always OFF.
- Power supply not connected, microscope on, battery charged: LED is GREEN.
- Power supply not connected, microscope on, battery almost low: LED is RED.

## 9.8 Use of the polarizer (optional)


1. Remove the specimen from the stage.
2. Looking inside the eyepieces, rotate the polarizer until the darkest position is achieved.
3. Once the dark is achieved (“extinction” or “Crossed Nicol” position) it is possible to begin the observation.

## 10. Use of the digital head

### 10.1 Operation procedure

#### 10.1.1 Connect the power adapter to the camera

1. Insert the power adapter into the power socket ① placed on the back side of the camera. (Fig. 16)
- After the connection, the red light will light up.


2. Press the power button ②. (Fig. 17)
  3. The indicator light will change from red to green, and the camera will start up.
- **Warning: disconnect the power supply if the equipment is not used for a long time.**

#### 10.1.2 Use of the rear buttons

1. Use the rear buttons to operate the camera functions and adjust the camera parameters.
2. After modifying the parameters, a dialog box appears allowing to save the modified parameters.
3. In the upper left corner of the screen, "Param. saved " appears. (Fig. 18)

#### 10.1.3 Snap

- The snap button ③ is on the right side of the camera, above the power button. (Fig. 17)
1. Press it to capture the current image on the screen and store it in the SD card.
  2. The screen shows "Snap Succeeded", which means that the picture has been taken successfully. (Fig. 19)

## 10.2 Menu and Functions

1. Press the MENU button to call out the menu.
2. The position of the current cursor (that is, the position of the highlighted icon) is the white balance function option. (Fig. 20)
3. Press  $\uparrow\downarrow$  for function selection, press  $\rightarrow$  to enter the sub-menu interface of the corresponding functions, press MENU to hide the interface, and save all parameters that have been modified.

Fig. 20

4. The specific functions of this menu are shown in Fig. 21.



Fig. 21


## 10.3 Operating Instructions

### 10.3.1 White Balance

After entering the White Balance menu, the default option is "Auto White Balance". (Fig. 22)


When the effect of automatic white balance is not ideal due to the difference of color temperature between different light sources, manual white balance can be used to adjust the parameters of color temperature, R, G and B respectively. (Fig. 23)


### 10.3.2 Exposure

After entering the exposure menu, the default option is "Auto exposure". (Fig. 24)

- Under automatic exposure, you can adjust the target "value" to adjust the degree of exposure.


- In manual exposure, you can also adjust the exposure by adjusting the value of "exposure time". (Fig. 25)

### 10.3.3 Line Measurement


This menu includes Cross Line, Calibration and Measurement. (Fig. 26)

#### a) Cross Line

- Four groups of cross lines are provided in red, blue, green and white colors. You can choose according to your requirements.
1. Enter the Cross Line menu. (Fig. 27)
- "ID" refers to the number of each group of cross lines.
  - "Display" informs whether the lines are displayed.
  - "X Position" and "Y Position" adjust the position of the center point of the lines.

You can select and press the "Turn off All Cross lines" to hide all cross lines.



#### b) Calibration

- Default calibration values are available for this system. However, due to the different objective standards of the microscope, the calibration value may have some discrepancies, so it is advisable to perform a new calibration.
- Calibration procedure requires a calibration micrometer.


The calibration procedure is described below.

1. Put the micrometer on the stage and adjust the microscope so that the micrometer scale is clearly displayed on the screen.
  2. In order to facilitate calibration, it is suggested to rotate the camera so that the micrometer is placed horizontally in the screen without being hidden by the menu.
  3. After entering the Calibration menu, the dialogue box in Fig. 28 appears.
  4. Adjust the positions of the calibration start and end points to make the calibration line coincide with the micrometer scale and try to select the maximum length that can be displayed, in order to make the most accurate measurement. (Fig. 28)
  5. The minimum range of the selected micrometer is 0.01 mm (10 micron). Fig.29 shows the image under a 10x objective lens. At this time, the "magnification" is set to "10X", the "unit" is marked as " $\mu\text{m}$ ", and the "length" is set to "40".
  6. Repeat the same procedure for all objective available on the microscope.
- After adjusting the parameters, exit the calibration interface: the calibration is completed.


### c) Measurement

- The system needs to be calibrated before measurement can be performed.
1. Enter the Measurement menu.
  2. Select the measurement magnification.
  3. Move the starting and ending point of X and Y position by pressing ← and →: the measurement length is displayed in real time. (Fig. 30)
- Changes in the percentage of the field of view do not affect measurements.


### 10.3.4 Color Adjustment

1. After entering the color adjustment menu (Fig. 31), Saturation, Hue, Brightness and Contrast can be adjusted to make the picture reach the required level.
  2. In order to facilitate the color adjustment, the “Color Value Reset” option is available in the menu.
- When selected and pressed, all the color values in the menu will be restored to the default value.

### 10.3.5 Record Video

1. Press M button to show menu in the screen (Fig. 32) and choose Record Video option by pressing ↑ or ↓.
  2. Press → button to start recording. The recording time is shown in the right part of the screen beside the dialog.
  3. If you want to quit recording, just press M button.
- **Before recording the video, you should check whether to insert an SD card with a FAT32 file system and free space. You cannot take a photo during the recording process.**

### 10.3.6 Other Functions

The menu includes functions of Monochrome, Flip vertical, Flip horizontal and FOV. (Fig. 33)

1. The numerical option can be adjusted by ← and →, and the switch options can be opened and closed by →.
2. After the adjustment is completed and the main menu is closed, the settings will be saved.
3. The function settings will be retained when the next boot is started.

#### a) Monochrome

The function produces visual images in varying tones of a single color (gray).


#### b) Flip

This function is divided into vertical flipping and horizontal flipping.

#### c) Field of View (FOV)

This function can adjust the range of observation field.

1. Use ← and → in the menu to adjust the size.
2. When the menu is not displayed, press ↑ and ↓ to adjust: the current percentage value appears in the upper left corner.



### 10.3.7 Other Settings

This menu contains Power Frequency, SD Management, Timed Snap, Language, factory Reset, Software Upgrade, Version. (Fig. 34)

#### a) Power Frequency


CMOS sensors have a rolling curtain effect that causes flicker problems, which can be solved by capturing a line of pixels as an integer (n) times the flicker period. Among them, 60Hz in North America and 50Hz in Europe. (Fig. 35)

- DC (DC): for DC (DC) light source, there is no light fluctuation, so there is no need to compensate the flashing light source.
- AC (50Hz): select AC (50Hz) to eliminate the dark strip of the lamp curtain caused by the 50Hz light source.
- AC (60Hz): select AC (60Hz) to eliminate the dark strip of the lamp curtain caused by 60Hz light source.

### b) SD Management


1. After inserting the SD card, the remaining space and total space of the SD card can be seen in "Available Space". (Fig. 36)
2. If "0.00 Gb / 0.00 Gb" is displayed (Fig. 37), the SD card was not successfully inserted; please try to reinsert it.

### c) Timed Snap

"Hours, Minutes and Seconds" refers to the time interval of timed snap, and "Counts" refers to the number of snaps to be done.

1. After setting the parameters, move the cursor to "Timed Snap Start" and press → to start the acquisition.
  2. At this point, the number jumps below. This is the number of photos that have been successfully taken so far. (Fig. 38)
- If the available space of SD card is insufficient during the process of timed snap, the process will be aborted.


### d) Language

The current version can switch between Chinese and English. (Fig. 39)

### e) Factory Reset

Press → to reset the menu settings to factory settings. (Fig. 40).

### f) Software Upgrade

After the system upgrade, the menu parameters will be restored to the factory settings. It is advisable to record the parameters before the upgrade, so they can be easily restored after the upgrade.

#### f.1) System normally upgraded

When the system upgrade files are published, the upgrade files can be saved into the SD card and the system upgrade can be carried out in this page.

The file name for the system upgrade is similar to this:


```
main_app_v1.0.bin    rootfs_uclibc_64k_v1.0.jffs2
```

1. Select “Yes” to upgrade, and the following message will appear: “UPGRADING...” while the system is being upgraded. (Fig 41)
- **WARNING: Please wait patiently for 2-3 minutes, during which time do not operate the device and keep power connected.**
- If “FILES ERROR” or “NO FILE” appears, please check whether the upgrade file is missing and the version number is corresponding. After checking, please upgrade the system again.

#### f.2) System upgrade failed


When a system upgrade fails, it goes into the tiny system, which is used for emergency updates.

1. After entering the tiny system, “*Upgrade Failed. Please try again according to the instructions.*” will be displayed on the screen.
2. Insert SD card with upgrade files, then the screen shows “*Files detected, Press Menu to Upgrade.*”
3. When you press the Menu button, the screen shows “*Upgrading... Please do not power off.*”
  - When the upgrade is completed, the system will be automatically restarted to complete the upgrade.
4. If the screen shows “*The version is illegal, please check the file.*” means the upgrade file is missing or the version number does not match.
5. Please check the file version number and update it again.



### g) Version

You can view the version information of the system. (Fig. 42)

## 10.4 Failure analysis and troubleshooting

1. When pressing the snap button, recording videos, taking photos at a fixed time, or entering the SD card management menu, "NO SD-CARD" will be displayed (Fig. 43). Please insert the SD card with file system FAT32 into the SD card slot on the right side of the camera, and then carry out corresponding operations.
2. When pressing the snap button, recording videos, and taking photos at a fixed time, "No Available Space" will be displayed (Fig. 44). Please free up some space on the SD card before inserting the card into the camera for use.
3. The camera's screen image shows horizontal lines. Open "Other Settings" - "Power Frequency" and select the appropriate power frequency.
4. The image is blurry and out of focus. Clean the objective or perform the focus again.
5. If the camera screen remains black when the power is turned on and restarts in about 15 seconds, the snap button or power switch may be locked. Return them to their normal positions.
6. If an unknown problem occurs and cannot be solved by yourself, please press the power button for 10 seconds to restart. If this operation doesn't solve the problem, please contact the manufacturer.

---

## 11. Maintenance

### Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

### To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

### Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

**For the best results, use the OPTIKA cleaning kit (see catalogue).**

If you need to send the microscope to Optika for maintenance, please use the original packaging.

## 12. Troubleshooting


Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
<b>I. Optical Section:</b>		
LED operates, but field of view remains dark	Power supply is unplugged	Connect
	Brightness is too low	Set brightness to a proper level
	Batteries are uncharged	Fully charge the batteries
Dirt or dust is visible in the field of view	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture diaphragm is stopped down too far	Open aperture diaphragm
Visibility is poor <ul style="list-style-type: none"> <li>• Image is not good</li> <li>• Contrast is poor</li> <li>• Details are indistinct</li> <li>• Image glares</li> </ul>	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture diaphragm is too closed or too open	Adjust aperture diaphragm
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly
One side of the image is out of focus	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The specimen is out of place (tilted)	Place the specimen flat on the stage
	The optical performance of the sample cover glass is poor	Use a cover glass of better quality
<b>II. Mechanical Section:</b>		
The coarse focus knob is hard to turn	The tension adjustment collar is too tight	Loosen the tension adjustment collar
The focus is unstable	The tension adjustment collar is too loose	Tighten the tension adjustment collar
<b>III. Electric section:</b>		
The LED doesn't turn on	No power supply	Check the power cord connection
	Batteries are uncharged	Fully charge the batteries
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord
<b>IV. Observation tube:</b>		
Field of view of one eye does not match that of the other	Interpupillary distance is incorrect	Adjust interpupillary distance
	Incorrect diopter adjustment	Adjust diopter
	Your view is not accustomed to microscope observation	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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**OPTIKA® S.r.l.**