B-380 Series

INSTRUCTION MANUAL

Model
B-383POL

Ver. 4.4 2023

Table of contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
4.	Unpacking	5
5.	Intended use	5
6.	Symbols and conventions	5
7.	Instrument description	6
8.	Assembling	8
9.	Summary of brightfield observation procedures	10
10.	Use of the microscope in brightfield	11
10.1	Light intensity adjustment	11
10.2	Adjust the interpupillary distance	11
10.3	Diopter adjustment	11
10.4	Coarse focus tension adjustment	11
10.5	Focus lock lever	12
10.6	Stage	12
10.7	Condenser centering	12
10.8	Aperture diaphragm	12
11.	Use of the microscope in polarized light	13
11.1	Centering the nosepiece	13
11.2	Checking the extinction of light	14
11.3	Use of Tint plates	14
11.4	Use of Bertrand lens	15
12.	Microphotography	16
12.1	Use of C-mount cameras	16
12.2	Use of reflex cameras	16
13.	Maintenance	17
14.	Troubleshooting	18
	Equipment disposal	20

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

- ① Frame
- ② Objectives
- ③ Observation head
- ④ Eyepieces
- ⑤ Bertrand lens
- ⑥ Analyzer
- ⑦ Tint plates
- ⑧ Tension adjustment tool
- ⑨ Allen wrenches
- ⑩ Dust cover
- ⑪ Power supply

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

7. Instrument description

Opposite side

8. Assembling

1. Insert the Bertrand lens ① in the frame and lock the locking screw ② with the provided Allen wrench. (Fig.1)

2. Insert the optical head on the Bertrand lens and lock the locking screw with the provided Allen wrench. (Fig. 2)

3. Insert both eyepieces into the tubes of the optical head. (Fig. 3)
 - **One of the two eyepieces is equipped with a crosshair for centering the entire optical system. It is advisable to insert the eyepiece with crosshair in the right eyepiece holder.**
 - The condenser is pre-installed before leaving the factory.
 - To remove the condenser use an Allen wrench diam. 1,5 and operate on the locking screw placed on the right side of the condenser holder.


4. Remove the dummy slider from Bertrand lens and insert the analyzer ③. (Fig. 4)

-
5. Insert the power supply jack in the connector placed at the rear side of the microscope. (Fig. 5)

9. Summary of brightfield observation procedures

10. Use of the microscope in brightfield

10.1 Light intensity adjustment

1. Operate on the light intensity dial ① to turn ON/OFF the microscope and to increase or decrease the illumination intensity. (Fig. 6)

10.2 Adjust the interpupillary distance

1. Hold the right and left parts of the observation head using both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. (Fig. 7)
- The graduation on the interpupillary distance indicator, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator’s eyes.
 - The range of the interpupillary distance is 48- 75 mm

10.3 Diopter adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
 2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the diopter adjustment ring ② to compensate. (Fig. 8)
- Highpoint eyepieces allow the use also to glass wearers.
 - NOTE: For optimal parfocality, we recommend using your glasses during normal microscope use

10.4 Coarse focus tension adjustment

- Adjust the tension using the provided tool.

The coarse knob tension is pre-set in the factory.

1. To modify the tension according to personal's needs, rotate the ring ③ using the provided tool (Fig. 9).
- Clockwise rotation increases the tension.
 - If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

10.5 Focus lock lever


The upper limit knob has two functions: prevent the contact between slide and objective and acts as “focus memory”.

1. After focusing the specimen, pull the lever ① toward the front of the microscope and lock it. (Fig. 9).
- In this way the focus upper limit is set.
2. Now one can lower the stage with coarse focus knob, replace the specimen and raise again the stage up to the upper limit: specimen will be in approximate focus and will need a fine adjustment to get the proper focus..
- **Fine focus movement is not affected by the coarse focus lock.**
 - **To unlock, move the lever in the opposite direction to the one used for the locking.**

10.6 Stage

- The rotating stage accepts specimens on slide.
 - It is possible to lock the specimen once placed on the stage using the stage clips ②. (Fig. 11)
1. After loosening the locking screw ③, the stage can be rotated by 360°

10.7 Condenser centering

- The condenser is installed and pre-centered in the factory.


Should a new centering is needed, operate in this way:

1. Insert 4x objective in the light path (in case 4x is not available use the lower magnification available).
2. Focus the specimen.
3. Close the aperture diaphragm using the ring ①, moving the ring to the value “4” related to the 4x objective. (Fig. 12)
4. Raise the condenser to the upper limit using the height adjustment knob ② placed on the left side of the condenser holder.
5. Center the condenser using the centering screws ③ until the field of view is evenly illuminated (in the field of view no dark and bright areas must be noticed).
6. Fully open the diaphragm.

10.8 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image. Move the diaphragm ring ① (Fig. 12) on the value corresponding to the objective in use. In this case the optimal setting of the condenser is achieved.
- With low contrast specimens set the numerical aperture to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's diaphragm in order to obtain an image like the one in Fig. 13.

11. Use of the microscope in polarized light


- The system allows observation in Orthoscopy (crossed Nicol) or in Conoscopy (crossed Nicol with the use of the Bertrand lens).
- For optimal performance in polarized light microscopy, accurate optical adjustments are essential before beginning the observation.

11.1 Centering the nosepiece

1. Loosen the stage rotation lock ① and rotate the stage until the graduated scale of the stage ② and the vernier scale ③ are aligned in the "0" position. (Fig. 14)
- **This operation serves to ensure a standard reference position for centering the nosepiece.**


2. Focus on a recognizable detail ④ in the field of view by placing it at the centre of the crosshair. (Fig. 15)
3. By turning the stage, the detail in focus will describe a circle ⑤. (Fig. 15)
4. Move again the stage on the "0" position and lock the locking screw ①.
5. Operating on the nosepiece centering screws ⑥ move the detail in diametrically opposite direction to the described circle. The movement must be about half the diameter of the circle described. (Fig. 16)
6. Manually move the specimen and put it again on the centre of the crosshair. Loosen again the stage rotation lock and rotate once more the stage.
7. If the centering has been carried out correctly, by turning the table the image of the focused detail does not move with respect to the centre of the crosshair. If not, repeat the operations described by 2. to 6. until the perfect coincidence of the centre of rotation of the stage with the centre of the crosshair for which the specimen remains in the centre of the crosshair turning the table.
8. Once centered with the 10x, rotate the nosepiece to insert in the optical path all the other objectives and verify the correct centering of the objectives, acting on the screws of centering the nosepiece ⑥, to make sure that all the objectives are perfectly centered with respect to the optical axis.

11.2 Checking the extinction of light

1. Remove the specimen from the light path and insert 10x objective.
 2. Loosen the locking screw of the polarizer ① and check that the scale is in the "0" position ②. (Fig. 17)
 3. Insert rotatable analyzer in the light path, loosen the locking screw ③ and put the vibration scale on 0° ④, then lock the locking screw ③. (Fig. 18)
 4. Loosen the polarizer locking screw ① and rotate the polarizer scale ② to obtain total extinction (total dark in the eyepieces or "Crossed Nicol position"). Tighten the screws ①. (Fig. 17)
- It may happen that the scale of the polarizer is not perfectly aligned on the reference mark but is shifted by one or two notches. This is not a defect but is due to the mechanical alignment of the polarizers during assembly.

11.3 Use of Tint plates


Three tint plates are supplied with the microscope:

- λ plate (1st order Red)
 - $\lambda/4$ plate
 - "Quartz wedge" plate (Q)
1. Put a specimen on the stage and focus.
 2. Insert in the right slot of the Bertrand lens ⑤ one of the tint plates ⑥. (Fig. 19)
 3. In polarized light, inserting one of the plates will have chromatic effects on the specimen.
- Using the λ plate (also called 1st order Red) the specimen will take a magenta tinge.
 - Using the $\lambda/4$ plate the specimen will take on a colour tending to pale yellow.
 - Using the Q plate the specimen will present a series of colored bands that will fade as the plate is inserted.

11.4 Use of Bertrand lens


- Bertrand lens allows observation in Orthoscopy and Conoscopy.
 - In the disengaged position (“O”) the lens allows observation in Orthoscopy, while in the inserted position (“B”) it is possible to make observations in Conoscopy.
1. Rotate the upper knurled ring of the Bertrand lens ① to engage the “B” position. (Fig. 20)
 2. Using one objective from 20x to 60x, focus the conoscopic image using the focus ring ②.
 3. If the conoscopic image is not perfectly centered with respect to the optical axis, centre the image using the centering screws ③.
- By turning the stage you will see black fringes that will appear and disappear depending on the rotation of the stage. These fringes are the crystallization axes of that specific crystal. (Fig. 21)

12. Microphotography

12.1 Use of C-mount cameras

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 22)

2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 23)

12.2 Use of reflex cameras

1. Insert the Reflex adapter ① into the relay tube to the microscope ②.
 2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
 3. Connect the Reflex camera ④ to the "T2" just installed (Fig. 24).
 4. Mount the other end of the relay tube ② into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 22)
- "T2" ring is not provided with the microscope, but is commercially available.
 - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
 - To calculate the magnification of the camera: objective magnification * camera magnification * lens magnification.
 - **If using an SLR camera, mirror movement may cause the camera to vibrate.**
 - **We suggest lifting the mirror, using long exposure times and a remote cord.**

13. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off- position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

14. Troubleshooting

Review the information in the table below to troubleshoot operating problems.


PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark.	Power supply is unplugged.	Connect
	Brightness is too low	Set brightness to a proper level
	Bertrand lens is in	Remove Bertrand lens from light path
	You are in a position of extinction	Disengage analyzer from the light path
Field of view is obscured or not evenly illuminated	Nosepiece is not correctly engaged.	Make sure that the nosepiece clicks properly into place.
	Tint plate or Bertrand lens are in a intermediate position.	Move to a click stop
Dirt or dust is visible in the field of view.	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture iris diaphragm is stopped down too far.	Open aperture iris diaphragm.
	The condenser is not well centered or it is in a wrong height	Set the condenser until a proper image is achieved
Visibility is poor. <ul style="list-style-type: none"> • Image is not sharp. • Contrast is poor • Details are indistinct • Image glares 	Nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture iris diaphragm is too closed or too open.	Adjust aperture iris diaphragm.
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly.
	For transmitted light observation, the coverglass thickness must not exceed 0.17mm	Use a coverglass with thickness 0.17mm
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position.
One side of the image is unfocused	Nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Slide is mounted not in a flat position (tilted)	Place the specimen in a flat position on the stage
	Poor quality of the glass slide	Use a glass slide with higher quality
Conoscopic image cannot be seen	Bertrand lens is not in the light path.	Insert the lens in the light path
Total extinction cannot be obtained	Analyzer is not in the light path	Insert analyzer in the light path.
II. Mechanical Section:		
Coarse focus knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
Focus is unstable	Tension adjustment ring is too loose	Tighten tension adjustment ring
III. Electrical Section		
LED doesn't turn on.	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection
IV. Observation tube		
Field of view of one eye does not match that of the other.	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope.

V. Microphotography

Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.