B-510 Series

INSTRUCTION MANUAL

Model
B-510BF4K

Ver. 1.0 2023

Table of contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
4.	Unpacking	5
5.	Intended use	5
6.	Symbols and conventions	5
7.	Instrument description	6
8.	Assembling	9
9.	Observation procedures	12
10.	Use of the microscope	13
10.1	Turning on the camera	13
10.2	Turning on the monitor	13
10.3	Turning on the microscope	13
10.4	Parfocality adjustment	14
10.5	Centering the condenser	14
10.6	Effects of the field diaphragm	15
10.7	Coarse focus tension adjustment	16
10.8	Coarse upper limit lever	16
10.9	Stage	16
10.10	Aperture diaphragm	16
10.11	Use of oil immersion objective	17
11.	Use of the camera	18
11.1	HDMI Mode	18
11.2	USB Mode	18
12.	C-HP4 User Interface (UI) and its Functions	19
12.1	The Camera Control Panel	20
12.2	The Measurement Toolbar	21
12.3	Camera calibration	23
12.4	Icons and Functions of the Camera Control Toolbar	25
12.4.1	Setting > Measurement	25
12.4.2	Setting > Magnification	26
12.4.3	Setting > Image Format	27
12.4.4	Setting > Video	27
12.4.5	Setting > Storage	28
12.4.6	Setting > Files	29
12.4.7	Setting > Language	29
12.4.8	Setting > Miscellaneous	30
13.	Micrometric Slide M-005	31
14.	Maintenance	32
15.	Troubleshooting	33
	Equipment disposal	34

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

- ① Microscope frame
- ② Objectives
- ③ Microscope power supply
- ④ Spacer
- ⑤ "C" mount
- ⑥ Camera with fixing plate
- ⑦ Camera power supply
- ⑧ Wireless mouse
- ⑨ HDMI monitor with fixing plate
- ⑩ Monitor power supply

- ⑪ Type "C" L-adapter
- ⑫ HDMI cable
- ⑬ HDMI-to-mini HDMI L-adapter
- ⑭ USB cable
- ⑮ SD card
- ⑯ Micrometric slide
- ⑰ Dust cover
- ⑱ Tension adjustment tool
- ⑲ Allen wrench
- ⑳ Immersion oil

NOTE: OPTIKA reserves the right to make corrections, modifications, enhancements, improvements and other changes to its products at any time without notice.

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION


This symbol indicates a potential risk and alerts you to proceed with caution.


ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

7. Instrument description

Opposite side


Camera control panel

- ① USB mouse connector
- ② USB2.0 connector
- ③ HDMI connector
- ④ ON-OFF button

- ⑤ SD card slot
- ⑥ Power supply plug
- ⑦ LED indicator

HDMI monitor

- ① HDMI camera connection
- ② Video Display / Audio / Touch Control


- ③ Power Supply / USB3.0 port
- ④ Keyboard / Mouse connection

NOTE: When using the monitor in this configuration, the Headphone, Type C ② and Micro USB ④ ports are not used.


8. Assembling

- The condenser is pre-installed in the factory. To remove the condenser use an Allen wrench 1,5 mm diam and operate on the locking screw placed on the right side of the condenser holder.

1. Insert the spacer ① above the stand and tighten the screw ②. (Fig. 1)

2. Using the supplied screws, join the fixing plate mounted on the camera ③ and the fixing plate mounted on the monitor ④. (Fig. 2)

3. Screw the "C" mount adapter ⑤ to the camera ⑥ (Fig. 3) and insert the round dovetail of the "C" mount into the empty hole of the spacer, then tighten the clamping screw ⑦. (Fig. 4)

- Screw each objective into the thread of the nosepiece, clockwise with increasing magnification. (Fig. 5)

- Insert the microscope power supply jack in the connector placed at the rear side of the microscope. (Fig. 6)
- Plug the power supply to the wall outlet.

- Plug the USB mouse receiver in the "USB MOUSE" port ① on the camera. (Fig. 7)

- Connect the cable with the HDMI side ① to the provided HD-MI-HDMI L-adapter ②. (Fig. 8)
- Connect the HDMI L-adapter to the "HDMI" port ② on the camera.

10. Plug mini HDMI cable to the mini HDMI port ③ placed on the right side of the monitor. (Fig. 9)


11. Insert the SD card into “SD” slot ④ for saving captured images and recorded videos. (Fig. 10)

12. Plug the monitor power supply connector ③ in the provided “L” adapter ④. (Fig. 11)

13. Plug the “L” adapter in the “Type C” port ⑤ placed on the right side of the monitor.

14. Plug the power supply to the wall outlet.


15. Insert the camera power supply jack in the “DC12V” socket ⑥. (Fig. 12)

16. Plug the power supply to the wall outlet.

9. Observation procedures

10. Use of the microscope

10.1 Turning on the camera

- The LED indicator ① (Fig. 13), when the power supply is connected but the camera is switched off, is red.
1. When the ON/OFF button is pressed, the LED ② starts flashing blue for a few seconds.
 2. Once the switching sequence is complete, the LED ② is steady blue.
 3. The camera is ready to work.

10.2 Turning on the monitor

1. Press the POWER button ③ placed on the left side of the monitor (Fig. 14) to turn on the monitor.
- The monitor takes a few seconds to become operational.

10.3 Turning on the microscope

1. Operate on the light intensity adjustment knob ④ to turn ON/OFF the microscope and to increase/decrease the illumination voltage. (Fig. 15)

10.4 Parfocality adjustment

- The parfocality must be adjusted by means of the focusable “C” mount adapter connected to the camera in order to achieve maximum performance.
1. Insert the lower magnification objective into the optical path.
 - **If the lower magnification objective is a 2x or 4x, the front lens of the condenser must be removed from the optical path.**
 2. While observing the image on the monitor, focus using the microscope’s focus knobs.
 3. Insert the higher magnification dry objective.
 4. While observing the image on the monitor, focus using the microscope’s focus knobs.
 5. Insert the lower magnification objective again.
 6. If the image is not perfectly in focus, turn the knurled ring ⑤ on the “C” mount (Fig. 16), until the image is perfectly in focus.
 - **Doing this operation DO NOT turn the microscope’s focus knobs.**
 7. The parfocality adjustment is completed.

Fig. 16

10.5 Centering the condenser

1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Insert the front lens of the swing-out condenser ①. (Fig. 17)
3. Move the mouse cursor to the bottom of the video window, the *Camera Control Toolbar* will pop up automatically. (Fig. 18)
4. From here, select the “Grid” command ⑤.
5. The “Grids” dialog box appears (Fig. 19). In the dialogue box select the value “1” in the “Row” field and the value “1” in the “Column” field.
6. A crosshair appears in the centre of the monitor. (Fig. 20)
7. Rotate the field diaphragm ring ② in counterclockwise direction, to fully close the diaphragm. (Fig. 17)
8. A bright spot appears in the screen. (Fig. 21)
9. Rotate the condenser height adjustment knob ③ (Fig. 17) to focus the edges of the diaphragm.
10. Rotate the two centering screws ④ (Fig. 17) to bring the bright spot in the center of the screen. (Fig. 22)
11. Gradually open the diaphragm until it is outside the screen.
12. At the end hide the crosshair by deselecting the check box “Visible” in the “Grids” dialog box.

Fig. 17





Fig. 18

10.6 Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image. Set the diaphragm according to the objective in use until it circumscribes the monitor size, in order to eliminate unnecessary light to camera.

10.7 Coarse focus tension adjustment

- **Adjust the tension using the provided tool.**
- The coarse knob tension is pre-set in the factory.
- 1. To modify the tension according to personal's needs, rotate the ring ① using the provided tool. (Fig. 23)
- Clockwise rotation increases the tension.
- If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

10.8 Coarse upper limit lever

The upper limit knob has two functions: prevent the contact between slide and objective and acts as "focus memory".

1. After focusing the specimen, pull the lever ② toward the front of the microscope and lock it. (Fig. 24)
 - In this way the focus upper limit is set.
 2. Now one can lower the stage with coarse focus knob, replace the specimen and raise again the stage up to the upper limit: specimen will be in approximate focus and will need a fine adjustment to get the proper focus.
- **Fine focus movement is not affected by the coarse focus lock.**
 - **To unlock, move the lever in the opposite direction to the one used for the locking.**


10.9 Stage

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverside 0,17 mm.


It is possible to place two slides side by side on the stage.

- **Open the spring arm of the slide holder ③ and place frontally the slides on the stage. (Fig. 25)**
- **Gently release the spring arm of the slide holder.**
- **A sudden release of the spring arm could cause the falling of the slide.**

10.10 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
- With low contrast specimens set the numerical aperture value ④ (printed on the condenser ring) to about 70%-80% of the objective's N.A. (Fig. 26)
- **Example: with objective PLAN 40x/0,65 set the scale to $0.65 \times 0.8 = 0,52$**

10.11 Use of oil immersion objective

1. Focus the specimen with a low power objective.
2. Lower the stage (remembering to lock the coarse upper limit lever).
3. Put a drop of oil (provided) on the area of the specimen to be observed. (Fig. 27)
 - **Make sure that there are no oil bubbles. Air bubbles in the oil damage the image quality.**
 - To check for bubbles: remove an eyepiece, fully open the aperture diaphragm and observe the objective exit pupil. (The pupil must be round and bright).
 - To remove the bubbles, gently move the nosepiece to the right and left to move the immersion objective a few times and allow the air bubbles to move.
4. Insert immersion objective.
5. Return the stage to the upper focusing point and obtain an optimal focus using the fine focus knob.
6. After use, gently remove the oil with a soft paper towel or a lightly moistened optic paper with a mixture of ethyl ether (70%) and absolute ethyl alcohol (30%).
 - **The immersion oil, if not immediately cleaned, could crystallize creating a glass-like layer. In this situation the observation of the specimen would be difficult if not impossible due to the presence of an additional thickness on the objective.**

11. Use of the camera

C-HP4 is a multiple interfaces (HDMI + USB2.0 + SD card) CMOS camera. HDMI and USB2.0 are used as the data transfer interface to HDMI display or computer.

The system has been optimized for use of the camera in HDMI 4K mode.

However, it is possible to use the system in USB mode by connecting the camera to an external PC.

- For HDMI output, the *Camera Control Panel + Measurement Toolbar* and *Camera Control Toolbar* are overlaid on the HDMI screen when the mouse moves to the related region. In this case, the USB mouse can be used to set the camera, browse and compare the captured image, play the video and perform the measurement.
- For USB output, connect camera to a computer with USB connection to transfer video in real time. From here the ProView software can control each camera function.

11.1 HDMI Mode

1. Plug the HDMI cable into the HDMI port (3) to connect the camera to HDMI display.
2. Plug the USB mouse into USB Mouse port (1) to get control of the camera by using built-in software HDMIPRO.
3. Plug 12V/1A power adapter into DC 12V slot (6) to supply power for the camera. The LED Indicator (7) will turn into red.
4. Insert SD card into SD Card slot (5) for saving captured images and recorded videos.
5. Press ON/OFF Button (4) to turn on the camera. The LED indicator will turn into blue.
6. Move mouse cursor to the left side of the video window, a *Camera Control Panel* will appear. It includes *Manual/ Automatic Exposure, White Balance, Sharpness* and other functions, please refer to 12.1 for details.
7. Move mouse cursor to the bottom of the video window and a *Camera Control Toolbar* will appear. Operations like *Zoom In, Zoom Out, Flip, Freeze, Cross Line, Comparison* and so on can be achieved. Please refer to 12.4 for details.
8. Move mouse cursor to the upper side of the video window, a *Measurement Toolbar* with calibration and other measurement tools will appear, please refer to 12.2 for details. The measurement data can be output with *.CSV format.

11.2 USB Mode

1. Plug the USB cable provided with the camera into USB2.0 (2) to connect the camera to the computer.
 2. Plug 12V/1A power adapter into DC 12V slot (6) to supply power for the camera. The LED Indicator (7) will turn into red.
 3. Press ON/OFF button (4) to turn on the camera. The LED indicator will turn into blue.
 4. Start the ProView software.
 5. Clicking the camera name in the camera list starts the live video.
- **NOTE: When the USB cable and the mouse are plugged into the camera at the same time, the USB cable has priority and the mouse is not available; when the USB cable is unplugged, the mouse can be used normally.**

12. C-HP4 User Interface (UI) and its Functions

The C-HP4 UI shown in Fig. 28 includes a *Camera Control Panel* on the left side of the video window, a *Measurement Toolbar* on the upper side of the video window, a *Camera Control Toolbar* on the bottom of the video window.


Fig. 28 - Camera Control User Interface


1. Move the mouse cursor to the left side of the video window: the *Camera Control Panel* will pop up automatically.
2. Move the mouse cursor to the upper side of the video window: a *Measurement Toolbar* will pop up for the calibration and measurement operations.
 - When user left-clicks the *Float/Fixed* button  on the Measurement Toolbar, the Measurement Toolbar will be fixed. In this case the Camera Control Panel will not pop up automatically even if users move mouse cursor to left side of the video windows.
 - Only when user left-clicks the  button on the *Measurement Toolbar* to exit from measuring procedure it will be possible to perform other operations on the *Camera Control Panel*, or *Camera Control Toolbar*.
 - During the measuring process, when a specific measuring object is selected an *Object Location & Attributes Control Bar*  will appear for changing location and properties of the selected objects.
3. Move the mouse cursor to the bottom of the video window, the *Camera Control Toolbar* will pop up automatically.

12.1 The Camera Control Panel

The *Camera Control Panel* controls the camera to achieve the best image quality according to the specific applications. It will pop up automatically when mouse cursor is moved to the left side of the video window (in measurement status, the Camera Control Panel will not pop up. Only when measurement process is terminated the Camera Control Panel will pop up by moving mouse cursor to the left side of the video window).

Left-click  button to achieve *Display/ Auto Hide* switch of the Camera Control Panel.

Control Panel	Function	Description
	Snap	Capture image from the current video window
	Record	Record video from the current video window
	Auto Exposure	When Auto Exposure is checked, the system will automatically adjust exposure time according to the value of exposure compensation
	Exposure Compensation	Available when Auto Exposure is checked. Slide to left or right to adjust Exposure Compensation according to the current video brightness to achieve proper brightness value
	Exposure Time	Available when Auto Exposure is unchecked. Slide to left or right to reduce or increase exposure time, adjusting brightness of the video
	Gain	Adjust Gain to reduce or increase brightness of video. The Noise will be reduced or increased accordingly
	Auto White Balance	White Balance adjustment according to the window video every time the button is clicked
	Manual White Balance	Slide Red or Blue cursors to manually set the video White Balance
	ROI White Balance	When the ROI region is moved, white Balance can be adjusted according to content inside the ROI region.
	Red	Slide to left or right to decrease or increase the proportion of Red in RGB on video
	Green	Slide to left or right to decrease or increase the proportion of Green in RGB on video
	Blue	Slide to left or right to decrease or increase the proportion of Blue in RGB on the video
	Sharpness	Adjust Sharpness level of the video window
	Denoise	Slide left or right to denoise the image
	Saturation	Adjust Saturation level of the video window
	Gamma	Adjust Gamma level of the video. Slide to the right side to increase gamma and to the left to decrease gamma
	Contrast	Adjust Contrast level of the video. Slide to the right side to increase contrast and to the left to decrease contrast
	DC	For DC illumination, there will be no fluctuation in light source so no need for compensating light flickering
	AC(50HZ)	Check AC(50HZ) to eliminate flickering "strap" caused by 50Hz illumination
	AC(60HZ)	Check AC(60HZ) to eliminate flickering "strap" caused by 60Hz illumination
Default	Restore all the settings in the Camera Control Panel to default values	

12.2 The Measurement Toolbar

The *Measurement Toolbar* will pop up when moving mouse cursor to any place near the upper side of the video window.


Fig. 29 - The Measurement Toolbar

Icon	Function
	<i>Float/ Fix</i> switch of the <i>Measurement Toolbar</i>
<input checked="" type="checkbox"/> Visible	Define measuring object in <i>Show up/ Hide</i> mode
Nanometer (nm)	Select the desired <i>Measurement Unit</i>
4X	Select <i>Magnification for Measurement after Calibration</i>
	<i>Object Select</i>
	<i>Angle</i>
	<i>4 Points Angle</i>
	<i>Point</i>
	<i>Arbitrary Line</i>
	<i>3 Points Line</i>
	<i>Horizontal Line</i>
	<i>Vertical Line</i>
	<i>3 Points Vertical line</i>
	<i>Parallel</i>
	<i>Rectangle</i>
	<i>Ellipse</i>
	<i>5 Points Ellipse</i>
	<i>Circle</i>
	<i>3 Points Circle</i>
	<i>Annulus</i>
	<i>Two Circles and Center Distance</i>
	<i>3 Points Two Circles and Center Distance</i>
	<i>Arc</i>
	<i>Text</i>

	<i>Polygon</i>
	<i>Curve</i>
	<i>Scale Bar</i>
	<i>Arrow</i>
	Execute <i>Calibration</i> to determine the corresponding relation between magnification and resolution, which will establish the corresponding relationship between measurement unit and the sensor pixel size. <i>Calibration</i> needs to be done with the help of a micrometer. For detailed steps of carrying out <i>Calibration</i> please refer to chapter 9.3.
	Export the <i>Measurement</i> to CSV file(*.csv)
	<i>Measurement Setup</i>
	<i>Delete All</i> the Measurement Objects
	<i>Exit</i> from <i>Measurement</i> mode
	When the measurement ends, left-click on a single measuring object and the <i>Object Location & Properties Control Bar</i> will show up. User could move the object by dragging the object with the mouse. But more accurate movement could be done with the control bar. The icons on the control bar mean <i>Move Left, Move Right, Move Up, Move Down, Color Adjustment</i> and <i>Delete</i> .

Note:


1. When user left-clicks *Display/Hide*  button on the *Measurement Toolbar*, it will be fixed. In this case the *Camera Control Panel* will not pop up automatically even if moving the mouse cursor to the left edge of the video window. Only when user left-click the  button on the *Measurement Toolbar* to exit from the measurement mode it will be possible to perform other operations on the *Camera Control Panel* or on the *Camera Control Toolbar*.
2. When a specific item is selected during the measurement process, the *Object Location & Attributes Control Bar*  will appear for changing the object location and properties of the selected objects.

12.3 Camera calibration


The camera must be calibrated before any measurements can be taken. To do this it is necessary to use the micrometer slide supplied with the camera.

Detailed procedure is the following:


1. Insert the objective with the lowest magnification available and put the stage micrometer on the stage.
2. Click on the “Calibrate”  button in the Measurement toolbar.
3. The following dialog box appears:

4. In the *Magnification* field enter the current magnification
 - Some magnification values (default ones 4x, 10x, 20x, 40x and 100x) are preset in the drop-down list.
 - If user must insert a different magnification value, double click in the *Magnification* field. The following dialog box appears:

- From here user can enter any desired magnification value.
5. Focus the ruler under the microscope.
 - Along with the previous dialog box a red line appears on the image. This line is the calibration line.

-
6. Align the two ends of the red line with the ruler scale (try to use as many notches as possible to keep the calibration accuracy).
 7. Set the desired unit (micrometer, millimeter or nanometer) in the *Calibration* dialog box.
 8. Insert the length of the drawn line in the *Actual Length* field (using the same procedure described above).
 - Micrometric slide has two rulers (see chapter 13). One is 1mm/100 (for biological microscopes) and the other is 10mm/100 (for stereomicroscopes).
 - When using the ruler 1mm/100 the interval between two long lines is 100 μm , while using the ruler 10mm/100 the interval between two long lines is 1000 μm .
 9. If everything is fine, click *OK* to end the calibration. The new magnification (the value you entered in the Magnification edit box) will be available in the *Magnification* drop-down list.
 10. Repeat the same procedure for each magnification available in your microscope.

12.4 Icons and Functions of the Camera Control Toolbar

Icon	Function	Icon	Function
	Zoom In the Video Window		Zoom Out the Video Window
	Horizontal Flip		Vertical Flip
	Color / Gray		Video Freeze
	Display Cross Line		Overlay
	Compare stored image with the current live image		Browse Images and Videos in the SD Card
	Settings		Check the HDMIPRO version

The *Setting* function is relatively more complicated than the other functions. Here are more info about it:

12.4.1 Setting > Measurement


Fig. 30 - Measurement Setting Page

<i>Global</i>	<i>Precision</i>	Used for setting digits behind the decimal point for measurement results.		
<i>Calibration</i>	<i>Line Width</i>	Used for defining width of the lines for calibration.		
	<i>Color</i>	Used for defining color of the lines for calibration.		
	<i>EndPoint</i>	<i>Type:</i>	Used for defining shape of the endpoints of lines for calibration.	
		<i>Null:</i>	Means no endpoints.	
<i>Rectangle:</i>		Means rectangle type of endpoints. It makes alignment more easily.		
<i>Point, Angle, Line, Horizontal Line, Vertical Line, Rectangle, Circle, Ellipse, Annulus, Two Circles, Polygon, Curve</i> Left-click the along with the Measurement command mentioned above will unfold the corresponding attribute settings to set the individual property of the Measurement Objects.				

12.4.2 Setting > Magnification


Fig. 31 - Magnification Settings Page

<i>Name</i>	The name of the <i>Magnification</i> , usually the Magnification of the microscope objective is used as the Magnification name when calibrating, such as <i>4X</i> , <i>10X</i> , <i>40X</i> , <i>100X</i> , etc. Besides, other user-defined information could be added into the Magnification name too, for example, microscope model, operator name, etc.
<i>Resolution</i>	Pixels per meter. Devices like microscopes have high resolution value.
<i>Clear All</i>	Click the <i>Clear All</i> button will clear the calibrated magnifications and resolutions.
<i>Delete</i>	Click <i>Delete</i> to delete the selected item for specific resolution.
<i>Up</i>	Click <i>Move Up</i> to move up the selected Magnification.
<i>Down</i>	Click <i>Move Down</i> to move down the selected Magnification.

12.4.3 Setting > Image Format


Fig. 32 - Image Format Settings Page

<i>Image Format</i>	<p>JPEG: The extension of JPEG file can get very high compression rate and display very rich and vivid images by removing redundant images and color data. In other words, it can get better image quality with a minimum least disk space. If measurement objects are available, the measurement objects will be burned into the image and the measurement cannot be edited.</p> <p>TIFF: Tag Image File Format (TIFF) is a flexible bitmap format that is mainly used to store images including photos and artistic images.</p>
<i>Measurement Object Save Method</i>	<p>Burn in Mode: The measurement objects are merged into the current image. User could not edit the measurement objects any more. This mode is not reversible.</p> <p>Layered Mode: The measurement objects are saved in different layer with current image data in the target file. User could edit the measurement objects in the target file with some software on the PC. This mode is reversible.</p>

12.4.4 Setting > Video


Fig. 33a - Video Settings Page - Playback


Fig. 33b - Video Settings Page - Video Encode

<i>Fast Forward / Reverse Interval</i>	The time interval of the playback of video files.
<i>Video Encode</i>	You can choose <i>H264</i> or <i>H265</i> encoding. H265 encoding can significantly reduce encoding bandwidth and save storage space under the same encoding quality.

12.4.5 Setting > Storage


Fig. 34 - SD Card Settings Page

<i>Storage Device</i>	<i>SD Card</i> : only SD Card is supported as the storage device.
<i>File System</i>	List the file system format of the current storage device.
<i>Format of the Storage Device</i>	<i>FAT32</i> : The file system of SD Card is FAT32. The maximum video file size of single file is 4GB. <i>exFAT</i> : The file system of SD Card is exFAT. The maximum video file size of single file is 4GB. <i>NTFS</i> : The file system of SD Card is NTFS. The maximum video file size of single file is 4GB. Use PC to format the SD Card and switch between FAT32, exFat and NTFS. <i>Unknown Status</i> : SD Card not detected or the file system is not identified.

12.4.6 Setting > Files


Fig. 35 - Files Settings Page

<i>Image File Name</i>	<i>Auto:</i> The image files will be saved automatically with the specified prefix. <i>Manual:</i> User has to specify the file name before image saving.
<i>Video File Name</i>	<i>Auto:</i> The video file will be saved automatically with the specified prefix. <i>Manual:</i> User has to specify the Video File Name before video recording.

12.4.7 Setting > Language


Fig. 36 - Language Selection Settings Page

<i>English</i>	Set language of the whole software into English
<i>Simplified Chinese</i>	Set language of the whole software into Simplified Chinese
<i>Traditional Chinese</i>	Set language of the whole software into Traditional Chinese
<i>Korean</i>	Set language of the whole software into Korean
<i>Thailand</i>	Set language of the whole software into Thai
<i>French</i>	Set language of the whole software into French
<i>German</i>	Set language of the whole software into German
<i>Japanese</i>	Set language of the whole software into Japanese
<i>Italian</i>	Set language of the whole software into Italian
<i>Russian</i>	Set language of the whole software into Russian

12.4.8 Setting > Miscellaneous



Fig. 37 - Miscellaneous Settings Page

<i>Ruler</i>	Select to display or hide the Ruler in the video window.
<i>Measurement</i>	Select to display the Measurement toolbar in the video window, otherwise, the Measurement toolbar will not be displayed.
<i>Overlay</i>	Select to support saving graphics Overlay information in fusion mode, or not to save graphics Overlay information in fusion mode.
<i>Auto Exposure</i>	The maximum exposure time during auto exposure process could be specified. Setting this item to a lower value could guarantee a faster frame rate during Auto Exposure.
<i>ROI Color</i>	Choose the <i>ROI</i> rectangle line color.
<i>Cursor</i>	Choice of <i>Cursor</i> size based on screen resolution or personal preference.
<i>Camera Parameters Import</i>	Import the <i>Camera Parameters</i> from the SD card to use the previously exported parameters.
<i>Camera Parameters Export</i>	Export the <i>Camera Parameters</i> to the SD card to use the previously imported parameters.
<i>Reset to factory defaults</i>	<i>Reset</i> camera parameters to factory status


13. Micrometric Slide M-005

Micrometric slide, 26x76mm, with 2 scales
(1mm/100div. for biological microscopes / 10mm/100div. for stereomicroscopes)

1 DIV=0.01mm

For biological microscopes calibration

1 DIV=0.1mm

For stereomicroscopes calibration

14. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region.
- The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.


15. Troubleshooting

Review the information in the table below to solve operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark.	Microscope power supply is unplugged	Connect
	Brightness setting is too low	Adjust brightness
	Camera power supply is unplugged	Connect
	Monitor power supply is unplugged	Connect
	Camera is turned off	Turn on camera
	Camera exposure parameters are not correctly set	Set camera with a proper exposure
	Switching objective the microscope is not parfocal	Perform parfocality adjustment again
Field of view is obscured or not evenly illuminated	Revolving nosepiece is not correctly engaged.	Make sure that the revolving nosepiece clicks properly into place.
Dirt or dust is visible in the field of view.	Dirt/dust on the specimen	Clean the specimen
Image looks double	Aperture iris diaphragm is stopped down too far.	Open aperture iris diaphragm.
	The condenser is not well centered or it is in a wrong height	Set the condenser according to Koehler settings.
Visibility is poor. <ul style="list-style-type: none"> • Image is poor. • Contrast is poor. • Details are indistinct. 	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture iris diaphragm is too closed or too open.	Adjust aperture iris diaphragm.
	Dust or dirt on lenses (condenser, objectives and slide)	Clean thoroughly
	For transmitted light observation, the coverglass thickness must not exceed 0.17mm	Use a coverglass with thickness 0.17mm
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position
One side of the image is unfocused	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Slide is mounted not in a flat position (tilted)	Place the specimen in a flat position on the stage
	Poor quality of the glass slide	Use a glass slide with higher quality
II. Mechanical Section:		
Coarse focus knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
Focus is unstable	Tension adjustment ring is too loose	Tighten tension adjustment ring
III. Electrical Section:		
LED doesn't turn on.	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
	Camera exposure time set incorrectly	Adjust exposure time
Light blinks	Power supply not well connected	Check for proper connection
IV. Microphotography:		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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10. Collegare il cavo mini HDMI alla porta mini HDMI © situata sul lato destro del monitor. (Fig. 9)
OPTIKA® S.r.l.