B-510 Series

INSTRUCTION MANUAL

Model
B-510MET
B-510METR

Ver. 1.5 2024

Table of Contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
3.1	B-510MET	4
3.2	B-510MET-R	5
4.	Unpacking	6
5.	Intended use	6
6.	Symbols and conventions	6
7.	Instrument description	7
7.1	B-510MET	7
7.2	B-510MET-R	9
8.	Assembling	11
9.	Summary of reflected light observation procedures	15
10.	Use of the microscope	16
10.1	Light intensity adjustment	16
10.2	Coarse focus tension adjustment	16
10.3	Coarse upper limit lever	16
10.4	Stage	16
10.4.1	B-510MET	16
10.4.2	B-510MET-R	17
10.5	Diopter adjustment	17
10.6	Adjusting the interpupillary distance	17
10.7	Use of eyeshields	17
10.8	Centering of reflected light diaphragms	18
10.8.1	Field diaphragm (FS)	18
10.8.2	Effects of the field diaphragm	18
10.8.3	Aperture diaphragm (AS)	19
10.9	Centering the condenser (B-510MET-R)	19
10.10	Aperture diaphragm	19
10.11	Use of filters	20
10.12	Oblique illumination	20
10.13	Use of polarized light	21
11.	Microfotography	22
11.1	Use of C-mount cameras	22
11.2	Use of reflex cameras	22
12.	Maintenance	23
13.	Troubleshooting	24
	Equipment disposal	26

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

3.1 B-510MET

① Microscope frame

② Eyepieces

③ Objectives

④ Observation head

⑤ Reflected light Illuminator

⑥ Allen wrench

⑦ Polarizer

⑧ Analyzer

⑨ Dust cover

⑩ Connection plate

⑪ Color filters

⑫ Power supply

⑬ Tension adjustment tool

3.2 B-510MET-R

- ① Microscope frame
- ② Eyepieces
- ③ Objectives
- ④ Observation head
- ⑤ Reflected light illuminator
- ⑥ Allen wrenches
- ⑦ Polarizer

- ⑧ Analyzer
- ⑨ Dust cover
- ⑩ Color filters
- ⑪ Connection plate
- ⑫ Condenser
- ⑬ Power supply
- ⑭ Tension adjustment tool

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


ELECTRICAL SHOCK


This symbol indicates a risk of electrical shock.

7. Instrument description


7.1 B-510MET


Opposite side

7.2 B-510MET-R


Opposite side

8. Assembling

1. Insert the reflected light illuminator ① on the frame and lock the locking knob ② with the provided Allen wrench. (Fig. 1)


2. Remove the cover ③ on the back side of the illuminator. (Fig. 2)

3. Insert the LED housing ④ and lock the locking screws ⑤. (Fig. 3-4)

4. Connect the illuminator cable ⑥ to the connector ⑦ placed on the back side of the frame. (Fig. 5)

5. Insert the optical head above the epi-illuminator and tighten the locking screw ⑧ with the provided Allen wrench. (Fig. 6)

 - **Hold the head with one hand during the locking in order to avoid that the head falls.**

(Only B-510MET-R)

- Condenser is pre-installed in the factory. To remove the condenser use an Allen wrench 1,5 mm diam and operate on the locking screw placed on the right side of the condenser holder.

6. Insert both eyepieces into the tubes of the optical head. (Fig. 7)

7. Screw each objective into the thread of the nosepiece, clockwise with increasing magnification. (Fig. 8)

8. Insert the polarizer "PO" with the inscription facing the user (before removing the empty slider). (Fig. 9)


9. Insert analyzer "AN" with the inscription facing up (first removing the white covering plates on the illuminator). (Fig.10-11)

- It is possible to combine the analyzer and the polarizer together (using the connection plate) so that the insertion of the two filters is simultaneous.


10. Screw the plate fixing screws to the holes on the end of the filters. (Fig. 12)

-
11. Insert the power supply jack in the connector placed on the back side of the frame. (Fig. 13)

9. Summary of reflected light observation procedures

10. Use of the microscope

10.1 Light intensity adjustment

Operate on the light intensity adjustment knob to turn ON / OFF the microscope and to increase / decrease the illumination voltage ①. (Fig. 14)

- **For B-510MET-R only.** There is a three-position switch on the rear of the stand: position “I” turns on the transmitted light, position “II” turns on the reflected light and position “O” turns off the microscope.

Fig. 14

10.2 Coarse focus tension adjustment

- **Adjust the tension using the provided tool.**

The tension of the coarse focusing knob is factory preset.

1. To modify the tension according to personal's needs, rotate the ring ②. (Fig. 15)
- Clockwise rotation increases the tension.
 - If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

Fig. 15

10.3 Coarse upper limit lever

The upper limit knob has two functions: prevent the contact between slide and objective and acts as “focus memory”.

1. After focusing the specimen, pull the lever ③ toward the front of the microscope and lock it. (Fig. 16)
- In this way the focus upper limit is set.
2. Now one can lower the stage with coarse focus knob, replace the specimen and raise again the stage up to the upper limit: specimen will be in approximate focus and will need a fine adjustment to get the proper focus.
- **Fine focus movement is not affected by the coarse focus lock.**
 - **To unlock, move the lever in the opposite direction to the one used for the locking.**

Fig. 16

10.4 Stage

10.4.1 B-510MET

Stage accepts metallurgical samples with thickness max 35 mm, placed on a standard support 26 x 76 mm.

- **Open the spring arm of the slide holder ④ and place frontally the slides on the stage. (Fig. 17)**
- **Gently release the spring arm of the slide holder.**
- **A sudden release of the spring arm could cause the falling of the slide.**


Fig. 17

10.4.2 B-510MET-R

Stage accepts metallurgical samples with thickness max 35 mm, placed on a standard support 26 x 76 mm, or standard slides 26x76 mm.(Fig. 18)

- **Open the spring arm of the slide holder ⑤ and place frontally the slides on the stage. (Fig. 18)**
- **Gently release the spring arm of the slide holder.**
- **A sudden release of the spring arm could cause the falling of the slide.**

10.5 Diopter adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
 2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the diopter adjustment ring ⑥ to compensate. (Fig. 19)
- **The adjustment range is ± 5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's diopter correction.**

10.6 Adjusting the interpupillary distance

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

- **The graduation on the interpupillary distance indicator ⑦, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator's eyes. (Fig. 34)**

The range of the interpupillary distance is 48-75 mm.

10.7 Use of eyeshields

- **Use with eyeglasses**
Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses. (Fig. 21)

- **Use without eyeglasses**

Raise eye shields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation. (Fig. 22)

10.8 Centering of reflected light diaphragms



10.8.1 Field diaphragm (FS)

1. Move the selector ① on the reflected light illuminator in the position corresponding to letter "R" (only B-510MET-R). (Fig. 23)
2. Place the specimen on the stage, insert 10x objective into the light path and focus.
3. Completely pull out the field diaphragm selector ② to fully close the diaphragm. (Fig. 24)
4. Using the provided Allen wrench use the two centering screws ③ to bring the bright spot in the center of the field of view.
5. Gradually open the diaphragm. The condenser is centered when the diaphragm image is symmetrical to the field of view. (Fig. 25)
6. In normal use, open the diaphragm until it circumscribes the field of view.



10.8.2 Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image. Set the diaphragm according to the objective in use until it circumscribes the field of view, in order to eliminate unnecessary light to eyepieces.

10.8.3 Aperture diaphragm (AS)

1. Move the selector ① on the reflected light illuminator in the position corresponding to letter "R" (**only B-510MET-R**). (Fig. 23)
2. Completely pull out the aperture diaphragm selector ④ to fully close the diaphragm. (Fig. 26)
3. Remove one eyepiece.
4. While looking into the empty eyepiece sleeve, use the provided Allen wrench ⑤ into the centering screws ⑤ to bring the bright spot in the center of the field of view. (Fig. 26)
5. The illuminator is centered when the diaphragm image is symmetrical to the field of view.
 - The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
 - With low contrast specimens move the diaphragm lever to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's ring in order to obtain an image like the one in Fig. 27.

10.9 Centering the condenser (B-510MET-R)

1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Insert the front lens of the swing-out condenser ①. (Fig. 28)
3. Rotate the field diaphragm ring ② in counterclockwise direction, to fully close the diaphragm.
4. Rotate the condenser height adjustment knob ③ to focus the edges of the diaphragm.
5. Rotate the two centering screws ④ to bring the bright spot in the center of the field of view.
6. Gradually open the diaphragm. The condenser is centered when the diaphragm image is symmetrical to the field of view.
7. In normal use, open the diaphragm until it circumscribes the field of view.

10.10 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
- With low contrast specimens set the numerical aperture value ① (printed on the condenser ring) to about 70%-80% of the objective's N.A. (Fig. 29) If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's ring in order to obtain an image like the one in fig. 27.

Example: with objective PLAN 40x / 0,65 set the scale to 0.65 x 0.8 = 0,52

10.11 Use of filters

1. Insert the desired filter for the observation in the slot placed in the left side of the illuminator (Fig. 30).
2. First click stop position leaves filter not inserted, while second position inserts filter into the light path.

Filter	Application
Blue	Converts the source color temperature to daylight
Green	Increases the contrast in black and white photography
Yellow	Contrast filter for the observation of wafers and semiconductors
Frosted white	Reduces lighting irregularities by creating a homogeneous field of view

Fig. 30

10.12 Oblique illumination

1. After focusing the specimen with the desired objective, gently move the oblique illumination sliders (Fig. 31).
 2. Stop sliding when a shadow begins to obscure the specimen.
- In this way an oblique illumination effect can be achieved, that highlights reliefs and surface defects of your specimen.


Fig. 31

10.13 Use of polarized light

1. Move the polarizer "PO" into the light path until the slider clicks into its position. (Fig. 32)

2. Move the rotatable analyzer "AN" into the light path until the slider clicks into its position. (Fig. 33)


3. Put a flat mirror on the stage and get into approximate focus the mirror surface.
4. While observing into the eyepieces, rotate the scale of the analyzer ① until you achieve the darkest image. (Fig. 34)
 - Now the light extinction (so-called "crossed Nicol position") is achieved.
5. Put a specimen on the stage, focus the specimen and begin observation.

11. Microphotography

11.1 Use of C-mount cameras

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 35)

2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 36)

11.2 Use of reflex cameras

1. Insert the Reflex adapter ② into the relay tube to the microscope ①.
2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
3. Connect the Reflex camera ④ to the "T2" ring just installed (Fig. 37).
4. Mount the other end of the relay tube ② into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 35)
 - "T2" ring is not provided along with the microscope, but is commercially available.
 - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
 - To calculate the magnification of the camera: objective magnification * camera magnification * lens magnification.
 - **If using an SLR camera, mirror movement may cause the camera to vibrate.**
 - **We suggest lifting the mirror, using long exposure times and a remote cord.**

12. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

13. Troubleshooting


Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark	Power supply is unplugged	Connect
	Brightness is too low	Set brightness to a proper level
	Color filter or oblique illumination slider are in a wrong position	Move to a click stop
Field of view is obscured or not evenly illuminated	Revolving nosepiece is not correctly engaged	Make sure that the revolving nosepiece clicks properly into place
	Oblique illumination slider is in a wrong position	Move the slider to a fully open position
Dirt or dust is visible in the field of view	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture iris diaphragm is stopped down too far	Open aperture iris diaphragm
	The field diaphragm (both for reflected light and transmitted light) is not well centered	Set the diaphragm according to Koehler settings
Visibility is poor: <ul style="list-style-type: none"> • Image is not good; • Contrast is poor; • Details are indistinct; • Image glares 	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture iris diaphragm is too closed or too open	Adjust aperture iris diaphragm
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly
	For reflected light observation, the specimen must not have any coverglass	Use specimens without coverglass
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position
One side of the image is unfocused	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Slide is mounted not in a flat position (tilted)	Place the specimen in a flat position on the stage
II. Mechanical Section:		
Coarse focus knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
Focus is unstable	Tension adjustment ring is too loose	Tighten tension adjustment ring
III. Electrical Section		
LED doesn't turn on	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection
IV. Observation tube		
Field of view of one eye does not match that of the other	Interpupillary distance is incorrect	Adjust interpupillary distance
	Incorrect diopter adjustment	Adjust diopter
	Your view is not accustomed to microscope observation	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope

V. Microphotography		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.