B-510 Series

INSTRUCTION MANUAL

Model
B-510POL
B-510POL-I

Ver. 2.3 2023

Table of Contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
3.1	B-510POL	4
3.2	B-510POL-I	5
4.	Unpacking	6
5.	Intended use	6
6.	Symbols and conventions	6
7.	Instrument description	7
7.1	B-510POL	7
7.2	B-510POL-I	8
8.	Assembling	9
8.1	B-510POL	9
8.2	B-510POL-I	11
9.	Transmitted light brightfield observation procedures	14
9.1	B-510POL	14
9.2	B-510POL-I	15
10.	Use of the microscope (transmitted light brightfield)	16
10.1	Light intensity adjustment	16
10.2	Coarse focus tension adjustment	16
10.3	Coarse upper limit lever	16
10.4	Stage	17
10.5	Diopter adjustment	17
10.6	Adjusting the interpupillary distance	17
10.7	Use of eye shields	17
10.8	Centering the condenser	18
10.9	Effects of the field diaphragm	18
10.10	Aperture diaphragm	19
11.	Transmitted light polarized observation procedures	20
11.1	B-510POL	20
11.2	B-510POL-I	21
12.	Reflected light polarized observation procedures	22
12.1	B-510POL-I	22
13.	Use of the microscope in polarized light	23
13.1	Centering the rotatable stage	23
13.2	Nosepiece centering	24
13.3	Checking the extinction of the light	24
13.3.1	B-510POL	24
13.3.2	B-510POL-I	25
13.4	Centering of reflected light diaphragms	26
13.4.1	Field diaphragm (FS)	26
13.4.2	Aperture diaphragm (AS)	26
13.5	Use of tint plates	27
13.6	Use of the Bertrand lens	27
13.7	Use of diffusing filter	27
14.	Microphotography	28
14.1	Use of C-mount cameras	28
14.2	Use of reflex cameras	28
15.	Maintenance	29
16.	Troubleshooting	30
	Equipment disposal	32

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

3.1 B-510POL

- ① Microscope frame
- ② Eyepieces
- ③ Objectives
- ④ Observation head
- ⑤ Bertrand lens
- ⑥ Analyzer
- ⑦ Nosepiece centering screws

- ⑧ Dust cover
- ⑨ Tint plates
- ⑩ Tension adjustment tool
- ⑪ Allen wrench
- ⑫ Dummy slider
- ⑬ Power supply

3.2 B-510POL-I

- ① Microscope frame
- ② Eyepieces
- ③ Objectives
- ④ Observation head
- ⑤ Bertrand lens
- ⑥ Analyzer
- ⑦ Nosepiece centering screws
- ⑧ Dust cover
- ⑨ Tint plates
- ⑩ Tension adjustment tool
- ⑪ Allen wrench
- ⑫ Reflected light illuminator
- ⑬ Reflected light polarizer
- ⑭ Power supply
- ⑮ Diffusing filter

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


ELECTRICAL SHOCK


This symbol indicates a risk of electrical shock.

7. Instrument description

7.1 B-510POL

7.2 B-510POL-I

8. Assembling

8.1 B-510POL

1. Insert the Bertrand lens ① on the frame and lock the locking screw ② with the provided Allen wrench. (Fig. 1)

2. Insert the optical head on the Bertrand lens and lock the locking screw with the provided Allen wrench. (Fig. 2)

- **Hold the head with one hand during the locking in order to avoid that the head falls.**

3. Insert both eyepieces into the tubes of the optical head. (Fig. 3)


- **One of the two eyepieces is equipped with a cross hair for centering the entire optical system. It is advisable to insert the eyepiece with cross hair in the right eyepiece sleeve.**
- Condenser is pre-installed in the factory. To remove the condenser use an Allen wrench 1,5mm diam and operate on the locking screw placed on the right side of the condenser holder.

4. Screw each objective into the thread of the nosepiece, clockwise with increasing magnification. (Fig. 4)

5. Remove the dummy slider from Bertrand lens ③ and insert the analyzer ④. (Fig. 5 - 6)

6. Insert the power supply jack in the connector placed at the rear side of the microscope. (Fig. 7)


8.2 B-510POL-I

1. Insert the reflected light illuminator ① on the frame and lock the locking screw ② with the provided Allen wrench. (Fig. 8)

2. Connect the illuminator cable jack ③ to the socket ④ placed in the back side of the frame. (Fig. 9)

3. Insert the Bertrand lens ⑤ on the reflected light illuminator and lock the locking screw ⑥ with the provided Allen wrench. (Fig. 9 - 10).

4. Insert the optical head on the Bertrand lens and lock the locking screw with the provided Allen wrench. (Fig. 11)
- **Hold the head with one hand during the locking in order to avoid that the head falls.**


5. Insert both eyepieces into the tubes of the optical head. (Fig. 12)
 - **One of the two eyepieces is equipped with a cross hair for centering the entire optical system. It is advisable to insert the eyepiece with cross hair in the right eyepiece sleeve.**
 - Condenser is pre-installed in the factory. To remove the condenser use an Allen wrench 1,5mm diam and operate on the locking screw placed on the right side of the condenser holder.


6. Screw each objective into the thread of the nosepiece, clockwise with increasing magnification. (Fig. 13)

7. Insert reflected light polarizer ⑦. (Fig. 14)

8. Remove the dummy slider from Bertrand lens ⑧ and insert the analyzer ⑨. (Fig. 15 - 16)


9. Insert the power supply jack in the connector placed at the rear side of the microscope. (Fig. 17).


9. Transmitted light brightfield observation procedures

9.1 B-510POL

9.2 B-510POL-I


10. Use of the microscope (transmitted light brightfield)

10.1 Light intensity adjustment

Operate on the light intensity adjustment knob to turn ON / OFF the microscope and to increase / decrease the illumination voltage ①. (Fig. 18)

- **For B-510POL-I model only.** There is a three-position switch on the rear of the stand: position “I” turns on the transmitted light, position “II” turns on the reflected light and position “O” turns off the microscope.

10.2 Coarse focus tension adjustment

- **Adjust the tension using the provided tool.**


The coarse knob tension is pre-set in the factory. To modify the tension according to personal's needs, rotate the ring ② using the provided tool. (Fig. 19) Clockwise rotation increases the tension. If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

10.3 Coarse upper limit lever


The upper limit knob has two functions: prevent the contact between slide and objective and acts as “focus memory”.

1. After focusing the specimen, pull the lever ③ toward the front of the microscope and lock it. (Fig. 20)
 - In this way the focus upper limit is set.
2. Now one can lower the stage with coarse focus knob, replace the specimen and raise again the stage up to the upper limit: specimen will be in approximate focus and will need a fine adjustment to get the proper focus.
 - **Fine focus movement is not affected by the coarse focus lock.**
 - **To unlock, move the lever in the opposite direction to the one used for the locking.**

10.4 Stage

1. The rotating stage accepts specimens on slide (B-510POL) or opaque specimens (B-510POL-I).
2. It is possible to lock the specimen once placed on the stage using the stage clips ①. (Fig. 21)
3. After loosening the locking screw ②, the stage can be rotated by 360°.

10.5 Diopter adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
 2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the diopter adjustment ring ③ to compensate. (Fig. 22)
- **The adjustment range is ± 5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's diopter correction.**

10.6 Adjusting the interpupillary distance

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

- **The graduation on the interpupillary distance indicator ④, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator's eyes. (Fig. 23)**

The range of the interpupillary distance is 48-75mm.

10.7 Use of eye shields

- **Use with eyeglasses**
Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses. (Fig. 24)

- **Use without eyeglasses**

Raise eye shields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation. (Fig. 25)

10.8 Centering the condenser


1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Insert the front lens of the swing-out condenser ①. (Fig. 26)
3. Rotate the field diaphragm ring ② in counterclockwise direction, to fully close the diaphragm.
4. Rotate the condenser height adjustment knob ③ to focus the edges of the diaphragm.
5. Rotate the two centering screws ④ to bring the bright spot in the center of the field of view.
6. Gradually open the diaphragm. The condenser is centered when the diaphragm image is symmetrical to the field of view.
7. In normal use, open the diaphragm until it circumscribes the field of view.

10.9 Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image.



Set the diaphragm according to the objective in use until it circumscribes the field of view, in order to eliminate unnecessary light to eyepieces. (Fig. 27)

10.10 Aperture diaphragm


- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
- With low contrast specimens set the numerical aperture value ① (printed on the condenser ring) to about 70%-80% of the objective's N.A. (Fig. 28) If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's ring in order to obtain an image like the one in fig. 29.

Example: with objective PLAN 40x / 0,65 set the scale to 0.65 x 0.8 = 0,52


11. Transmitted light polarized observation procedures

11.1 B-510POL


11.2 B-510POL-I

12. Reflected light polarized observation procedures

12.1 B-510POL-I

13. Use of the microscope in polarized light


- The system allows observation in Orthoscopy (crossed Nicol) or in Conoscopy (crossed Nicol with the use of the Bertrand lens).
- For optimal performance in polarized light microscopy, accurate optical adjustments are essential before beginning the observation.

13.1 Centering the rotatable stage


1. Loosen the stage rotation locking screw ① and rotate the stage until the graduated scale of the stage ② and the vernier scale are aligned ③ in the "0" position. (Fig. 30)
- This operation serves to ensure a standard reference position for centering the rotating stage.

2. Focus on a recognizable detail ④ in the field of view by placing it at the center of the eyepiece cross hair. (Fig. 31)
3. By turning the stage, the detail in focus will draw a circle ⑤. (Fig. 31)
4. Move again the stage on the "0" position and lock the locking screw ①. (Fig. 30)

5. Using the centering screws of the stage ⑥ move the detail in diametrically opposite direction to the described circle. The movement must be about half the diameter of the circle drawn. (Fig. 32)
6. Manually move the specimen and bring it back to the center of the cross hair. Loosen the stage locking screw again and rotate the stage again.
7. If the centering has been carried out correctly, by turning the stage, the image of the focused detail does not move with respect to the center of the cross hair. If not, repeat the operations described by (1) to (6) until the perfect coincidence of the center of rotation of the stage with the center of the cross hair for which the specimen remains in the center of the cross hair when turning the stage.

13.2 Nosepiece centering

1. After centering the stage with the 10X objective, return the recognizable detail used for centering to the centre of the cross hair.
2. Rotate the revolver by inserting all the other objectives in the optical path and check that the detail is always in the centre of the cross hair.
3. If this is not the case, use the centering screws on the ① nosepiece to ensure that all the objectives are perfectly centered in relation to the optical axis. (Fig. 33)

13.3 Checking the extinction of the light

13.3.1 B-510POL


1. Remove the specimen from the light path and insert 10x.
 2. Loosen the locking screw of the polarizer ① and check that the scale is in the "0" position ②. (Fig. 34)
 3. Insert rotatable analyzer in the light path, loosen the locking screw ③ and put the vibration scale on 0° ④, then lock the locking screw ③. (Fig. 35)
 4. Rotate the polarizer scale ② to obtain total extinction (total dark in the eyepieces). Then tighten the screw ①. (Fig. 34)
- It may happen that the scale of the polarizer is not perfectly aligned on the reference mark but is shifted by one or two notches. This is not a defect but it is due to the mechanical alignment of the polarizers during assembly.


13.3.2 B-510POL-I

Extinction in reflected light

1. Move the selector ① on the reflected light illuminator in the fully inserted position, corresponding to letter "R". (Fig. 36)

2. Insert the reflected light polarizer ②. (Fig. 37)
3. Place on the stage a flat mirror and insert 10x objective.

4. Insert rotatable analyzer in the light path, loosen the locking screw ③ and put the vibration scale on 0° ④, then lock the locking screw ③. (Fig. 38)
- It may happen that the scale of the polarizer is not perfectly aligned on the reference mark but is shifted by one or two notches. This is not a defect but it is due to the mechanical alignment of the polarizers during assembly.



Extinction in transmitted light

1. Move the selector ① on the reflected light illuminator in the fully removed position, corresponding to letter "T". (Fig. 36)
2. Repeat procedure described in steps 1. to 4. for the B-510POL.

13.4 Centering of reflected light diaphragms


13.4.1 Field diaphragm (FS)

1. Move the selector ① on the reflected light illuminator in the fully inserted position, corresponding to letter "R". (Fig. 36)
2. Place a specimen on the stage, insert 10x objective and focus the specimen.
3. Rotate the knurled ring of the field diaphragm ① in the direction shown by the arrow to fully close the diaphragm. (Fig. 39)
4. Using the provided Allen wrench use the two centering screws ② to bring the bright spot in the center of the field of view.
5. Gradually open the diaphragm. The illuminator is centered when the diaphragm image is symmetrical to the field of view. (Fig. 40)
6. In normal use, open the diaphragm until it circumscribes the field of view.

13.4.2 Aperture diaphragm (AS)

1. Rotate the knurled ring of the aperture diaphragm ③ in the direction shown by the arrow to fully close the diaphragm.
2. Remove one eyepiece.
3. While observing in the empty eyepiece holder, use the provided Allen wrench and use the centering screws ④ to bring the bright spot in the center of the field of view. (Fig. 41)
4. The illuminator is centered when the diaphragm image is symmetrical to the field of view.
 - The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
 - With low contrast specimens rotate the ring on the illuminator to set the numerical aperture value to about 70%-80% of the objective's N.A. If necessary, remove one eyepiece and, looking into empty sleeve, adjust the ring of the illuminator in order to obtain an image like the one in fig. 29.

13.5 Use of tint plates

Three tint plates are supplied with the microscope:



- Lambda (λ) plate (1st order Red)
 - Lambda/4 ($\lambda/4$) plate
 - "Quartz wedge" (Q) plate
1. Insert in the right slot of the Bertrand lens ① one of the tint plates ②. (Fig. 42)
 2. In polarized light, inserting one of the plates will have chromatic effects on the specimen.
 - Using the λ plate (also called 1st order Red) the specimen will take a magenta tinge.
 - Using the $\lambda/4$ the specimen will take on a color tending to pale yellow.
 - Using the Q plate the specimen will show a series of colored bands that will fade as the plate is inserted.

13.6 Use of the Bertrand lens

Bertrand lens allows observation in Orthoscopy and Conoscopy. In the disengaged position ("O") the lens allows observation in Orthoscopy, while in the engaged position ("B") it is possible to make observations in Conoscopy.

1. Rotate the upper knurled ring of the Bertrand lens ① to engage the "B" position. (Fig. 43)
 2. Using one objective from 20x to 60x, focus the conoscopic image using the focus ring ②.
 3. If the conoscopic image is not perfectly centered with respect to the optical axis, center the image using the centering screws.
- By turning the stage you will see black fringes that will appear and disappear depending on the rotation of the stage. These fringes are the crystallization axes of that specific crystal. (Fig. 44)

13.7 Use of diffusing filter

Depending on the type of sample to be observed, it may be useful to remove or insert the diffusing filter ④ that is present in the back of the illuminator.


1. Insert the slide into the illuminator until it reaches the end of its stroke to insert the diffuser filter in the optical path. (Fig. 45)
2. Pull out one click (until the first "click") the slide to remove the filter from the optical path, but always leaving the slide in place.
3. If you intend to completely remove the slide from the illuminator, remove it completely from its housing.

14. Microphotography

14.1 Use of C-mount cameras

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 46)

2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 47)

14.2 Use of reflex cameras

1. Insert the Reflex adapter ① into the relay tube to the microscope ②.
 2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
 3. Connect the Reflex camera ④ to the "T2" just installed. (Fig. 48)
 4. Mount the other end of the relay tube ② into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 46)
- "T2" ring is not provided with the microscope, but is commercially available.
 - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
 - To calculate the magnification of the camera: objective magnification * camera magnification * lens magnification.
 - **If using an SLR camera, mirror movement may cause the camera to vibrate.**
 - **We suggest lifting the mirror, using long exposure times and a remote cord.**

15. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

16. Troubleshooting


Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark.	Power supply is unplugged.	Connect
	Brightness is too low	Set brightness to a proper level
	Bertrand lens is in	Remove Bertrand lens from light path
	You are in a position of extinction	Disengage analyzer from the light path
Field of view is obscured or not evenly illuminated	Revolving nosepiece is not correctly engaged.	Make sure that the revolving nosepiece clicks properly into place.
	Tint plate, a filter or Bertrand lens are in an intermediate position.	Move to a click stop
Dirt or dust is visible in the field of view.	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture iris diaphragm is stopped down too far.	Open aperture iris diaphragm.
	The condenser is not well centered or it is in a wrong height	Set the condenser according to Kohler settings.
Visibility is poor. <ul style="list-style-type: none"> • Image is not sharp. • Contrast is poor. • Details are indistinct. 	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture iris diaphragm is too closed or too open.	Adjust aperture iris diaphragm.
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly.
	For transmitted light observation, the coverglass thickness must not exceed 0.17mm	Use a coverglass with thickness 0.17mm
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position.
One side of the image is unfocused	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Slide is mounted not in a flat position (tilted)	Place the specimen in a flat position on the stage
	Poor quality of the glass slide	Use a glass slide with higher quality
Conoscopic image cannot be seen	Swing-out lens of the condenser is not in the light path	Insert the lens in the light path
	Bertrand lens is not in the light path.	Insert the lens in the light path.
Total extinction cannot be obtained	Analyzer is not in the light path	Insert analyzer in the light path.
II. Mechanical Section:		
Coarse focus knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
Focus is unstable	Tension adjustment ring is too loose	Tighten tension adjustment ring
III. Electrical Section:		
LED doesn't turn on.	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection

IV. Observation tube:		
Field of view of one eye does not match that of the other.	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope.
V. Microphotography:		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.