HDMI Cameras – H Series

INSTRUCTION MANUAL

Model
C-HP
C-HPSC
C-HE
C-HESC

Ver. 3.0 2019

Table of Contents

1. Warning	3
2. Symbols and conventions	3
3. Safety Information	3
4. Intended use	3
5. C-HP Basic features	4
5.1 Quick Instructions for C-HP camera	5
5.2 USB Video Mode	5
5.3 HDMI Mode	5
6. C-HP USER INTERFACE (UI) AND ITS FUNCTIONS	6
6.1 The Camera Control Panel	7
6.2 Icons and Functions of the Camera Control Toolbar	8
6.3 The Measurement Toolbar	11
7. C-HE BASIC FEATURES	13
7.1 Quick Instructions for C-HE camera	14
7.1 HDMI Mode	14
8. C-HE USER INTERFACE (UI) AND ITS FUNCTIONS	15
8.1 The Camera Control Panel	16
8.2 Icons and Functions of the Camera Control Toolbar	17
9. C-HPSC / C-HESC	18
9.1 Overview	18
9.2 Components	18
9.3 Assembling	19
10. MICROMETRIC SLIDE M-005	20
Equipment disposal	21

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

3. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

4. Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

5 C-HP Basic features

C-HP is a multiple interfaces (HDMI + USB2.0 + SD card) CMOS camera. HDMI + USB2.0 are used as the data transfer interface to HDMI display or computer.

For HDMI output, the Camera Control Panel+ Measurement Toolbar and Camera Control Toolbar are overlaid on the HDMI screen when the mouse move to the related region.

In this case, the USB mouse can be used to set the camera, browse and compare the captured image, play the video and perform the measurement.

For USB Video camera mode, plug in the micro USB host cable to the camera USB video port and computer USB port, then the video stream can be transfer to computer with the advanced software ProView. With ProView, you can control the camera, process the video and image as OPTIKA's other USB series camera.


The C-HP's basic characteristics are as follows:

For HDMI output:

- All in 1(HDMI + USB + SD card) C-mount camera with Sony high sensitivity CMOS sensor;
- Simultaneous HDMI & USB output;
- Built-in mouse control;
- Built-in image capture & video record to SD card;
- Built-in camera control panel, including exposure(manual/auto)/gain, white balance(lockable), color adjustment, sharpness control;
- Built-in video and image measurement;
- Built-in toolbar including zoom, mirror, comparison, freeze, cross, browser functions;
- Built-in image & video browsing, display & play;
- Real time clock (RTC)

For USB Video output:

- Ultra-Fine color engine with perfect color reproduction capability(USB);
- With advanced video & image processing application ProView, which including professional image processing such as 2D measurement, HDR, image stitching, EDF(Extended Depth of Focus), image segmentation & count, image stacking, color composite and denoising (USB);
- Support standard UVC protocol for Windows (USB);

C-HP - SW UI for Mouse Control

5.1 Quick Instructions for C-HP camera

Figure 1 - The Layout of C-HP Rear Cover

- | | |
|-----------------------|---------------------|
| ① HDMI connector | ④ SD card slot |
| ② USB mouse connector | ⑤ Power supply plug |
| ③ USB PC connector | ⑥ ON-OFF button |

Before starting the camera, please connect the standard C-HP C-mount camera to the camera adapter and insert it into the microscope's photo tube which will relay microscope object's middle image to the camera sensor.

5.2 USB Video Mode

1. Plug 12V/ 1A power cable into DC 12V slot to supply power for the camera. The LED Indicator will turn into red;
2. Press ON/OFF Button to start the camera and the LED indicator will turn into blue;
3. Plug the USB cable which comes with the camera into USB Video to connect the video to the computer;
4. Open ProView software, start C-HP by clicking the camera model name listed in Camera List. For more details please refer to the ProView help manual;

5.3 HDMI Mode

1. Plug the HDMI cable into the HDMI Port to connect the C-HP camera to HDMI display;
2. Plug a USB mouse into USB Mouse to get control of the camera by using built-in software ProView;
3. Plug 12V/ 1A power adapter into Power Interface to supply power for the camera. The LED Indicator will turn into red;
4. Insert SD card into SD Card Slot for saving captured images and recorded videos;
5. Press ON/ OFF Button to start the camera. The LED Indicator will turn into blue;
6. Move mouse cursor to the left side of the video window, a Camera Control Panel will appear. It includes Manual/ Automatic Exposure, White Balance, Sharpness and other functions, please refer to 2.1 for details;
7. Move mouse cursor to the bottom of the video window and a Camera Control Toolbar will appear. Operations like Zoom In, Zoom Out, Flip, Freeze, Cross Line, Comparison and etc. can be realized. Please refer to 2.2 for details;
8. Move mouse cursor to the upper side of the video window, a Measurement Toolbar with calibration and other measurement tools will appear, please refer to 2.3 for details. The measurement data can be output with *.CSV format.

6 C-HP User Interface (UI) and Its Functions

The C-HP UI shown in Fig. 2 includes a *Camera Control Panel* on the left side of the video window, a *Measurement Toolbar* on the upper side of the video window, a *Camera Control Toolbar* on the bottom of the video window.



Figure 2 - The C-HP Camera Control UI

Notes

1. When users move mouse cursor to the left side of the video window, the *Camera Control Panel* will pop up automatically.
2. When users move mouse cursor to the bottom of the video window, the *Camera Control Toolbar* will pop up automatically;
3. When user moves mouse cursor to the bottom of the video window, the *Camera Control Toolbar* will pop up automatically.
4. Move the mouse cursor to the upper side of the video window: a *Measurement Toolbar* will pop up for the calibration and measurement operations. When user left-clicks the *Float/Fixed* button  on the Measurement Toolbar, the Measurement Toolbar will be fixed. In this case the Camera Control Panel will not pop up automatically even if users move mouse cursor to left side of the video windows. Only when user left-clicks the  button on the Measurement Toolbar to exit from measuring procedure will they be able to do other operations on the *Camera Control Panel*, or *Synthesis Camera Control Toolbar*. During the measuring process, when a specific measuring object is selected an *Object Location & Attributes Control Bar*  will appear for changing location and properties of the selected objects.

6.1 The Camera Control Panel

The Camera Control Panel controls the camera to achieve the best image quality according to the specific applications. It will pop up automatically when mouse cursor is moved to the left side of the video window (in measurement status, the Camera Control Panel will not pop up. Only when measurement process is terminated will the Camera Control Panel pop up by moving mouse cursor to the left side of the video window). Left-clicking button  to achieve *Display/ Auto Hide* switch of the Camera Control Panel.

Camera Control Panel	Function	Description	
	Snap	Capture image from the current video window	
	Record	Record video from the current video window	
	Auto Exposure	When Auto Exposure is checked, the system will automatically adjust exposure time according to the value of exposure compensation	
	Target	Available when Auto Exposure is checked. Slide to left or right to adjust Target according to the current video brightness to achieve proper brightness value	
	Exposure Time	Available when Auto Exposure is unchecked. Slide to left or right to reduce or increase exposure time, adjusting brightness of the video	
	Gain	Adjust Gain to reduce or increase brightness of video. The Noise will be reduced or increased accordingly	
	Red	Slide to left or right to decrease or increase the proportion of Red in RGB on video	
	Green	Green is base for reference and cannot be adjusted	
	Blue	Slide to left or right to decrease or increase the proportion of Blue in RGB on the video	
	White Balance	Auto:	White Balance adjustment according to the window video
		Manual:	Slide the Red or Blue to manually set the video White Balance
		ROI:	Set the White Balance according to the ROI. The ROI can be resized and moved
	Sharpness	Adjust Sharpness level of the video window	
	Saturation	Adjust Saturation level of the video window	
	Gamma	Adjust Gamma level of the video. Slide to the right side to increase gamma and to the left to decrease gamma	
	Contrast	Adjust Contrast level of the video. Slide to the right side to increase contrast and to the left to decrease contrast	
	DC	For DC illumination, there will be no fluctuation in light source so no need for compensating light flickering	
	AC(50HZ)	Check AC(50HZ) to eliminate flickering “strap” caused by 50Hz illumination	
	AC(60HZ)	Check AC(60HZ) to eliminate flickering “strap” caused by 60Hz illumination	
	Default	Set all the settings in the Camera Control Panel to default values	

6.2 Icons and Functions of the Camera Control Toolbar

Icon	Function	Icon	Function
	Zoom In the Video Window		Zoom Out the Video Window
	Reset Zoom to Original 1X		Horizontal Flip
	Vertical Flip		Video Freeze
	Compare Image with the current video		Display Cross Line
	Browse Images and Videos in the SD Card		Settings
	Check the Version of ProView		

The  Setting function is relatively more complicated than the other functions. Here are more info about it:


Figure 3 - Measurement Settings Page

Global: Used for setting digits behind the decimal point for measurement results;

Calibration--> Line Width: Used for defining width of the lines for calibration;

Color: Used for defining color of the lines for calibration;

EndPoint: Type: Used for defining shape of the endpoints of lines for calibration.

Null: means no endpoints, *Rectangle* means rectangle type of endpoints. It makes alignment more easily;

Point, Angle, Line, Horizontal Line, Vertical Line, Rectangle, Circle, Ellipse, Annulus, Two Circles, Polygon, Curve:

Left-click the  beside the measuring patterns mentioned above will unfold the corresponding attribute settings to set the individual property of the measuring objects.


Figure 4 - Magnification Calibration Management Settings Page

- Name:** Names such as 10X, 40X, 100X are based on magnification of the microscopes. For continuous zoom microscopes, ensure that the selected magnification coincides with the scale alignment line on the microscope zoom knob;
- Resolution:** Pixels per meter. Devices like microscopes have high resolution value;
- Clear All:** Click the Clear All button will clear the calibrated magnifications and resolutions;
- Delete:** Click Delete to delete the selected item for specific resolution;


Figure 5 - Setting of SD Card Setting Page

- Current File System:** The maximum file FAT32 can store is of 4G Bytes; for EXFAT, it's 2048G Bytes. Suggest converting FAT32 file into EXFAT format on a PC; Unknown Status: SD card not detected or the file system is not identified.

Figure 6 - Setting of Language Selection Setting Page

- English:* Set language of the whole software into English;
- Simplified Chinese:* Set language of the whole software into Simplified Chinese;
- Traditional Chinese:* Set language of the whole software into Traditional Chinese;
- Korean:* Set language of the whole software into Korean;
- Thailand:* Set language of the whole software into Thailand;


Figure 7 - The Real Time Clock Setting

- Year:* Current Year
- Month:* Current Month
- Day:* Current Day
- Hour:* Current Hour
- Minute:* Current Minute
- Second:* Current Second

6.3 The Measurement Toolbar

The Measurement Toolbar will pop up when moving mouse cursor to any place near the upper side of the video window.

Figure 8 - The Measurement Toolbar Button on the Upper Side of the Video window

Icon	Function
	Float/ Fix switch of the Measurement Toolbar
<input checked="" type="checkbox"/> Visible	Define measuring object in Show up/ Hide mode
Nanometer (nm)	Select the desired Measurement Unit
4X	Choose the same Magnification as the microscope to ensure accuracy of measurement result when measurement unit is not in Pixel units
	Object Select
	Point
	Angle
	Arbitrary Line
	Parallel
	Horizontal Line
	Vertical Line
	Rectangle
	Circle
	Ellipse
	Annulus
	Two Circles and Center Distance
	Arc
	Polygon
	Curve
	Make Calibration to determine the corresponding relation between magnification and resolution, this will establish the corresponding relationship between measurement unit and the sensor pixel size. Calibration needs to be done with the help of a micrometer. For detailed steps of carrying out calibration please refer to ProView help manual
	Export the measurement information to CSV file(*.csv)
	Delete All the Measurement Objects

	Setting
	Exit from Current Measurement Mode
	When the measurement ends, left-click on a single measuring object and the Object Location & Properties Control Bar will show up. The icons on the control bar mean Move Left, Move Right, Move Up, Move Down, Color Adjustment and Delete

Note:

1. When user left-clicks Display/Hide button  on the Measurement Toolbar, the Measurement Toolbar will be fixed. In this case the Camera Control Panel will not pop up automatically even if moving mouse cursor to the left side of the video window. Only when users left click the  button on the Measurement Toolbar to exit from the measurement mode they will be able to doing other operations in the Camera Control Panel, the Auto Focus Control Panel or the Synthesis Camera Control Toolbar.
2. When a specific measuring object is selected during the measuring process, the Object Location & Attributes Control Bar  will appear for changing the object location and properties of the selected objects.

7 C-HE Basic features

C-HE is a HDMI + SD card CMOS camera.
HDMI is used as the data transfer interface to HDMI display.


For HDMI output, the Camera Control Panel+ Measurement Toolbar and Camera Control Toolbar are overlaid on the HDMI screen when the mouse move to the related region.

The USB mouse is used to set the camera.

The C-HE's basic characteristics are as follows:

For HDMI output:

- All in 1(HDMI + SD card) C-mount camera with Aptina high sensitivity CMOS sensor;
- Built-in mouse control;
- Built-in image capture to SD card;
- Built-in camera control panel, including exposure (manual/auto)/gain, color adjustment, sharpness control;
- Built-in toolbar including zoom, mirror, freeze, cross;

C-HE - SW UI for Mouse Control

7.1 Quick Instructions for C-HE camera

Figure 8 - The Layout of C-HE Rear Cover

- | | |
|-----------------------|---------------------|
| ① HDMI connector | ③ Power supply plug |
| ② USB mouse connector | ④ SD card slot |

Before starting the camera, please connect the standard C-HE C-mount camera to the camera adapter and insert it into the microscope's photo tube which will relay microscope object's middle image to the camera sensor.

7.1 HDMI Mode

1. Plug the HDMI cable into the HDMI Port to connect the C-HE camera to HDMI display;
2. Plug a USB mouse into USB Mouse to get control of the camera by using built-in software LiteView;
3. Insert SD card into SD Card Slot for saving captured images;
4. Plug 12V/ 1A power adapter into Power Interface to supply power for the camera. The LED Indicator will turn into blue;
5. Move mouse cursor to the left side of the video window, a Camera Control Panel will appear. It includes Manual/ Automatic Exposure, White Balance, Sharpness and other functions, please refer to 8.1 for details;
6. Move mouse cursor to the bottom of the video window and a Camera Control Toolbar will appear. Operations like Zoom In, Zoom Out, Flip, Freeze, Cross Line, and etc. can be realized. Please refer to 8.2 for details;

8 C-He User Interface (UI) and Its Functions

The C-He UI shown in Fig. 2 includes a *Camera Control Panel* on the left side of the video window and a *Camera Control Toolbar* on the bottom of the video window.



Figure 9 - The C-He Camera Control UI

Notes


1. When users click on the  button, the *Camera Control Panel* will pop up automatically.
2. When users move mouse cursor to the bottom of the video window, the *Camera Control Toolbar* will pop up automatically;

8.1 The Camera Control Panel

The Camera Control Panel controls the camera to achieve the best image quality according to the specific applications. It will pop up automatically when mouse cursor is moved to the left side of the video window.

Camera Control Panel	Function	Description
	Auto Exposure	When Auto Exposure is checked, the system will automatically adjust exposure time according to the value of exposure compensation
	Exposure Time	Available when Auto Exposure is unchecked. Slide to left or right to reduce or increase exposure time, adjusting brightness of the video
	Gain	Adjust Gain to reduce or increase brightness of video. The Noise will be reduced or increased accordingly
	Red	Slide to left or right to decrease or increase the proportion of Red in RGB on video
	Green	Green is base for reference and cannot be adjusted
	Blue	Slide to left or right to decrease or increase the proportion of Blue in RGB on the video
	White Balance	Search a zone on the slide wher only the glass is visible and click White Balance
	Sharpness	Adjust Sharpness level of the video window
	Denoise	Adjust Noise level of the video window
	Saturation	Adjust Saturation level of the video window
	Gamma	Adjust Gamma level of the video. Slide to the right side to increase gamma and to the left to decrease gamma
	Contrast	Adjust Contrast level of the video. Slide to the right side to increase contrast and to the left to decrease contrast
	DC	For DC illumination, there will be no fluctuation in light source so no need for compensating light flickering
	AC(50HZ)	Check AC(50HZ) to eliminate flickering "strap" caused by 50Hz illumination
	AC(60HZ)	Check AC(60HZ) to eliminate flickering "strap" caused by 60Hz illumination
Default	Set all the settings in the Camera Control Panel to default values	

8.2 Icons and Functions of the Camera Control Toolbar

Icon	Function	Icon	Function
	Snap the image and save it to SD card		Horizontal Flip
	Zoom In the Video Window		Video Freeze
	Zoom Out the Video Window		Display Cross Line
	Vertical Flip		Toggle ON/OFF Camera Control Panel

9. C-HPSC / C-HESC

9.1 Overview


9.2 Components

①

②

④

③

- ① C-HP / C-HE Camera
- ② HDMI Monitor

- ③ Mounting bracket
- ④ Screws

9.3 Assembling

1. Screw the four provided allen screws on the back side of the monitor. (Fig. 1)


2. Screw the two provided cross screws on the back side of the camera. (Fig. 2)

3. Install the "C" mount and connect the HDMI and power cables to monitor and camera.
4. System is ready.


10. Micrometric Slide M-005

**Micrometric slide, 26x76mm, with 2 scales
(1mm/100div. for biological microscopes / 10mm/100div. for stereomicroscopes)**

1 DIV=0.01mm

For biological microscopes calibration

1 DIV=0.1mm

For stereo microscopes calibration

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.