HDMI Cameras – H Series

INSTRUCTION MANUAL

Model
C-HUB4K

Ver. 1.3 2023

Table of Contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
4.	Unpacking	5
5.	Intended use	5
6.	Symbols and conventions	5
7.	C-HUB4K Basic features	6
7.1	Quick Instructions for C-HUB4K camera	6
7.1.1	HDMI Mode	6
7.1.2	USB Mode	7
7.1.3	WLAN Mode (AP mode)	7
7.1.4	LAN Ethernet Mode	8
7.1.5	Connecting multi-cameras to a router via LAN port/WLAN STA mode for network application	11
8.	C-HUB4K User Interface (UI) and its Functions	14
8.1	The Camera Control Panel	15
8.2	The Measurement Toolbar	16
8.3	Icons and Functions of the Camera Control Toolbar	17
8.3.1	Setting > Network > General	18
8.3.2	Setting > Network > LAN	18
8.3.3	Setting > Network > WLAN	20
8.3.4	Setting > Measurement	20
8.3.5	Setting > Magnification	21
8.3.6	Setting > Image Format	21
8.3.7	Setting > Video	22
8.3.8	Setting > Storage	23
8.3.9	Setting > Files	23
8.3.10	Setting > Time	24
8.3.11	Setting > Language	24
8.3.12	Setting > Miscellaneous	25
8.4	Camera calibration	26
9.	Micrometric Slide M-005	28
	Equipment disposal	29

1. Warning

This device is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

- ① Camera
- ② HDMI cable
- ③ WLAN 5G adapter
- ④ Wireless mouse
- ⑤ Micrometric slide
- ⑥ SD card 32GB
- ⑦ Power supply
- ⑧ USB2.0 cable

4. Unpacking

The system is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items fall out and get damaged.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

7. C-HUB4K Basic features

C-HUB4K is a multiple interfaces (HDMI + USB2.0 + LAN + WiFi + SD card) CMOS camera. HDMI, USB2.0 and LAN are used as the data transfer interface to HDMI display or computer.

For HDMI output, the *Camera Control Panel + Measurement Toolbar* and *Camera Control Toolbar* are overlaid on the HDMI screen when the mouse moves to the related region.

In this case, the USB mouse can be used to set the camera, browse and compare the captured image, play the video and perform the measurement.

For USB output, there are three possibilities:

1. Connect USB flash drive to save pictures and videos.
2. Connect 5G WLAN modules to transfer image wirelessly in real time (AP/STA)
3. Connect computer with USB connection to transfer image in real time.

For LAN output, the camera can be connected with an Ethernet cable directly to a PC or to a router. From here the ProView/LiteView software can control each camera function.

7.1 Quick Instructions for C-HUB4K camera

Before turning on the camera, please connect the standard C-mount camera to the camera adapter and insert it into the microscope's trinocular port.

Fig. 1 - C-HUB4K Rear Panel

- | | |
|-----------------------|---------------------|
| ① USB mouse connector | ⑥ SD card slot |
| ② USB2.0 connector | ⑦ ON-OFF button |
| ③ USB Video connector | ⑧ LED indicator |
| ④ HDMI connector | ⑨ Power supply plug |
| ⑤ LAN connector | |

You can use the C-HUB4K camera in 5 different ways. Each application requires different hardware environment.

7.1.1 HDMI Mode

1. Plug the HDMI cable into the "HDMI" port (4) to connect the C-HUB4K camera to HDMI display.
2. Plug the USB mouse into "USB Mouse" port (1) to get control of the camera by using built-in software HDMIPRO.
3. Plug 12V/1A power adapter into "DC 12V" slot (9) to supply power for the camera.
 - The LED Indicator (8) will turn into red.
4. Insert SD card into "SD Card" slot (6) or a USB flash disk into "USB2.0" port (2) for saving captured images and recorded videos.
5. Press "ON/OFF" Button (7) to turn on the camera.
 - The LED indicator will turn into blue.
6. Move mouse cursor to the left side of the video window, a *Camera Control Panel* will appear. It includes Manual/ Automatic Exposure, White Balance, Sharpness and other functions, please refer to 8.1 for details.
7. Move mouse cursor to the bottom of the video window and a *Camera Control Toolbar* will appear. Operations like Zoom In, Zoom Out, Flip, Freeze, Cross Line, Comparison and so on can be achieved. Please refer to 8.3 for details.
8. Move mouse cursor to the upper side of the video window, a *Measurement Toolbar* with calibration and other measurement tools will appear, please refer to 8.2 for details. The measurement data can be output with *.CSV format.

7.1.2 USB Mode

- For Windows user (Windows 7/8/10 [32/64 bit]), use ProView.
 - For Mac OS and Linux user (Mac OS 10.10 or above or Linux distributions with kernel 2.6.27 or higher), use LiteView.
1. Plug the USB cable provided with the camera into “USB Video” ③ to connect the C-HUB4K camera to the computer.
 - **Do not connect the USB cable into the “USB Mouse” port. Connecting the camera to the “USB Mouse” port, it will not be connected to the PC and it will not display any image on the screen.**
 2. Plug 12V/1A power cable into “DC 12V” slot ⑨ to supply power for the camera.
 - The LED Indicator ⑧ will turn into red.
 3. Press “ON/OFF” Button ⑦ to turn on the camera.
 - The LED indicator will turn into blue.
 4. Start the software.
 5. Clicking the camera name in the *Camera List* starts the live image.
- **When the USB cable is connected, the mouse will not work. If you want to use the mouse, please unplug the USB cable and restart the camera.**

7.1.3 WLAN Mode (AP mode)

In WLAN (WiFi) mode, the camera can be used either via PC or Android mobile device (tablet or phone).

In case of use via PC:

- For Windows user (Windows 7/8/10 [32/64 bit]), use ProView.
 - For Mac OS and Linux user (Mac OS 10.10 or above or Linux distributions with kernel 2.6.27 or higher), use LiteView.
- The PC should be a WLAN enabled one.

In the case of use via mobile device (tablet or phone), the Optika Mobile View App is used, which can be downloaded from the Google Play store.

Use with PC

1. Plug the HDMI cable into the “HDMI” port ④ to connect the C-HUB4K camera to HDMI display.
2. Plug the USB mouse into “USB Mouse” port ① to get control of the camera by using built-in software HDMIPRO.
3. Plug 12V/1A power adapter into “DC12V” slot ⑨ to supply power for the camera.
 - The LED Indicator ⑧ will turn into red.
4. Press “ON/OFF” Button ⑦ to turn on the camera.
 - The LED indicator will turn into blue.
5. Move the mouse to the bottom of the UI and click the  button on the *Camera Control Toolbar* at the bottom of the video window. A small window called *Settings* will pop up as shown below.



Fig. 2 - WLAN configuration Page


6. Click *Network>WLAN* property page and choose the *AP* option in the *WiFi Mode* edit box (the factory default configuration is *AP* mode).
7. Install ProView or LiteView on your PC.
8. Plug the USB WLAN adapter into the “USB2.0” port ② of the camera.

9. Connect the PC with the WLAN AP point that the camera provides.
 - The network name (SSID) and the WLAN password (default one is 12345678) can be found on the camera's *Setting>Network>WLAN* page in AP mode.
10. Start ProView / LiteView software and check the configuration. Normally, active cameras are automatically recognized. The live image of each camera is displayed.
11. For the display, the *Camera List* tool window is used in the ProView / LiteView software.

Use with mobile device

1. Plug 12V/1A power adapter into “DC12V” slot ⑨ to supply power for the camera.
 - The LED Indicator ⑧ will turn into red.
2. Press “ON/OFF” Button ⑦ to turn on the camera.
 - The LED indicator will turn into blue.
3. Plug the USB WLAN adapter into the “USB2.0” port ② of the camera.

4. In the “Settings” of your device select the WiFi device **C-HUB4K-xxxxx**.
 - The password for connecting the camera is “12345678”.
5. Launch the Optika Mobile View App.
6. Manage the camera using the App to control all function.

7.1.4 LAN Ethernet Mode

This application uses the camera as the network camera. User must configure the IP of the camera and PC manually and ensure their IP addresses are in the same net. The subnet mask and gateway of the camera and PC must be the same.

1. Plug the HDMI cable into the “HDMI” port ④ to connect the C-HUB4K camera to HDMI display.
2. Plug the USB mouse into “USB Mouse” port ① to get control of the camera by using built-in software HDMIPRO.
3. Plug 12V/1A power adapter into “DC 12V” slot ⑨ to supply power for the camera.
 - The LED Indicator ⑧ will turn into red.
4. Press “ON/OFF” Button ⑦ to turn on the camera.
 - The LED indicator will turn into blue.
5. Move the mouse to the bottom of the UI and clicking the  button on the *Camera Control Toolbar* at the bottom of the video window. A small window called *Settings* will pop up as shown below.
6. Click the *LAN* property page, uncheck the *DHCP* item. Input the *IP Address*, *Subnet Mask* and *Default Gateway* for the camera.


Fig. 3 - Camera IP configuration Page

- Designate the Internet Protocol Version 4 (TCP/IPv4) Settings page's IP address on the PC with similar configuration as shown below but with different IP address.




Fig. 4 - PC'S IP configuration Page


- After the above configurations are finished, user can connect the camera to the computer through the Ethernet adapter as shown below:
 - Connect the LAN port with the Ethernet cable to the PC's network port.

- Insert the SD card (supplied) or USB flash disk (not supplied) into the camera SD card slot / USB2.0 port.


9. Start ProView / LiteView software and check the configuration.
10. Clicking the camera name in the *Camera List* starts the live video.

7.1.5 Connecting multi-cameras to a router via LAN port/WLAN STA mode for network application

In LAN/WLAN STA mode, the camera connects to router LAN/WLAN by LAN port / WLAN STA mode. If a router with LAN/WLAN capability is used, users could connect the router with Ethernet cable/WLAN and control the camera.

1. When connecting via LAN the connection and configuration are just the same as in Sec. 7.1.1. or 7.1.4, but here, users need to check *DHCP*.
 - If *Multicast* is disabled or is not supported, users should only select *Unicast*.
 - If *Multicast* is supported by the network, users could select *Multicast* to achieve a better performance, especially in the case that multi-users connect to the same camera.
 - In addition, please guarantee that the broadcasting function is enabled in the network.


Fig. 5 - LAN configuration

- When connecting via WLAN the connection and configuration are just the same as in Sec. 7.1.3 but here, users need to check STA.



Fig. 6 - WLAN configuration


- Install ProView or LiteView on your PC.
 - Plug the Ethernet cable into the camera's LAN port and the other end to the PC (for those connected to router with WLAN STA mode).

- Or plug the USB WLAN adapter into the camera's USB2.0 port (for those connected to router with WLAN STA mode).

- Finally as shown below, 2 HDMI cameras are connected to the router with LAN cable and 2 HDMI cameras are connected to the same router with WLAN STA mode (the number of the cameras, the connection mode (LAN or WLAN STA) connected to the router are determined by the router performance).

5. Make sure that your PC is connected to the LAN or WLAN of the router.
6. Start ProView / LiteView software and check the configuration. Normally, active cameras are automatically recognized. The live image of each camera is displayed.
7. For the display, the *Camera List* tool window is used in the ProView/LiteView software.


Fig. 7 - PROVIEW and C-HUB4K Cameras connected in LAN port/WLAN STA mode

- **Note on data security**

The data transfer of the HDMI 4K camera in LAN or WLAN is not encrypted. Anyone who is connected to the network and has installed the Optika software, can see the live image of all active cameras. Operate the camera with the ProView software, if you want to make sure that nobody in the network can see the camera's live image.

- **About the routers/switches**

It is suggested that routers/switches supporting 802.11ac 5G segment should be selected to achieve better wireless connection experience.

8. C-HUB4K User Interface (UI) and its Functions

The C-HUB4K UI shown in Fig. 8 includes a *Camera Control Panel* on the left side of the video window, a *Measurement Toolbar* on the upper side of the video window, a *Camera Control Toolbar* on the bottom of the video window.


Fig. 8 - Camera Control User Interface

1. Move the mouse cursor to the left side of the video window: the *Camera Control Panel* will pop up automatically.
2. Move the mouse cursor to the upper side of the video window: a *Measurement Toolbar* will pop up for the calibration and measurement operations.
 - When user left-clicks the *Float/Fixed* button  on the Measurement Toolbar, the Measurement Toolbar will be fixed. In this case the Camera Control Panel will not pop up automatically even if users move mouse cursor to left side of the video windows.
 - Only when user left-clicks the  button on the *Measurement Toolbar* to exit from measuring procedure it will be possible to perform other operations on the *Camera Control Panel*, or *Camera Control Toolbar*.
 - During the measuring process, when a specific measuring object is selected an *Object Location & Attributes Control Bar*  will appear for changing location and properties of the selected objects.
3. Move the mouse cursor to the bottom of the video window, the *Camera Control Toolbar* will pop up automatically.

8.1 The Camera Control Panel

The *Camera Control Panel* controls the camera to achieve the best image quality according to the specific applications. It will pop up automatically when mouse cursor is moved to the left side of the video window (in measurement status, the Camera Control Panel will not pop up. Only when measurement process is terminated will the Camera Control Panel pop up by moving mouse cursor to the left side of the video window). Left-clicking button  to achieve *Display/ Auto Hide* switch of the Camera Control Panel.

Control Panel	Function	Description
 <p>The screenshot shows the Camera Control Panel with the following settings: Snap and Record buttons; Auto Exposure checked; Exposure Compensation: 71; Exposure Time: 8ms; Gain: 0; Red: 101; Green: 102; Blue: 75; Auto selected; Sharpness: 0; Denoise: 0; Saturation: 50; Gamma: 6; Contrast: 60; Brightness: 50; AC(60Hz) selected; and a Default button.</p>	Snap	<i>Capture</i> image and save it to the SD card
	Record	<i>Record</i> video and save it to the SD card
	Auto Exposure	When <i>Auto Exposure</i> is checked, the system will automatically adjust exposure time according to the value of exposure compensation
	Exposure Compensation	Available when <i>Auto Exposure</i> is checked. Slide to left or right to adjust <i>Exposure Compensation</i> according to the current video brightness to achieve proper brightness value
	Exposure Time	Available when <i>Auto Exposure</i> is unchecked. Slide to left or right to reduce or increase exposure time, adjusting brightness of the video
	Gain	Adjust <i>Gain</i> to reduce or increase brightness of video. The Noise will be reduced or increased accordingly
	Red	Slide to left or right to decrease or increase the proportion of <i>Red</i> in RGB on video
	Green	Slide to left or right to decrease or increase the proportion of <i>Green</i> in RGB on video
	Blue	Slide to left or right to decrease or increase the proportion of <i>Blue</i> in RGB on the video
	Auto	<i>White Balance</i> adjustment according to the window video every time the button is clicked
	Manual	Slide the <i>Red</i> or <i>Blue</i> to manually set the video <i>White Balance</i>
	ROI (Region of Interest)	Check the <i>ROI</i> item will display a red <i>ROI</i> rectangle on the video window, drag it to the interested area will perform the <i>White Balance</i> according to the area video data
	Sharpness	Adjust <i>Sharpness</i> level of the video window
	Denoise	Slide left or right to <i>denoise</i> the image
	Saturation	Adjust <i>Saturation</i> level of the video window
	Gamma	Adjust <i>Gamma</i> level of the video. Slide to the right side to increase Gamma and to the left to decrease Gamma
	Contrast	Adjust <i>Contrast</i> level of the video. Slide to the right side to increase contrast and to the left to decrease Contrast
	Brightness	Adjust <i>Brightness</i> level of the video. Slide to the right side to increase and to the left to decrease Brightness
	DC	For <i>DC</i> illumination, there will be no fluctuation in light source so no need for compensating light flickering
	AC(50HZ)	Check <i>AC(50HZ)</i> to eliminate flickering “strap” caused by 50Hz illumination
AC(60HZ)	Check <i>AC(60HZ)</i> to eliminate flickering “strap” caused by 60Hz illumination	
Default	Restore all the settings in the <i>Camera Control Panel</i> to default values	

8.2 The Measurement Toolbar

The *Measurement Toolbar* will pop up when moving mouse cursor to any place near the upper side of the video window.

Fig. 9 - The Measurement Toolbar

Icon	Function
	<i>Float/ Fix</i> switch of the <i>Measurement Toolbar</i>
<input checked="" type="checkbox"/> Visible	Define measuring object in <i>Show up/ Hide</i> mode
Nanometer (nm)	Select the desired <i>Measurement Unit</i>
4X	Select <i>Magnification for Measurement after Calibration</i>
	<i>Object Select</i>
	<i>Angle</i>
	<i>4 Points Angle</i>
	<i>Point</i>
	<i>Arbitrary Line</i>
	<i>3 Points Line</i>
	<i>Horizontal Line</i>
	<i>Vertical Line</i>
	<i>3 Points Vertical line</i>
	<i>Parallel</i>
	<i>Rectangle</i>
	<i>Ellipse</i>
	<i>Circle</i>
	<i>3 Points Circle</i>
	<i>Annulus</i>
	<i>Two Circles and Center Distance</i>
	<i>3 Points Two Circles and Center Distance</i>
	<i>Arc</i>
	<i>Text</i>
	<i>Polygon</i>

	Curve
	Scale Bar
	Arrow
	Execute <i>Calibration</i> to determine the corresponding relation between magnification and resolution, which will establish the corresponding relationship between measurement unit and the sensor pixel size. <i>Calibration</i> needs to be done with the help of a micrometer.
	Export the <i>Measurement</i> to CSV file(*.csv)
	<i>Measurement Setup</i>
	Delete All the Measurement Objects
	Exit from <i>Measurement</i> mode
	When the measurement ends, left-click on a single measuring object and the <i>Object Location & Properties Control Bar</i> will show up. User could move the object by dragging the object with the mouse. But more accurate movement could be done with the control bar. The icons on the control bar mean <i>Move Left, Move Right, Move Up, Move Down, Color Adjustment</i> and <i>Delete</i> .

Note:

1. When user left-clicks *Display/Hide* button  on the *Measurement Toolbar*, it will be fixed. In this case the *Camera Control Panel* will not pop up automatically even if moving the mouse cursor to the left edge of the video window. Only when user left-click the  button on the *Measurement Toolbar* to exit from the measurement mode it will be possible to perform other operations on the *Camera Control Panel* or on the *Camera Control Toolbar*.
2. When a specific item is selected during the measurement process, the *Object Location & Attributes Control Bar*  will appear for changing the object location and properties of the selected objects.

8.3 Icons and Functions of the Camera Control Toolbar

Icon	Function	Icon	Function
	Zoom In the Video Window		Zoom Out the Video Window
	Horizontal Flip		Vertical Flip
	Color / Gray		Video Freeze
	Display Cross Line		Image Overlay
	Compare Image with the current video		Browse Images and Videos in the SD Card
	Settings		Check the HDMIPRO Version

The  *Setting* function is relatively more complicated than the other functions. Here are more info about it:

8.3.1 Setting > Network > General


Fig. 10 - General Network Settings Page

<i>Name</i>	The current camera name recognized as the network name
-------------	--

8.3.2 Setting > Network > LAN


Fig. 11 - LAN Network Settings Page

<i>DHCP</i>	Dynamic host control protocol allows DHCP server to automatically assign IP information to the camera. Only in Sec. 7.1.4 LAN networking this item should be checked, so that cameras can automatically get IP information from routers/switches to facilitate networking operation
<i>Unicast/Multicast</i>	By default, <i>Unicast</i> function is used. Only in Sec. 7.1.5 networking environment, when the router/switch has <i>Multicast</i> function, camera can switch to multicast mode, which can save the network bandwidth consumed by the camera and facilitate the connection of more cameras in the same network
<i>IP Address</i>	<p>Every machine on a network has a unique identifier. Just as you would address a letter to send in the mail, computers use the unique identifier to send data to specific computers on a network. Most networks today, including all computers on the Internet, use the TCP/IP protocol as the standard for how to communicate on the network. In the TCP/IP protocol, the unique identifier for a computer is called IP address.</p> <p>There are two standards for IP address: IP Version 4 (IPv4) and IP Version 6 (IPv6). All computers with IP addresses have an IPv4 address, and many are starting to use the new IPv6 address system as well.</p> <p>Users must manually conFig. their IP addresses on the camera side and computer side. The IP addresses set on the camera side and computer side should be in the same network segment. The specific settings are shown</p> <p>Fig. 12: It's usually a private address. Private address is a non-registered address used exclusively within an organization. The internal private addresses retained are listed below: Class A 10.0.0-10.255.255; Class B 172.16.0-172.31.255.255; Class C 192.168.0-192.168.255.255. The suggested <i>IP address</i> is Class C</p>

<i>Subnet Mask</i>	Subnet Mask is used to distinguish network domain from host domain in 32-bit IP address
<i>Default Gateway</i>	A default gateway allows computers on a network to communicate with computers on another network. Without it, the network is isolated from the outside. Basically, computers send data that is bound for other networks (one that does not belong to its local IP range) through the default gateway; Network administrators configure the computer's routing capability with an IP range's starting address as the default gateway and point all clients to that IP address

Uncheck the *DHCP* and select the *Unicast* item, user still need to set the *IP address*, *Subnet mask* and *Default Gateway* as shown below:


Fig. 12 - Manual DHCP in Unicast Settings Page

Uncheck the *DHCP* and select the *Multicast* item, user still need to set the *IP address*, *Subnet Mask* and *Default Gateway* as shown below:


Fig. 13 - Manual DHCP in Multicast Settings Page

8.3.3 Setting > Network > WLAN


Fig. 14 - WLAN Network Settings Page

<i>WiFi Mode</i>	AP/STA mode to select
<i>SSID</i>	Channel for the AP mode and SSID for the STA mode. Here, the SSID is the router's SSID
<i>Password</i>	Camera Password for the AP mode. Router Password for the STA mode

8.3.4 Setting > Measurement

This page is used for the definition of the *Measurement Object* properties


Fig. 15 - Measurement Settings Page

<i>Global</i>	Used for setting digits behind the decimal point for measurement results		
<i>Calibration</i>	<i>Line Width</i>	Used for defining width of the lines for calibration	
	<i>Color</i>	Used for defining color of the lines for calibration	
	<i>EndPoint</i>	<i>Type:</i>	Used for defining shape of the endpoints of lines for calibration
		<i>Null:</i>	Means no endpoints
		<i>Rectangle:</i>	Means rectangle type of endpoints. It makes alignment more easily
<i>Point, Angle, Line, Horizontal Line, Vertical Line, Rectangle, Circle, Ellipse, Annulus, Two Circles, Polygon, Curve:</i>			
Left-click the  beside the measuring patterns mentioned above will unfold the corresponding attribute settings to set the individual property of the measuring objects			

8.3.5 Setting > Magnification

This page's items are formed by the *Measurement Toolbar* --> *Calibration* command.


Fig. 16 - Magnification Settings Page

<i>Name</i>	Names such as 10X, 40X, 100X are based on magnification of the microscopes. For continuous zoom microscopes, ensure that the selected magnification coincides with the scale alignment line on the microscope zoom knob
<i>Resolution</i>	Pixels per meter. Devices like microscopes have high resolution value
<i>Clear All</i>	Click the Clear All button will clear the calibrated magnifications and resolutions
<i>Delete</i>	Click Delete to delete the selected item for specific resolution

8.3.6 Setting > Image Format


Fig. 17 - Image Format Settings Page

<i>Image Format</i>	<i>JPEG</i>	A JPEG file can get very high compression rate and display very rich and vivid images by removing redundant images and color data. In other words, it can get better image quality with the least disk space. If measurement objects are available, the measurement objects will be burned into the image and the measurement cannot be edited
	<i>TIFF</i>	TIFF is a flexible bitmap format mainly used to store images including photos and artistic images
<i>Measurement Object Save Method</i>	<i>Burn in Mode</i>	The measurement objects are merged into the current image. User could not edit the measurement objects anymore. This mode is not reversible
	<i>Layered Mode</i>	The measurement objects are saved in different layer with current image data in the target file. User could edit the measurement objects in the target file with some software on the PC. This mode is reversible

8.3.7 Setting > Video


Fig. 18 - Video Settings Page - Playback


Fig. 19 - Video Settings Page - Video Encode

<i>Playback</i>	Fast Forward/Reverse interval in second units for Video Playback
<i>Video Encode</i>	Select the Video Encode format. It can be H264 or H265. Compared with H264, H265 has a higher H265 compression ratio which is primarily used to further reduce the design flow rate, in order to lower the cost of storage and transmission

8.3.8 Setting > Storage


Fig. 20 - SD Card / USB Flash Disk Settings Page

<i>File System Format of the Storage Device</i>	List the file system format of the current storage device	
	<i>FAT32</i>	The file system is FAT32. The maximum file size of single file is 4GB
	<i>exFAT</i>	The file system is exFAT. The maximum file size of single file is 16EB
	<i>NTFS</i>	The file system is NTFS. The maximum file size of single file is 2TB.
	<i>Unknown Status</i>	SD card /USB Flash Disk not detected or the file system is not identified

- **Note:** For USB Flash Drive, USB 3.0 interface is preferred.

8.3.9 Setting > Files


Fig. 21 - File Name Settings Page

<i>Image or Video File Name Paradigm</i>	Provide Auto or Manual naming paradigm for Image or Video file
<i>Auto</i>	With specified name as the Prefix, HDMIPRO will add digits after the Prefix for the Image or Video file
<i>Manual</i>	A file dialog will pop up to enter the Image or Video file name for the captured Image or Video

8.3.10 Setting > Time


Fig. 22 - Time Settings Page

<i>Time</i>	User can set <i>Year, Month, Day, Hour, Minute</i> and <i>Second</i> in this page
-------------	---

8.3.11 Setting > Language

Fig. 23 - Language Settings Page

<i>English</i>	Set language of the whole software into English
<i>Simplified Chinese</i>	Set language of the whole software into Simplified Chinese
<i>Traditional Chinese</i>	Set language of the whole software into Traditional Chinese
<i>Korean</i>	Set language of the whole software into Korean
<i>Thailand</i>	Set language of the whole software into Thai
<i>French</i>	Set language of the whole software into French
<i>German</i>	Set language of the whole software into German
<i>Japanese</i>	Set language of the whole software into Japanese
<i>Italian</i>	Set language of the whole software into Italian
<i>Russian</i>	Set language of the whole software into Russian

8.3.12 Setting > Miscellaneous



Fig. 24 - Miscellaneous Settings Page

<i>Clarity Factor</i>	Check this to show the <i>Clarity Factor</i> on the video window screen and inform if the camera is correctly focused or not
<i>Ruler</i>	Select to display the ruler in the video window, otherwise not to display the ruler
<i>Overlay</i>	Select to support saving graphics overlay information in fusion mode, otherwise it will not support
<i>Grids</i>	Select to support saving grids information in fusion mode, otherwise not to support
<i>USB video output switch back to mouse operation</i>	Select automatic restart or manual restart to switch from USB video output to mouse operation
<i>ROI Color</i>	Choosing the <i>ROI</i> rectangle line color
<i>Cursor</i>	Choosing the <i>Cursor</i> size according to the screen resolution or personal preference
<i>Auto Exposure</i>	Define the maximum automatic exposure time
<i>Auto Exposure Region</i>	Select the AE (Auto Exposure) reference area
<i>Camera Parameters Import</i>	Import the Camera Parameters from the SD Card or USB flash drive to use the previously exported Camera Parameters
<i>Camera Parameters Export</i>	Export the Camera Parameters to the SD Card or USB flash drive to use the previously exported Camera Parameters
<i>Reset to factory defaults</i>	Restore camera parameters to its factory status


8.4 Camera calibration

The camera must be calibrated before any measurements can be taken. To do this it is necessary to use the micrometer slide supplied with the camera.

1. Insert the objective with the lowest magnification available and put the stage micrometer on the stage.
2. Click on the “*Calibrate*” button in the Measurement toolbar.
 - The following dialog box appears:


3. In the *Magnification* field enter the current magnification
 - Some magnification values (default ones 4x, 10x, 20x, 40x and 100x) are preset in the drop-down list.
 - If user must insert a different magnification value, double click in the *Magnification* field. The following dialog box appears:

- From here user can enter any desired magnification value.


4. Focus the ruler under the microscope.
 - Along with the previous dialog box a red line appears on the image. This line is the calibration line.

5. Align the two ends of the red line with the ruler scale (try to use as many notches as possible to keep the calibration accuracy).
6. Set the desired unit (micrometer, millimeter or nanometer) in the *Calibration* dialog box.
7. Insert the length of the drawn line in the *Actual Length* field (using the same procedure described above).
 - Micrometric slide has two rulers (see chapter 9). One is 1 mm/100 (for biological microscopes) and the other is 10mm/100 (for stereomicroscopes).
 - When using the ruler 1mm/100 the interval between two long lines is 100 μm , while using the ruler 10mm/100 the interval between two long lines is 1000 μm .
8. If everything is fine, click *OK* to end the calibration. The new magnification (the number you entered in the Magnification edit box) will be available in the *Magnification* drop-down list.
9. Repeat the same procedure for each magnification available in your microscope.


9. Micrometric Slide M-005

Micrometric slide, 26x76mm, with 2 scales
(1mm/100div. for biological microscopes / 10mm/100div. for stereomicroscopes)

1 DIV=0.01mm

For biological microscopes calibration


1 DIV=0.1mm

For stereomicroscopes calibration

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. “According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal.”

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.