Camera series C-WF

# INSTRUCTION MANUAL

| Model |
|-------|
| C-WF  |
| C-WFR |

Ver. 1.5 2023

---

## Table of contents

|       |  |    |
|-------|--|----|
| 1.    | Warning  | 3  |
| 2.    | Safety Information                               | 3  |
| 3.    | Package content                                  | 4  |
| 3.1   | C-WF   | 4  |
| 3.2   | C-WFR  | 4  |
| 4.    | Unpacking  | 5  |
| 5.    | Intended use                                     | 5  |
| 6.    | Symbols and conventions                          | 5  |
| 7.    | Instrument description                           | 6  |
| 7.1   | C-WF   | 6  |
| 7.2   | C-WFR  | 6  |
| 8.    | Assembling                                       | 7  |
| 8.1   | Assembling on microscope                         | 7  |
| 8.1.1 | Trinocular port with focusable “C” mount adapter | 7  |
| 8.1.2 | Trinocular port with projection lens             | 7  |
| 8.1.3 | Binocular or monocular head                      | 8  |
| 8.2   | Battery installing/replacement (C-WFR)           | 10 |
| 9.    | Use of camera                                    | 11 |
| 9.1   | WiFi camera connection                           | 11 |
| 9.1.1 | Windows operating system                         | 11 |
| 9.1.2 | Android operating system                         | 11 |
| 9.2   | Use of the camera (C-WF) using PC connection     | 11 |
| 9.3   | Use of the camera (C-WFR) using PC connection    | 12 |
| 9.4   | Parfocality adjustment                           | 12 |
| 9.4.1 | Biological microscope                            | 12 |
| 9.4.2 | Stereo microscope with “C” mount adapter         | 12 |
| 9.4.3 | Stereo microscope with projection lens           | 13 |
| 10.   | Micrometric Slide M-005                          | 14 |
|       | Equipment disposal                               | 15 |

---

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

## 2. Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

### 3. Package content

#### 3.1 C-WF

- ① Camera
- ② Projection lens
- ③ Reduction ring 30.5 mm
- ④ Reduction ring 30 mm
- ⑤ Micrometric slide
- ⑥ Power supply with USB2.0 socket

#### 3.2 C-WFR

- ① Camera
- ② Projection lens
- ③ Reduction ring 30.5 mm
- ④ Reduction ring 30 mm
- ⑤ Micrometric slide
- ⑥ Power supply with USB2.0 socket

---

## 4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 5. Intended use

### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

## 6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


### ELECTRICAL SHOCK



This symbol indicates a risk of electrical shock.

## 7. Instrument description

### 7.1 C-WF

### 7.2 C-WFR

## 8. Assembling

### 8.1 Assembling on microscope

#### 8.1.1 Trinocular port with focusable "C" mount adapter

1. Remove the dust cover from camera and from "C" mount adapter.
2. Screw the "C" mount adapter to the camera's thread. (Fig. 1)

Fig. 1

3. Loosen the clamping screw on the trinocular port of the microscope, remove the dust cover and install the camera "C" mount unit previously assembled. (Fig. 2)


Fig. 2

4. Tighten the clamping screw firmly. (Fig. 3)

Fig. 3

#### 8.1.2 Trinocular port with projection lens

1. Remove the dust cover from camera and from projection lens.
2. Screw the projection lens to the camera's thread. (Fig. 4)


Fig. 4

- **Installing on biological microscopes**

3. Insert the final part of the projection lens into the empty hole of the photo tube. (Fig. 5)

- **Installing on stereo microscopes**


4. Insert the final part of the projection lens into the empty hole of the photo tube. (Fig. 6)

### 8.1.3 Binocular or monocular head

- Any OPTIKA monocular or binocular microscope/stereo microscope (eyepiece sleeve with diameter 23 mm).
- Any non-OPTIKA monocular or binocular microscope/stereo microscope (eyepiece sleeve with diameter 23 mm).

1. Remove the dust cover from camera and from projection lens.
2. Screw the projection lens to the camera's thread. (Fig. 7)

3. Remove one eyepiece from the observation head. (Fig. 8)

4. Insert the final part of the projection lens into the empty eyepiece sleeve. (Fig. 9)


5. In case of eyepiece sleeve with diameter 30 or 30.5 mm use the adapter rings. (Fig. 10)

## 8.2 Battery installing/replacement (C-WFR)

- **Batteries are not provided. It is necessary to use 4x AA NiMH 1.2V batteries.**

1. Unplug the power supply.
2. Unscrew the top cover ① of the battery compartment. (Fig. 11)

3. Remove the battery holder and insert the 4 batteries (respecting the polarities). (Fig. 12)
- In case of battery replacement, remove the exhausted batteries and insert new ones.

4. Connect the battery pack to the connector in the battery compartment. (Fig. 13)
5. Insert the battery pack into the compartment and screw in the top cover.

---

## 9. Use of camera

The camera is driven by a software.

Can be used with both Windows and Android operating systems.

1. When used with Windows operating system, two levels of software are available: PROVIEW and LITEVIEW.

Inside the package is enclosed a Function Table showing the several software functions.

It will be the user's responsibility to decide which level of software best meets his needs.

The software can be downloaded from the site:

<http://www.optikamicroscopes.com/optikamicroscopes/optika-support/download-drivers-softwares/>

Once the file has been downloaded, you will have to run the setup.exe file.

At the end of the installation it is possible to start the software.

- **NOTE: no driver installation is required for the cameras. The software setup procedure automatically installs all the needed drivers for the correct operation of the camera.**

The software's User Manual is available in PDF format within the application itself and can be opened using the "F1" function key.

You must have Acrobat Reader installed to view the manual.

The manual contains all the operating instructions for using the camera and for the various functions of the software.

2. When used with Android operating system, one App is available: Optika Mobile View.

The App can be downloaded from the Google Play store.

### 9.1 WiFi camera connection

The camera is powered via USB cable (C-WF version) or via USB cable and/or rechargeable batteries (C-WFR version).

#### 9.1.1 Windows operating system

- **The camera connects to PC only and exclusively with a WiFi connection.**
- **There is no possibility to connect the camera to PC with a USB cable.**

1. Connect the USB power cable to the camera.
2. In the "Settings" of your PC select the WiFi device **C-WFR-xxxxx**.
  - The password for connecting the camera is "12345678".
3. Launch the PROVIEW or LITEVIEW software.
4. Manage the camera using the software to control all function.

#### 9.1.2 Android operating system


1. Connect the USB power cable to the camera.
2. In the "Settings" of your device select the WiFi device **C-WFR-xxxxx**.
  - The password for connecting the camera is "12345678".
3. Launch the Optika Mobile View App.
4. Manage the camera using the App to control all function.

### 9.2 Use of the camera (C-WF) using PC connection

1. Connect the USB power cable to the camera.
2. Launch the PROVIEW or LITEVIEW software.
3. From "Camera List" panel in the software, click on the camera name (C-WFR).
4. Live image appears.
5. Manage the camera using the software to control all function.

### 9.3 Use of the camera (C-WFR) using PC connection

1. Connect the USB power cable to the camera.
    - The camera is equipped with an ON-OFF switch ① (Fig. 14) which allows you to recharge the batteries (when the selector switch is set to “OFF”) or recharge the batteries and turn the camera on (when the selector switch is set to “ON”).
    - If the selector switch is set to “OFF” the operating LED ② (Fig. 15) is off. When the selector switch is set to “ON” the LED lights up. This indicates that the camera is active.
    - The camera also works when the power cable is disconnected. In this case the camera is only powered with batteries.
  2. Launch the PROVIEW or LITEVIEW software
  3. From “Camera List” panel in the software, click on the camera name (C-WFR).
  4. Live image appears.
  5. Manage the camera using the software to control all function.
- **NOTE: if the batteries are not inserted, the camera does not work even if the power cable is connected.**

### 9.4 Parfocality adjustment

In order to have the same focus when observing the specimen through the eyepieces and on the screen/monitor, please check that the microscope is properly installed and follow the instructions below.

#### 9.4.1 Biological microscope

1. Use a low power objective and focus the specimen.
2. Switch to the highest dry objective available on the microscope (40x or 60x) and focus the specimen again.
3. Turn on the live-view on the camera, without changing the focus on the microscope.
4. While observing the image on the screen/monitor, adjust the focus by turning the knurled knob on the C-mount adapter. (Fig. 16)

#### 9.4.2 Stereo microscope with “C” mount adapter

1. Use a low power magnification and focus the specimen.
2. Reach the highest magnification available using the zoom knob and then focus the specimen again.
3. Turn on the live-view on the camera, without changing the focus on the microscope.
4. While observing the image on the screen/monitor, adjust the focus by turning the knurled knob on the C-mount adapter. (Fig. 16)

### 9.4.3 Stereo microscope with projection lens

1. Use a low power magnification and focus the specimen.
2. Reach the highest magnification available using the zoom knob and then focus the specimen again.
3. Turn on the live-view on the camera, without changing the focus on the microscope.
4. While observing the image on the screen/monitor, unscrew the fixing knob ① of the trinocular port (Fig. 17) and move up or down the movable part of the adapter ② until the focus is achieved on the camera. (Fig. 18)
5. Lock the fixing knob ①. (Fig. 17)


The proper parfocality adjustment is obtained when the same focus is reached when looking into the eyepieces and on the screen/monitor.

---


## 10. Micrometric Slide M-005

**Micrometric slide, 26x76mm, with 2 scales**  
(1mm/100div. for biological microscopes / 10mm/100div. for stereo microscopes)

1 DIV=0.01mm

**For biological microscopes calibration**

1 DIV=0.1mm

**For stereo microscopes calibration**

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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**OPTIKA® S.r.l.**