Camera series C-WF

INSTRUCTION MANUAL

| Model |
|-------|
| C-WF2 |

Ver. 1.0 2024

Table of contents

| | | |
|-------|--|----|
| 1. | Warning | 3 |
| 2. | Safety Information | 3 |
| 3. | Package content | 4 |
| 3.1 | C-WF2 | 4 |
| 4. | Unpacking | 5 |
| 5. | Intended use | 5 |
| 6. | Symbols and conventions | 5 |
| 7. | Instrument description | 6 |
| 7.1 | C-WF2 | 6 |
| 8. | Assembling | 7 |
| 8.1 | Assembling on microscope | 7 |
| 8.1.1 | Trinocular port with focusable “C” mount adapter | 7 |
| 8.1.2 | Trinocular port with projection lens | 7 |
| 8.1.3 | Binocular or monocular head | 8 |
| 9. | Use of camera | 10 |
| 9.1 | PC camera connection | 10 |
| 9.2 | LAN camera connection | 10 |
| 9.3 | WiFi camera connection (AP mode) | 11 |
| 9.3.1 | Windows operating system | 11 |
| 9.3.2 | IOS / Android operating system | 11 |
| 9.4 | WiFi camera connection (STA mode) | 12 |
| 9.5 | Parfocality adjustment | 13 |
| 9.5.1 | Biological microscope | 13 |
| 9.5.2 | Stereo microscope with “C” mount adapter | 13 |
| 9.5.3 | Stereo microscope with projection lens | 13 |
| 10. | Configure the Camera’s WiFi STA Mode’s SSID and Password | 14 |
| 11. | Micrometric Slide M-005 | 15 |
| | Equipment disposal | 16 |

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

3.1 C-WF2

① Camera

② Projection lens

③ Reduction ring 30.5 mm

④ Reduction ring 30 mm

⑤ Micrometric slide

⑥ USB to micro USB cable

4. Unpacking

The microscope is housed in a molded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

7. Instrument description

7.1 C-WF2

8. Assembling

8.1 Assembling on microscope

8.1.1 Trinocular port with focusable "C" mount adapter

1. Remove the dust cover from camera and from "C" mount adapter.
2. Screw the "C" mount adapter to the camera's thread. (Fig. 1)

** These images represent the installation of a generic camera. The C-WF2 is installed in exactly the same way.

Fig. 1

3. Loosen the clamping screw on the trinocular port of the microscope, remove the dust cover and install the camera "C" mount unit previously assembled. (Fig. 2)


Fig. 2

4. Tighten the clamping screw firmly. (Fig. 3)

Fig. 3

8.1.2 Trinocular port with projection lens

1. Remove the dust cover from camera and from projection lens.
2. Screw the projection lens to the camera's thread. (Fig. 4)


Fig. 4

- **Installing on biological microscopes**

3. Insert the final part of the projection lens into the empty hole of the photo tube. (Fig. 5)


- **Installing on stereo microscopes**

4. Insert the final part of the projection lens into the empty hole of the photo tube. (Fig. 6)

8.1.3 Binocular or monocular head

- Any OPTIKA monocular or binocular microscope/stereo microscope (eyepiece sleeve with diameter 23 mm).
 - Any non-OPTIKA monocular or binocular microscope/stereo microscope (eyepiece sleeve with diameter 23 mm).
1. Remove the dust cover from camera and from projection lens.
 2. Screw the projection lens to the camera's thread. (Fig. 7)

3. Remove one eyepiece from the observation head. (Fig. 8)

4. Insert the final part of the projection lens into the empty eyepiece sleeve. (Fig. 9)

5. In case of eyepiece sleeve with diameter 30 or 30.5 mm use the adapter rings. (Fig. 10)

9. Use of camera

The camera is driven by a software.
It can be used with Windows, IOS and Android operating systems.

1. When used with Windows operating system, two levels of software are available: PROVIEW and LITEVIEW.

Inside the package is enclosed a Function Table showing the several software functions.
It will be the user's responsibility to decide which level of software best meets his needs.

The software can be downloaded from the site:
<http://www.optikamicroscopes.com/optikamicroscopes/optika-support/download-drivers-softwares/>

Once the file has been downloaded, you will have to run the setup.exe file.
At the end of the installation it is possible to start the software.

- **NOTE: no driver installation is required for the cameras. The software setup procedure automatically installs all the needed drivers for the correct operation of the camera.**

The software's User Manual is available in PDF format within the application itself and can be opened using the "F1" function key.
You must have Acrobat Reader installed to view the manual.

The manual contains all the operating instructions for using the camera and for the various functions of the software.

2. When used with IOS or Android operating system, one App is available: Optika Mobile View.

The App can be downloaded from the App Store or Google Play Store.

- | | |
|-----------------|------------------|
| ① WiFi antenna | ④ LAN port |
| ② AP/STA button | ⑤ Micro USB port |
| ③ AP/STA LED | ⑥ Power LED |

9.1 PC camera connection

The camera is powered via micro USB cable.

1. Connect the micro USB side of the provided cable to the micro USB camera port (5).
2. Connect the USB side of the provided cable to one of the USB port in your PC.
3. Launch the Proview or Liteview software.
4. Manage the camera using the software to control all function.


9.2 LAN camera connection

The camera is powered via mini USB cable.
In LAN mode, the camera connects to a router by Ethernet cable.

You can connect one or more cameras to the router by Ethernet cable for network applications.

1. Install the Proview or Liteview software on your PC. Or, install the free Optika Mobile View App on the mobile device.
2. Connect the micro USB side of the provided cable to the micro USB camera port (5).

3. Plug the USB side of the cable to a USB power adapter (not supplied).
4. Plug the power adapter to a wall outlet.
5. Connect the camera and the router by plugging the Ethernet cable into the LAN port ④: if the connection is successful, the network indication LED will be off and the LAN indicator will be light up.
6. Make sure that your PC or your mobile device is connected to the router (by LAN or WiFi).
7. Start the Proview or Liteview software or Optika Mobile View App and check the configuration.
8. Normally, active C-WF2 cameras are automatically recognized. The connected cameras will be listed in the *Camera List* group of the Camera Control Bar of Proview or Liteview on PC.
9. On mobile device, connected cameras will be listed in the Camera Thumbnail page of Optika Mobile View App.
10. Click the corresponding camera to start live image.

9.3 WiFi camera connection (AP mode)

The camera is powered via micro USB cable.

- A camera can be connected to a maximum of three devices at the same time.
- The WiFi signal generated by the camera may conflict with the signal generated by a second camera in the same room. In order to avoid interference between the signals (slowing down of the live image or black screen), it is recommended not to install more than 6 cameras in a room of approximately 60 sq. mt. Furthermore, the mobile device connected to the camera should be close enough to the camera itself.
- If more than 6 cameras are to be installed in the same room, it is necessary to proceed in STA mode.

9.3.1 Windows operating system

- The PC must be equipped with a WiFi receiver to handle the WiFi signal generated by the camera.
1. Connect the micro USB side of the provided cable to the micro USB camera port ⑤.
 2. Plug the USB side of the cable to a USB power adapter (not supplied).
 3. Plug the power adapter to a wall outlet.
 4. Press the AP/STA button ②. The operating LED ③ will light up in green when AP mode is selected.
 5. In the "Network Settings and Internet" of your PC select the WiFi device **C-WF2-xxxxx**.
 - The password for connecting the camera is "12345678".
 6. Launch the Proview or Liteview software.
 7. Manage the camera using the software to control all function.


9.3.2 IOS / Android operating system

1. Connect the micro USB side of the provided cable to the micro USB camera port ⑤.
2. Plug the USB side of the cable to a USB power adapter (not supplied).
3. Plug the power adapter to a wall outlet.
4. Press the AP/STA button ②. The operating LED ③ will light up in green when AP mode is selected.
5. In the "Settings" of your device select the WiFi device **C-WF2-xxxxx**.
 - The password for connecting the camera is "12345678".
6. Launch the Optika Mobile View App.
7. Manage the camera using the App to control all function.

9.4 WiFi camera connection (STA mode)

The camera is powered via micro USB cable.

- In STA mode, the camera is supposed to connect to a router.
- A WiFi router can usually provide 10 signals, which means 10 cameras can be connected to the router.
- Assuming 4 WiFi routers are installed in a room, we ensure that 40 cameras can coexist in a room without interference from WiFi signals.

1. Connect the micro USB side of the provided cable to the micro USB camera port ⑤.
2. Plug the USB side of the cable to a USB power adapter (not supplied).
3. Plug the power adapter to a wall outlet.
4. Press the AP/STA button ②. The operating LED ③ will light up in blue when STA mode is selected.
 - This means the camera is in STA mode and has connected to the router (See Chapter 10 on how to set the router's name and password in the camera).
5. Make sure that your PC or your mobile device is connected to the router (by LAN or WiFi).
6. Start the Proview or Liteview software or Optika Mobile View App and check the configuration.
7. Normally, active C-WF2 cameras are automatically recognized. The connected cameras will be listed in the *Camera List* group of the Camera Control Bar of Proview or Liteview on PC.
8. On mobile device, connected cameras will be listed in the Camera Thumbnail page of Optika Mobile View App.
9. Click the corresponding camera to start the stream.

• Note on data security

The data transfer of the C-WF2 camera in LAN or WiFi is not encrypted. Anyone who is connected to the network and has installed the Optika software or Optika App, can see the live image of all active C-WF2 cameras.

• About the routers/switches

It is suggested that routers/switches supporting 802.11ac 5G segment should be selected to achieve better wireless connection experience.

9.5 Parfocality adjustment

In order to have the same focus when observing the specimen through the eyepieces and on the screen/monitor, please check that the microscope is properly installed and follow the instructions below.

9.5.1 Biological microscope

1. Use a low power objective and focus the specimen.
2. Switch to the highest dry objective available on the microscope (40x or 60x) and focus the specimen again.
3. Turn on the live-view on the camera, without changing the focus on the microscope.
4. While observing the image on the screen/monitor, adjust the focus by turning the knurled knob on the C-mount adapter. (Fig. 11)

9.5.2 Stereo microscope with “C” mount adapter

1. Use a low power magnification and focus the specimen.
2. Reach the highest magnification available using the zoom knob and then focus the specimen again.
3. Turn on the live-view on the camera, without changing the focus on the microscope.
4. While observing the image on the screen/monitor, adjust the focus by turning the knurled knob on the C-mount adapter. (Fig. 11)

** These images represent the set-up of a generic camera. The C-WF2 is set in exactly the same way.

9.5.3 Stereo microscope with projection lens

1. Use a low power magnification and focus the specimen.
2. Reach the highest magnification available using the zoom knob and then focus the specimen again.
3. Turn on the live-view on the camera, without changing the focus on the microscope.
4. While observing the image on the screen/monitor, unscrew the fixing knob ① of the trinocular port (Fig. 12) and move up or down the movable part of the adapter ② until the focus is achieved on the camera. (Fig. 13)
5. Lock the fixing knob ①. (Fig. 12)

The proper parfocality adjustment is obtained when the same focus is reached when looking into the eyepieces and on the screen/monitor.

Fig. 11

Fig. 12


Fig. 13


10. Configure the Camera's WiFi STA Mode's SSID and Password

Referring to Section 9.4, the detailed steps are as follows:


1. Connect the C-WF2 camera working in WiFi AP mode using iOS or Android devices.
2. Launch the Optika Mobile View App.


3. Press the configuration button  in the camera thumbnail above to pop up the interface for selecting WiFi SSID.
4. A new dialog box will appear:

5. Tap on the WiFi icon  on the left side of the screen.
6. A dialog box containing all the available WiFi router (SSID) will appear.


7. Select the SSID to be connected (for example OPTIKA_1) and input the Password (12345678).
8. Press OK to finish the SSID and Password setup process.

9. If the pop-up dialog box does not have the desired SSID, tap the Add WLAN item on the bottom and input SSID name and Password manually (the current SSID is OPTIKA_1 and the password is 12345678).
10. After completion, press the button to switch to STA mode. The camera WiFi will be connected to the designated router in STA mode for operation, and the camera network LED indicator will turn blue.

11. Micrometric Slide M-005

Micrometric slide, 26x76mm, with 2 scales
(1mm/100div. for biological microscopes / 10mm/100div. for stereo microscopes)

1 DIV=0.01mm

For biological microscopes calibration


1 DIV=0.1mm

For stereo microscopes calibration

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.