**IM-3 Series**

# INSTRUCTION MANUAL

<b>Model</b>
IM-3MET

Ver. 3.0 2019

---

## Table of Contents

<b>1. Warning</b>	<b>3</b>
<b>2. Symbols and conventions</b>	<b>3</b>
<b>3. Safety Information</b>	<b>3</b>
<b>4. Intended use</b>	<b>3</b>
<b>5. Overview</b>	<b>4</b>
<b>6. Unpacking</b>	<b>5</b>
<b>7. Assembling</b>	<b>5</b>
7.1 Microscope assembling	6
7.1.1 Installing the objectives	6
7.1.2 Installing the stage insert	6
7.1.3 Installing the eyepieces	7
7.1.4 Installing the lamp housing	7
7.1.5 Installing the filter holder	7
7.1.6 Installing the bulb	8
7.1.7 Connecting the power cord	9
7.1.8 Replacing the fuse	9
7.1.9 Installing mechanical stage or stage extension	9
7.1.10 Installing polarizer and analyzer	10
<b>8. Brightfield observation procedures</b>	<b>11</b>
<b>9. Use of the microscope</b>	<b>12</b>
9.1 Microscope switch on	12
9.2 Light intensity adjustment	12
9.3 Centering the bulb	12
9.4 Coarse focus tension adjustment	13
9.5 Diopter adjustment	13
9.6 Adjusting the interpupillary distance	13
9.7 Use of eye shields	13
9.8 Centering the aperture diaphragm	14
9.9 Effects of the field diaphragm	14
9.10 Aperture diaphragm	14
9.11 Light path selection	15
9.12 Use of color filters	15
9.13 Use of polarized light	16
<b>10. Microphotography</b>	<b>17</b>
10.1 Use of C-mount cameras	17
10.2 Use of Reflex cameras	17
<b>11. Maintenance</b>	<b>18</b>
<b>12. Troubleshooting</b>	<b>19</b>
<b>Equipment disposal</b>	<b>20</b>

---

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

## 2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

### ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

## 3. Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

## 4. Intended use


### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

## 5. Overview

## 6. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 7. Assembling

Once opened the box, the microscope parts are the following:



- |                                 |                                 |
|---------------------------------|---------------------------------|
| ① Microscope body               | ⑦ Analyzer                      |
| ② Eyepieces                     | ⑧ Polarizer                     |
| ③ Objectives                    | ⑨ Filter holder                 |
| ④ Stage metal insert            | ⑩ Lamp housing connection cable |
| ⑤ Lamp housing                  | ⑪ Power cord                    |
| ⑥ Color filters (LBD and IF550) |                                 |

## 7.1 Microscope assembling

### 7.1.1 Installing the objectives

1. Rotate the coarse focusing knob ① until the nose-piece reaches its lowest position. (Fig. 1)
  - **For a safe transport, the nosepiece is placed in the lowest position before shipment from the factory.**
2. Screw the lowest magnification objective on to the turret from the right side, then turn the nose-piece clockwise. Mount the other objectives in the same way, following the sequence from low to high magnification. (Fig. 2)
  - **Clean the objectives regularly. In inverted microscopes, the objectives are very sensitive to dust.**
  - **To prevent dust and dirt from entering the microscope, cover all the unused holes with dust caps ②. (Fig.3)**

### 7.1.2 Installing the stage insert

1. Install the stage insert in the stage opening. (Fig. 4)


### 7.1.3 Installing the eyepieces

1. Insert both eyepieces into the tubes of the optical head. (Fig. 5)


### 7.1.4 Installing the lamp housing

1. Insert the lamp house and tighten the fixing screws ①. (Fig. 6)


### 7.1.5 Installing the filter holder


- The microscope is equipped with a filter holder located on the back side of the microscope. (Fig. 7).
- The filter holder has three positions: ② filter holder for color filters, ③ frosted filter, ④ filter holder for color filters.

1. Remove one of the knobs ⑤ from the filter holder. (Fig. 8)

2. Insert the filter holder into the slot in front of the lamp housing ⑥. (Fig. 9)
3. When the filter holder is in place, reinstall the filter holder knob.


### 7.1.6 Installing the bulb


1. Open the lamp housing cover. (Fig. 10)

2. Insert the halogen bulb into its holder ⑦. (Fig. 11)
  - **Do NOT touch the bulb with the bare hands, this could reduce the bulb efficiency and lifetime.**

3. Plug the power cable from the lamp housing to the microscope body. (Fig. 12)

### 7.1.7 Connecting the power cord

1. Put the main switch ① to “O” (OFF) before connecting the power cord. (Fig.13)


Fig. 13

2. Insert the cable into the power socket of the microscope. (Fig.14)
3. Plug the power cord into the mains socket. Check for a safe connection.
  - **Please use the provided power cord.**
  - **If lost or damaged, please refer to qualified service.**
  - **Connect the power cord to a grounded power supply only.**


Fig. 14

### 7.1.8 Replacing the fuse

- Before replacing the fuse, turn the main switch to “O” (OFF) and unplug the power cord.
1. Rotate the fuse support out of the holder using a straight screwdriver. Insert a new fuse in the support, then rotate the support back into the holder. (Fig.15)
- **Fuse rating: see back of the microscope.**

Fig. 15

### 7.1.9 Installing mechanical stage or stage extension

- **Stage extension and mechanical stage are optional accessories.**
  - Stage extension can be installed on either side of the stage to enlarge the working surface.
  - Mechanical stage can only be installed on the right side of the stage.
1. Screw the bolts in the fixing holes of the attachments, then mount the unit from below the stage. (Fig. 16)

Fig. 16

### 7.1.10 Installing polarizer and analyzer


1. Insert the polarizer “S-PO” into the slot ① placed on the left side of the lamp housing support. (Fig. 17)

2. Insert the analyzer “S-ANR” in the slot ② placed on the right side of the microscope body. (Fig. 18)


## 8. Brightfield observation procedures

## 9. Use of the microscope


### 9.1 Microscope switch on

Put the main switch ① to "I" (ON) to turn on the microscope. (Fig. 19)



### 9.2 Light intensity adjustment

Rotate the light intensity adjustment knob ② to increase / decrease the illumination voltage. (Fig. 20)

### 9.3 Centering the bulb


1. Remove one objective from the nosepiece and put the empty hole in the optical path.
2. Put a piece of paper on the stage, slightly unscrew the knurled ring of the focus adjustment lever ③ on the lamp housing and rotate it counterclockwise until the bulb filament is in focus on the paper. (Fig. 21)
  - If the bulb filament is not in the center of the light path, the resulting image will be not evenly illuminated. In this case a re-centering is needed.
3. Rotate the centering screws ④ on the right side of the lamp housing to bring the bulb filament in the center of the field of view.
4. At the end fully rotate the focus adjustment lever ③ clockwise and screw the knurled ring to fix the lever into its position.
5. Reinstall the objective.

## 9.4 Coarse focus tension adjustment

The coarse knob tension is pre-set in the factory.

1. To modify the tension according to personal's needs, rotate the ring ①. (Fig. 23)
- Clockwise rotation increases the tension.
  - If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

## 9.5 Diopter adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
  2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptic adjustment ring ② to compensate. (Fig. 24)
- **The adjustment range is  $\pm 5$  diopter. The number shown on the diopter ring should correspond to the operator's dioptic correction.**


## 9.6 Adjusting the interpupillary distance

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

- **The graduation on the interpupillary distance indicator ③, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator's eyes. (Fig. 25)**

The range of the interpupillary distance is 48-75mm.

## 9.7 Use of eye shields

### • Use with eyeglasses

Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses. (Fig. 26)

- **Use without eyeglasses**

Raise eye shields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation. (Fig. 28)

Fig. 27

### 9.8 Centering the aperture diaphragm

1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Pull the field diaphragm lever ① to fully close the field diaphragm. (Fig. 28)
3. Rotate the two centering screws ② to bring the bright spot in the center of the field of view.(Fig. 28)
4. Gradually open the diaphragm by inserting the lever ①. The diaphragm is centered when the diaphragm image is symmetrical to the field of view.(Fig. 29)
5. In normal use, open the diaphragm until it circoscribes the field of view.


Fig. 28

### 9.9 Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image.

Set the diaphragm according to the objective in use until it circoscribes the field of view, in order to eliminate unnecessary light to eyepieces.


Fig. 29

### 9.10 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
- With low contrast specimens pull the lever ③ (Fig. 30) to slightly close the aperture diaphragm. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the diaphragm's lever in order to obtain an image like the one in Fig. 31.



Fig. 30

### 9.11 Light path selection


Pull the light path selector lever ① to the left or to the right using your thumb, selecting the light path you need. (Fig. 32)

LIGHT PATH SELECTOR LEVER	LIGHT SPLIT	APPLICATION
In	20% binocular and 80% for PHOTO/TV port	Binocular and camera observation can be operated simultaneously
Out	100% used for binocular observation	Binocular observation

### 9.12 Use of color filters

1. Insert the filter in one of the filter holder. (Fig. 33).

2. Move the slider to the left or to the right in order to insert the desired one. (Fig. 34)

FILTER	USE
Blue	Converts the color temperature of the bulb into daylight
Green	Increases the contrast in black and white photography
Frosted	Reduces lighting irregularities by creating a homogeneous illumination

Fig. 34

### 9.13 Use of polarized light

1. Remove the specimen from the stage.
2. Put a mirror on the stage and get into approximate focus on the mirror surface.
3. Insert polarizer and analyzer into the light path.
4. Rotate the knurled ring of the analyzer ① and, looking into the eyepiece, reach the darkest position available. (Fig. 35)
5. Once the dark is achieved ("extinction" or "Crossed Nicol" position) it is possible to begin the observation.
6. Remove the mirror from the stage, put a specimen and focus.


Fig. 35

## 10. Microphotography

### 10.1 Use of C-mount cameras

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 36)

2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 37)

### 10.2 Use of Reflex cameras

1. Insert the Reflex adapter ① into the relay tube to the microscope ②.
  2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
  3. Connect the Reflex camera ④ to the "T2" just installed (Fig. 38).
  4. Mount the other end of the relay tube ② into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 36)
- "T2" ring is not provided with the microscope, but is commercially available.
  - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
  - To calculate the magnification of the camera: objective magnification \* camera magnification \* lens magnification.
  - **If using an SLR camera, mirror movement may cause the camera to vibrate.**
  - **We suggest lifting the mirror, using long exposure times and a remote cord.**

---

## 11. Maintenance

### Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

### To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

### Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

**For the best results, use the OPTIKA cleaning kit (see catalogue).**

If you need to send the microscope to Optika for maintenance, please use the original packaging.

## 12. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
<b>I. Optical Section:</b>		
Bulb operates, but field of view remains dark.	Power supply is unplugged.	Connect
	Brightness is too low	Set brightness to a proper level
Field of view is obscured or not evenly illuminated	Revolving nosepiece is not correctly engaged.	Make sure that the revolving nosepiece clicks properly into place.
	Filter holder, polarizer or analyzer are partially inserted	Insert or remove the filter slider
Dirt or dust is visible in the field of view.	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture iris diaphragm is stopped down too far.	Open aperture iris diaphragm.
	The field diaphragm is not well centered	Set the diaphragm according to Koehler settings.
Visibility is poor. <ul style="list-style-type: none"> <li>• Image is not poor.</li> <li>• Contrast is poor.</li> <li>• Details are indistinct.</li> </ul>	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture iris diaphragm is too closed or too open.	Adjust aperture iris diaphragm.
	Dust or dirt on lenses (objectives, eyepieces)	Clean thoroughly.
	Focus is not even	Specimen holder is not flat. Move the specimen to a flat position.
One side of the image is unfocused	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Specimen is not in a flat position	Place the specimen in a flat position on the stage
<b>II. Mechanical Section:</b>		
Coarse focus knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
Focus is unstable	Tension adjustment ring is too loose	Tighten tension adjustment ring
<b>III. Electrical Section:</b>		
Bulb doesn't turn on.	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection
<b>IV. Observation tube:</b>		
Field of view of one eye does not match that of the other.	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find helpful to look up and into distance for a moment before looking back into microscope.
<b>V. Microphotography:</b>		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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**OPTIKA® S.r.l.**