Solar Battery Pack

# INSTRUCTION MANUAL

<b>Model</b>
M-069

Ver. 7.3    2023

---

## Table of contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
4.	Unpacking	5
5.	Intended use	5
6.	Symbols and conventions	5
7.	Instrument description	6
8.	Connection guide	8
9.	Recharging the battery	8
10.	Use of the solar panel	9
10.1	Powering the microscope	9
10.2	Using the LED torch light	9
11.	Troubleshooting	10
	Equipment disposal	11

---

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

## 2. Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

### 3. Package content

① Solar panel

② USB-Jack 2.1 mm cable

③ USB-mini USB cable

④ Power supply

---

## 4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 5. Intended use

### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

## 6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION


This symbol indicates a potential risk and alerts you to proceed with caution.

### ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

## 7. Instrument description

## 8. Connection guide

Plug your microscope to the solar panel using the provided USB/Jack cable:

- **PLEASE NOTE:** Although the solar panel has two USB outputs, you can only connect one microscope. Optimum operation of the microscope is not guaranteed if two devices are connected to the solar panel at the same time.

## 9. Recharging the battery

- **Fully recharge the battery via USB cable and power supply when using it for the first time.**

Solar panel can be recharged in three ways:

1. Expose the front solar panel to direct light (sunlight or ambient light).
2. The green LED indicator will turn on. (Fig. 1)

3. Connect via USB/micro-USB cable to the external USB power supply.
4. Connect via USB/micro-USB cable to a PC USB port.
  - The blue LED indicators will turn on (according to current the battery charge status of the solar panel). (Fig. 2)
  - If the battery is not fully charged one of the blue light will blink indicating the recharging.
  - When the battery is fully charged all the four blue LED will be lit.

## 10. Use of the solar panel

### 10.1 Powering the microscope

1. Connect the USB/micro-USB cable from the solar panel to the microscope.
2. The blue LED will immediately turn on. (Fig. 3)
  - This means that the solar battery is giving power to the microscope.
3. The battery charge level decreases when the microscope is used. Blue LED will turn off as the charge level decreases.

Fig. 3

### 10.2 Using the LED torch light

1. Press twice the switching button ①. (Fig. 4)

Fig. 4

2. The LED torch light ② will turn on. (Fig. 5)
3. Press the button ① to turn off the torch light.

Fig. 5

## 11. Troubleshooting


Review the information in the table below to solve operating problems.

PROBLEM	CAUSE	SOLUTION
After charging for long time, indicators don't show full capacity	Charging time is not enough	Charge it 8 hours continuously
	Power supply or battery pack is faulty	Refer to technical assistance
Microscope light is not working	No power left	Recharge the battery
	Something wrong with cables and connections	Check that USB is well inserted into the output connector of the battery. Check that the jack connector is well inserted into the input of the microscope
	Wrong settings of microscope switches	Check that the main power switch of the microscope is in ON position. Check that the brightness adjustment knob of the microscope is set to the desired value

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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**OPTIKA® S.r.l.**