**ACCESSORIES Series**

# INSTRUCTION MANUAL

	<b>Model</b>
	M-1157

Ver. 1.0 2019

---

## Summary

<b>1. Warning</b>	<b>3</b>
<b>2. Intended use</b>	<b>3</b>
<b>3. Unpacking</b>	<b>3</b>
<b>4. Condenser overview</b>	<b>4</b>
4.1 Condenser	5
4.2 Top lenses	5
4.3 Optical inserts	6
<b>5. Assembling</b>	<b>7</b>
5.1 Assembling the top lens	7
5.2 Installing optical inserts	7
5.3 Use of indicator markers	9
5.4 Installing the condenser	9
5.5 Centering the condenser	9
<b>6. Observation</b>	<b>10</b>
6.1 Brightfield observation (BF)	10
6.2 Darkfield observation (DF)	11
6.3 Phase Contrast observation (PH)	12
6.4 Simple Polarized Light observation (PO)	13
6.5 Differential Interference Contrast Nomarski (DIC)	14
<b>7. Observation modes and usable objectives</b>	<b>16</b>
<b>8. Troubleshooting</b>	<b>17</b>

---

## 1. Warning

This system is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

## 2. Intended use

### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.



## 3. Unpacking

The device is housed in a proper packaging. Take some care to avoid that the optical items fall out and get damaged. Lift the system from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 4. Condenser overview

The universal condenser M-1157 allows observation with different methods: Brightfield, Darkfield, Phase Contrast, Simple Polarized Light, Differential Interference Contrast Nomarski (DIC).


The configuration of the condenser varies according to the chosen observation methods.

Once opened the box, the system parts are the following:

#### 4.1 Condenser

① Condenser

② Magnetick stickers

#### 4.2 Top lenses

① Top lens 0.2  
② Top lens 0.9

③ Top lens 1.4

### 4.3 Optical inserts

See paragraph 7 for usable objectives.

① Phase Contrast inserts

- SL1: PH ring 10X/20X
- SL2: PH ring 40X/60
- SL3: PH ring 100X
- SL4: PH ring 4X

② Darkfield inserts

- DF: darkfield ring "dry"
- DF OIL: darkfield ring "oil"

③ DIC inserts

- U-DIC10: DIC prism 10X
- U-DIC20: DIC prism 20X
- U-DIC40/60: DIC prism 40X / 60X

## 5. Assembling

### 5.1 Assembling the top lens

- Depending on the observation needs, use a top lens dedicated to the observation method or objectives used.
- Fully screw the desired top lens onto the housing. (Fig. 1)
- Avoid excessively forcing when screwing in.**

Top lens	Usable objectives
0.2	2X
0.9	4X - 100X (*)
1.4 (**)	100X

\*) The A.N. with the objectives 60X and 100X may be slightly insufficient without affecting the observation.  
(\*\*) Only for observation with oil.

Fig. 1

### 5.2 Installing optical inserts

- Using the Allen key supplied with the microscope, unscrew the fixing screw ① of the condenser cover. (Fig. 2)
  - Remove the condenser cover to access the housings of the optical inserts.
  - Insert the optical insert into one of the available housings.
- To mount the inserts see the following table.


Fig. 2

#### Only for DIC inserts

- Each condenser housing has a small alignment notch ②. (Fig. 3)
- Align the foot on the underside of the DIC prism ③ (Fig. 4) with the notch on the condenser ②.
- Do not touch the DIC prism during assembly.**



Fig. 3

2. Place the insert. (Fig. 5)
3. Replace the cover and tighten the fixing screw ①. (Fig. 2)

Observation Mode	Insert	Condenser position
PH	SL1, SL2, SL3, SL4	1 - 2 - 3
DF	DF, DF OIL	4, 8
DIC	U-DIC10, U-DIC20, U-DIC 40/60	5 - 6 - 7 (*)

(\*) DIC prisms can be installed also in the positions 4 and 8.


### 5.3 Use of indicator markers

- Once the inserts are installed, use the magnetic markers to quickly identify the position of each single inserto.
1. Align the positioning number of the optical insert installed in the condenser ② with the reference number ① on the cover (Fig. 7) and then apply the magnetic marker.
  - When one of the eight available capacitor positions is inserted, a dot “○” appears inside the corresponding window on the condenser cover. This indicates that the optical insert installed in that position is currently present in the optical path.


### 5.4 Installing the condenser

- Refer to the instruction manual of the microscope.
1. Raise the stage with the coarse focusing knob, avoiding contact with the objectives.
  2. Lower the condenser holder using the condenser height adjustment knob.
  3. If necessary, lower the top lens.
  4. Loosen the condenser locking screw ③. (Fig. 8)
  5. Insert the condenser into the condenser holder and place the alignment pin in the groove on the back of the condenser holder.
  6. Tighten the condenser fixing screw.

### 5.5 Centering the condenser

- Refer to the instruction manual of the microscope.
1. Place the specimen on the stage, insert 10x objective into the light path and focus.
  2. Insert the front lens of the swing-out condenser ①. (Fig. 9)
  3. Rotate the field diaphragm ring ② in counterclockwise direction, to fully close the diaphragm.
  4. Rotate the condenser height adjustment knob ③ to focus the edges of the diaphragm.
  5. Rotate the two centering screws ④ to bring the bright spot in the center of the field of view.
  6. Gradually open the diaphragm. The condenser is centered when the diaphragm image is symmetrical to the field of view. (Fig. 10)
  7. In normal use, open the diaphragm until it circumscribes the field of view.

## 6. Observation

### 6.1 Brightfield observation (BF)

- Usable objectives with different top lenses

Objective	2X	4X	10X	20X	40X	60X	100X
Top lens							
0.2	○	-	-	-	-	-	-
Top lens position	IN		-				
0.9	○	○	○	○	○	○*	○*
Top lens position	OUT		IN				
1.4	-	-	-	-	-	○	○
Top lens position	OUT					IN	

- Remove all filters from the optical path.

- Rotate the condenser turret to enter the “BF” position (no optical insert in the optical path).
  - Pull the lever of the polarizer ① to remove it from the optical path. (Fig. 11)
  - Follow the indications in the table “Usable objectives with different top lenses” to optimize the observation.
- When the top lens is removed, the field diaphragm functions as an aperture diaphragm.
- Place a specimen on the stage.
  - Focus the specimen.
  - Begin the observation.
  - Adjust the aperture diaphragm according to the type of sample to optimize image quality.


Fig. 11

- Using the oil top lens, read the upper scale of the condenser (indicated by “OIL”), while using the dry top lenses, read the lower scale. (Fig. 12)

Fig. 12

## 6.2 Darkfield observation (DF)

- Usable objectives with different inserts

Objective	4X	10X	20X	40X	60X	100X
Insert						
DF	-	○	○	○	-	-
Top lens position	OUT	IN			-	-
DF OIL	-	-	-	-	○	○(*)
Top lens position	OUT	-			IN	

(\*) For observation with a 100X objective in darkfield oil we recommend the use of a 100X oil objective with iris.

- Remove all filters from the optical path.

1. Rotate the condenser turret to enter the "DF" position.
2. Pull the lever of the polarizer ① to remove it from the optical path. (Fig. 11)
3. Place a specimen on the stage and focus the specimen.
4. Observing in the eyepieces, lower or raise the condenser until a homogeneous illumination of the sample is obtained and therefore an optimal dark field effect is achieved.
5. Using the centering screws on the condenser ① (Fig. 13) centre the dark field ring until the field of view is evenly illuminated.

- Darkfield requires a huge amount of light. Switching from darkfield to brightfield, one could be dazzled. Do not keep your eyes on the eyepieces when moving the condenser turret from DF to BF.

- "Dry" darkfield observation, that is, without the use of oil, is only possible with objectives with N.A. lower than 0,7.

- Observing in darkfield, it may be necessary to slightly raise the condenser from the normal position to obtain a more homogeneous illumination. This is not a defect.


Fig. 13




### 6.3 Phase Contrast observation (PH)

- Usable objectives with different phase rings

Objective	4X	10X	20X	40X	60X	100X
Ring						
SL1	-	○	○	-	-	-
Top lens position	OUT	IN		-	-	-
SL2	-	-	-	○	○	-
Top lens position	OUT	-		IN		-
SL3	-	-	-	-	-	○
Top lens position	OUT	-	-	-	-	IN
SL4	○	-	-	-	-	-
Top lens position	IN	-	-	-	-	-

- Remove all filters from the optical path.

- Pull the lever of the polarizer ① to remove it from the optical path. (Fig. 11)
  - Rotate the condenser turret to enter the position relative to the needed phase ring (SL1, SL2, SL3, SL4).
  - Follow the indications in the table “Usable objectives with different phase rings” for objective/ring compatibility.
  - Insert the objective corresponding to the ring into the optical path.
  - Place a specimen on the stage and focus the specimen.
  - Remove one eyepiece and insert the centering telescope. (Fig. 14)
  - Rotate the upper part of the centering telescope until the two phase rings (one dark and one bright) visible in the telescope are in focus. (Fig. 14-15)
  - Using centering screws on the condenser ① (Fig. 13), center the phase rings to make the bright ring ② be concentric to the dark ring ③. (Fig. 16)
  - Repeat the same operation with other objectives to check the ring centering.
  - At the end remove the centering telescope, reinstall the eyepiece and begin observation.
- With 40x and 100x objectives it may be useful to slightly raise the condenser, to obtain a better projection of the phase rings. This is not a defect.

## 6.4 Simple Polarized Light observation (PO)

### • Usable objectives

Objective	2X	4X	10X	20X	40X	60X	100X
0.9	△	○	○	○	○	○	○
Top lens position	OUT		IN				

- For simple polarized light observation, a transmitted light analyzer is required.
1. Rotate the condenser turret to enter the “BF” position (no optical insert in the optical path).
  2. Remove the dummy slider from the nosepiece and insert the analyzer into the dummy slider, then insert the ① assembly into the slot ②. (Fig. 17)
  3. Push the lever of the polarizer ③ to insert it in the optical path. (Fig. 18)
  4. Remove the specimen from the stage.
  5. Loosen the polarizer rotation fixing screw ④. (Fig. 18)
  6. Turn the polarizer wheel ⑤ under the condenser (Fig. 19) to achieve maximum darkening of the eyepieces, and then tighten the polarizer rotation locking screw ④.
  7. Place the sample on the stage and focus.
  8. Begin observation.


Fig. 17


Fig. 18


Fig. 19

## 6.5 Differential Interference Contrast Nomarski Observation (DIC)


### • Usable objectives with different DIC inserts

Objective	10X	20X	40X	60X
Insert				
DIC10	○	-	-	-
DIC20	-	○	-	-
DIC40/60	-	-	○	○


The Nomarski DIC observation in transmitted light requires the kit consisting of the following accessories: Universal condenser ① (containing the dedicated DIC prisms according to the objectives in use), transmitted light Analyzer ②, DIC slider ③. (Fig. 20)

1. Rotate the condenser turret to enter the “BF” position (no optical insert in the optical path).
2. Using the knob ④, insert the built-in polarizer into the condenser and loosen the polarizer rotation fixing screw ⑤. (Fig. 21)

3. Remove the dummy slider from the nosepiece and insert the analyzer into the dummy slider, then insert the ⑥ assembly into the slot ⑦. (Fig. 22)

4. Remove the specimen from the stage.
5. Turn the polarizer knob ⑧ under the condenser to achieve maximum darkening of the eyepieces, and then tighten the polarizer locking screw ⑤. (Fig. 23)

Fig. 23

6. Once the maximum darkening is achieved, remove the slider from the nosepiece, remove the analyzer from the dummy slider and insert it into the DIC prism. Now insert the DIC slider ⑨ into the slot ⑦. (Fig. 24)


Fig. 24

7. Rotate the condenser turret ⑩ to insert the DIC prism matching the objective in use. (Fig. 25)


Fig. 25

8. Put a specimen on the stage and focus.
9. Begin the observation by turning the knob on the DIC slider ⑪ to obtain a three-dimensional effect of the sample. (Fig. 26)

Fig. 26

## 7. Observation modes and usable objectives

Observation mode	Insert	Usable objectives
BF	--	IOS W-PLAN, IOS W-PLAN F, IOS U-PLAN F, IOS U-PLAN F PH, IOS U-PLAN F APO
DF	DF / DF OIL(*)	IOS W-PLAN, IOS W-PLAN F, IOS U-PLAN F, IOS U-PLAN F APO
PH	SL1, SL2, SL3, SL4	IOS U-PLAN F PH
PO	--	IOS W-PLAN, IOS W-PLAN F, IOS U-PLAN F, IOS U-PLAN F PH, IOS U-PLAN F APO
DIC	DIC10, DIC20, DIC40/60	IOS U-PLAN F, IOS U-PLAN F APO

(\*) For observation with a 100X objective in darkfield oil we recommend the use of a 100X oil objective with iris.

## 8. Troubleshooting

Review the information in the table below to solve operating problems.

PROBLEM	CAUSE	SOLUTION
When using 10x to 100x objectives, the field diaphragm image does not appear.	The slide is too thick	Use slides with a thickness not exceeding 1.2 mm
	Top lens is not in	Put the top lens in the light path
In bright field the image is too bright and with low resolution	Aperture diaphragm too closed	Open aperture diaphragm in a proper way
	Top lens is not in	Put the top lens in the light path
Condenser phase ring does not overlap to the objective phase ring	An incorrect insert is inserted in the optical path	Rotate the condenser turret to insert the appropriate optical insert for your objective
	The objective does not match the phase ring in use	Use the objective corresponding to the phase ring
Contrast in darkfield is poor	Top lens is not in	Put the top lens in the light path
	Aperture diaphragm is closed	Fully open the aperture diaphragm
	An incorrect insert is inserted in the optical path	Rotate the condenser turret to insert the appropriate optical insert for your objective
	You're using an inappropriate objective	See paragraph 7. "Observation modes and usable objectives".
	Darkfield ring not well centered	Center the ring correctly
Image quality when observing in polarized light is not high	Polarizer is not inserted in the optical path	Insert the polarizer in the optical path
	Analyzer is not inserted in the optical path	Insert the analyzer in the optical path
	One insert is present in the optical path	Remove the insert from the optical path
While observing in DIC you can not see the typical interference colors	Polarizer is not inserted in the optical path	Insert the polarizer in the optical path
	Analyzer is not inserted in the optical path	Insert the analyzer in the optical path
	DIC prism is not inserted in the optical path	Insert the DIC prism in the optical path
	DIC slider is not inserted in the optical path	Insert the DIC slider in the optical path
	Polarizer and analyzer are not in the total extinction position (crossed Nicol)	Reposition the polarizer
While observing in DIC the typical interference colors are not homogeneous	The condenser is not well centered or is at an incorrect height	Adjust the condenser according to Koehler's setting
	An incorrect insert is inserted in the optical path	Rotate the condenser turret to insert the appropriate optical insert for your objective
	You're using an inappropriate objective	See paragraph 7. "Observation modes and usable objectives".

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**OPTIKA® S.r.l.**