ACCESSORIES Series

# INSTRUCTION MANUAL

Model
M-185

---

## Table of content

1.	Warning	3
2.	Package content	3
3.	Unpacking	3
4.	Intended use	3
5.	Assembling	4
6.	Use of the condenser	5

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

## 2. Package content

① Darkfield Dry condenser

## 3. Unpacking

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 4. Intended use

### Standard models


For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### IVD Models


Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

## 5. Assembling

1. Using the height adjustment knob ① placed on the left side of the condenser holder (please refer to the microscope instruction manual), lower the condenser up to the lowest position achievable. (Fig. 1)



2. Using a 1.5mm Allen wrench, loosen the condenser fixing screw ② placed on the right side of the condenser holder. (Fig. 2)
3. Remove the existing condenser and place the darkfield condenser.
4. Lock the fixing screw ② to secure the condenser in its place.
5. Raise the condenser to the upper position using the height adjustment knob ①.

## 6. Use of the condenser

- “Dry” darkfield observation, that is, without the use of oil, is only possible with objectives with N.A. lower than 0,7 (40x objective).
  - The darkfield observation with 100x oil objective is not possible with this condenser.
1. Insert 10x objective.
  2. Place a specimen on the stage and focus.
  3. Observing into eyepieces adjust the height of the condenser until a homogeneous illumination of the specimen can be achieved, thus obtaining a proper darkfield effect. (Fig. 3)
  4. Using the condenser centering screws ① placed on the condenser holder (Fig. 4) center the condenser so that the illumination is as even as possible over the entire field of view.
  5. Once the setting has been performed, begin the observation.

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**OPTIKA® S.r.l.**