**IM Series**

# INSTRUCTION MANUAL

	<b>Model</b>
	M-797
	M-798

Ver. 1.2    2023

---

## Table of Contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
3.1	M-797	4
3.2	M-798	4
4.	Unpacking	5
5.	Intended use	5
6.	Symbols and conventions	5
7.	Instrument description	6
7.1	M-797	6
7.2	M-798	6
8.	Assembling	7
8.1	Installing the filter holder	7
8.2	Installing the fluorescence filters	7
8.2.1	M-797	7
8.2.1	M-798	8
8.3	Installing the filter covers	8
8.4	Installing the filter selection lever (M-798)	8
8.5	Installing the epi-illuminator	8
8.6	Installing the lamp housing	9
8.7	Installing the bulb	9
9.	Use of the microscope in fluorescence (reflected light)	11
9.1	Centering the mercury bulb	11
9.2	Centering the field diaphragm	13
9.3	Effects of the field diaphragm	13
9.4	Available fluorescence filter cubes	13
10.	Simultaneous observation in Phase Contrast + Fluorescence	14
	Equipment disposal	15

---

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

## 2. Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

### 3. Package content

#### 3.1 M-797

- ① Fluorescence power supply
- ② Lamp housing
- ③ Epi-illuminator
- ④ 2 position filter holder
- ⑤ Fluorescence filters (B-G)

- ⑥ Slider covers
- ⑦ HBO bulb
- ⑧ UV shield
- ⑨ Fluorescence connection cable
- ⑩ Power cord

#### 3.2 M-798

- ① Fluorescence power supply
- ② Lamp housing
- ③ Epi-illuminator
- ④ 4 position filter holder with filters included (B-G)
- ⑤ Filters selector lever

- ⑥ Slider covers
- ⑦ HBO bulb
- ⑧ UV shield
- ⑨ Fluorescence connection cable
- ⑩ Power cord

---

## 4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

## 5. Intended use

### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

## 6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


### ELECTRICAL SHOCK


This symbol indicates a risk of electrical shock.

## 7. Instrument description

### 7.1 M-797


### 7.2 M-798


## 8. Assembling

### 8.1 Installing the filter holder

1. Unscrew the screw ① on both side of the microscope to remove the plastic covers ②. (Fig. 1)

2. Unscrew the three screws ③ holding the plate on the microscope body. (Fig. 2)

3. Insert the filter holder (either 2 positions or 4 positions) in the microscope, with the "L" facing the back side of the microscope. (Fig. 3)
4. Fix the filter holder by screwing the three screws ③, removed before from the plate.


### 8.2 Installing the fluorescence filters

#### 8.2.1 M-797

Insert the filter cubes by sliding them from the sides of the microscope: B from the left and G from the right. (Fig. 4)


- Pay attention to install the filter cubes with the emission filter facing down.

### 8.2.1 M-798

Insert the filter holder by sliding from the side of the microscope. The insertion side is not relevant. (Fig. 5)

- Pay attention to install the filter holder with the emission filter facing down.

### 8.3 Installing the filter covers

1. Once the fluorescence filters are in place, cover the filter holder with the provided covers.
- For M-797: the two cover plates have the same shape and can be mounted either on the left or on the right.
  - For M-798: install the cover with the printed encryption on the left side of the microscope. (Fig. 6)

### 8.4 Installing the filter selection lever (M-798)

Screw the filter selection lever inserting the lever into the hole of the left cover. (Fig. 7)

### 8.5 Installing the epi-illuminator

1. Remove the black plastic cover from the microscope rear.
2. Insert the fluorescence illuminator from the back. In order to ease the insertion, just tilt the assembly at about 45° and move it forward.
3. Fix it using the 3 provided Allen screws. (Fig. 8)


## 8.6 Installing the lamp housing

Insert the lamp house and fix it with the Allen screw ①. (Fig. 9)


## 8.7 Installing the bulb

1. Open the lamp housing using the door lock screw ② and remove the lamp holder. (Fig. 10)


2. Remove the plastic block ③ from the lamp holder (or the exhausted lamp in case of replacement) by loosening the two locking screws ④. (Fig. 11)

3. Insert the mercury bulb ⑤ (respect the polarity of the bulb), tighten the locking screws and refit the lamp holder inside the lamp housing. (Fig. 12)

- **Do not touch the bulb with bare hands.**

4. Connect the cable from the external power supply to the lamp housing. (Fig. 13)

5. Connect the power cable to the external power supply. (Fig. 14)


- The input voltage for fluorescence power supply is 110-240Vac.
- Please use the standard power cable provided. Select suitable one when missing or damaged.
- Connect the power supply correctly, be sure to have a good earth connection.
- Before connecting the power cord, secure the lamp housing cable to the power supply.
- If the power cord is connected before, there may be a risk of electrical shock.

- Disconnect all electrical cables before installing or replacing the bulb.
- The bulb has an anode and a cathode of different sizes. Respect the polarity during assembly.

6. In order to prevent possible damages from UV radiation, mount the orange protection screen as shown. (Fig. 15)


## 9. Use of the microscope in fluorescence (reflected light)

### 9.1 Centering the mercury bulb


- **Wait around 5 minutes before proceeding with this operation to allow the bulb to properly warm up.**
1. Turn on the HBO mercury bulb by operating the power supply switch ①. (Fig. 16)
  2. Turn the nosepiece into an empty position (without objectives) and remove the protective cap, or remove an objective from the nosepiece.


3. Place a piece of white paper on the stage and insert the filter cube "B" into the optical path. (Fig. 17)



4. Acting on the focus screw of the collector lens ② and on the centering screws ③ try to obtain the light spot of the bulb's arc. (Fig. 18-19)


5. Using the focus screw of the collector lens ②, put the image of the arc projected onto the paper. The light spot must be brighter and sharper as possible. (Fig. 19)

6. Using the centering screws ③ on the side of the lamp housing, center the image of the arc. (Fig. 20-21)

7. Using the focusing screw of the collector lens ② enlarge the image until a homogeneous illumination is achieved. (Fig. 22). At this point, insert an objective into the optical path and, looking into the eyepieces, optimize the illumination always using the screws ② and ③.

8. After replacing the exhausted bulb, reset the time counter on the power supply by pressing the "Reset" button ①. (Fig. 23)

## 9.2 Centering the field diaphragm

1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Rotate the field diaphragm lever ①, to fully close the diaphragm. (Fig. 24)
3. Rotate the two centering screws ② to bring the bright spot in the center of the field of view.
4. Gradually open the diaphragm. The illuminator is centered when the diaphragm image is symmetrical to the field of view.
5. In normal use, open the diaphragm until it circumscribes the field of view.


Fig. 24

## 9.3 Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image.

Set the diaphragm according to the objective in use until it circumscribes the field of view, in order to eliminate unnecessary light to eyepieces. (Fig. 25)


Fig. 25

## 9.4 Available fluorescence filter cubes

### • M-797

FILTER NAME	EXCITATION FILTER	DICHROIC MIRROR	BARRIER FILTER	APPLICATIONS
B	460-490 nm	495 nm	520LP nm	<ul style="list-style-type: none"> <li>• FITC: fluorescent antibodies</li> <li>• Acridine orange: DNA, RNA</li> <li>• Auramine</li> </ul>
G	540-580 nm	585 nm	590LP nm	<ul style="list-style-type: none"> <li>• Rhodamine, TRITC: fluorescent antibodies</li> <li>• Propidium iodide: DNA, RNA</li> <li>• RFP</li> </ul>

### • M-798

FILTER NAME	EXCITATION FILTER	DICHROIC MIRROR	BARRIER FILTER	APPLICATIONS
B	460-490 nm	500 nm	520LP nm	<ul style="list-style-type: none"> <li>• FITC: fluorescent antibodies</li> <li>• Acridine orange: DNA, RNA</li> <li>• Auramine</li> </ul>
G	527-553 nm	565 nm	575LP nm	<ul style="list-style-type: none"> <li>• Rhodamine, TRITC: fluorescent antibodies</li> <li>• Propidium iodide: DNA, RNA</li> <li>• RFP</li> </ul>
UV	325-375 nm	415 nm	435LP nm	<ul style="list-style-type: none"> <li>• Nuclear counterstaining</li> </ul>
V	390-420 nm	440 nm	455LP nm	<ul style="list-style-type: none"> <li>• Acridine orange: DNA, RNA</li> </ul>

## 10. Simultaneous observation in Phase Contrast

---


## + Fluorescence

- **Fluorescence models allow observation in transmitted light Phase contrast in combination with reflected light Fluorescence. Samples with rapid decay must first be observed in Fluorescence and then in Phase Contrast. The combined observation allows you to easily identify some areas of the sample that emit fluorescence.**
1. Turn on the power supply for the HBO fluorescent bulb and wait 5 minutes before the arc stabilizes.
  2. Move the filter selector to an empty position or, if the filter holder is full, to the position containing the UV filter.
  3. Insert the desired PH objective and move the phase contrast slider to the position containing the corresponding phase ring.
  4. Focus the sample.
  5. Adjust the light intensity of the transmitted light.
  6. Move the fluorescence filter selector to the desired position.
  7. To obtain the proper observation of the sample, adjust the light intensity of the transmitted light to modulate the intensity of the fluorescence with the one of the phase contrast.

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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**OPTIKA® S.r.l.**