SFX/STEREO Series

INSTRUCTION MANUAL

Model
SFX-85V

Ver. 1.0 2024

Table of Contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
4.	Unpacking	5
5.	Intended use	5
6.	Symbols and conventions	5
7.	Instrument description	6
8.	Assembling	7
9.	Use of the microscope	8
9.1	Switching on the microscope	8
9.2	Turning on the light	8
9.4	Placing the specimen	8
9.4	Focus	8
9.5	Magnification	9
9.6	Contrast disc	9
10.	Use of the software	10
10.1	Function menu	10
10.1.1	Camera control	10
10.1.2	Measurement	13
10.1.3	Calibration	14
11.	Micrometric Slide M-005	16
12.	Maintenance	17
13.	Troubleshooting	18
	Equipment disposal	19

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

- ① Microscope body
- ② Black/white plastic stage plate
- ③ Glass stage plate
- ④ Micro SD card
- ⑤ Wireless mouse
- ⑥ Dust cover
- ⑦ Micrometric slide
- ⑧ Power supply

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

7. Instrument description

8. Assembling

No installation procedure is required for this model.

1. Insert the SD card in the slot ① until it clicks in place. (Fig. 1)

2. Plug the mouse receiver in the USB port ② placed in the bottom side of the monitor. (Fig. 2)


3. Connect the power supply to the microscope and to the wall socket.
 - **The connector on the monitor has no function in this configuration. The monitor turns on automatically when the microscope is turned on.**

9. Use of the microscope

9.1 Switching on the microscope

1. Put the main switch ① in the ON position ("1"). (Fig. 3)
2. This turns on both the microscope's illumination system and the internal software.
3. The monitor turns on. The start-up of the internal software takes a few seconds.
4. At the end of the start-up procedure the image that the microscope is displaying appears.

9.2 Turning on the light

1. Repeatedly press one of the touch button ② in order to change the light intensity.
 - "INC" button controls the incident illumination, while "TRA" button controls the transmitted illumination.
2. It is possible to use both illuminations at the same time.
 - LED ③ indicates the power status of the microscope. (Fig. 4)
 - LED ON (green): one of the illumination (INC or TRA) is on.
 - LED OFF (no color): no illumination (INC or TRA) is on.

9.3 Placing the specimen

1. Place the specimen on the microscope stage and lock it using the slide-clamp if it is necessary.
 - Ensure that the specimen is centered over the stage opening. (Fig. 5)

9.4 Focus

1. Rotate the focusing knob ④ to bring the sample in focus. (Fig. 6)

9.5 Magnification

1. Select the desired magnification by rotating the zoom changer ⑤. (Fig. 7)

9.6 Contrast disc

- Two discs are supplied: a frosted glass disc and a white/black plastic disc.
 - For transmitted light applications the glass disc is used, while for incident light applications the black/white disc is used.
1. If observing dark objects place the disk with the white side facing up.
 2. If observing bright objects place the disk with the black side facing up. (Fig. 8)

10. Use of the software

10.1 Function menu

1. Move the mouse to the top and bottom of the display screen, the function menu interface will appear.
- **Note: the key switch on the right side of the monitor is used to turn the display on or off, not the microscope on or off.**
 - **Note: Before using a new screen or function, you must first close the previous screen or function.**

10.1.1 Camera control

Move the mouse to the lower part of the screen to pop up a main interface menu as shown below.

Icon	Function
	Toggle On/Off Parameter adjustment
	Camera parameters
	Image adjustment
	View storage media and software version
	Turn off the monitor
	Snap an image
	Record a video
	Preview
	Zoom in
	Zoom out
	Flip horizontal
	Flip vertical
	Color / black and white image
	Freeze
	Crosshair
	Comparison of two images
	Comparison of four images
	HDR (High Dynamic Range)
	White balance
	Exposure setting
	Factory data reset
	Language selection
	Date / Time setting


Clicking on the  button, user enters in the Camera Parameters adjustment panel.

Here it is possible to set the several parameters for the live image.

Parameter	Description
Auto Exposure	When “Auto” is checked, the system will automatically adjust exposure time according to the value of Exposure Target
Exposure Target	Available when “Auto” is checked. Slide to left or right to adjust Exposure Compensation according to the current video brightness to achieve proper brightness value
Exposure (ms)	Available when “Auto” is unchecked. Slide to left or right to reduce or increase exposure time, adjusting brightness of the video
Auto WB	When the ‘Auto’ button is selected, the camera continuously performs an automatic white balance Click on the ‘one push’ button and the camera needs a reference white sample to perform an optimal white balance
R_GAIN	Slide to left or right to decrease or increase the proportion of Red in RGB on video
G_GAIN	Slide to left or right to decrease or increase the proportion of Green in RGB on video
B_GAIN	Slide to left or right to decrease or increase the proportion of Blue in RGB on video
Temp Color	Move the slider to change the colour temperature of the light
Power Frequency	Check “50Hz” or “60Hz” to eliminate flickering “strap” caused by 50Hz or 60Hz illumination

Clicking on the  button, user enters in the Image adjustment panel

Parameter	Description
Contrast	Adjust Contrast level of the image. Slide to the right side to increase contrast and to the left to decrease contrast
Saturation	Adjust Saturation level of the image. Slide to the right side to increase saturation and to the left to decrease saturation
Gamma	Adjust Gamma level of the image. Slide to the right side to increase gamma and to the left to decrease gamma
Sharpness	Adjust Sharpness level of the image. Slide to the right side to increase gamma and to the left to decrease sharpness
Denoise	Adjust Noise level of the image Slide left or right to denoise the image
Flags Reset	A series of grids are pre-set and will appear in overlay. Click to activate or modify them according to your needs

10.1.2 Measurement

Move the mouse to the top of the screen to pop up a main interface menu as shown below.

Icon	Function	Explanation
	<i>Point drawing</i>	Draw points and display the coordinates of points. Click to confirm.
	<i>Free line</i>	Click to determine the starting point, then drag the mouse to the end point to be measured. Click a second time to determine the end point. Click a third time to confirm.
	<i>Horizontal line</i>	
	<i>Vertical line</i>	
	<i>Arrow</i>	Click to determine the starting point. Drag the mouse to determine the length, inclination and direction of the arrow. Click a second time to determine the end point. Click a third time to confirm the position of the arrow.
	<i>Point spacing</i>	Click to determine the starting point, then drag the mouse to the end point to be measured. Click a second time to determine the end point. Click a third time to confirm.
	<i>Parallel line</i>	Click to draw a first line and drag the mouse. Click a second time to finish drawing the first line. Move the mouse to the bottom of the line and click a third time to draw the second line. Finish drawing the second line with a fourth click. Click a fifth time to confirm.
	<i>Perpendicular line</i>	Click to draw a first line and drag the mouse. Click a second time to finish drawing the first line. Now move the mouse and you will get the line perpendicular to the first drawn line. Click a third time to draw the end of the perpendicular line. Click a fourth time to confirm
	<i>Angle</i>	Click to draw a first line and drag the mouse. Click a second time to finish drawing the first line. Now move the mouse and you will get the oblique line that determines the angle with the vertex starting from the end of the first drawn line. Click a third time to draw the angle. Click a fourth time to confirm
	<i>Rectangle</i>	Click to determine the first vertex of the rectangle. Drag the mouse and click a second time to finish drawing the rectangle. Click a third time to confirm.
	<i>Polygon</i>	Click as many times as there are sides of the polygon to be measured. Click one last time to close the polygon. Click again to confirm.
	<i>Arc</i>	Click to draw a first line and drag the mouse. Click a second time to finish drawing the first line. Now move the mouse to determine the arc to be measured. Click a third time to finish drawing the arc. Click a fourth time to confirm
	<i>Circle (Radius)</i>	Click in the centre of the circle to be measured and drag the mouse to trace the circle. Click a second time to finish drawing. Click a third time to confirm
	<i>Circle (Two points)</i>	Click on the circumference of the circle to be measured and drag the mouse to trace the circle. Click a second time to finish drawing. Click a third time to confirm
	<i>Circle (Three points)</i>	Click three times on the circumference of the circle to be measured and drag the mouse to trace the circle. Click a fourth time to finish drawing. Click a fifth time to confirm.
	<i>Distance between line and circle</i>	Click three times on the circumference of the circle to be measured and drag the mouse to trace the circle. Click a fourth time to finish drawing. Now drag the mouse to the line whose distance is to be measured and click a fifth time. Click a sixth time to confirm


	<i>Concentric circles</i>	Click three times on the circumference of the circle to be measured and drag the mouse to draw the circle. Click a fourth time to finish drawing. Now drag the mouse to draw a second circle concentric to the first and click a fifth time. Click a sixth time to confirm
	<i>Center distance</i>	Click three times on the circumference of the first circle to be measured and drag the mouse to trace the circle. Click a fourth time to finish plotting. Move the mouse to the second circle and make three clicks on the circumference. Click to finish drawing the second circle. Click a final time to confirm
	<i>Text</i>	Click to draw the text box. Click inside the text box to activate the “virtual keyboard” and click with the mouse to type the desired text. Click Enter and the text is displayed
	<i>Color</i>	Select the color of the line. There are 8 colors to choose
	<i>Thickness</i>	Select the thickness of the line. There are 4 sizes to choose
	<i>Clear screen</i>	Click to delete all existing overlays
	<i>Delete all</i>	Click to delete all items on the screen
	<i>Screen capture</i>	Save the screen capture to the storage space
	<i>File output</i>	Save the file to storage media
	<i>Calibration</i>	Enter the calibration menu

10.1.3 Calibration

In order to make measurements, it is necessary to calibrate the system using a micrometer slide.

Note: Calibration is required for each available magnification on the microscope


Once the calibration button  is selected, the calibration dialog box appears:

1. Click on the  icon to add a calibration name, and then click on the box beside the “Name” field to assign a name to the the calibration item, usually named using the current magnification (e.g. 0.7x).
2. Click the blank box to pop up the “virtual keyboard”, as shown below.


3. When assigning the calibration name, click the dialog box to fill in the required content.
4. After confirming, press "Enter" and then press "Close" to close the dialog box.
5. Click the blank box beside "Pixel" field, and press the left mouse button on the screen image to draw the calibration line on the micrometer scale.
 - The calibration line should be an integer as far as possible to facilitate calculation. (In principle, within the scope of visibility and convenient operation, the larger the span of the calibration line is, the better, in order to eliminate the errors generated in the operation).
6. The following figure shows the calibration line selected during 4X objective observation.
7. Click one point to confirm the starting point and click again to confirm the end.

8. When the calibration line is 458.00px, release the mouse and this parameter will be automatically filled in the blank box beside the "Pixel" field.
9. The actual length of the calibration line is 1000 μm , so enter 1000 in the blank beside "Length".
10. The unit is set to " μm ".
11. Now, the blank field beside "Scale" will be calibrated to automatically calculate and display the proportion based on the actual length and pixel length.

11. Micrometric Slide M-005

Micrometric slide, 26x76mm, with 2 scales
(1mm/100div. for biological microscopes / 10mm/100div. for stereo microscopes)

1 DIV=0.01mm

For biological microscopes calibration

1 DIV=0.1mm

For stereo microscopes calibration

12. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in OFF position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.


13. Troubleshooting

Review the information in the table below to solve operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
The illumination is ON, but the field of view is dark.	The plug is not connected to the illumination	Connect the cable
	The brightness is too low	Adjust to a proper setting
The edge of the field of view is vignetted or the brightness is asymmetric.	The incident illuminator is not correctly oriented	Change the angle of the incident illuminator
Dust and stains can be seen in the field of view.	There are stains and dust on the specimen	Clean the specimen
	There are stains and dust on the eyepiece	Clean the eyepiece
Poor image quality: <ul style="list-style-type: none">• The image is not sharp• The contrast is not high• The details are not clear• Image glares	The lenses (additional lens, objective) are dirt	Thoroughly clean all the optical system
One side of the image is out of focus.	The specimen is out of place (tilted)	Place the specimen flat on the stage.
II. Mechanical Section:		
The focus knob is hard to turn.	The tension adjustment collar is too tight	Loosen the tension adjustment collar
The focus is unstable.	The tension adjustment collar is too loose	Tighten the tension adjustment collar
III. Electric Section:		
The LED doesn't turn on.	No power supply	Check the power cord connection
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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
The screenshot shows a 'Calibration' dialog box with the following elements:

- Title:** Calibration
- Buttons:** A '+' button, a trash icon, a 'Show' button, a 'Hide' button, and a 'Confirm' button.
- Input Fields:**
 - Name:** A text input field.
 - Length:** A text input field followed by a dropdown arrow.
 - Pixel:** A text input field with the label 'Px' to its right.
 - Scale:** A text input field.
 - Pos:** A text input field followed by a dropdown arrow.

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OPTIKA® S.r.l.