**SZP Series**

# INSTRUCTION MANUAL

Model
SZP-6
SZP-6E
SZP-8
SZP-8E
SZP-10
SZP-10E
ST-172

Ver. 2.0 2020

---

## Table of Contents

<b>1. Warning</b>	<b>3</b>
<b>2. Symbols and conventions</b>	<b>3</b>
<b>3. Safety Information</b>	<b>3</b>
<b>4. Intended use</b>	<b>3</b>
<b>5. Instrument description</b>	<b>4</b>
5.1 SZP-6 / SZP-8 / SZP-10	4
5.2 SZP-6E / SZP-8E / SZP-10E	5
5.3 ST-172	6
<b>6. Unpacking</b>	<b>7</b>
<b>7. Assembling</b>	<b>7</b>
7.1 SZP-6 / SZP-8 / SZP-10	7
7.2 SZP-6E / SZP-8E / SZP-10E	8
7.3 ST-172	8
7.4 Assembling procedure	9
<b>8. Using the microscope</b>	<b>11</b>
8.1 Adjusting interpupillary distance	11
8.2 Dioptic compensation	11
8.3 Magnification	11
8.4 Use of additional lens	12
8.5 Use of aperture diaphragm (ST-172)	12
<b>9. Maintenance</b>	<b>14</b>
<b>10. Troubleshooting</b>	<b>15</b>
<b>Equipment disposal</b>	<b>16</b>

---

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

## 2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### **CAUTION**

This symbol indicates a potential risk and alerts you to proceed with caution.

### **ELECTRICAL SHOCK**

This symbol indicates a risk of electrical shock.

## 3. Safety Information

### **Avoiding Electrical Shock**

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

## 4. Intended use

### **Standard models**


For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### **IVD Models**


Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

## 5. Instrument description

### 5.1 SZP-6 / SZP-8 / SZP-10

## 5.2 SZP-6E / SZP-8E / SZP-10E

5.3 ST-172

## 6. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 7. Assembling

Once opened the box, the microscope parts are the following:

### 7.1 SZP-6 / SZP-8 / SZP-10

- ① Microscope body
- ② Objective
- ③ Observation head

- ④ Eyepieces
- ⑤ Allen wrench
- ⑥ Dust cover

## 7.2 SZP-6E / SZP-8E / SZP-10E

- ① Microscope body
- ② Objective
- ③ Observation head

- ④ Eyepieces
- ⑤ Allen wrench
- ⑥ Dust cover

## 7.3 ST-172

- ① Aperture diaphragm

## 7.4 Assembling procedure

1. Insert the microscope body in the round hole of the stand.
2. Screw the objective in the bottom part of the microscope body. (Fig. 1)

3. Install the aperture diaphragm attachment ST-172 (if provided) by aligning the positioning pin ② and the positioning groove ③.
4. Lock the fixing screw with the provided Allen wrench. (Fig. 2-3)

5. Install the head (Fig. 3) by aligning the positioning pin ② and the positioning groove ③. (Fig. 2)
6. Lock the fixing screw. (Fig. 4)

---

7. Insert eyepieces into empty eyepieces sleeves. (Fig. 5)

## 8. Using the microscope


### 8.1 Adjusting interpupillary distance

1. Hold the right and left eyepiece tube with both hands and adjust the interpupillary distance by moving the two parts until one circle of light can be seen. (Fig. 6)
2. If two circles appear, the interpupillary distance is too big, and if two overlapped circles appear, the interpupillary distance is too small. (Fig. 6)


### 8.2 Dioptic compensation

- **This compensation makes possible for people wearing glasses to adjust the microscope to their eyes and use the microscope without glasses.**
1. Put the zoom down to the lowest magnification ① and focus the specimen with the focusing knobs ②. (Fig. 7)
  2. Put the zoom to the maximum magnification and repeat the focusing.
  3. Return to the lowest magnification: the specimen will be out of focus.
  4. Adjust the diopter compensation ring of the right eyepiece ③ (Fig. 8) until the image of the right eyepiece is clear and sharp. Repeat the procedure for the left eyepiece ④.
  5. Check the focus of the image for the whole zoom range. It should be perfectly parfocal (focus is maintained during the change of magnification).

### 8.3 Magnification

1. Select the desired magnification by adjusting the zoom knob ⑤. (Fig. 9)
- Change the eyepieces and/or add an appropriate additional lens if a different magnification range is needed.
  - The microscope body is equipped with a “click-stop” function that allows to obtain a precise setting of the desired magnification.

## 8.4 Use of additional lens

1. Unscrew the 1X objective and screw the desired additional lens on the microscope body. (Fig. 1)

Total magnification used can be calculated as:  
 Eyepiece magnification \* Zoom magnification \* Objective lens magnification.

### SZP-6 / SZP-6E

Eyepiece		10x		15x		20x	
Field number (mm)		24		15		10	
Objective	Total Mag	F.O.V (mm)	Total Mag	F.O.V (mm)	Total Mag	F.O.V (mm)	
0.3X	2.4X - 15X	100 - 16	3,6x - 22,5X	62,5 - 10	4,8x - 30X	41,7 - 6,7	
0.5X	4X - 25X	60 - 9,6	6x - 37,5X	37,5 - 6	8x - 50X	25 - 4	
1X	8X - 50X	30 - 4,8	12x - 75X	18,8 - 3	16x - 100X	12,5 - 2	
2X	16X - 100X	15 - 2,4	24x - 150X	9,4 - 1,5	32x - 200X	6,25 - 1	

### SZP-8 / SZP-8E

Eyepiece		10x		15x		20x	
Field number (mm)		24		15		10	
Objective	Total Mag	F.O.V (mm)	Total Mag	F.O.V (mm)	Total Mag	F.O.V (mm)	
0.3X	2,4X - 19,2X	100 - 12,5	3,6X - 28,8X	62,5 - 7,8	4,8X - 38,4X	41,7 - 5,2	
0.5X	4X - 32X	60 - 7,5	6X - 48X	37,5 - 4,7	8X - 64X	25 - 3,1	
1X	8X - 64X	30 - 3,75	12X - 96X	18,8 - 2,3	16X - 128X	12,5 - 1,5	
2X	16X - 128X	15 - 1,9	24X - 192X	9,4 - 1,2	32X - 256X	6,25 - 0,8	

### SZP-10 / SZP-10E

Eyepiece		10x		15x		20x	
Field number (mm)		24		15		10	
Objective	Total Mag	F.O.V (mm)	Total Mag	F.O.V (mm)	Total Mag	F.O.V (mm)	
0.3X	2,4x - 24x	100 - 10	3,6x - 36x	62,5 - 6,25	4,8x - 48x	41,7 - 4,2	
0.5X	4x - 40x	60 - 6	6x - 60x	37,5 - 3,75	8x - 80x	25 - 2,5	
1X	8x - 80x	30 - 3	12x - 120x	18,8 - 1,87	16x - 160x	12,5 - 1,25	
2X	16x - 160x	15 - 1,5	24x - 240x	9,45 - 0,94	32x - 320x	6,25 - 0,625	

## 8.5 Use of aperture diaphragm (ST-172)

By varying the iris aperture, the observer can adjust the brightness of the image for better contrast and greater depth of field.

- Over-closing the iris causes the image quality to deteriorate.
1. Move the aperture diaphragm lever ① to the right or left. Rotation towards the symbol ○ opens the diaphragm while rotation towards the ⊗ symbol closes it.
  2. Adjust the diaphragm during observation to optimize the effect on contrast and depth of field.
- **Do not close the diaphragm too tightly so as not to have a poor image quality and not to have an excessively dark image.**
3. The scale on the diaphragm serves as a reference for future use.

---

## 9. Maintenance

### Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

### To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

### Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

**For the best results, use the OPTIKA cleaning kit (see catalogue).**

If you need to send the microscope to Optika for maintenance, please use the original packaging.

## 10. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
<b>I. Optical Section:</b>		
The edge of the field of view is vignetted or the brightness is asymmetric	The incident illuminator is not correctly oriented	Change the angle of the incident illuminator
Dust and stains can be seen in the field of view	Stain and dust on the specimen	Clean the specimen
	Stain and dust on the eyepiece	Clean the eyepiece
Poor image quality: <ul style="list-style-type: none"> <li>The image is not sharp</li> <li>The contrast is not high</li> <li>The details are not clear</li> <li>Image glares</li> </ul>	The lenses (objective, eyepieces) are dirt	Thoroughly clean all the optical system
One side of the image is out of focus	The specimen is out of place (tilted)	Place the specimen flat on the stage
<b>II. Observation tube:</b>		
The field of view of the two eyes is different	The interpupillar distance is not correct	Adjust the interpupillar distance
	The dioptic correction is not right	Adjust the dioptic correction
	The viewing technique is not correct, and the operator is straining the eyesight	When look into the eyepieces, do not stare at the specimen but look at the whole field of view. Periodically, move the eyes away to look at a distant object, then back into the eyepieces

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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**OPTIKA® S.r.l.**