IM-5 Series

INSTRUCTION MANUAL

Model
IM-5MET

Ver. 1.5 2025

Table of contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
4.	Unpacking	5
5.	Intended use	5
6.	Symbols and conventions	5
7.	Instrument description	6
8.	Assembling	8
8.1	Installing the objectives	8
8.2	Use of nosepiece adapter rings	8
8.3	Removing the stage lock	9
8.4	Installing the stage insert	9
8.5	Installing the eyepieces	10
8.6	Installing the lamp housing	10
8.7	Installing the bulb	11
8.8	Connecting the power cord	11
9.	Brightfield observation procedures	13
10.	Use of the microscope in Brightfield (BF)	14
10.1	General switch on	14
10.2	Brightness adjustment	14
10.3	Switching Brightfield observation	14
10.4	Adjust the interpupillary distance	15
10.5	Diopter adjustment	15
10.6	Use of eyeshields	15
10.7	Light path selection	16
10.8	Coarse focus tension adjustment	16
10.9	ECO button	16
10.10	Centering the field diaphragm (FS)	17
10.11	Effects of the field diaphragm	17
10.12	Centering the Aperture Diaphragm (AS)	17
10.13	Use of color filters	18
11.	Use of the microscope in Darkfield (DF)	18
11.1	Switching Darkfield observation	18
12.	Use of the microscope in Polarized light (PO)	19
13.	Microphotography	20
13.1	Use of C-mount cameras	20
13.2	Use of Reflex cameras	20
14.	Maintenance	21
15.	Troubleshooting	22
	Equipment disposal	24

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.


2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

- ① Microscope body
- ② Eyepieces
- ③ Stage
- ④ Objectives
- ⑤ Lamp housing
- ⑥ Daylight filter (LBD)
- ⑦ Analyzer
- ⑧ Polarizer
- ⑨ Filter holder

- ⑩ DIC dummy slider
- ⑪ Power supply
- ⑫ Power cord
- ⑬ Revolving nosepiece
- ⑭ Nosepiece adapter rings
- ⑮ Dust cover
- ⑯ Stage clips
- ⑰ Allen wrench
- ⑱ Metal insert for stage

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

7. Instrument description


Opposite side

8. Assembling

8.1 Installing the objectives

1. Turn the coarse focusing knob ① until the nosepiece reaches its lowest position. (Fig. 1)
- **For a safe transport, the nosepiece is placed in the lowest position and the tension adjustment collar ② is adjusted to the appropriate tension when the microscope leaves the factory. (Fig. 1)**
2. Screw the lowest magnification objective on to the nosepiece from the right side, then turn the turret clockwise. Install the other objectives in the same way, following the sequence from low to high magnification.
- **Note: the objectives can also be installed through the stage opening. (Fig. 2)**
- Clean the objectives regularly. In inverted microscopes, the objectives are very sensitive to dust.
- To prevent dust and contamination from entering the microscope, cover all the unused holes with dust caps ③. (Fig. 3)
- When operating, use the low magnification objective (5x or 10X) to search and focus the specimen, then switch to higher magnifications.

8.2 Use of nosepiece adapter rings

The nosepiece can allocate objectives for brightfield/darkfield observation (BD). However, it is possible to install objectives for brightfield using the adapter rings.

1. Screw the adapter ring into the empty hole of the nosepiece. (Fig. 4 - 5)

Fig. 5

2. Screw the objective into the adapter ring. (Fig. 6)
3. Perform the observation as usual.

- The brightfield objective cannot be used for darkfield, therefore switching the mirror cube selector to “DF” could bring to a completely dark image.
- **Using the adapter rings it could happen an imperfect parfocality between brightfield and brightfield/darkfield objectives. This is not a defect.**

Fig. 6

8.3 Removing the stage lock

1. Remove the 4 screws of the metal retaining plate ①, using the supplied Allen screwdriver. (Fig. 7)


Fig. 7

8.4 Installing the stage insert

- Make sure that the plate is horizontal.
1. Install the stage insert in the stage opening. (Fig. 8)

Fig. 8

8.5 Installing the eyepieces

Insert both eyepieces into the tubes of the optical head. (Fig. 9)

8.6 Installing the lamp housing

1. Insert the lamp house and tighten the provided Allen wrench ①. (Fig. 10)

2. Slide the filter holder into the tube in front of the lamp housing. (Fig.11)

3. Connect the lamp housing connector to the socket placed in the back of the microscope. (Fig. 12)


- Check for a safe connection.

8.7 Installing the bulb

- Before installing the bulb, set the main switch to “OFF” and disconnect all electrical cables.

1. Using the provided Allen wrench, release the lamp housing cover ②. (Fig. 13)


2. Press the spring clips ③ of the lamp holder and insert the halogen bulb into the holes ④. (Fig. 14)

3. Release the spring clips to lock the bulb in place

- **Do NOT touch the bulb with the bare hands: this could reduce the lamp efficiency and lifetime.**


4. Once installed the bulb, put the lamp housing cover into its original position and lock the screw ②.

8.8 Connecting the power cord

1. Turn the main switch ① to “OFF” before connecting the power cord. (Fig. 15)


2. Insert the plug of the provided power supply into the power socket of the microscope. (Fig. 16)

3. Plug the power cord into the socket of the power supply. (Fig. 17)
4. Plug the power cord into the mains socket.
 - Check for a safe connection.
 - **Please use the provided power cord.**
 - **If lost or damaged, please refer to qualified service.**
 - **Connect the power cord to a grounded power supply only.**


9. Brightfield observation procedures

10. Use of the microscope in Brightfield (BF)


10.1 General switch on

To activate the reflected light illuminator, put the main switch ①, located on the right side of the stand, to the position "ON". (Fig. 18)

10.2 Brightness adjustment

1. Turn the brightness adjustment knob ② to increase and decrease the brightness. (Fig. 19)


2. While rotating the brightness adjustment knob, the LED indicator on the front panel ③ will turn ON or OFF some segments. (Fig. 20)

10.3 Switching Brightfield observation

Move the lever ④ in "BF" position for brightfield illumination. (Fig. 21)

10.4 Adjust the interpupillary distance

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view match.

- **The graduation on the interpupillary distance indicator ①, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator’s eyes. (Fig. 22)**

The range of the interpupillary distance is 48-75 mm.

10.5 Diopter adjustment

- **Both eyepieces are provided with diopter adjustment ring. Put each ring in the “0” position before proceeding with the adjustment.**
 1. Look into the right eyepiece with your right eye only, and focus on the specimen.
 2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the diopter adjustment ring ② to compensate. (Fig. 23)
- **The adjustment range is ± 5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator’s diopter correction.**

10.6 Use of eyeshields

- **Use with eyeglasses**

Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses. (Fig. 24)

- **Use without eyeglasses**

Raise eyeshields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation. (Fig. 25)

10.7 Light path selection

- The observation head is equipped with an optical path selector that allows the light to be distributed to the eyepieces and to the photo / TV port.
1. Move the selector ① to the left (In) or to the right (Out) to distribute the light. (Fig. 26)

POSITION	LIGHT
IN	100% TV
OUT	50% EYEPIECES - 50% TV

Fig. 26

10.8 Coarse focus tension adjustment

The tension of the coarse focusing knob is factory preset.

1. To modify the tension according to personal's needs, rotate the ring ②. (Fig. 27)
- Clockwise rotation increases the tension.
 - If the tension is too loose, the revolver could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

Fig. 27

10.9 ECO button


1. Press the ECO button (Fig. 28) to activate the "ECO" function.
 - Once activated, the system will automatically shut-off after 20 minutes from the activation.
2. To deactivate the function, press the ECO button again.

Fig. 28

10.10 Centering the field diaphragm (FS)


1. Move the selector on the reflected light illuminator in the position "BF". (Fig. 21)
2. Place the specimen on the stage, insert 10x objective into the light path and focus.
3. Completely pull out the field diaphragm lever ① to fully close the diaphragm. (Fig. 29)
4. Using the provided Allen wrench use the two centering screws ② to bring the bright spot in the centre of the field of view. (Fig. 29)
5. Gradually open the diaphragm. The diaphragm is centered when the diaphragm image is symmetrical to the field of view. (Fig. 30)
6. In normal use, open the diaphragm until it circumscribes the field of view.

10.11 Effects of the field diaphragm



Field diaphragm adjusts the illuminated area to obtain a high contrast image.

Set the diaphragm according to the objective in use until it circumscribes the field of view, in order to eliminate unnecessary light to eyepieces. (Fig. 30)

10.12 Centering the Aperture Diaphragm (AS)

1. Move the selector on the reflected light illuminator in the position "BF". (Fig. 21)
2. Completely pull out the aperture diaphragm lever ③ to fully close the diaphragm. (Fig. 31)
3. Remove one eyepiece.
4. While looking into the empty eyepiece sleeve, use the provided Allen wrench into the centering screws ④ to bring the bright spot in the centre of the field of view. (Fig. 31)
5. The illuminator is centered when the diaphragm image is symmetrical to the field of view.
 - The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
 - With low contrast specimens move the diaphragm lever to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the AS lever in order to obtain an image like the one in Fig. 32.

10.13 Use of color filters

1. Insert the color filter in the slot of the filter holder. (Fig. 33)
2. Move the slider into one of the three available position until it clicks into the right position.

FILTER	APPLICATION
Blue	Converts the source colour temperature to daylight
Frosted white	Reduces lighting irregularities by creating a homogeneous field of view

11. Use of the microscope in Darkfield (DF)

11.1 Switching Darkfield observation

- The microscope allows to work in brightfield (BF) or darkfield (DF).
 - For observation in darkfield it is necessary to use specific objectives called “BD”.
1. Move the lever in “DF” position for darkfield illumination. (Fig. 21)
- Darkfield objectives also allow observation in brightfield.

12. Use of the microscope in Polarized light (PO)

The performance of polarizer may deteriorate when it has been exposed to light for a long period (about continuous 2000 hours). If this happens, replace the polarizer.

1. Move the lever in “BF” position for brightfield illumination. (Fig. 21)
- **Polarized light observation cannot be performed when the BF/DF selector is in “DF” position.**
2. Insert the polarizer “PO” with the inscription facing the user.
3. Move the polarizer “PO” into the light path until the slider clicks into its position. (Fig. 34)

Fig. 34

4. Insert analyzer “AN” with the inscription facing up.
5. Move the rotatable analyzer “AN” into the light path until the slider clicks into its position. (Fig. 35)

Fig. 35

6. Put a flat mirror on the stage and get into approximate focus the mirror surface.
7. While observing into the eyepieces, rotate the scale of the analyzer ① until you achieve the darkest image. (Fig. 36)
- Now the light extinction (so-called “crossed Nicol position”) is achieved.
8. Put a specimen on the stage, focus the specimen and begin observation.



Fig. 36

13. Microphotography

13.1 Use of C-mount cameras

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 37)

2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 38)

13.2 Use of Reflex cameras

1. Insert the Reflex adapter ① into the relay tube to the microscope ②.
 2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
 3. Connect the Reflex camera ④ to the "T2" ring just installed (Fig. 39).
 4. Mount the other end of the relay tube ② into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 37)
- "T2" ring is not provided along with the microscope, but is commercially available.
 - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
 - To calculate the magnification of the camera: objective magnification * camera magnification * lens magnification.
 - **If using an SLR camera, mirror movement may cause the camera to vibrate.**
 - **We suggest lifting the mirror, using long exposure times and a remote cord.**

14. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.
- The equipment is approved according to EC safety regulations. However, users bear full responsibility for the safe use of the microscope.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

15. Troubleshooting


Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
Bulb does not light	Power cord/supply is unplugged	Connect into the power outlet
Bulb operates, but field of view remains dark	Brightness is too low	Set brightness to a proper level
	DF selector on the reflected light illuminator is inserted but no DF objectives are available	Move the selector to BF position
	Colour filter slider is in a wrong position	Move to a click stop
Field of view is obscured or not evenly illuminated	Light path selector knob is set to the camera position	Move the selector to the eye position
	The nosepiece is not in the correct position	Turn the nosepiece to a click stop
	The colour filter is partially inserted	Insert the filter to full depth
	Field diaphragm is stopped down too far	Open the field diaphragm until it circumscribes the field of view
	DF selector on the reflected light illuminator is inserted but no DF objectives are available	Move the selector to BF position
Dust and stains can be seen in the field of view	Dirt/dust on the specimen	Clean thoroughly
	Dirt/dust on the eyepieces	
Image looks double	Aperture diaphragm is stopped down too far	Open aperture diaphragm
	The field diaphragm is not well centered	Set the diaphragm according to Kohler settings
Poor image quality: <ul style="list-style-type: none"> • The image is not sharp • The contrast is not high • The details are not clear • Image glares 	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	Aperture diaphragm is too closed or too open	Adjust aperture diaphragm
	The lenses (objective, eyepieces) are dirty	Thoroughly clean all the optical system
	For reflected light observation, the specimen must not have any coverglass	Use specimens without coverglass
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position
One side of the image is out of focus.	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The specimen is out of place (tilted)	Place the specimen flat on the stage
Image appears to waver	Revolving nosepiece is not corrected mounted	Push slide dovetail all the way until it is stopped
	Objective is not correctly engaged in light path	Make sure that revolving nosepiece clicks into place correctly
	Field diaphragm not properly centered	Center the field diaphragm properly
Field of view becomes only slightly brighter when the voltage is raised	Field diaphragm not properly centered	Center the field diaphragm properly

II. Mechanical Section:		
The coarse focus knob is hard to turn	The tension adjustment collar is too tight	Loosen the tension adjustment collar
The revolver goes down by itself during observation	The tension adjustment collar is too loose	Tighten the tension adjustment collar
III. Electrical Section:		
Bulb does not light	Power cord/supply is unplugged	Connect into the power outlet
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection
IV. Observation tube:		
Field of view of one eye does not match that of the other	Interpupillary distance is incorrect	Adjust interpupillary distance
	Incorrect diopter adjustment	Adjust diopter
	Your view is not accustomed to microscope observation	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope
V. Microphotography:		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.