High-performance Full HD HDMI camera

INSTRUCTION MANUAL


Model	Code
OPTIKAM HDMI PRO	4083.13

v 1.0 2017

M-005

Micrometer slide for software calibration, 1mm/10um, 10mm/100um.

1 DIV=0.1mm


For stereo microscopes calibration

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1 DIV=0.01mm


For biological microscopes calibration

For biological microscopes calibration

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Ø 0.15mm

For checking the calibration done by drawing a line or a circle

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Table of Contents

Package Contents.

Installation of camera drivers and software.

Some notes.

How to test your Camera.

Select the right driver.

Getting started with Optika Vision Lite, a very easy measuring and documentation software.

Getting started with optika IsView, a complete software for measuring, adjusting, capturing images and video clip.

HDMI user guide.

Equipment disposal.

Package Contents

DESCRIPTION	QUANTITY
HDMI Camera	1
12V 2000mA power supply	1
SD card	1
USB cable	1
HDMI cable	1
Calibration slide	1
CD	1
Mouse	1

Installation of camera drivers and software

1. Operating system requirements: all the Windows O.S.
2. Please insert your CD and run:
\\optika vision lite\ setup.exe
3. Now the Optika Vision Lite icon appears on your desktop.
4. Please Run
\\optika IsView\ setup.exe
5. Now the Optika IsView icon appear on your desktop
6. You don't need to install any driver.
7. Plug the Camera in your PC: the camera will be identified automatically.

Some notes

Often a PC has some ports on the front and some on the back: the ports on the back are usually faster. We suggest to use these ones.

If you connect the camera into another USB port, Windows needs to recognize its Driver again even if it has already done this operation for the first USB port.

It may happen sometimes that Windows doesn't recognize a Camera using one of your PC USB ports, while gives no problem using another port.

How to test your Camera

In order to check if your Camera works correctly, it's not necessary to put it into the microscope: it's enough to connect it into the USB port, run Vision Lite using the icon on your desktop, select the right driver and check whether the shown image is not totally black.

Select the right driver

Please select: USB CAMERA.

With Optika Vision Lite you have to tell the Software which is the Driver by selecting the microscope-shaped icon. With Optika IsView the live image coming from your camera will automatically appear.

Getting started with Optika Vision Lite, a very easy measuring and documentation software.

Run Optika Vision Lite by clicking on its icon on your desktop.

Select the right Driver for capturing the image. Note that after choosing the driver and before capturing the image you can set the brightness, saturation, gamma ecc..
Then you can capture the image.

The image is now ready for being manipulated using Optika Vision Lite.

Getting started with optika IsView, a complete software for measuring, adjusting, capturing images and video clip.

Run Optika IsView by clicking on your desktop.


The live image with easy buttons for set your camera live view will automatically appear.

HDMI user guide

1. USB interface.

A: Connect a mouse to the USB port. Use the mouse to control the camera directly.

B: Connect the USB port to PC to make the camera work as a Driver-Free camera.

2. HDMI interface.

Use the HDMI cable to connect the camera to the HDMI monitor.

3. Power interface.

Please use the provided 12V2A power supply. When power is plugged in the red light is on. When switch on the camera blue light is on.

4. ON/OFF KEY

Press and hold the ON/OFF key until the blue light is ON, press the ON/OFF key again to turn the camera off.


5. SD CARD

To get faster and more stable data transfer recommend to use "Class10" SD card we supply a 8 giga SD card.

6. C MOUNT

Standard C-mount optical port.

HDMI Operation


Step 1. Connect the camera.


1. Plug in the 12V2A power supply. Use HDMI cable to connect the camera to the monitor. Press and hold ON/OFF key until blue light is on.
2. Connect the mouse to the USB port. Move the cursor to get the settings on the screen.
3. Insert the SD card.

Step 2. Move the cursor to the left of the screen.

When move the cursor to the left of the screen, 'Capture' and 'Setting' icons will appear (See image on the left side). Click 'Setting' to get parameter setting menu (See image on the right side).


Step 3. Move the cursor to the right of the screen.


Step 4. Check the images and videos you have captured on the SD card.


1. Check the album and delete the images.

2. Check the video and delete the video. To delete the video successfully, be sure this video is not in use.


Step 5. Compare preview with the captured images.



Connect the camera to the PC

1. Use the special USB cable for connecting the camera to the PC.
2. Plug in the 12V 2A Power Supply.
3. Press and hold the on/off key until the blue light is on.
4. Run the Optika IsView software, please open the on-line help to discover all the software features.


Step 6. How to do calibration and measurement in HDMI mode:

Basic steps:

1. Click on the calibration icon  to obtain the calibration table.
2. Move the cursor out of the calibration table to start the calibration process.
 1. The live image should correspond to the calibration slide or to any known-sized specimen.
 2. Draw a line to obtain a reference length, and enter the length value into the calibration table.
3. Five options are available for the length unit : MIL, CM,MM,UM and INCH.
3. Select the objective magnification in  it is possible to enter an extra value for the magnification, if needed.
 1. These data should be entered only if the user intends to use the same camera to obtain images by means of different objectives.
In this case, there is no need to create a calibration file for each image at any objective setting; just change the “Mag” option to obtain the corresponding calibration file.
 2. f only one objective is needed for your application, you should select “N/A”
 3. If you cannot find a suitable objective magnification in the list , you can use the keyboard to directly enter a new mag. group.
4. Enter a name for the newly created calibration file 
You can create up to 10 calibration groups.
5. Click  to complete the calibration settings.
6. Click  to go to the measurement window.

Click on the live image using your right key to get to the Measurement menu.

1	Scale bar	On/off the scale bar on the picture.
2	Coordinate	To get the coordinates of the selected points.
3	Line	Arbitrary direction, horizontal and vertical line measurement.
4	Freehand line	To get a freehand line length.
5	Rectangle	Rectangular measure of perimeter and area. Selectable to achieve the data of perimeter or area.
6	Circle	Center Point: It uses both center point and point on the circle to draw a circle. Two Point: To draw a circle according to the diameter. Three Point: To use 3 points on the circle to draw a circle. Selectable to obtain data on radius, perimeter or area. For a pre-set setting, diameter, perimeter and area are selected.
7	Polygon	To measure a polygon perimeter and area. Selectable to obtain data of perimeter or area.
8	Perpendicular	It measures the perpendicular length.
9	Concentric circles	It gets the diameters of two concentric circles.
10	Bicircle	It gets the distance between the central points of two circles.
11	Angle	To measure the angle.
12	Setting	1) To set the measurement line width & colour and the font colour. 2) Magnifier On/off: To switch the magnifier on/off. When the magnifier is on, the area in which the cursor is located will be zoomed and placed in the corner to help an accurate location of the measurement point.
13	Save results	Export the measurement results to an Excel or Text file. The exported results will be saved in the / MEASURE folder on the SD card.
14	Set scale	It allows to create 10 groups of calibration files, among which the user can chose the required measurements.
15	Exit	It allows the exit from the measurement.


Note:

All the measurements on the live image will be removed when exiting from there.

The measurement result can be saved on the captured images provided these are saved.

Equipment disposal

Art.13 Dlsg 25 july 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA S.r.l.