OPTIStar

By **OPTIKA**

OPTISTAR By OPTIKA Series

INSTRUCTION MANUAL

Model
OS-STX
OS-MS2
OS-ST30FX

Ver. 1.0 2021

Table of Contents

1.	Warning	3
2.	Safety Information	3
3.	Unpacking	3
4.	Intended use	3
5.	Symbols and conventions	3
6.	Instrument description	4
6.1	OS-STX	4
6.2	OS-MS2	4
6.3	OS-ST30FX	5
7.	Assembling	5
8.	Use of the microscope	6
8.1	Turning on the light	6
8.1.1	OS-STX	6
8.1.2	OS-ST30FX	6
8.2	Place the specimen on the stage	6
8.3	Adjust interpupillary distance	6
8.4	Focus	7
8.5	Diopter compensation	7
8.5.1	OS-ST30FX	7
8.6	Magnification	7
8.6.1	OS-ST30FX	7
8.7	Contrast disc	7
8.8	Inserting the batteries	8
9.	Maintenance	9
10.	Troubleshooting guide	10
	Equipment disposal	11

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

4. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

5. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


ELECTRICAL SHOCK


This symbol indicates a risk of electrical shock.

6. Instrument description


6.1 OS-STX

6.2 OS-MS2

6.3 OS-ST30FX

7. Assembling

No installation procedure is required for these models.

- **Only for OS-ST30FX:** The microscope comes with an electrical illuminator. Insert the power supply jack in the socket, turn on the switch on the main body and select your light source.

8. Use of the microscope

8.1 Turning on the light

8.1.1 OS-STX

Press the button ① to turn on or off the light. (Fig. 1)

- **Light intensity cannot be changed.**

8.1.2 OS-ST30FX

1. Press the button ② to turn on or off the transmitted light. (Fig. 2)
 2. Press the button ③ to turn on or off the incident light. (Fig. 2)
- **Light intensity cannot be changed.**

8.2 Place the specimen on the stage


Place the specimen on the microscope stage and lock it using the slide-clamp if it is necessary. Ensure that the specimen is centered over the stage opening. (Fig. 3)

8.3 Adjust interpupillary distance

Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by moving the two parts until one circle of light can be seen.

- If two circles appear, the interpupillary distance is too big, and if two overlapped circles appear, the interpupillary distance is too small. (Fig. 4)
- **This operation is not needed for OS-MS2.**

8.4 Focus

Rotate the focus knob ① to bring the specimen into focus. (Fig. 5)


8.5 Diopter compensation

This compensation makes it possible for people with glasses to adjust the microscope to their eyes and use the microscope without glasses.

8.5.1 OS-ST30FX

1. Focus the specimen using the focus knob ① while observing with the right eyepiece. (Fig. 5)
2. Adjust the diopter compensation ring ② of the left eyepiece tube until the image is clear and sharp. (Fig. 6)

8.6 Magnification

8.6.1 OS-ST30FX

Select the desired magnification by rotating the objective. (Fig. 7)

- Total magnification used can be calculated as:
Eyepiece magnification x Objective lens magnification.

8.7 Contrast disc

You can use the black/white disc in order to enhance the contrast of the image when using the incident illumination.

1. If you want to observe bright opaque objects, place the disc with the black side facing up. (Fig. 8)
2. If you want to observe dark opaque objects, place the disc with the white side facing up.


8.8 Inserting the batteries

(OS-STX only)

Before you start using the instrument, you must insert the batteries.

1. Open the battery cover ① (located on the bottom of the instrument) and insert the batteries, taking care to respect the polarities. (Fig. 9)
2. When finished, close the battery cover again.

9. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the dedicated cleaning kit.

10. Troubleshooting guide

Review the information in the table below to solve operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
The illumination is ON, but the field of view is dark	The plug is not connected to the illumination	Connect the cable
	The brightness is too low	Adjust to a proper setting
The edge of the field of view is vignetted or the brightness is asymmetric	The incident illuminator is not correctly oriented	Change the angle of the incident illuminator
Dust and stains can be seen in the field of view	There are stains and dust on the specimen	Clean the specimen
	There are stains and dust on the eyepiece	Clean the eyepiece
Poor image quality: <ul style="list-style-type: none"> The image is not sharp The contrast is not high The details are not clear 	The lenses (additional lens, objective, eyepieces) are dirt	Thoroughly clean all the optical system
One side of the image is out of focus	The specimen is out of place (tilted)	Place the specimen flat on the stage
II. Electric Section:		
The LED doesn't turn on	No power supply	Check the power cord connection
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord
III. Viewing tube assembly:		
The field of view of the two eyes is different	The interpupillar distance is not correct	Adjust the interpupillar distance
	The diopter correction is not right	Adjust the diopter correction
	The viewing technique is not correct, and the operator is straining the eyesight	When look into the eyepieces, do not stare at the specimen but look at the whole field of view. Periodically, move the eyes away to look at a distant object, then back into the eyepieces

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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