**OPTIStar**

By **OPTIKA**

---

OPTISTAR By OPTIKA Series

# INSTRUCTION MANUAL

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Ver. 1.0 2021

---

## Table of contents

1.	Warning	3
2.	Safety Information	3
3.	Unpacking	3
4.	Intended use	3
5.	Symbols and conventions	3
6.	Instrument description	4
6.1	OS-SZ1 / OS-SZ2	4
6.2	OS-SZ3 / OS-SZ4	4
6.3	OS-SZ5 / OS-SZ6	5
6.4	Ringlight Illuminator	5
7.	Assembling	6
7.1	OS-SZ1 / OS-SZ2	6
7.1.1	Installing focus system	6
7.1.2	Connecting incident light cable	6
7.1.3	Installing the head	6
7.1.4	Installing the eyepieces	6
7.1.5	Installing the phototube (only trinocular model)	7
7.1.6	Connecting the power supply	7
7.2	OS-SZ3 / OS-SZ4	7
7.3	OS-SZ5 / OS-SZ6	10
7.4	Installing Ringlight illuminator (OS-SZ3 / OS-SZ6)	12
8.	Use of the microscope	13
8.1	Adjusting interpupillary distance	13
8.2	Focusing	13
8.3	Adjusting the tension of the focus knob	13
8.4	Diopter compensation	13
8.5	Magnification	14
8.6	Use of additional lens	14
8.7	Black/white stage plate (OS-SZ1 / OS-SZ2)	15
8.8	Use of illumination (OS-SZ1 / OS-SZ2)	15
8.9	Use of overhanging stand	16
8.9.1	OS-SZ3 / OS-SZ4	16
8.9.2	OS-SZ5 / OS-SZ6	17
8.10	Use of Ringlight Illuminator	18
9.	Microphotography	19
10.	Maintenance	20
11.	Troubleshooting	21
	Equipment disposal	22

---

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

## 2. Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

## 3. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 4. Intended use

### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

## 5. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


### ELECTRICAL SHOCK


This symbol indicates a risk of electrical shock.

## 6. Instrument description


### 6.1 OS-SZ1 / OS-SZ2

### 6.2 OS-SZ3 / OS-SZ4

**6.3 OS-SZ5 / OS-SZ6**

**6.4 Ringlight Illuminator**

## 7. Assembling

### 7.1 OS-SZ1 / OS-SZ2


#### 7.1.1 Installing focus system

Place the focusing system on the pillar and, once the desired height is reached, tighten the locking knob ① located on the back of the focusing system. (Fig. 1)

#### 7.1.2 Connecting incident light cable

Insert the cable plug into the socket at the top of the pillar ②. (Fig. 2)

#### 7.1.3 Installing the head

Install the head and screw the fixing knob ③. (Fig. 3)

#### 7.1.4 Installing the eyepieces

Remove the dust caps from the eyepiece sleeves and insert the eyepieces. (Fig. 4)

### 7.1.5 Installing the phototube (only OS-SZ2)

1. Remove the dust cover by unlocking the fixing screw ①. (Fig. 5)
2. Insert the ring of the photo port ② making sure to leave the fixing screw of the photo adapter ③ facing right.
3. Lock the fixing screws ①.

### 7.1.6 Connecting the power supply

Connect the power supply plug to the socket on the back of the microscope base. (Fig. 6)

### 7.2 OS-SZ3 / OS-SZ4

1. Screw the pillar on the base. (Fig. 7)

2. Tighten the screw to lock the pillar. (Fig. 8)

3. Insert the drop preventing ring and fix it at the desired height by screwing the fixing knob. (Fig. 9)

4. Insert the horizontal arm and secure it with the fixing screw ①. (Fig. 10-11)

5. Install the head holder. Unscrew the locking knob ③ and insert the arm for the focus support from above ② into the hole of the horizontal arm. (Fig. 12-13)

6. Once fully inserted, tighten the fixing screw ④ (Fig. 13)

7. Insert from below the focusing system, tighten the fixing screw ⑤ and re-tighten the locking knob ③ from below. (Fig. 14-15)
8. Repeat steps from 7.1.3 to 7.1.5 for the installation of the head.

### 7.3 OS-SZ5 / OS-SZ6

1. Screw the pillar on the base. (Fig. 16)

2. Tighten the screw to lock the pillar. (Fig. 17)


3. Insert the drop preventing ring and fix it at the desired height by screwing the fixing knob. (Fig. 18)



4. Insert the horizontal arm and secure it with the fixing screw ①. (Fig. 19)

5. For a bigger safety, lock the fixing screw ② with the provided Allen wrench. (Fig. 20)

6. Install the head holder. Unscrew the prevention knob ③ and insert the head holder from below in the pillar. Lock the fixing knob ④. At the end screw again the prevention knob ③. (Fig. 21-22)
7. Repeat steps from 7.1.3 to 7.1.5 for the installation of the head.

#### 7.4 Installing Ringlight illuminator (OS-SZ3 / OS-SZ6)

1. Screw the fixing ring ① on the bottom of the Stereo head. (Fig. 23)

2. Place the illuminator on the fixing ring, using the three locking screws ②. (Fig. 24)

3. Connect the power supply jack to the connector on the illuminator. (Fig. 25)
4. Connect the power supply to the wall socket.

## 8. Use of the microscope

### 8.1 Adjusting interpupillary distance

1. Hold the right and left eyepiece tube with both hands and adjust the interpupillary distance by moving the two parts until one circle of light can be seen. (Fig. 26)
- If two circles appear, the interpupillary distance is too big, and if two overlapped circles appear, the interpupillary distance is too small.

### 8.2 Focusing


Put the sample to be observed on the stage plate and focus the sample using the focusing knobs ①.

- If necessary adjust the height of the microscope head along the vertical stand.
- Remember to lock the fixing knob ② after setting the height of the microscope. (Fig. 27)


### 8.3 Adjusting the tension of the focus knob

- **This adjustment allows to increase or decrease the tension of the knob by avoiding an involuntary descent of the microscope body under its own weight. Adjust the tension just above the point where the focus is stable.**
1. Grab the knobs ③ both hands and, while holding the left knob, turn the right knob.
- The tension increases or decreases depending on the direction of rotation of the right knob. (Fig. 28)

### 8.4 Diopter compensation

- **This compensation makes possible for people wearing glasses to adjust the microscope to their eyes and use the microscope without glasses.**
1. Put the zoom down to the lowest magnification ④ and focus the specimen with the focusing knobs ⑤. (Fig. 29)
  2. Put the zoom to the maximum magnification and repeat the focusing.
  3. Return to the lowest magnification: the specimen will be out of focus.


4. Adjust the diopter compensation ring of the right eyepiece ⑥ (Fig. 30) until the image of the right eyepiece is clear and sharp. Repeat the procedure for the left eyepiece.
5. Check the focus of the image for the whole zoom range. It should be perfectly parfocal (focus is maintained during the change of magnification).

### 8.5 Magnification

1. Select the desired magnification by adjusting the zoom knob ①. (Fig. 31)
- Change the eyepieces and/or add an appropriate additional lens if needed.

### 8.6 Use of additional lens

1. Screw the desired additional lens on the microscope body. (Fig. 32)
- The stroke of the focus adapter could not compensate the different working distances of the several additional lenses.
- In the event that the focus support fails to focus the sample, the entire body of the microscope must be raised or lowered.

Total magnification used can be calculated as:  
 Eyepiece magnification \* Zoom magnification \* Objective lens magnification.

Eyepiece	10x		20x	
Field Number (mm)	22		9	
Objective	Total Mag.	Field Of View (mm)	Total Mag.	Field Of View (mm)
0.5X	3.5X-25X	62.86-8.80	7X-50X	25.71-3.60
1X	7X-50X	31.43-4.40	14X-100X	12.86-1.80
2X	14X-100X	15.71-2.20	28X-200X	6.43-0.90

### 8.7 Black/white stage plate (OS-SZ1 / OS-SZ2)

- You can use a black/white contrast disc to increase the contrast of the image when working in incident light. (Fig. 33)
1. If you are looking at bright samples, place the disc with the black part facing up.
  2. If you are looking at dark samples, place the disc with the white part facing up.

### 8.8 Use of illumination (OS-SZ1 / OS-SZ2)

1. Move the main switch ① to turn ON the power supply. (Fig. 34)
  2. Turn the transmitted light knob ② to change the intensity of the transmitted light LED.
  3. Turn the incident light knob ③ to change the intensity of the incident light LED.
- It is possible to use both illuminations at the same time.

## 8.9 Use of overhanging stand

### 8.9.1 OS-SZ3 / OS-SZ4

#### Moving the horizontal arm

1. Unlock the knob on the right side of the horizontal arm ①. (Fig. 35)

2. The arm can be extended or shortened according to specific needs. (Fig. 36)


#### Swivel the head

1. Unlock the fixing knob ② and rotate the head to the desired swivel angle (left or right), then tighten the knob again. (Fig. 37)

#### Rotating the horizontal arm


1. Loosen the horizontal arm fixing knob ③ and rotate the arm, then tighten again the fixing knob. (Fig. 38)
- **NOTE: 180° rotation of the microscope with respect to the base could cause a rollover of the entire system.**

## 8.9.2 OS-SZ5 / OS-SZ6

### Moving the horizontal arm

1. Unlock the knob on the left side of the horizontal arm ①. (Fig. 39)
2. The arm can be extended or shortened according to specific needs. (Fig. 40)

### Tilting the head holder

1. Loosen the knob ② (Fig. 41) on the right side of the head holder. The tilting pillar can be moved. Once the desired angle is achieved, tighten the knob again.

### Swivel the head

1. Unlock the fixing knob ③ and rotate the head to the desired swivel angle (left or right), then tighten the knob again. (Fig. 42-43)


Fig. 43

### Rotating the horizontal arm

1. Loosen the horizontal arm fixing knob ④ and rotate the arm, then tighten again the knob. (Fig. 44)
- **NOTE: 180° rotation of the microscope with respect to the base could cause a rollover of the entire system.**

Fig. 44

### 8.10 Use of Ringlight Illuminator

1. Use ON/OFF switch ① to turn the LEDs on. (Fig. 45)
2. Use the knob ② to adjust the light intensity.

Fig. 45

---

## 9. Microphotography

- Refer to specific camera Instruction manual for detailed explanation about microphotography procedure.
- Instruction manual can be downloaded from this link: <https://www.optikamicroscopes.com/optikamicroscopes/download/optistar-by-optika-digital-cameras/>

---

## 10. Maintenance

### Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

### To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

### Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

**For the best results, use the dedicated cleaning kit.**

## 11. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
<b>I. Optical Section:</b>		
The illumination is ON, but the field of view is dark.	The plug is not connected to the illumination	Connect the cable
	The brightness is too low	Adjust to a proper setting
The edge of the field of view is vignetted or the brightness is asymmetric.	The incident illuminator is not correctly oriented	Change the angle of the incident illuminator
Dust and stains can be seen in the field of view.	There are stains and dust on the specimen	Clean the specimen
	There are stains and dust on the eyepiece	Clean the eyepiece
Poor image quality: <ul style="list-style-type: none"> <li>The image is not sharp</li> <li>The contrast is not high</li> <li>The details are not clear</li> <li>Image glares</li> </ul>	The lenses (additional lens, objective, eyepieces) are dirt	Thoroughly clean all the optical system
One side of the image is out of focus.	The specimen is out of place (tilted)	Place the specimen flat on the stage.
<b>II. Mechanical Section:</b>		
The focus knob is hard to turn.	The tension adjustment collar is too tight	Loosen the tension adjustment collar
The focus is unstable.	The tension adjustment collar is too loose	Tighten the tension adjustment collar
<b>III. Electric Section:</b>		
The LED doesn't turn on.	No power supply	Check the power cord connection
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord
<b>IV. Viewing tube assembly:</b>		
The field of view of the two eyes is different	The interpupillary distance is not correct	Adjust the interpupillary distance
	The diopter correction is not right	Adjust the diopter correction
	The viewing technique is not correct, and the operator is straining the eyesight	When look into the eyepieces, do not stare at the specimen but look at the whole field of view. Periodically, move the eyes away to look at a distant object, then back into the eyepieces
<b>V. Microphotography and video:</b>		
The image is unfocused	Incorrect focusing	Adjusting the focus system as in the present manual
Bright patches appear on the image	Stray light is entering the microscope through the eyepieces and through the camera viewfinder	Cover the eyepieces and the viewfinder with a dark cloth

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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