## **ST-155 / ST-156** **Microscope stand with** **LED illumination**

 **OPERATION MANUAL**


 **GUIDA UTENTE**

 **MANUAL DE INSTRUCCIONES**

 **MANUEL D'INSTRUCTIONS**

<b>1.0 DESCRIPTION</b>	page 3
<b>2.0 INTRODUCTION</b>	page 4
<b>3.0 UNPACKING AND ASSEMBLY</b>	page 4
<b>4.0 USING THE STAND</b>	page 5
<b>5.0 MAINTENANCE</b>	page 6
<b>6.0 TECHNICAL SPECIFICATIONS</b>	page 7
<b>7.0 RECOVERY AND RECYCLING</b>	page 8

## 2.0 INTRODUCTION

The OPTIKA ST-155 and ST-156 are high quality microscope stands with built-in LED illuminators (transmitted and incident), both with intensity control.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

There are, at the moment, 2 models in the series:

**ST-155:** Modern stand with large base. LED trans-illuminator and LED incident illuminator, both with intensity control.

**ST-156:** Same as ST-155 but with coaxial coarse and fine focusing system.

## 3.0 UNPACKING AND ASSEMBLY

The components for the ST-155 and ST-156 series are shipped in a Styrofoam box for protection. Open the box with care to avoid that any components fall out and get damaged.

If the instrument has been damaged during the transport, please contact both the carrier and your supplier immediately.

Unpack the stand and put it on a good stable table.

Put the stereo head (sold separately) onto the focusing carrier of the stand and lock it by the lock screw on the front of the focusing carrier. During the operation of the microscope, this screw should always be locked to maintain maximum stability.

### 4.1 Illumination system

The ST-155 and ST-156 come with an integrated illumination system with incident and transmitted light.

Before starting, read the section 5.3 about electrical safety precautions. Then, insert the power cable into the plug on the rear and press the ON/OFF switch.

Use the brightness adjustment knobs to control the light intensity (separately for incident and transmitted illuminators).

The incident LED illuminator can be oriented in order to obtain the best condition of light.

The base is equipped with two clips that can be used to block the specimen.

### 4.2 Focus

Put the sample to be observed on the stage plate and adjust the zoom to the desired value on the optical head.

Try to focus the sample using the focusing knobs of the stand. On the ST-156 model you can also fine-adjust the focus through the coaxial knob.

## 5.0 MAINTENANCE

### 5.1 Microscopy environment

This illuminator stand is recommended to be used in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

### 5.2 To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope or stand yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

### 5.3 Electrical safety precautions


- Before plugging in the power cable with the supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the ON/OFF switch is in OFF position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

## 6.0 TECHNICAL SPECIFICATIONS

INPUT VOLTAGE	110 – 240VAC 50HZ
EYEPIECE	100MM DIAMETER PLATE. LED MATRIX ILLUMINATOR. COLOR: COOL WHITE 6300K. BRIGHTNESS ADJUSTABLE.
INCIDENT ILLUMINATOR	LED MATRIX ILLUMINATOR. COLOR: COOL WHITE 6300K. BRIGHTNESS ADJUSTABLE. ANGLE ADJUSTABLE.
STAND DIMENSIONS	W X L X H: 330MM X 285MM X 305MM
FOCUS MOVEMENT	RACK & PINION MECHANISM. ST-155: COARSE FOCUS. ST-156: COARSE AND FINE COAXIAL FOCUS.

Art.13 Dlsg 25 July 2005 N°151. “According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal.”

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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