OPTIKA

M I C R O S C O P E S
I T A L Y

Ver. 3.0.0

SZM-D

 **OPERATION MANUAL**

 **GUIDA UTENTE**

 **MANUEL D'INSTRUCTIONS**

 **GUIA USUARIO**

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OPTIKA MICROSCOPES - ITALY

www.optikamicroscopes.com - info@optikamicroscopes.com

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

Safety guidelines

This manual contains important information and warnings regarding safety about installation, use and maintenance of the microscope. Please read this manual carefully before using the equipment. To ensure safe use, the user must read and follow all instructions in this manual. OPTIKA products are designed for safe use in normal operating conditions. The equipment and accessories described in the manual are manufactured and tested according to industry standards for safety instrumentation laboratory. Misuse can cause personal injury or damage to the instrument. Keep this manual at hand close to the instrument, for an easy consultation.

Electrical safety

Before connecting the power cord to wall outlet, ensure that your mains voltage for your region corresponds to the voltage supply of the instrument, and that the illuminator's switch is in position OFF. The user must observe the safety regulations in force in his region. The instrument is equipped with CE safety marking, in any case the user has full responsibility concerning the safe use of that instrument.

Warning/Caution symbols used in this manual

The user should be aware of safety aspects when using the instrument. Warning or hazard symbols are shown below. These symbols are used in this manual.

DANGER

The instructions on this symbol to avoid possible severe personal injuries.

WARNING

Warning of use; the incorrect operation on the instrument can cause damages to the person or instrument.

WARNING

Possibility of electric shock.

HOT!

Attention: high temperature surfaces. Avoid direct contact.

NOTE

Technical notes or usage tips.

**Connect the mains plug into the socket at the base**

Make sure, before you turn the illumination on, that the voltage selector is set to the mains voltage for your region.

The power cord should be used only on network sockets equipped with adequate grounding. Contact a technician to check the state of your electrical system. If there is no need to install additional accessories, the instrument is now ready for use. Once positioned and installed with the necessary components, the microscope is ready to be used. Your microscope is a laboratory instrument designed to last. Handle it always carefully and avoid abrupt vibrations or shocks. Always disconnect the power cable from the microscope when not in use for long time, while you clean it or when you perform any maintenance.

**AVOID DISASSEMBLING THE INSTRUMENT**

Do not disassemble the instrument. This entails the cancellation of the warranty and may cause malfunction.

SZM-D - OPERATION MANUAL**1.0 INSTALLATION OF CAMERA DRIVERS AND SOFTWARE**

- 1.1 Operating system requirements: Windows XP, Windows Vista or Windows 7. If your PC is Windows XP: before the installation of your Camera please check if Service Pack2 is installed. If not, please update your XP from the Microsoft web site. If your PC has a USB port less than 2.0, there may be some problems in capturing image from the camera due to the slow transmission.
- 1.2 Please insert your CD and run:
\\ vision lite \ setup.exe
- 1.3 Now Optika Vision Lite icon appears on your desktop.
- 1.4 Choose from the folder “\drivers” the model of camera you want to use and select the file “setup.exe”.
- 1.5 The installer may show warnings about missing Windows certification: ignore such warnings and select “Continue”. Note : do NOT cancel the installation.
- 1.6 At the end of the installation process restart the PC if requested.
- 1.7 Plug the Camera in your PC: the camera will be identified and Windows will recognize the driver which has just been installed from the CD.
Then ignore any message about Windows compatibility and select “Continue” .
Note : do NOT cancel the installation.
- 1.8 At the end of the identification process, disconnect your Camera and try to plug it again, in order to

check if the Driver was installed correctly (no Windows error messages should appear). Should any error message appear on your PC, you will find a “\Troubleshooting” folder on your CD which shows how to solve any problem which Windows could give.

2.0 SOME NOTES

You should know that your Windows PC needs to install a Driver before capturing images from the camera: after installing the driver from the CD, Windows needs to recognize the camera before using it.

Often a PC has some ports on the front and some on the back: the ports on the back are usually faster. We suggest to use these ones.

All our cameras, even though they may have an external connection, do not need any external power supply because it comes from the USB cable of your PC.

If you connect the camera into another USB port, Windows needs to recognize its Driver again even if it has already done this operation for the first USB port.

Generally, every time you use a peripheral device and after that you connect your Camera, then Windows needs to recognize the Driver again.

It may happen sometimes that Windows doesn't recognize a Camera using one of your PC USB ports, while gives no problem using another port.

To prevent any installation problem we have put on the CD all the software and drivers in specific folders with an appropriate names.

3.0 HOW TO TEST YOUR CAMERA

In order to check if your Camera works correctly, it's not necessary to put it into the microscope: it's enough to connect it into the USB port, run Vision Lite using the icon on your desktop, select the right driver and check whether the shown image is not totally black.

4.0 SELECT THE RIGHT DRIVERS

Select the driver :

OPTIKA TCA 1.31

To use Optika Vision Lite you have to tell the Software which is the Driver by selecting the microscope-shaped icon.

5.0 GETTING STARTED WITH OPTIKA VISION LITE, A VERY EASY MEASURING AND DOCUMENTATION SOFTWARE

Run Optika Vision Lite by clicking on its icon on your desktop.

Select the right Driver for capturing the image. Note that after choosing the driver and before capturing the image you can set the brightness, saturation, gamma ecc..


Then you can capture the image.

The image is now ready for being manipulated using Optika Vision Lite.

Optika vision Lite also allows to make simple video clips.

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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