 **OPERATION MANUAL**

 **MANUEL D'INSTRUCTIONS**

 **GUIDA UTENTE**

 **MANUAL DE INSTRUCCIONES**

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<b>1.0 DESCRIPTION</b>	page 4
<b>2.0 INTRODUCTION</b>	page 5
<b>3.0 UNPACKING AND ASSEMBLY</b>	page 5
<b>4.0 USING THE MICROSCOPE</b>	page 6
<b>5.0 MAINTENANCE</b>	page 7
<b>6.0 ELECTRICS</b>	page 8
<b>7.0 OPTIONAL ACCESSORIES</b>	page 8
<b>8.0 RECOVERY AND RECYCLING</b>	page 10

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

### Safety guidelines

This manual contains important information and warnings regarding safety about installation, use and maintenance of the microscope. Please read this manual carefully before using the equipment. To ensure safe use, the user must read and follow all instructions in this manual. OPTIKA products are designed for safe use in normal operating conditions. The equipment and accessories described in the manual are manufactured and tested according to industry standards for safety instrumentation laboratory. Misuse can cause personal injury or damage to the instrument. Keep this manual at hand close to the instrument, for an easy consultation.

### Electrical safety

Before connecting the power cord to wall outlet, ensure that your mains voltage for your region corresponds to the voltage supply of the instrument, and that the illuminator's switch is in position OFF. The user must observe the safety regulations in force in his region. The instrument is equipped with CE safety marking, in any case the user has full responsibility concerning the safe use of that instrument.

### Warning/Caution symbols used in this manual

The user should be aware of safety aspects when using the instrument. Warning or hazard symbols are shown below. These symbols are used in this manual.

**DANGER**

The instructions on this symbol to avoid possible severe personal injuries.

**WARNING**

Warning of use; the incorrect operation on the instrument can cause damages to the person or instrument.

**WARNING**


Possibility of electric shock.

**HOT!**

Attention: high temperature surfaces. Avoid direct contact.

**NOTE**

Technical notes or usage tips.

### **SZM-1**

### **SZM-4**

#### TRINOCULAR HEAD

PHOTO/VIDEO CAPTURING TUBE

The OPTIKA zoom stereo microscopes SZM series are high performance Greenough microscopes with continuous zoom system 0.7:4.5.

These instruments are made for examining any three dimensional object in industrial, biological and educational applications. They are equipped with two separate optical systems that give an excellent field depth to the image, a stereoscopic effect.

## 3.0 UNPACKING AND ASSEMBLY

The components for the SZM series are shipped in a Styrofoam box for protection. Open the box with care to avoid that any components fall out and get damaged.

If the instrument has been damaged during the transport, please contact both the carrier and your supplier immediately.

Unpack the stand and put it on a good stable table. Loosen the lock screw of the focusing carrier, adjust the height of the focusing carrier and lock the lock screw again. For the SZM-3/4, make sure that the support collar is secured firmly below the focusing carrier along the vertical post. This is important as this collar functions to avoid the accidental falling of the microscope head along the stand.

Put the stereo head onto the focusing carrier and lock it by the lock screw on the side of the focusing carrier. During the operation of the microscope, this screw should always be locked to maintain maximum stability.

Unpack the eyepieces and other optical parts carefully and place the eyepieces into the eyepiece tubes. When handling the optical parts, avoid touching any lens surface with naked hand or fingers. Any fingerprints or grease stains will negatively affect the image quality.

Connect the mains plug into the socket at the base

Make sure, before you turn the illumination on, that the voltage selector is set to the mains voltage for your region.

The power cord should be used only on network sockets equipped with adequate grounding. Contact a technician to check the state of your electrical system. If there is no need to install additional accessories, the instrument is now ready for use. Once positioned and installed with the necessary components, the microscope is ready to be used. Your microscope is a laboratory instrument designed to last. Handle it always carefully and avoid abrupt vibrations or shocks. Always disconnect the power cable from the microscope when not in use for long time, while you clean it or when you perform any maintenance.

### AVOID DISASSEMBLING THE INSTRUMENT

Do not disassemble the instrument. This entails the cancellation of the warranty and may cause malfunction.

### 4.1 Illumination system

The SZM-1 and SZM-2 come with an integrated illumination system with incident and transmitted light. For the microscopes without integrated illumination systems it is possible to use an external cold light source. By using the correct light it is possible to obtain a better image of your sample.

Before starting, read the section 5.3 about electrical safety precautions. Then, insert the power cord and use the brightness controls to turn on the lamps and control the light intensity.

The base is equipped with two clips that can be used to block the specimen.

### 4.2 Adjust interpupillary distance

Hold the right and left parts of the observation head by both hands and adjust the interpupillary distance by moving the two parts until one circle of light can be seen. If two circles appear, the interpupillary distance is too big, and if two overlapped circles appear, the interpupillary distance is too small.

### 4.3 Focus and regulation of the focusing tension

Put the sample to be observed on the plate and turn the zoom to the highest magnification using the bilateral zoom knob. Try to focus the sample using the focusing knobs.

### 4.4 Dioptic compensation

This compensation makes it possible for people with glasses to adjust the microscope to their eyes and use the microscope without glasses. Turn the zoom down to the lowest magnification. Adjust the diopter compensation ring of the right eyepiece tube until the image of the right eyepiece is clear and sharp. Repeat the procedure for the left eyepiece. Then, check the focus of the image for the whole zoom range. It should now be perfectly parfocal (focus is always maintained during the change of magnification).

### 4.5 Magnification and working distance

The total magnification varies from 1.75x to 180x depending on the eyepieces and auxiliary objectives used. With the bilateral zoom knob, the user can change the magnification in a factor of 6.4 with perfect parfocality and the image is always centred.

Select the desired magnification by adjusting the bilateral zoom knob. Change the eyepieces and/or add an appropriate magnifying objective lens if necessary. The total magnification used can be calculated by the following equation:

$$\text{Total magnification} = \text{Eyepiece magnification} \times \text{Zoom magnification} \times \text{Objective lens magnification}$$

Normal working distance for the standard configuration (1x objective lens) is 95 mm.

### 4.6 Video capturing (optional)

Two kinds of observation heads, binocular or trinocular, are available for the SZM series. They can all be connected to cameras via an adaptor, for digital and analogical photo and video capturing.

Before taking a picture or filming video, pull out the light path selector lever so the light will be deflected into the photo tube. At the same time, no light will enter into the right eyepiece tube for observation. Please refer to the adaptor and camera manuals for further details.

### 5.1 Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 0-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

### 5.2 To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

### 5.3 Electrical safety precautions

- Before plugging in the power cord with the supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- For models with halogen lamp, do not turn the power on and off, off and on immediately as this will shorten the life span of the bulb and may cause damage to the electrical system.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### 5.4 Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.  
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

### 5.5 Changing the lamp [only for models with halogen lamp]

Before changing the bulb, make sure that the power switch is off and the power cord has been disconnected from the main supply.

For incident light: unscrew the lamp collector piece and pull it out, remove the old light bulb from the socket, and carefully plug the new bulb in. Screw the collector piece back after changing the bulb.

For transmitted light: remove the glass plate by pushing down its rear side, remove the old bulb from the socket and carefully plug the new bulb in.

Never touch the glass surface of the bulb with naked hand, any grease stain brought onto the bulb by the hand will negatively affect the heat dissipation, and thus, greatly shorten the life span of the bulb. Clean the bulb surface with alcohol and tissue if the user has touched the bulb surface accidentally.

<b>Universal Power:</b>	230Vac [optional 110Vac], 50/60 Hz
<b>Transmitted light:</b>	12V/15W Halogen bulb LED model: high power diffusing disc, 6300K
<b>Incident light:</b>	12V/15W Halogen bulb with dichroic mirror LED model: high power single chip LED, 6300K

## 7.0 OPTIONAL ACCESSORIES

### 7.1 Eyepieces and auxiliary objectives

With the standard equipment is included a pair of 10x eyepieces. There are also 15x and 20x eyepieces available as optional accessories.

To change the eyepieces, remove the original eyepieces and replace with the new pair.

The standard objective is 1x and as optional accessories there are 0,5x, 1,5x, and 2x additional objectives. An additional objective can be needed to obtain a specific magnification or a different working distance.

The additional objective is mounted simply by screwing it onto the stereo-head. The height of the microscope has to be re-adjusted as the working distance is changed when additional objective is used.

### 7.2 Dark field condenser

The dark field condenser has to be used with transmitted light system. By putting this accessory onto the transmitted light outlet (removing the frosted glass or stage plate), it creates the dark field effect. It is especially useful in jewellery or gem study and special techniques in biomedical applications.

### 7.3 Photo/Video adapter

This accessory, mounted on the trinocular port, allows a reflex camera of any brand to be used on the microscope in order to capture photos. T2 ring, specific to each model of camera, must be inserted into the adapter before mounting it on the camera. T2 ring can be purchased in any photcamera store.

### 7.4 "C-mount" adapter


This accessory, mounted on the trinocular port, allows a digital or analog camera to be used on the microscope in order to capture photos. The camera must have a C-mount or CS-mount connection type.

## 7.5 List of available accessories

<b>ST-081</b>	Eyepieces (pair) WF10x/20 mm
<b>ST-082</b>	Eyepieces (pair) WF15x/15 mm
<b>ST-083</b>	Eyepieces (pair) WF20x/10 mm
<b>ST-084</b>	Micrometric eyepiece WF10x/20 mm
<b>ST-085</b>	Additional lens 0,5x
<b>ST-091</b>	Additional lens 0,75x
<b>ST-086</b>	Additional lens 1,5x
<b>ST-087</b>	Additional lens 2x
<b>ST-088</b>	Polarising set (filters and rotating stage)
<b>ST-040</b>	Darkfield condenser
<b>ST-041</b>	Sample clip
<b>ST-100</b>	Hand moving stage
<b>ST-089</b>	Reflex camera adapter
<b>ST-090</b>	CCD camera adapter (for 1/3" sensor)
<b>ST-090.1</b>	CCD camera adapter (for 1/2" sensor)
<b>ST-036</b>	Eye cups (pair) type 2
<b>ST-012</b>	White/black object-plate, type 2 dia. 95 mm
<b>ST-014</b>	Glass stage, type 2, dia. 95 mm
<b>ST-038</b>	Halogen bulb, 12V/15W
<b>ST-037</b>	Halogen bulb, 12V/15W, with dichroic mirror
<b>ST-033</b>	Dust cover type 13
<b>ST-666</b>	Heating stage, with digital temperature controller
<b>ST-173</b>	APS-C Reflex camera adapter

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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