## XC-100L

 **OPERATION MANUAL**

 **GUIDA UTENTE**

 **MANUAL DE INSTRUCCIONES**

 **MANUEL D'INSTRUCTIONS**

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<b>1.0 DESCRIPTION</b>	page 3
<b>2.0 INTRODUCTION</b>	page 4
<b>3.0 UNPACKING AND ASSEMBLY</b>	page 4
<b>4.0 USING THE MICROSCOPE</b>	page 4
<b>5.0 MAINTENANCE</b>	page 5
<b>6.0 RECOVERY AND RECYCLING</b>	page 6

**XC-100L**

## 2.0 INTRODUCTION

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily laboratory use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

## 3.0 UNPACKING AND ASSEMBLY

Unpack the microscope from the box and put it on a stable surface. Unscrew the bottom part of the penlight and insert two AAA batteries. Screw the lamp holder onto the microscope and fixate the lamp to the lamp holder with the screws. The microscope is now ready to use.

## 4.0 USING THE MICROSCOPE

- 4.1 Put the microscope on top of the object or surface you want to examine.
- 4.2 Switch on the light on the bottom of the lamp.
- 4.3 Rotate the focusing knob to bring the object into focus.
- 4.4 To avoid that the objective hits against the surface that are to be examined and in order to always use the same focusing distance, the minimum height of the objective can be set. When the image is in focus, lower the focus-stop ring on the microscope tube as much as is possible and fixate it using the screw.

### 5.1 Always think about

- The following environment is required: Indoor temperature: 0-40°C, Maximum relative humidity: 85 % (non condensing)
- Keep the microscope away from dust and shocks while in use.
- Turn off the light immediately after use.
- Use a soft lens tissue to clean the optics after use.
- Only if needed, use a cloth moistened with water and a mild detergent, rinsing with water and drying immediately with a lint-free cloth.
- After use, cover the microscope with the included dust-cover, and keep it in a dry and clean place.


### 5.2 Do not!

- Wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Use solvents, neither on the microscope, nor on the optics.
- Disassemble objective or eyepieces to attempt to clean them.
- Mishandle or impose unnecessary force on the microscope.
- Clean the unit with volatile solvents or abrasive cleaners.
- Attempt to service the microscope yourself.

**5.3** If you need to send the microscope to Optika for maintenance, please use the original packaging.

The appliance reports the symbol:

This symbol means that the appliance can be a precious source of raw materials. Therefore, it must not be disposed of as waste, but separately collected for the recycling and the recovery of the materials it contained in it. Such materials, if improperly dispersed into the environment, can be harmful to the environment and to human health.

The producer of the equipment, Optika Microscopes, recovers, re-uses and recycles the raw materials contained in the equipment. Such recovery, however, needs your help.

When, at the end of its operating life, you shall decide to dispose of the apparatus, do not try to open it, nor to use parts of it in ways other than reported in this Manual, but bring it back to the Distributor you bought it from, or to other Optika Microscopes distributors. The Distributor shall collect the apparatus free of charge.

The recovery of the raw materials shall then take place in accordance with the European Directive 2002 / 96 / EC and all other relevant Directives. Never disassemble, nor dispose of as waste, apparatus reporting the “crossed bin” symbol indicated above.

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