OPTIKA

M I C R O S C O P E S

I T A L Y

Ver. 1.1.0

 **OPERATION MANUAL**


 **MANUAL DE INSTRUCCIONES**

OPTIKA MICROSCOPES - ITALY

www.optikamicroscopes.com - info@optikamicroscopes.com

1.0 DESCRIPTION	page 3
2.0 INTRODUCTION	page 4
3.0 UNPACKING AND ASSEMBLY	page 4
4.0 USING THE MICROSCOPE	page 7
5.0 MAINTENANCE	page 10
6.0 ELECTRICAL SPECIFICATIONS	page 14
7.0 RECOVERY AND RECYCLING	page 15

2.0 INTRODUCTION

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily laboratory use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

3.0 UNPACKING AND ASSEMBLY

Components of the microscope

Unpack the components of the microscope from the box and put them on a stable surface.

MICROSCOPE AS YOU FIND IT IN THE BOX

The following components are provided:

- Microscope stand with focus and optics already mounted.
- 10X micrometric eyepiece.
- C-mount and eyepiece adaptors. The C-mount is mounted on the microscope when you open the box.

EYEPiece
ADAPTOR


C-MOUNT
ADAPTOR

- Black/white opaque disc.
- One spare fuse F500mA 250V.
- Dust cover.
- Power cable.

Assembling the microscope

1. Put the microscope stand with the mounted optical zoom on a stable plane.
2. If you want to use the eyepiece, put on the top of the microscope the eyepiece adaptor and tighten the screw, in order to fix it in position.

TIGHTEN THIS SCREW


3. Insert the eyepiece.


4. Insert the power cable on the rear (see section 5 and 6 for electrical maintenance and specifications).

1. Adjust the height of the microscope holder so that the distance between the front lens of the objective and the glass plate on the stand is about 10cm. You can move the holder by loosening the big screw on the rear. Later you can fine-adjust the focus using the focus knobs.


LOOSEN THIS
SCREW TO MOVE
THE HOLDER


2. Adjust the safety ring 3-4cm under the holder. This prevents the objective from falling on the plate and getting damaged.

3. Put the sample under observation on the glass plate just below the zoom objective.
4. Put the ON/OFF switch in I position, then press the selector switch to choose between incident illumination (I position) or transmitted illumination (II position). Adjust the brightness of the light with the proper knob on the side. Note that you can also change the orientation angle of the incident illuminator.
5. Rotate the zoom until you reach the maximum magnification. Then rotate the focus knobs until you get a sharp image of your sample. Now you should obtain perfect focussed images at any level of zoom (parfocal zoom).


6. If you want to change the transparent glass plate with the black/white opaque plate:

UNSCREW THE
FIXING NUT

PULL OUT THE GLASS
PLATE AND REPLACE IT
WITH THE OPAQUE PLATE

7. If you want to use a digital camera with C-mount standard, screw on the camera's C-mount the provided adaptor and then replace the eyepiece. Refer to the digital camera user manual for more information.

5.1 Electrical safety precautions

- Before plugging in the power cable with the supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in OFF-position.
- Do not turn the power on and off, off and on immediately as this will shorten the life span of the bulb and may cause damage to the electrical system.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

5.2 Always think about

- The following environment is required: Indoor temperature: 0-40°C, Maximum relative humidity: 85 % (non condensing).
- Keep the microscope away from dust and shocks while in use.
- Turn off the light immediately after use.
- Use a soft lens tissue to clean the optics after use.
- Only if needed, use a cloth moistened with water and a mild detergent, rinsing with water and drying immediately with a lint-free cloth.
- After use, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

5.3 Do not!

- Wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Use solvents, neither on the microscope, nor on the optics.
- Disassemble objective or eyepieces to attempt to clean them.
- Mishandle or impose unnecessary force on the microscope.
- Clean the unit with volatile solvents or abrasive cleaners.
- Attempt to service the microscope yourself.

5.4 Lamp replacement


Before replacing a lamp, be sure that it is completely cooled.

Do not touch an halogen lamp with your bare hands, this could shorten the lifetime of the bulb.

- To replace the **transmitted** light bulb, lift the microscope then open the lamp door. Extract the lamp from the holder and replace it with the new one.

LAMP DOOR


- To replace the **incident** light bulb, first remove the microscope tube from the holder:


Then unscrew the blue filter in order to access the lamp:

Extract the dichroic lamp and insert the new one in the holder:

LAMP HOLDER [BOTTOM VIEW]

Screw the blue filter in front of the lamp, and replace the microscope tube in its position.

5.5 Fuse replacement

If both incident and transmitted illuminators don't work, check that the 2 fuses are not broken, unscrewing the fuse holders on the bottom of the microscope. Be sure to disconnect the power cable before this operation.

FUSE HOLDERS

5.6 If you need to send the microscope to OPTIKA for maintenance, please use the original packaging.


6.0 ELECTRICAL SPECIFICATIONS

Voltage Supply:	230Vac, 50Hz
Power consumption:	10W (at maximum brightness)
Lamps:	halogen 12V – 10W
Fuses:	F500mA 250V (x2)

Art.13 Dlsg 25 July 2005 N°151. “According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal.”

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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