



Model:

CL-31

Typology:

ILLUMINATOR

Description:

Double flexible arm (L = 60 cm) with double LED 2 W.

Chromaticity coordinates (x,y):

0.320, 0.320 (pure white).

Optical power:

Luminous flux: 160 lumen.

Illuminance (at 10 cm distance): 50,000 lux.

Aperture angle:

35°.

Power supply:

Brightness control potentiometer.

External power supply provided.

Input: 110/240 Vac 50/60 Hz; output: 12 Vdc 1500 mA.