



1080p full HD stand-alone Real Time Autofocus Camera

Camera specifications	Stand alone use	Yes
	HDMI Camera resolution (MP)	2
	Camera resolution (n° of pixels: W x H)	1920x1080
	HDMI Signal output	Yes
	Trigger Signal output	Yes
	Color / Monochrome	Color
	Sensor Size	1/2.8"
	Sensor Technology	CMOS
	Sensor Type	SONY EXMOR
	Optical Zoom	1x-14x
	Digital Zoom	1x-6x
	Rolling Shutter	Yes
	Autofocus	Yes
	Image Format	16/9
	Pixel Size (µm)	2.9x2.9
	Frame rate full resolution (fps)	60@1920x1080
	G sensitivity	1300mV at 1/30sec
	Dark Signal	0.34mV at 1/30sec
	ADC conversion	8 Bit
	Color Depth	8 Bit
	Exposure time	Auto
	Binning	1x1
	IR filter	Yes
Cooling Temperature	None	
Camera Power	5V 2A	
Dimensions (mm)	88 x 72 x 125	
Weight (Kg)	0.57	
Camera accessories	HDMI cable (2m)	Yes
	Wireless Mouse	Yes
	16GB SD card	Yes
Environment specifications	Working Temperature	-10°C up to +50°C
	Storage Temperature	-20°C up to +60°C
	Relative Humidity	30% to 80%
	Storage Humidity	10% to 60%

Warranty

Years

2