

Model:
Eyepieces

Description:
Widefield eyepieces, high eyepoint type



M-1001



M-1002



M-1003



M-781

Magnification	10x (M-1001, M-1002, M-781) 15x (M-1003)
Field Number	Models available: - 22 mm (M-1001 and M-781) - 24 mm (M-1002) - 16 mm (M-1003)
Optical system	High eyepoint type (for observation with eyeglasses).
Tube diameter	30 mm.
Notes	With eyepiece cups (except M-1003). M-781: with micrometer 10 mm, 0.1 mm divisions.