

Model:

Observation heads

Description:

Optical observation heads for B-810 and B-1000



*Trinocular head
M-1188 for B-810*



*Trinocular head
M-1010 for B-810
M-1011 for B-1000*



Binocular ERGO head M-1012



Trinocular ERGO head M-1013

Output ports	Models available: - Trinocular head with fixed photo port (ratio 50/50) (M-1188) - Trinocular head with 2 positions (ratios 100/0, 50/50) (M-1010) - Trinocular head with 3 positions (ratios 100/0, 50/50, 0/100) (M-1011) - Binocular ergonomic head (M-1012) - Trinocular ergonomic head (ratios 100/0, 50/50) (M-1013)
Inclination	Trinocular head: 30° Binocular ERGO head: 30° - 60° Trinocular ERGO head: 5° - 35°
Rotation	Heads can be rotated by 360°
Interpupillary adjustment	50 mm - 75 mm
Diopter adjustment	+ - 5 diopter adjustment
Notes	Photo tube (internal diameter 23mm) included in Trinocular ERGO head. Eyepieces NOT included with heads.