

**Model:**

Fluorescence attachment with LED source

**Description:***LED fluorescence with 4 filter positions*

<b>Source</b>	Blue and Green LED. Power: 3.6 W. Lifetime: about 50,000 h.  LEDs are factory-precentered.
<b>Filters</b>	4-position sliding filter holder.  Standard mounted filtersets: - B: EX 450-490, DM 495, EM 520LP; - G: EX 500-540, DM 565, EM 575LP;
<b>Power supply</b>	Included. Voltage: 100/240 Vac, 50/60 Hz, Output 6 V 2.5 A. Max power required: 8 W.
<b>Dimensions</b>	230 mm overall length.