



Model:
Metallographic attachment

Description:
Metallographic attachment for incident brightfield and darkfield observation.

Source	Halogen lamp. Power: 100 W. Lifetime: about 3,000 hours.
Filters	Slots for analyzer and polarizer filters (filters included). Slider for ND neutral density filter insertion (filter included).
Diaphragms	Aperture and field diaphragms.
Darkfield	Front slider for incident darkfield insertion.
Power supply	Included. Voltage: 100/240 Vac, 50/60 Hz Fuse: T3,15A 250 V. Max power required: 130 W.
Dimensions	600 mm overall length.