

**Model:**

Metallographic attachment

**Description:**

*Metallographic attachment for incident brightfield and darkfield observation, with LED illuminator.*

<b>Illumination</b>	High efficiency 18 W LED. Pre-centered. Manual brightness adjustment. LED average life time approx. 50,000 hours. External Power Supply: Input 110/240 Vac, 50/60 Hz, 1 A. Max required power: 24 W.
<b>Filters</b>	Slots for analyzer and polarizer filters (filters included). Slider for ND neutral density filter insertion.
<b>Diaphragms</b>	Aperture and field diaphragms, centrable.
<b>Darkfield</b>	Front slider for incident darkfield insertion.
<b>Dimensions</b>	430 mm overall length.