



Model:
Camera adapter

Description:
CCD camera adapter 0.5x

Optical magnification	0.5x
Connection to the camera	C-mount standard
CS-mount	Allowed, but the C-mount ring (optional) is suggested
Connection to the microscope	23 mm eyepiece connection
Allowed connections	30 mm & 30.5mm eyepiece connection (with M-113.1 & M-113.2)
Sensor size compatibility	1/2" 1/2.3" 1/2.5"
External fast identification	Black
Length	55 mm
Color	Black
Allowed	The camera is free to move 360 degree