

**Model:**

MPC (Mineral Solid Surface) rackless mechanical stage

Description:

Double-layer mechanical stage with coaxial translation knobs, with belt-drive movement. Innovative mineral/plastic material, anti-scratch and acid-resistant.

Dimensions	220x150 mm.
Translation	X: 80 mm. Y: 50 mm. Coaxial drive control length: 137 mm.
Slide holder	For 2 standard slides.
Graduated axis	On both X and Y: 0.1 mm vernier scale.
Notes	Internal belt-drive system on X axis (no protruding parts during movement). Innovative mineral/plastic material, anti-scratch and acid-resistant.