

**Model:**

MPC (Mineral Solid Surface) Belt Drive Mechanical Stage

**Description:**

*Double-layer mechanical stage with coaxial translation knobs, with belt-drive movement. Innovative mineral/plastic material, anti-scratch and acid-resistant.*

<b>Dimensions</b>	243x160 mm.
<b>Translation</b>	X: 100 mm. Y: 60 mm. Coaxial drive control length: 110 mm (extendable to 125 mm).
<b>Slide holder</b>	For 2 standard slides.
<b>Graduated axis</b>	On both X and Y; 0.1 mm vernier scale.
<b>Notes</b>	Internal belt-drive system (no protruding parts during movement). Innovative mineral/plastic material, anti-scratch and acid-resistant.