

**Model:**

Rotating mechanical stage + Attachable translation stage

Description:

Rotating stage for polarizing microscopes, with X-Y translation mechanism

Dimensions	Diameter 165 mm.
Translation	Attachable Mechanical Stage M-1146: X: 41 mm Y: 32 mm Low-profile translation knobs (stage can be rotated without interference with the objectives).
Slide holder	For 1 standard slide.
Graduated axis	On both X and Y; 0.1 mm vernier scale. 360° graduation with 0.1° vernier scale.
Notes	Click-stop every 45° during rotation. Fixing screw to lock the position.