



*8-position Universal Condenser  
for B-1000*

<b>Material</b>	Type	Aluminum / Aluminum die-cast / Plastic
<b>Top lens</b>	Type	Swing-out
	Top lens <b>M-1205</b>	0.2
	Top lens <b>M-1206</b>	0.9
	Top lens <b>M-1207</b>	1.4
<b>Positions</b>	Brightfield	Yes
	Darkfield	Yes
	Phase Contrast	Yes (3 positions)
	DIC	Yes (3 positions)
<b>Darkfield stops</b>	Darkfield stop (dry) <b>M-1216</b>	Yes (for objectives with N.A. < 0.7)
	Darkfield stop (oil) <b>M-1217</b>	Yes
<b>Phase stops</b>	4x phase ring <b>M-1215</b>	Yes
	10x/20x phase ring <b>M-1212</b>	Yes
	40x/60x phase ring <b>M-1213</b>	Yes
	100x phase ring <b>M-1214</b>	Yes
<b>DIC prisms</b>	10x DIC prism <b>M-1208</b>	Yes
	20x DIC prism <b>M-1209</b>	Yes
	40x/60x DIC prism <b>M-1210</b>	Yes
<b>Diaphragm</b>	Type	Iris aperture
<b>Polarizing filter</b>	Type	Insertable
	Rotatable	Yes, 360°
<b>Product Weight</b>	(kg)	0.65