



**Model:**  
DIC set

**Description:**  
*Koehler DIC set for transmitted light*

<b>Observation Mode</b>	Differential Interference Contrast based on Koehler technique (best performance with condenser apertures smaller than around 0.4 times the objective aperture).
<b>DIC slider</b>	Nomarski prism slider with prism control knob for contrast adjustment.
<b>Filters</b>	Kit includes the following filters: <ul style="list-style-type: none"> <li>• Polarizer filter</li> <li>• Analyzer filter</li> <li>• Interference Green filter (IF550)</li> </ul>
<b>Notes</b>	DIC requires objectives M-1075, M-1076, M-1077, M-1078 and M-1079 in order to achieve the best performance. <i>Please contact OPTIKA for further information.</i>